

**QUANTITATIVE CONJUGATE IMAGING OF IODINE-123 AND
TECHNETIUM-99m LABELED BRAIN AGENTS
IN THE BASAL GANGLIA**

A Thesis
Presented to
The Academic Faculty

By

Desiree N. Jangha

In Partial Fulfillment
of the Requirements for the Degree
Doctor of Philosophy in the
George W. Woodruff School
of Mechanical Engineering

Georgia Institute of Technology

August 2006

**QUANTITATIVE CONJUGATE IMAGING OF IODINE-123 AND
TECHNETIUM-99m LABELED BRAIN AGENTS
IN THE BASAL GANGLIA**

Approved by:

Dr. C K Chris Wang, Co-Advisor
Mechanical Engineering
Georgia Institute of Technology

Dr. John N. Aarsvold, Co-Advisor
Nuclear Medicine
Emory University

Dr. Nolan Hertel
Mechanical Engineering
Georgia Institute of Technology

Dr. Farzad Rahnema
Mechanical Engineering
Georgia Institute of Technology

Dr. Paul J. Benkeser
Biomedical Engineering
Georgia Institute of Technology

Date Approved: March 29, 2006

This body of work is dedicated to my grandmother, Luella F. Foster,
whom instilled in me the importance of education.

ACKNOWLEDGEMENTS

The completion of my dissertation would not have been possible without a large number of supportive individuals. I first like to give thanks to my Creator. I can do all things through Christ who strengthens me, Philippians 4:13.

I have been blessed with extraordinary advisors. Dr. John Aarsvold has been a true blessing. I thank him for his outstanding guidance, encouragement, and many late nights. I thank Dr. Chris Wang for his time and frank advisement. I give special thanks to Mr. Robert Mintzer for sharing his time (late nights as well) and extensive knowledge of computer programming and science. I thank Dr. John Valentine for introducing me to Dr. Aarsvold and Mr. Mintzer. I thank my reading committee, Dr. Paul J. Benkeser, Dr. Nolan Hertel, and Dr. Farzad Rahnema for their advisement and review of my research. I acknowledge Ms. Arlene Smith, Ms. Cosetta Williams, Ms. Jeannie Baggett, Ms. Glenda Johnson, and Ms. Terri Keita for their tremendous support and encouraging words. I thank Dr. Wepfer for his time and advisement. I thank the Nuclear and Radiological Engineering/Health Physics department for affording me a wonderful educational opportunity. I thank the David and Lucille Packard Foundation and Georgia Institute of Technology for funding my graduate studies.

My family deserves a large thank you. I thank my mother, Joyce R. Dixon, for her encouragement, active support and patience. (I will be forever grateful to her for saving my paper at the midnight hour.) I thank my sisters, Jasmine Marie Batth and Erin Leigh

Batth for their support and comedic breaks. I thank my grandparents Ms. Luella F. Foster and Mr. and Mrs. William R. Batth.

During this journey I met and married my wonderful husband, Sundiata Kwesi Jangha. He has been most supportive of my education endeavors. I am extremely grateful for my children, Naimah Afi and Amari Sundiata, who were constantly asking me, “Are you a doctor, yet?” They have encouraged me to finish my thesis as soon as possible. I hope through watching me they have learned tenacity. I thank my mother-in-law, Ms. Gerry Jangha, my grandmother-in-law, Ms Georgia Whitted, my sister-in-law, Awa Georgia Jangha, and my brother-in-law, Mr. Kweku Jangha and his family.

I thank the following individuals for pushing me and sometimes pulling me to this point: Ms. Stephanie Alexander, Ms. Tanya LaShawn Bridges, Dr. Aisha (Fields) Buchanan, Dr. Chekesha Clingman, Mr. Roosevelt Elivert, Ms. Carmen Greene, Ms. Likisha Griffin, Cleotha and John Griffith, Jerolyn and Ron Kight, Ms. Asha Robinson, Mrs. Xin Zhang, and Ms. Anett Young. I specially thank Ms. Gena Woodruff for all her help and advice. She has been a blessing to me on numerous occasions. I am truly grateful.

TABLE OF CONTENTS

ACKNOWLEDGEMENTS	IV
LIST OF TABLES	X
LIST OF FIGURES	XXIII
SUMMARY	XXVI
CHAPTER 1: INTRODUCTION.....	1
CHAPTER 2: BACKGROUND	4
2.1 SINGLE-PHOTON EMISSION IMAGING	5
2.1.1 Photon Interactions	6
2.1.1.1 Photoelectric Absorption	6
2.1.1.2 Compton Scattering	7
2.1.2 Radiopharmaceuticals	9
2.1.3 Gamma Camera	10
2.1.3.1 Components	10
2.1.3.2 Gamma Camera Function	13
2.1.3.3 Gamma Camera Performance	16
2.1.4 Attenuation and Scatter	18
2.1.5 Qualitative Single-photon Emission Imaging	20
2.1.6 Planar Emission Imaging	22
2.1.7 SPECT	23
2.1.8 Quantitative Single-photon Emission Imaging	25
2.2 CLASSIC CONJUGATE IMAGING	30
2.2.1 The Planar Model	32
2.2.2 The Volume Model	37
2.2.3 Sum of Two Sources	41
2.2.4 Applications	43
2.2.5 Extension of Classic Conjugate Imaging	45
2.2.5.1 Non-Uniformly Attenuating Medium	46
2.2.5.2 Background Activity	48
2.2.5.3 Multiple Sources	50
2.3 STRIATAL DOPAMINE AND PARKINSONISM AND SCHIZOPHRENIA	53
2.3.1 Striata	54

2.3.2	Parkinsonism.....	58
2.3.3	Schizophrenia.....	62
2.4	IMAGING THE STRIATA.....	66
2.4.1	Non-specific Binding.....	69
2.4.2	Pre- and Post-Synaptic Radiopharmaceuticals	70
2.4.2.1	Biodistribution of ^{123}I β -CIT.....	72
2.4.2.2	Biodistribution of ^{123}I FP-CIT	73
2.4.2.3	Biodistribution of $^{99\text{m}}\text{Tc}$ TRODAT.....	75
2.4.2.4	Biodistribution of ^{123}I IBZM.....	77
CHAPTER 3: CONJUGATE IMAGING OF TWO ALIGNED SOURCES.....		79
3.1	GEOMETRIC MODELS	80
3.1.1	A Model of Striata and Head	80
3.1.2	Simplified Model of Striata and Head	80
3.2	PARAMETERS REQUIRED.....	82
3.3	PIXEL-BASED COMPUTATIONS AND ANALYSES.....	83
3.4	CONJUGATE IMAGING OF TWO ALIGNED SOURCES.....	84
3.5	ERROR ANALYSIS	93
3.5.1	Uncertainty of Bone Attenuation Coefficient μ_b	94
3.5.2	Uncertainty of A_1 -to- A_2 Ratio	95
3.5.3	Uncertainty of Brain Tissue Attenuation Coefficient μ	96
3.5.4	Uncertainty of Position	98
CHAPTER 4: PHANTOM EXPERIMENTS		103
4.1	METHODS	104
4.1.1	Conjugate Imaging.....	104
4.1.1.1	Equipment.....	104
4.1.1.2	Activity Concentrations	105
4.1.1.3	Data Acquisition	106
4.1.1.4	Data Processing.....	107
4.1.2	Quantitative SPECT.....	110
4.1.2.1	Equipment.....	112
4.1.2.2	Activity Concentrations	112
4.1.2.3	Data Acquisition	112
4.1.2.4	Data Processing.....	113
4.2	RESULTS	114
4.2.1	Conjugate Imaging.....	114
4.2.2	Quantitative SPECT Study	117
4.3	DISCUSSION	119
CHAPTER 5: SHIELDING CONSIDERATIONS FOR A NOVEL DEDICATED SYSTEM.....		123
5.1	MONTE CARLO SIMULATIONS.....	125
5.1.1	Primate Model.....	128

5.1.2	Camera Model.....	128
5.2	SHIELDING DETERMINATION.....	132
5.2.1	Data Analysis.....	133
5.2.2	Results.....	136
5.2.3	Discussion.....	141
5.3	SHIELDING CHARACTERIZATION.....	143
5.3.1	Data Analysis.....	143
5.3.2	Results.....	147
5.3.3	Discussion.....	153
CHAPTER 6: EVALUATION AND CHARACTERIZATION OF THE MODIFIED CONJUGATE IMAGING TECHNIQUE		155
6.1	MONTE CARLO SIMULATIONS.....	156
6.1.1	Data Acquisition	156
6.1.2	Image Generation.....	156
6.1.3	Conjugate Imaging.....	162
6.1.4	ROI Study	164
6.2	RESULTS	166
6.2.1	Absolute Activity Estimation.....	166
6.2.1.1	13 pixels x 13 pixels ROI.....	166
6.2.1.2	19 pixels x 19 pixels ROI.....	169
6.2.2	Relative Activity Estimation.....	171
6.3	DISCUSSION	173
6.4	CLINICAL RELEVANCE	177
CHAPTER 7: SCATTER CORRECTIONS		178
7.1	TRIPLE ENERGY WINDOW.....	179
7.2	MONTE CARLO SIMULATIONS.....	182
7.2.1	TEW Scatter Correction.....	183
7.3	RESULTS	188
7.3.1	¹²³ I Radiopharmaceutical Studies	188
7.3.1.1	Scatter-Corrected Emission Only Study	188
7.3.1.2	Scatter-Corrected Emission and Camera Sensitivity Study.....	194
7.3.2	^{99m} Tc Radiopharmaceutical Study	199
7.3.2.1	Scatter-Corrected Emission Only Study	199
7.3.2.2	Scatter-Corrected Emission and Camera Sensitivity Study.....	201
7.4	DISCUSSION	203
7.5	CLINICAL RELEVANCE	206
CHAPTER 8: CONCLUSIONS AND FUTURE WORK		207
8.1	CONCLUSIONS.....	207
8.2	FUTURE WORK	209

APPENDIX A	211
APPENDIX B	233
REFERENCES.....	343

LIST OF TABLES

Table 2.1. Radiopharmaceuticals used in single-photon emission imaging studies of dopaminergic neurotransmission in humans.	67
Table 3.1 Positioning Errors	101
Table 3.2 Summary of Errors	102
Table 4.1. Effective Attenuation Coefficient and Background Activity Concentration	115
Table 4.2. Estimated Above Background Activity in Left and Right Spheres.....	115
Table 4.3. Sum of Estimated Above Background Activities and Estimated Total Activity	115
Table 4.4. No Background Activity Study: Estimates of Effective Attenuation Coefficient and Background Activity Concentration	116
Table 4.5. No Background Activity Study Estimated Above Background Activity in Left and Right Spheres.....	117
Table 4.6. No Background Activity Study Sum of Estimated Above Background Activities and Estimated Total Activity	117
Table 4.7. Percent Errors for quantitative SPECT Studies I	118
Table 4.8. Percent Errors for quantitative SPECT Studies II.....	119
Table 5.1 Percent Uptake of Initial Dose in Critical Organs	135
Table 6.1 Percent Uptake of Radiopharmaceuticals in Brain, Striata, and Critical Organs.....	157
Table 6.2. Percent Errors for Background Activity Estimates for Activity Uptakes	168
Table 6.3. Percent Errors for Activity Estimates for Narrow-beam Attenuation Study.....	168

Table 6.4. Percent Errors for Activity Estimates for Broad-beam Attenuation Study	168
Table 6.5. Percent Errors for Background Activity Estimates for Activity Uptakes	170
Table 6.6. Percent Errors for Activity Estimates for Narrow-beam Attenuation Study.....	170
Table 6.7 Percent Errors for Activity Estimates for Broad-beam Attenuation Study	171
Table 7.1. Scatter energy windows for ^{123}I studies.....	184
Table 7.2. Scatter energy windows for $^{99\text{m}}\text{Tc}$ studies.....	188
Table 7.3. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 175-180.....	189
Table 7.4. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 185-190.....	189
Table 7.5. Percent Error in Striatal Activity Estimates; Lower:137-142, Upper:175-180.....	189
Table 7.6. Percent Error in Striatal Activity Estimates; Lower:139-144, Upper: 173-178.....	190
Table 7.7. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 000-000.....	190
Table 7.8. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper:175-180.....	192
Table 7.9. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 185-190.....	192
Table 7.10. Percent Error in Striatal Activity Estimates; Lower:137-142, Upper:175-180.....	192
Table 7.11. Percent Error in Striatal Activity Estimates; Lower:139-144, Upper: 173-178	193
Table 7.12. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 000-000	193

Table 7.13. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 175-180	195
Table 7.14. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 185-190	195
Table 7.15. Percent Error in Striatal Activity Estimates Lower:137-142, Upper: 175- 180.....	195
Table 7.16. Percent Error in Striatal Activity Estimates; Lower:139-144, Upper: 173-178	196
Table 7.17. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 000-000	196
Table 7.18. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 175-180	197
Table 7.19. Percent Error of Striatal Uptakes; Lower:130-142, Upper: 185-190...	197
Table 7.20. Percent Error in Striatal Activity Estimates; Lower:137-142, Upper: 175-180	197
Table 7.21. Percent Error in Striatal Activity Estimates; Lower:139-144, Upper: 173-178	198
Table 7.22. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 000-000	198
Table 7.23. ^{99m} Tc TRODAT Percent Error in Striatal Activity Estimates.....	200
Table 7.24. ^{99m} Tc TRODAT Percent Error in Striatal Activity Estimates.....	200
Table 7.25. ^{99m} Tc TRODAT Percent Error in Striatal Activity Estimates.....	202
Table 7.26. ^{99m} Tc TRODAT Percent Error in Striatal Activity Estimates.....	202
Table A.1. ¹²³ I β-CIT 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions	212
Table A.2. ¹²³ I β-CIT 1:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	213
Table A.3. ¹²³ I β-CIT 2:1 Activity Uptake Ratio, Primary and High-Energy Emissions	214
Table A.4. ¹²³ I β-CIT 2:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	215

Table A.5. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Primary and High-Energy Emissions	216
Table A.6. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	217
Table A.7. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions	218
Table A.8. ^{123}I FP -CIT 1:1 Activity Uptake Ratio, Primary 159 keV Emissions...	219
Table A.9. ^{123}I FP -CIT 2:1 Activity Uptake Ratio Primary and High-Energy Emissions	220
Table A.10. ^{123}I FP -CIT 2:1 Activity Uptake Ratio, Primary 159 keV Emissions.	221
Table A.11. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Primary and High-Energy Emissions	222
Table A.12. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Primary 159 keV Emissions.	223
Table A.13. ^{123}I IBZM 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions	224
Table A.14. ^{123}I IBZM 1:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	225
Table A.15. ^{123}I IBZM 2:1 Activity Uptake Ratio, Primary and High-Energy Emissions	226
Table A.16. ^{123}I IBZM 2:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	227
Table A.17. ^{123}I IBZM 5:1 Activity Uptake Ratio Primary and High-Energy Emissions	228
Table A.18. ^{123}I IBZM 5:1 Activity Uptake Ratio, Primary 159 keV Emissions.....	229
Table A.19. $^{99\text{m}}\text{Tc}$ TRODAT 1:1 Activity Uptake Ratio.....	230
Table A.20. $^{99\text{m}}\text{Tc}$ TRODAT 2:1 Activity Uptake Ratio.....	231
Table A.21. $^{99\text{m}}\text{Tc}$ TRODAT 5:1 Activity Uptake Ratio.....	232
Table B.1. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	234
Table B.2. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	235

Table B.3. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	236
Table B.4. ^{123}I FP -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	237
Table B.5. ^{123}I FP -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	238
Table B.6. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	239
Table B.7. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	240
Table B.8. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	241
Table B.9. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	242
Table B.10. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	243
Table B.11. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	244
Table B.12. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	245
Table B.13. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	246
Table B.14. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	247
Table B.15. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	248
Table B.16. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	249
Table B.17. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	250

Table B.18. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	251
Table B.19. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	252
Table B.20. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	253
Table B.21. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	254
Table B.22. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	255
Table B.23. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	256
Table B.24. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	257
Table B.25. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	258
Table B.26. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	259
Table B.27. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	260
Table B.28. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	261
Table B.29. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	262
Table B.30. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	263
Table B.31. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	264
Table B.32. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	265

Table B.33. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	266
Table B.34. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	267
Table B.35. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	268
Table B.36. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	269
Table B.37. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	270
Table B.38. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	271
Table B.39. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	272
Table B.40. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	273
Table B.41. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	274
Table B.42. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	275
Table B.43. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	276
Table B.44. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	277
Table B.45. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	278
Table B.46. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	279
Table B.47. $^{99\text{m}}\text{Tc}$ TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV	280

Table B.48. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV	281
Table B.49. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV	282
Table B.50. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV	283
Table B.51. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV	284
Table B.52. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV	285
Table B.53. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV	286
Table B.54. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV	287
Table B.55. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV	288
Table B.56. ^{123}I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	289
Table B.57. ^{123}I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	290
Table B.58. ^{123}I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	291
Table B.59. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	292
Table B.60. ^{123}I FP -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	293

Table B.61. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	294
Table B.62. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	295
Table B.63. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	296
Table B.64. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	297
Table B.65. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	298
Table B.66. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	299
Table B.67. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	300
Table B.68. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	301
Table B.69. ^{123}I FP -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	302
Table B.70. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	303
Table B.71. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	304

Table B.72. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV	305
Table B.73. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV	306
Table B.74. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	307
Table B.75. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	308
Table B.76. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	309
Table B.77. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	310
Table B.78. ^{123}I FP -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	311
Table B.79. ^{123}I FP -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	312
Table B.80. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	313
Table B.81. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	314
Table B.82. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV	315

Table B.83. ^{123}I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	316
Table B.84. ^{123}I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	317
Table B.85. ^{123}I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	318
Table B.86. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	319
Table B.87. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	320
Table B.88. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	321
Table B.89. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	322
Table B.90. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	323
Table B.91. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV	324
Table B.92. ^{123}I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	325
Table B.93. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	326

Table B.94. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	327
Table B.95. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	328
Table B.96. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	329
Table B.97. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	330
Table B.98. ^{123}I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	331
Table B.99. ^{123}I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	332
Table B.100. ^{123}I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV	333
Table B.101. $^{99\text{m}}\text{Tc}$ TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV	334
Table B.102. $^{99\text{m}}\text{Tc}$ TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV	335
Table B.103. $^{99\text{m}}\text{Tc}$ TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV	336
Table B.104. $^{99\text{m}}\text{Tc}$ TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV	337

Table B.105.^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV	338
Table B.106.^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV	339
Table B.107.^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV	340
Table B.108.^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV	341
Table B.109.^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV	342

LIST OF FIGURES

Figure 2.1. Photoelectric absorption.	7
Figure 2.2. Compton scattering.	8
Figure 2.3. Basic components of a gamma camera system.	11
Figure 2.4. Collimators and cross section of a hexagonal collimator.	12
Figure 2.5. Acquisition of emission data.	21
Figure 2.6. SPECT projection acquisition.	24
Figure 2.7. Acquisition of transmission data.	27
Figure 2.8. Conjugate imaging: Acquisition of emission data.	33
Figure 2.9. Conjugate Imaging: Acquisition of flood data.	34
Figure 2.10. Conjugate Imaging: Acquisition of transmission data.	36
Figure 2.11. Volume source with thickness t	38
Figure 2.12. A source with thickness t embedded in an object with thickness T	39
Figure 2.13. Two sources separated by a distance s in a uniformly attenuating object.	42
Figure 2.14. Two activity distributions embedded in a non-uniformly attenuating medium.	48
Figure 2.15. A primary lesion surrounded by regions of background activity in a non-uniformly attenuating medium.	49
Figure 2.16. Illustration of the basal ganglia.	56
Figure 2.17. Striatal neurons.	57
Figure 2.18. Percentage yields of ^{123}I decays.	68

Figure 2.19. ^{123}I β -CIT and ^{123}I FP-CIT striatal activity uptake.....	74
Figure 2.20. $^{99\text{m}}\text{Tc}$ TRODAT-1 striatal activity uptake.....	76
Figure 2.21. ^{123}I IBZM striatal activity uptake.	78
Figure 3.1. Seven-volume model of the head.	81
Figure 3.2. Simplified seven-region model.....	82
Figure 4.1. Experimental set-up..	108
Figure 4.2. The transmission data acquisition set-up.	109
Figure 4.3. Examples of flood, transmission/emission, and transmission activity count maps for heads II and III.....	111
Figure 5.1. Schematic of gantry of mobile dedicated conjugate imaging system....	124
Figure 5.2. Photograph of conjugate imaging gantry.....	124
Figure 5.3. Configuration of idealized primate and camera used in the simulations.	127
Figure 5.4. Detailed diagram of idealized camera components.....	129
Figure 5.5. Small gamma camera for use in a novel conjugate imaging system.....	131
Figure 5.6. Cross-sectional schematic of small gamma camera for use in a novel conjugate imaging system.....	131
Figure 5.7. Pulse height spectrums as detected by ideal gamma camera.	137
Figure 5.8. Pulse height spectrums as detected by gamma camera with finite resolution.	138
Figure 5.9. Comparison of percentage increase for ^{123}I β -CIT activity uptake in the striata of controls, early-stage and late-stage Parkinson's disease patients....	139
Figure 5.10. Comparison of percentage increase for ^{123}I FP-CIT activity uptake in the striata of controls, early-stage and late-stage Parkinson's disease patients.	140
Figure 5.11. Comparison of percentage increase for ^{123}I IBZM activity uptake in the striata of controls, early-stage and late-stage Parkinson's disease patients....	141

Figure 5.12. Count distribution across the face of the scintillation crystal from activity uptake in the GI tract and liver for a 159 keV energy-windowed image..	148
Figure 5.13. Count distribution across the face of the scintillation crystal from activity uptake in the lungs for a 159 keV energy-windowed image.....	149
Figure 5.14. Comparison of percentage increase in striatal counts in a 4 cm ² ROI due to activity uptake of ¹²³ I β-CIT, ¹²³ I FP-CIT, and ¹²³ I IBZM..	150
Figure 5.15. Comparison of percentage increase in striatal counts in a 9 cm ² ROI due to activity uptake of ¹²³ I β-CIT, ¹²³ I FP-CIT, and ¹²³ I IBZM..	151
Figure 5.16. Comparison of percentage increase in striatal counts in a 16 cm ² ROI due to activity uptake of ¹²³ I β-CIT, ¹²³ I FP-CIT, and ¹²³ I IBZM..	152
Figure 6.1. Images of a calibration source in the 159 keV 20% energy window as seen by the left gamma camera.....	159
Figure 6.2. Images of brain activity in the 159 keV 20% energy window as seen as by the left gamma camera..	159
Figure 6.3. Images of activity in the right striatum in the 159 keV 20% energy window as seen by the left gamma camera.....	160
Figure 6.4. Images of activity in the left striatum in the 159 keV 20% energy window as seen by the left gamma camera.	160
Figure 6.5. Images of activity in the GI tract (left) and lungs (right) in the 159 keV 20% energy window as seen by the left gamma camera.	161
Figure 6.6. ¹²³ I β-CIT images 1:1 simulated data.....	161
Figure 6.7. Relative uptake ratio for 13 pixels x 13 pixels ROI.....	172
Figure 6.8. Relative uptake ratio for 19 pixels x 19 pixels ROI.....	173
Figure 7.1. Locations of the energy windows used in standard TEW.....	208
Figure 7.2. ¹²³ I emission spectrum in water.....	214
Figure 7.3. ^{99m} Tc photon emissions in water.....	215

SUMMARY

In the research reported in this dissertation, the concept of classic conjugate imaging, a non-tomographic nuclear medicine technique, is modified such that activity of a radiopharmaceutical distribution in the striata can be estimated. A mathematical model is developed that extended the application of classic conjugate imaging to estimation of two distinct and aligned activity distributions. Error analysis of the mathematical model is performed to characterize the accuracy of the model and to benchmark the limitations of the model. Phantom experiments are performed to demonstrate the practical application of the model and to evaluate its' accuracy. A Monte Carlo simulation model of conjugate imaging of activity uptake in the striata of a primate is developed to evaluate the accuracy of the modified conjugate imaging technique as applied in the use of a dedicated conjugate imaging system. In addition, the simulation model is used to determine and characterize the shielding design of the small field of view gamma cameras comprising the dedicated conjugate imaging system. The application of scatter correction is investigated to address the downscatter of high-energy photon emissions into the photopeak window and the inclusion of scattered primary photons in the photopeak window.

In this dissertation, it is shown that the modified conjugate imaging technique developed can be used to estimate accurately activity uptake in each of two distinct and aligned activity distributions. The accuracy of the technique was shown to be comparable to that of clinical quantitative SPECT. The modified conjugate imaging technique used

with the dedicated conjugate imaging system may, therefore, be a viable quantitative nuclear medicine technique for activity estimation of radiopharmaceutical uptake in the striata of Parkinsonian and schizophrenic patients. The portability and low cost relative to SPECT systems make the dedicated conjugate imaging system advantageous for clinics with Parkinsonian and schizophrenic patients, who are unable to travel due to physical or mental limitation.

Chapter 1

INTRODUCTION

Parkinsonism is a term for progressive neurodegenerative disorders, which describes a variety of movement disorders pathologically characterized by the degeneration of dopaminergic neurons in the striatal dopaminergic system. Schizophrenia is a chronic neuropsychiatric disorder associated with the hyper-production of striatal dopamine. During the last two decades nuclear medicine has been in the forefront of investigations of dopamine loss characteristic of Parkinsonism and dopamine hyper-production associated with schizophrenia. The striata, part of the basal ganglia, a collection of nuclei located deep in the cerebrum, are greatly responsible for controlling movements and maintaining postures by producing and transmitting dopamine. Studies have suggested that nuclear medicine imaging can play a vital role in the clinical management of Parkinsonism and schizophrenia, specifically in the initial diagnosis, intermittent follow-ups, and drug efficacy investigations for clinical management of these diseases.

The primary goal of the research presented in this dissertation is the development and characterization of conjugate imaging of two aligned activity distributions. The pursuit of this research is motivated by the study and evaluation of activity uptake in the striata of individuals with Parkinsonian or schizophrenic symptoms. Of particular interest is the influence of radiopharmaceuticals. Most of the radiopharmaceuticals

presently available for striatal imaging are labeled with iodine-123 (^{123}I) and technetium-99m ($^{99\text{m}}\text{Tc}$). Differences in biodistribution and emission spectrums should be addressed to ensure accurate assessment of activity uptake in the striata.

There are two main objectives for this study. The first objective is the development and assessment of a model describing conjugate imaging to the estimation of activities in each of two aligned sources. The second objective is the evaluation of the modified conjugate imaging technique as it applies to the estimation of activity uptake in the striata of a primate using a dedicated conjugate imaging system. Six specific tasks are completed to achieve the objectives. The first three tasks are related to the first main objective. The last three tasks are related to the second main objective.

In task one, a mathematical model that extends the application of classic conjugate imaging to include the estimation of activity above background in each of two aligned source distributions is developed. In task two, error analysis of the model is performed. In task three, phantom experiments are performed to evaluate the accuracy of the mathematical model and to demonstrate the practical applicability of the technique.

Tasks four, five, and six utilize a Monte Carlo simulation model developed to simulate conjugate imaging of activity uptake in the striata of a primate using a dedicated conjugate imaging system. The Monte Carlo simulation model was developed using the Monte Carlo radiation transport code MCNP-4C. In task four, the Monte Carlo simulation model is used to design and characterize the shielding design of two small field of view (FOV) gamma cameras comprising the dedicated conjugate imaging system. In task five, the Monte Carlo simulation model is used to assess the accuracy of the estimations of activity uptake in the striata using the modified conjugate imaging

technique and dedicated conjugate imaging system. In task six, scatter correction is investigated. The application of the triple energy window (TEW) scatter correction technique to data generated from the Monte Carlo simulations using ^{123}I and $^{99\text{m}}\text{Tc}$ radiopharmaceuticals is assessed.

This dissertation has eight chapters. Chapter 2 provides background relevant to this research. The work and results of this research are presented in Chapters 3-7. Chapter 8 summarizes the research and provides suggestions for future work.

Chapter 2

BACKGROUND

During the last two decades, nuclear medicine has been used in investigations at the forefront of research of the abnormal dopamine function characteristic of Parkinsonian syndromes and schizophrenia (Booi, Tissingh et al. 1999). In these investigations the dopaminergic system in the striata is assessed in vivo through qualitative and semi-quantitative analyses. Qualitative imaging techniques have limitations, as radiopharmaceutical uptake assessment is based mainly on human eye discernment and physician experience. Quantitative nuclear medicine techniques rely less on the physician's eye and more on the comparison of numerical values representing radiopharmaceutical uptake. Quantitative studies often result in more accurate assessment of disease development and drug therapy treatment. However, quantitative studies performed to date have focused on the relative activity ratio between the left and right striatum. The modified conjugate imaging technique introduced in this dissertation determines absolute activity estimates, in addition to relative activity estimates. The modified conjugate imaging technique is a single-photon emission technique that utilizes the geometry of the striata to estimate activity uptake in the striata.

This chapter reviews the background topics relevant to the techniques adopted in this dissertation. The subjects include single-photon emission imaging with specific

focus on classic conjugate imaging, striatal dopamine, Parkinsonism and Schizophrenia, and finally conjugate imaging of striata.

2.1 Single-photon Emission Imaging

This section provides a brief discussion on single-photon emission nuclear medicine including a description of standard techniques and equipment. The focus is single-photon emission imaging, qualitative and quantitative. An understanding of single-photon emission imaging provides a foundation for understanding conjugate imaging. While planar emission imaging is the backbone of this work, single-photon emission computed tomography (SPECT) is also presented, as it is a standard nuclear medicine technique used in brain imaging, the primary clinical application relevant to this dissertation. Measurements from the modified conjugate imaging technique are compared to measurements from quantitative SPECT studies.

Single-photon emission imaging is the form of nuclear medicine imaging that uses single photon emitting radiopharmaceuticals. There are two primary forms of single-photon emission imaging, planar single-photon emission imaging and SPECT. Planar single-photon emission imaging generates two-dimensional (2-D) projections (images) of a distribution of a radiopharmaceutical. SPECT generates three-dimensional (3-D) volume representations of a radiopharmaceutical distribution. Such 3-D volumes are generated through tomographic reconstruction of 2-D projections acquired from multiple angular positions about a body being imaged. Single-photon emission imaging uses gamma cameras to detect the photon emissions.

This section first discusses three aspects of single-photon emission image generation: photon interactions, radiopharmaceuticals, and gamma cameras. It then describes the difference between qualitative and quantitative images, and the mathematical basis of qualitative and quantitative planar and SPECT imaging. The section concludes with a brief discussion of the task of quantitative assessment of striata via the single-photon emission technique of conjugate imaging.

2.1.1 Photon Interactions

Single-photon emission imaging exploits photon interactions in matter to generate in vivo images of the body and its physiological functions. The energy of photons used to generate single-photon emission images range from 80 keV to 800 keV. These photons interact in matter primarily through photoelectric absorption and Compton scattering. The probability of each interaction occurring is dependent on the energy of photons and the material it is transversing.

2.1.1.1 Photoelectric Absorption

In photoelectric absorption, an incident photon is absorbed by an atom. The energy of the photon is transferred to an inner shell electron, and the electron is ejected from the atom with kinetic energy equal to the incident photon energy minus the binding energy of the electron. Figure 2.1 is a sketch of the photoelectric process. Photoelectric absorption is a dominant process for low energy photons with energy 100 keV and below. The cross section of photoelectric absorption is proportional to $Z^3/h\nu^3$ where $h\nu$ is the photon energy. Consequently, the probability that a photoelectric interaction will occur




Figure 2.1. Photoelectric absorption. An incident photon of energy $h\nu$ is absorbed by the atom and an inner shell electron is ejected. The electron is ejected at angle θ and has energy T , which is equal to the energy of the incident photon minus the electron's binding energy.

increases significantly in materials with high atomic numbers and decreases with an increase in photon energy.

2.1.1.2 Compton Scattering

In Compton scattering, an incident photon collides with a loosely bound electron. The photon scatters off the electron at angle ϕ transferring some of its energy to the electron. The electron is scattered at angle θ with kinetic energy equal to the energy difference between the incident photon and the scattered photon. Figure 2.2 is a sketch of the Compton scatter process.




Figure 2.2. Compton scattering. An incident photon scatters off a loosely bound electron. The scattered photon scatters off at angle ϕ with energy of $h\nu'$. The electron is ejected at angle θ and has energy T , which is equal to the difference between the energy of the incident photon $h\nu$ and the energy of the scattered photon $h\nu'$.

The energy of the scattered photon is expressed as

$$h\nu' = \frac{h\nu}{1 + \left(\frac{h\nu}{m_0 c^2} \right) (1 - \cos \phi)}, \quad (2.1)$$

where $h\nu'$ is the energy of the scattered photon; $h\nu$ is the energy of the incident photon; $m_0 c^2$ is the rest mass of an electron; and ϕ is the scattering angle of the scattered photon.

The kinetic energy of the scattered electron is

$$T = h\nu - h\nu'. \quad (2.2)$$

Compton scattering is the dominant form of interaction for photons with energy ranging from about 100 keV to about 10 MeV making Compton scattering the dominant form of interaction in single-photon emission imaging.

2.1.2 Radiopharmaceuticals

Radiopharmaceuticals are pharmaceutical agents labeled with radioisotopes. They are designed to accumulate in a specific organ or volume of tissue through participation in a physiological function that is specific to that organ or volume of tissue. In this manner, a radiopharmaceutical transports a radioisotope to the target organ or volume. The radiopharmaceuticals used in imaging the striata participate in the process of neurotransmission of dopamine.

Radiopharmaceuticals suitable for routine clinical single-photon imaging have several basic properties in common. The production of single-photon emission radiopharmaceutical is relatively inexpensive. The chemical combination of the pharmaceutical agent and the radioisotope does not change the chemical reactivity of the pharmaceutical agent as regards to its participation in a designated physiological function. In other words, the radiopharmaceutical participates in the function of interest without disrupting that function. A suitable clinical pharmaceutical has a high affinity for the target volume, and conversely, a low affinity for non-target organs or tissue, especially non-target organs in close proximity to the target volume. Affinity refers to the degree to which the radiopharmaceutical binds to the target. In addition to having high affinity for a target, the radiopharmaceutical has a low dissociation rate. The dissociation rate is the rate at which the radiopharmaceutical unbinds from the target site. Before

imaging can take place, the kinetics of the radiopharmaceutical must reach equilibrium. A pharmaceutical with a low dissociation rate remains bound to the target site long enough for equilibrium to be reached and the imaging study completed. Thus, it is important that the rate of dissociation be relatively small. In addition, the effective half-life of the radioisotope must be long enough to allow the completion of an imaging procedure and short enough to keep the radiation dose to the patient low.

Once the radiopharmaceutical is administered, photons emitted from the radiopharmaceutical travel through the body. Many interact in the body through photoelectric absorption and Compton scattering. Some photons that exit the body, unscattered or scattered, are detected by a gamma camera. This detection is the first step in the image generation process.

2.1.3 Gamma Camera

2.1.3.1 Components

The gamma camera, also called the Anger Camera after its inventor, Hal Anger, is the most commonly used device in nuclear medicine (Anger 1958). The gamma camera consists of a collimator, scintillation crystal, lightguide, photomultiplier tubes (PMTs), preamplifiers, amplifiers, positioning circuitry, and pulse-height analyzers (PHAs). Figure 2.3 is an illustration of a gamma camera system, which includes a computer for image formation and display. The components of the gamma camera are encased in housing made of materials with high atomic numbers, for example lead or tungsten. This housing shields the camera from ambient radiation. The collimator is an array of holes cast in a material with a high atomic number Z and high-density ρ , similar to the




Figure 2.3. Basic components of a gamma camera system. The photon emitted from the radioactivity distribution within the patient enters the gamma camera through the collimator. The photons interact in the scintillation crystal, which results in the emission of light photons. These light photons enter PMTs through the light guide. The PMTs convert light signal into an electric signal. The energy discrimination and positioning electronics use the electric signal to determine energy and position to generate images. The images are displayed on the computer. Figure courtesy of Larry Zeng, Ph.D., at the University of Utah. (Zeng 2003)




Figure 2.4. Collimators and cross section of a hexagonal collimator. The septa are made of a high atomic number and high-density material. The cross section resembles a honeycomb. Photons are emitted isotropically from activity distribution within the patient. Not all the photons are used to generate the image. A small percentage of photons that encounter the collimator are used to generate the image. Photons traveling through the collimator are usable. Photons traveling at an angle through the collimator have a high probability of being absorbed in the septa. These photons do not contribute to image formation. Figure courtesy of Larry Zeng, Ph.D., at the University of Utah. (Zeng 2003)

shielding. The walls of the holes are called septa. Alloys of lead are the most common materials used as septa. Figure 2.4 includes an example of a cross section of a parallel-hole collimator with hexagonal-shaped holes. The collimator allows photons traveling in the direction parallel to the collimator holes to pass through to the scintillation crystal and absorbs photons traveling in other directions. Photons passing through the collimator encounter the scintillation crystal. When a photon interacts with the atoms in the crystal, the crystal emits visible light photons approximately proportional to the amount of energy absorbed due to the photon interacting in the crystal. A light guide is an optical coupling compound that aids in the transmission of light photons from the scintillation crystal to the camera's PMTs by minimizing reflection loss. A PMT is a vacuum tube that contains an anode, photocathode, and dynodes. A PMT generates an electric current (output signal) approximately proportional to the number of visible light photons that enter the PMT. The preamplifiers and amplifiers are electronics used to propagate and amplify PMT output signals prior to digitization, discrimination, and storage. The positioning circuitry estimates, from the output signals of the camera's PMTs, the position of the gamma photon interaction in the crystal. The electronics of the PHA reject detected photons with energies not within a specified range. Data collected from the positioning circuitry and PHA are used to generate images. These images are digitally displayed on computer monitors.

2.1.3.2 Gamma Camera Function

In a typical single-photon emission study, the radiopharmaceutical is administered to a patient orally or intravenously. The radiopharmaceutical is given time

to localize in the desired organ or volume in the patient. Gamma photons are emitted isotropically from the activity distribution. Many of these photons go undetected, as they never come in contact with the gamma camera. Photons with the potential to be detected are those that enter the gamma camera through its collimator or aperture. The collimator mechanically selects photons for possible detection. Photons that travel through the holes reach the scintillation crystal. Some interact with the crystal; some pass through it. Ideally, photons that impinge on the collimator at an angle are absorbed in the walls of the honeycomb-like structure of the septa. Figure 2.4 includes an illustration of collimation. Collimated photons that interact in the scintillation crystal transfer energy through Compton or photoelectric interactions to electrons in the scintillation crystal. These electrons, in turn, lose their energy by interacting with the crystal atoms through excitation and ionization. The excited atoms return to ground state by emission of visible photons. The number of visible photons emitted is proportional to the amount of energy transferred to the atom indirectly by the incident gamma photon. The light output of the crystal is transferred to the PMTs via the lightguide. A gamma camera usually has several PMTs. The exact number depends on the size of the scintillation crystal and the size of the PMTs. In a typical gamma camera, there are 37 to 93 PMTs in a hexagonal array. The lightguide distributes the light output from photon interactions in the crystal across several PMTs. The photocathode of each PMT converts the visible light photons into electrons. About four to six light photons must strike the photocathode to generate one electron. The electrons produced at the photocathode are then multiplied as they propagate through the dynodes. Dynodes are located along the PMT. Each dynode is held at a higher positive voltage than the previous dynode. When an electron strikes a dynode,

two or four new electrons are generated. These new electrons add to the number of electrons traveling through the tube to the anode. Although each PMT outputs electrons multiplicatively greater in number than the electrons released at its photocathode, the number of electrons needs to be further amplified. A preamplifier, attached to the PMT, enhances the electrical signal so the signal is sufficiently strong that it can pass through a cable to the main amplifier. The amplifier increases the signal by about three orders of magnitude. The positioning circuitry linked to the PMTs determines the position of the energy deposition of a photon in the crystal by summing weighted signals from the PMTs. The amount of light detected from each PMT is in direct relation to the PMT's proximity to the location of the gamma-ray interaction in the scintillation crystal, i.e. the closer the PMT is to the event, the larger the PMT signal. As a result, the PMT signals are weighted in relation to their location to the event. The sum of the PMT output signals, known as the Z pulse, is proportional to the energy of the initial photon interacting in the crystal. The PHA accepts Z pulses within a specified energy range and rejects Z pulses outside the energy range. This range is referred to as an energy window. The detected photons within the photopeak window are the energy-window events used to generate images. The size of the energy window is dependent on the energy resolution of the gamma camera. The poorer the energy resolution, the wider the window. If the Z pulse is within a specified energy range, the Z pulse and the corresponding location are digitized and stored.

Typical nuclear medicine images are stored as $n \times n$ image matrices, where each image element is called a pixel. Each pixel corresponds to a small continuous range of (x, y) locations in real object space. The image matrix is similar to a counter; a count is

added to a pixel if the Z pulse occurred within the real object coordinate range represented by the pixel. Nuclear medicine images are usually stored as 64 x 64, 128 x 128, 256 x 256, or 512 x 512 matrices. While the physical image size remains the same; the pixel size varies according to the number of total pixels. The pixels in a 64 x 64 image matrix represent a larger range of spatial points in real object space than the pixels in a 512 x 512 image matrix. As a result, a 512 x 512 image records better spatial resolution, if the device hardware can accurately record such resolution.

2.1.3.3 Gamma Camera Performance

Not all gamma cameras are designed the same. As a result the performance of gamma cameras varies. The design of a gamma camera is largely dependent on the energy of photons to be imaged. The collimator and scintillation crystal thickness are selected to correspond with the energy of the photon to be imaged. The collimator and the scintillation crystal have great influence on the performance of a gamma camera.

Collimators play a significant role in the performance of gamma cameras as they are the entry point into the gamma camera. While the performance of the gamma camera is influenced by several collimator characteristics, such as hole diameter and bore, septal wall thickness is of interest in this research. The selection of septal wall thickness is also dependent on the energy of the photon to be imaged. Imaging high-energy photons requires thick septal walls to minimize septal penetration – photons permeating the septal walls. Thinner septal walls are acceptable for imaging lower energy photons. However, the selection of an adequate collimator becomes more complicated, if the energy spectrum of the radiopharmaceutical has energy peaks higher than the principle energy

peak used for image generation. This is because higher energy gammas may downscatter into the photopeak of the spectrum. This is the case with ^{123}I labeled radiopharmaceuticals. In addition to the primary emission photon of 159 keV, 3% of ^{123}I gamma photons are in the energy range of 248 keV - 784 keV. For ^{123}I quantitative single-photon emission imaging studies of striatal dopamine neurotransmission, in which spatial resolution is important, a high-resolution collimator is needed. Typically, a low energy high-resolution parallel-hole collimator is used for such studies as imaging focuses on ^{123}I 's primary gamma emission of 159 keV. However, the septal walls of such a collimator are too thin to absorb the energy of the high-energy gamma emissions of ^{123}I . Thus, septal penetration occurs and the high-energy gamma emissions may interact in the scintillation crystal increasing image noise.

Thickness of the scintillation crystal also influences the performance of a gamma camera, as it is the essential component for detection. A thicker crystal has higher detection efficiency than a thinner one of the same material. More photons are absorbed in a thicker crystal because photons that would have escaped a thinner crystal are now interacting in the crystal. A thicker crystal also decreases the spatial resolution of the gamma camera. Although more photons are absorbed in a thicker crystal, light photons leaving the crystal and entering the PMTs are spread out a larger area. In addition, many photons transfer energy through several Compton scatter events before escaping the crystal or eventually being absorbed by photoelectric interactions at distance from the location of the initial interaction.

2.1.4 Attenuation and Scatter

Ideally, images generated from gamma cameras would be a true spatial and quantitative representations of the activity distributed throughout the target volume. However, in practice the images generated from gamma cameras are less-accurate representations due to attenuation and scatter. That is, the grayscale values of images generated from such data do not reflect absolutely the values of radioisotope uptakes in particular volumes of interest.

Primary photons are photons emitted from a source of radioactivity that have not participated in an interaction. Primary photons that are absorbed in the scintillation crystal via the photoelectric effect provide useful spatial information about the activity distribution. Unfortunately, a large number of primary photons do not reach the camera due to Compton scattering or photoelectric absorption in the body. The loss of these photons is called attenuation. However, many photons are scattered in the body and are then detected by the gamma camera. Some photons exiting the body, unscattered or scattered, scatter in the camera before detection. The Z pulses of many of these scattered photons are often below the specified threshold and are therefore rejected by a PHA. Due to the finite energy resolution of the gamma camera, some Z pulses of scattered photons are not rejected. This is problematic because scattered photons are generally located a distance away from the Compton scatter event. The inclusion of scattered photons with approximately 159 keV energy degrades image quality as scattered photons do not provide accurate spatial information of the radiopharmaceutical distribution. In the case of radiopharmaceuticals that emit photons with energies higher than those of the primary photons, contributions from scattered higher-energy photons contaminant data from

primary photons. Such higher-energy photons are energetic enough to penetrate collimator septa and to deposit sufficient energy in the scintillation crystal that they contribute to the image. ^{123}I is an example of a radioisotope that emits photons with a range of energies above the primary photon emission of 159 keV. ^{123}I is an integral part of this dissertation and is discussed in more detail in section 2.4.

Image data that lacks data from primary photons due to attenuation and includes scattered photon data are qualitative images. Even though qualitative images do not reflect fully the amount of activity distributed in volumes of interest, they do have clinical value. Qualitative single-photon emission images are a significant component of the management of patients in present nuclear medicine clinics. Although qualitative images can sometimes be sufficiently informative, quantitative images are at times necessary as they provide more definitive information, such as absolute quantities or ratios. For example, a longitudinal study of a disease documented only with sequential qualitative images may require that significant changes in the images exist before clinicians are able to discern any change. Use of only qualitative images proves to be problematic if a disease progresses slowly. Small changes are lost in qualitative images. Quantitative imaging provides a means of recognizing change and of measuring its amount. Quantitative imaging involves compensating for the effects of attenuation and scatter. Parkinsonism is an example of a disease that progresses slowly. Studies have suggested that quantitative imaging can be used to monitor successfully Parkinsonism and related drug therapies (Ichise and Ballinger 1996; Bourguignon, Pauwel et al. 1997; Messa, Volonte et al. 1998; Booij, Tissingh et al. 1999; Benamer, Patterson et al. 2000; Booij, Speelman et al. 2001; Costa 2001; Grosset 2001; Tatsch 2001; Tzen, Chin-Song et

al. 2001; Frey 2002; Pinborg, Videbaek et al. 2002; Poewe and Wenning 2002; Abella 2003; Burn and O'Brien 2003; Catafau 2003; Jersild and Trevino 2003; Leenders 2003; Morrish 2003; Oertel, Gerstner et al. 2003; Piccini 2003; Pirker 2003; Seibyl 2003).

Thus, of interest, are quantitative approaches to striatal uptake of relevant radiopharmaceuticals.

2.1.5 Qualitative Single-photon Emission Imaging

Single-photon emission imaging, planar or tomographic, is based fundamentally on the attenuated Radon integral. The standard geometry and coordinate systems of patient and gamma camera set-up for single-photon emission imaging are illustrated in figure 2.5. Two coordinate systems, one based on the frame of reference of the patient and the other based on the frame of reference of the gamma camera, are depicted. The patient's frame of reference is specified by the spatial coordinates $\mathbf{x} = (x, y, z)$ and is assumed to be stationary. The z coordinate is out of the page. The gamma camera's frame of reference, specified by the spatial coordinates (x_r, y_r, z) , rotates about the z -axis with a rotation angle ϕ . The attenuated Radon integral has the following form:

$$E(x_r, z) = \int_L f(x) e^{-\int_{L(x)} \mu(x) dy_r'} dy_r, \quad (2.3)$$

where $E(x_r, z)$ is the emission measurement, $f(x)$ is the spatial distribution of the radiopharmaceutical, $\mu(x)$ is the distribution of attenuation coefficients, and dy_r is the differential length along the line L in the y_r direction. Equation 2.3 consists of two




Figure 2.5. Acquisition of emission data. The geometry and coordinate systems of patient and gamma camera set-up used to acquire an emission measurement. The patient is represented by the object. The detector in this illustration is acquiring an emission measurement from point x . The detector detects photons that are attenuated as they pass through the body from x to the detector. $\mu(x)$ is the distribution of attenuation coefficients within the object along the line segment $L(x)$. The emission measurement is described by the attenuated Radon integral. The Radon integral sums the product of the radiotracer

distribution $f(x)$ and the fraction of transmitted photons from x along L , $e^{-\int_{L(x)} \mu(x) dy_r}$. Figure courtesy of Larry Zeng, Ph.D., at the University of Utah. (Zeng 2003)

integrals. The exponential integral is the summation of the attenuation coefficients along the line segment $L'(\mathbf{x})$ as the tissue along $L'(\mathbf{x})$ may not be homogenous. The integral of $f(\mathbf{x})$ and the exponential term (the attenuation term) are the total radioactivity along the line segment $L(\mathbf{x})$ as seen by the gamma camera. Depending on the image study the gamma camera or gamma cameras can be rotated about the z_r axis through angle ϕ . At each angular location ϕ an emission measurement, an attenuated Radon integral, is acquired for each pixel of the image. Each measurement acquired by the camera in a specific angular location is a measurement of radiation traveling in the y_r direction.

2.1.6 Planar Emission Imaging

Planar emission imaging is the process by which images of single views of a radiopharmaceutical distribution are generated. An image of a single view of a 3D radiopharmaceutical distribution is called a planar projection. The gamma camera is placed at a fixed angular position ϕ and emission data are collected. The emission measurement in each pixel of a planar projection (planar image) is described approximately by the attenuated Radon integral, representing the sum of the activity along the line extending from a point within the distribution to the face of the gamma camera. Due to the integration of the activity along the line, the 3D distribution of the activity in the patient is depicted as a 2D image. One planar projection of the patient comprises the collection of attenuated Radon integrals in the y_r direction.

Clinically, planar emission imaging is used to observe structure and function through static studies and dynamic studies. A static study is the acquisition of a single projection; a dynamic study, in this context, is the collection of a series of planar images

acquired intermittently over a specific time interval. A thyroid study is a common static study used to assess the function of a thyroid. Generally, three planar images are acquired – anterior, left anterior oblique, right anterior oblique. A common dynamic study is a renal study. A renal study is used to evaluate the function of a kidney by observing the kidney extracting the administered radiopharmaceutical from the blood and excreting it to the bladder. Usually, several sequential planar images are required for accurate assessment, as time-activity curves of total activity in the kidney or kidneys are needed.

2.1.7 SPECT

Single-photon emission computed tomography is the process of generating, from a collection of planar projections, cross-sectional views of a radioactivity distribution. These views can be combined to form a 3D representation of the distribution. SPECT is a two-step process: 1) image data acquisition and 2) image data reconstruction. The complete process of SPECT requires one or multiple gamma cameras, a means to display images, a gantry for camera rotation about the patient, and software to perform image reconstruction. For step one, a gamma camera or multiple gamma cameras are rotated about the patient (z-axis) through angle ϕ to acquire various projection views. The angular range is generally 180 or 360 degrees. The range depends on the organ imaged. Most commercial SPECT systems have two gamma cameras on a rotating gantry. Single and triple gamma camera systems are also available. Figure 2.6 is an illustration of a dual head gamma camera system acquiring projection data from different angles. Once the projections have been collected, the data they contain are manipulated to generate 2D




Figure 2.6. SPECT projection acquisition. Dual-head gamma camera system acquiring eight projections at different angles. A projection acquired for each gamma camera. Figure courtesy of Larry Zeng, Ph.D., at the University of Utah. (Zeng 2003)

transaxial cross-sectional images. In addition to transaxial slices, coronal and sagittal slices are often generated for SPECT studies by reslicing the volume of transaxial slices.

The process of generating transaxial slices is called image reconstruction. Image reconstruction is an inverse problem that uses the measured data, the projections, to generate the transaxial slices of the activity distribution. Filtered backprojection is the most common reconstruction method applied in clinical settings (King and Farncombe, 2003; Rosenthal, 1995; Tsui, 1996; Zeng, 2001). Iterative image reconstruction schemes are used in most research labs and their use is becoming more popular for clinical studies (Vandenberghe, D'Asseler et al. 2001; Zeng 2001).

As stated previously, in section 2.2.4, a planar projection is described mathematically as a collection of attenuated Radon integrals. However, the attenuated Radon integral has proved too complex for practical use in clinical image reconstruction. Consequently, the attenuated Radon integral is often replaced with a Radon integral in the reconstruction process. If the attenuation coefficient distribution $\mu(\mathbf{x})$ in equation 2.3 is assumed to be zero, then equation 2.3 is the Radon integral of $f(\mathbf{x})$. That is, equation 2.3 becomes

$$E(x_r, z) = \int_L f(x) dy_r, \quad (2.4)$$

where $f(\mathbf{x})$ is the spatial distribution of the activity and $E(x_r, z)$ is the emission measurement along the single line segment L measured at the single point (x_r, z) .

Although the reconstruction calculation is less rigorous and less accurate with the use of the Radon integral than with the use of the attenuated Radon integral, the resulting SPECT images are clinically usable representations of the true activity distribution in the body. They are not quantitatively accurate and are sometimes misleading in the clinic setting, but nonetheless they are often diagnostically useful.

SPECT studies are used routinely in nuclear medicine departments. Often times, SPECT studies are performed when planar studies reveal abnormalities that require a better assessment. Bone SPECT studies are examples of studies performed in addition to corresponding planar bone studies. SPECT studies are also performed as primary diagnostic procedures. Such is the case for myocardial perfusion and brain perfusion assessment. Myocardial perfusion SPECT studies are the most common SPECT studies performed. They are used to assess a patient's heart at rest and under stress. Brain perfusion SPECT is performed to identify clinical indications of cerebrovascular disease, dementia, and seizure.

2.1.8 Quantitative Single-photon Emission Imaging

Quantitative single-photon emission imaging incorporates correction techniques to estimate and correct for attenuation and scatter. Attenuation correction uses an attenuation map, usually now generated from transmission imaging, to account for counts

lost due to attenuation. Scatter correction estimates the fraction of scattered photons detected and subtracts that fraction from the total number of photons detected.

Transmission imaging generates estimates of attenuation coefficients by using an external source. The external source is placed in such a position that the emissions of the source must pass through the patient to reach the detector. This configuration is illustrated in figure 2.7. The transmission measurement, $T(x_r, z)$, at a specific point (x_r, z) is described as

$$T(x_r, z) = \int_L \mu(\mathbf{x}) dy_r, \quad (2.5)$$

where $\mu(\mathbf{x})$ is the attenuation coefficient at \mathbf{x} . Equation 2.5 is the Radon integral of $\mu(\mathbf{x})$ along the single line L measured at the single point (x_r, z) . The collection of Radon integrals in the y_r direction is a transmission projection. Transmission tomography uses the collection of transmission projections in all y_r directions for a fixed z to generate a transmission image of a transaxial slice of the body. A single transmission projection is used to correct attenuation in planar imaging studies. Transmission tomography is used to correct attenuation in SPECT studies.

Classic conjugate imaging is a quantitative planar imaging technique that uses data acquired by two opposed gamma cameras to estimate the activity of a radiopharmaceutical distribution at an unknown depth (Barnes 1996; King and Farncombe 2003). Attenuation is corrected by using a measured effective attenuation map determined by the transmission data. Classic conjugate imaging will be discussed in




Figure 2.7. Acquisition of transmission data. Geometry and coordinate systems for patient and detector associated with a transmission projection. The patient is represented by the object with an attenuation coefficient $\mu(\mathbf{x})$. The patient is placed between the detector and source. The photons emitted from the source travel through the body. Some photons exiting the body are detected by the detector. The transmission measurement along L is described by the Radon integral in equation 2.5. Figure courtesy of Larry Zeng, Ph.D., at the University of Utah. (Zeng 2003)

more detail in section 2.2 as it is the foundation for the modified conjugate imaging technique presented in this dissertation. Although transmission maps are used in SPECT studies, an attenuation constant has often been uniformly applied for volumes where the distribution of attenuation coefficients is relatively uniform (Tsui 1996; Catafau 2003). This is done, for example, in the reconstruction of brain studies. For volumes with a heterogeneous distribution of attenuation coefficients, an attenuation map is applied.

Scatter contributes significantly to image data. Jaszczak estimated that about 40% of acquired image data are scattered radiation when ^{99m}Tc radiopharmaceuticals are used (Jaszczak, Greer et al. 1984; Rosenthal, Cullom et al. 1995). This percentage can be larger for radiopharmaceuticals that also emit photons with higher energies than its primary photon emission, such as ^{123}I (Gilland, Jaszczak et al. 1994; Dobbeleir, Hambye et al. 1999). These high-energy photons can downscatter into the photopeak window increasing the scatter fraction within that window. Although scatter-correction schemes have been developed and applied to SPECT studies, to date, scatter-correction techniques have yet to be implemented routinely in clinical settings.

King and Farncombe classify scatter correction techniques into two main categories - scatter-estimate and reconstruction-based scatter compensation (RBSC) techniques (King and Farncombe 2003). Techniques that fall into the scatter-estimation group estimate scatter contribution using acquired emission data. RBSC techniques incorporate a modeled scatter response function into an iterative reconstruction scheme. RBSC techniques use principles of scatter interactions to estimate the scatter contribution. Since RBSC techniques are used with SPECT studies, these techniques will not be discussed further.

The scatter-estimation group can be sub-divided into energy-distribution and spatial-distribution methods. Energy-distribution methods estimate scatter contribution in the photopeak window using emission data in other energy-windows. This method is a pixel-based method, i.e. it is applied to each image pixel. Several energy-distribution methods have been used, but the Compton scatter window and triple-energy window (TEW) techniques are two of the most popular. Jaszczak et al. introduced the Compton scatter window technique in which emission data from an energy window below the photopeak window are used to correct for scatter (Jaszczak, Greer et al. 1984). Ogawa et al. introduced the TEW scatter technique that uses emission data acquired in two energy windows, one abutting the lower boundary of the photopeak window and the other abutting the upper boundary of the photopeak window (Ogawa, Harata et al. 1991). This technique, detailed in this work, is often applied to single-photon emission imaging studies using radiopharmaceuticals with high-energy photon emissions, in addition to their primary photon emission. In the spatial-distribution method, emission data acquired in the photopeak window is used to estimate the source distribution and to model the source contribution to the scatter distribution. In the convolution-subtraction method, an example of a spatial-distribution method, a model of the scatter distribution function is convolved with acquired emission data to obtain an estimate of the scatter distribution. This scatter distribution is subtracted from the emission data.

Scatter contributions are not specifically addressed in the methodology of classic conjugate imaging. To account for scatter in classic conjugate imaging, one can modify the technique to include build-up factors or the projection data must be corrected using scatter-correction techniques (Wu and Siegel 1984; Barnes 1996; King and Farncombe

2003). The build-up factors can be calculated from empirical measurements or Monte Carlo computer simulations. However, the inclusion of build-up factors in the classic conjugate imaging technique increases the complexity of conjugate imaging. Although most scatter-correction techniques are applied primarily in SPECT studies, some can be modified for application to planar projection data.

The image data generated in quantitative imaging more accurately reflect the true activity distribution than do data of qualitative techniques. The primary outcome sought in use of quantitative imaging, however, is not an image, but an estimate of the radioactivity in a radiopharmaceutical distribution. The estimation of activity can be presented as an absolute activity measurement or as one of a variety of other parameters such as binding potential (Karesh 1996). The binding potential is the ratio between specific binding and non-specific binding. It is the most common measurement made in quantitative imaging.

2.2 Classic Conjugate Imaging

Classic conjugate imaging is a non-tomographic nuclear medicine technique that uses gamma camera projection data from two diametrically opposed views to estimate the activity in a source embedded at an unknown depth in an attenuating medium. The technique can also be used to estimate the sum of activities in multiple sources embedded at unknown depths in an attenuating medium. The projection data can be acquired simultaneously by two gamma cameras or sequentially by one gamma camera, if the tracer dissociation rate is sufficiently slow. For the remainder of this dissertation, it is assumed two gamma cameras are used. Three sets of projection data – emission,

transmission, and flood - are acquired for each view. The flood data are used to determine the camera sensitivities for each camera. The transmission data are used to estimate the transmission factor for the imaged object. The emission data are geometrically averaged and used together with the calculated camera sensitivities and effective attenuation of the medium to estimate the activity of the source.

At the roots of conjugate imaging is conjugate counting introduced in 1937 by Robley Evans (Evans 1937). Evans used anterior and posterior count data collected by a Geiger counter to estimate the total activity of radium deposits in the bones of radium dial painters. For each measurement taken, the body of a radium dial painter was curved in an arc of 1 meter radius of curvature and the Geiger counter was positioned at the radius of curvature. The advent of imaging, rectilinear and planar, introduced new imaging geometries. Sorenson evaluated the influence of source-to-detector distance on the estimation of activity using conjugate imaging for different imaging set-up geometries including Evans arc geometry (Sorenson 1974). The geometry investigated by Sorenson that is most relevant to this work is a planar imaging geometry, i.e. parallel-hole collimator imaging.

The unique aspect of conjugate imaging is the one-dimensional treatment used to determine quantitative measurements. In ideal parallel-hole collimator imaging in which the detection of scattered photons and depth dependent blurring are considered negligible, the area circumscribed by each bore in the collimator can be treated as an individual detector. As such the counts in an image pixel can be determined by one-dimensional analysis in which the depth of activity is the only variable. The derivation of conjugate imaging presented in this chapter and my discussion of an extension in Chapter VI

assume an ideal parallel-hole collimator imaging geometry, which results a one-dimensional treatment of the activity estimation problem and in pixel-by-pixel calculations.

This section discusses the methodology and mathematics of classic conjugate imaging. The discussion includes conjugate imaging of planar and volume sources, applications of classic conjugate imaging, and modified conjugate imaging techniques based on classic conjugate imaging.

2.2.1 The Planar Model

We assume a planar source of radioactivity embedded in a medium with a thickness T as diagrammed in figure 2.8. The medium is uniformly attenuating with an effective linear attenuation coefficient μ . (That is the transmission factor at depth d is $e^{-\mu d}$.) The planar source A is at an unknown depth d . Emission data are acquired by placing the object between two opposed gamma cameras as shown in figure 2.8. The cameras in figure 2.8 are denoted as “Camera L” (left) and “Camera R” (right). This notation will be used throughout the remainder of the dissertation.




Figure 2.8. Conjugate Imaging: Acquisition of emission data. Planar source with activity A is located at a depth d within a medium of thickness T. The medium is uniformly attenuating with an effective linear attenuation coefficient μ . The object is placed between two gamma cameras, camera L and camera R for emission data acquisition.

The emission data can be expressed as

$$E_L = C_L (Ae^{-\mu d}) \quad (2.6)$$

$$E_R = C_R (Ae^{-\mu(T-d)}), \quad (2.7)$$

where E_L and E_R are the emission data from the left and right cameras, respectively; C_L and C_R are the camera sensitivities of the left and right cameras, respectively; A is the source activity; T is the total object thickness; and d is the depth of the source as measured from the right surface of the medium. The units for the emission data are counts per second (cps). The camera sensitivity is expressed as counts per second divided by disintegrations per second (cps/dps). The activity is expressed as dps. Equations 2.6 and 2.7 can be used to solve for activity A. The result is

$$A = \left(\frac{E_L E_R}{C_L C_R e^{-\mu T}} \right)^{1/2}. \quad (2.8)$$

Since equation 2.8 is not dependent on the depth parameter d , conjugate imaging proves useful when a quantitative measurement of a source is desired but the depth of the source is unknown. Equation 2.8 cannot be solved without knowledge or measurements of the camera sensitivities, C_L and C_R , and the transmission factor, $e^{-\mu T}$. Knowledge of these requires two additional sets of data, flood and transmission data.

Flood data are acquired by placing a uniform flood source between and equidistant from the two cameras as seen in figure 2.9. The flood source has a known




Figure 2.9. Conjugate Imaging: Acquisition of flood data. A flood source with activity A_F uniformly distributed is placed between camera L and camera R. The data acquired are used to determine the camera sensitivities.

activity A_F . The flood data can be expressed as

$$F_L = C_L A_F \quad (2.9)$$

$$F_R = C_R A_F . \quad (2.10)$$

The camera sensitivities, C_L and C_R , can then be obtained as

$$C_L = \frac{A_F}{F_L} \quad (2.11)$$

$$C_R = \frac{A_F}{F_R} . \quad (2.12)$$

Transmission data are needed to compute the transmission factor of the medium, $e^{-\mu T}$.

Transmission data are obtained by each camera, as the data acquired by each camera are most likely different as the cameras likely have different sensitivities. When transmission data are collected by Camera L, the flood source is placed between Camera R and the object. When collected by Camera R, the flood source is placed between Camera L and the object. This is depicted in figure 2.10. In the absence of radioactivity in the object, the transmission data can be expressed as

$$T_L = C_L A_F e^{-\mu T} \quad (2.13)$$

$$T_R = C_R A_F e^{-\mu T} . \quad (2.14)$$

The transmission factor $e^{-\mu T}$ can be determined from the above equations by computing




Figure 2.10. Conjugate Imaging: Acquisition of transmission data. The object with planar source located at depth d is placed between camera L and camera R. The flood source is placed between camera L and the object to acquire transmission data for camera R. Transmission data from camera L are acquired with the transmission source between camera R and the object.

the geometric mean of the transmission data, T_L and T_R , and manipulating the resulting formula. When this is done, one obtains

$$e^{-\mu T} = \frac{1}{A_F} \sqrt{\frac{T_L T_R}{C_L C_R}}. \quad (2.15)$$

Although ideally, transmission data would be acquired without the presence of radioactivity in the object, clinically this is usually not feasible. Accurate registration of transmission and emission data is necessary and accurate repositioning of the patient is difficult. As a result, the transmission data measured in most cases contain data from a source or sources together with data from the transmission source. The measured transmission data are expressed as

$$T_L' = C_L A_F e^{-\mu l} + E_L \quad (2.16)$$

$$T_R' = C_R A_F e^{-\mu l} + E_R, \quad (2.17)$$

where the superscript ' indicates a measured quantity. The emission data must then be subtracted from the measured transmission data to obtain the true transmission data used in Equations 2.13-2.14.

2.2.2 The Volume Model

In this section, the extension of the planar source model to the volume source model is discussed. Activity injected into an object will not accumulate in a planar or point distribution. The activity will gather in an organ or tumor with some thickness. Self attenuation is important in sources with a measureable thickness, as potential counts may be lost due to gamma interactions within the source.

Consider activity A uniformly distributed throughout an attenuating medium of thickness t as seen in figure 2.11. The activity measured or seen by a gamma camera is referred to as effective activity, A_{eff} . An infinitesimal slice of the source as seen by an ideal detector can be expressed as

$$dA_{\text{eff}} = dA e^{-\mu x} \quad 0 \leq x \leq t, \quad (2.18)$$

where
$$dA = \left(\frac{A}{t} \right) dx. \quad (2.19)$$

In the above, dA_{eff} is the infinitesimal amount of effective activity in an infinitesimal slice dx; dA is the infinitesimal amount of activity in an infinitesimal slice dx; $e^{-\mu x}$ is the




Figure 2.11. Volume source with thickness t . Activity is distributed uniformly throughout an attenuating medium with attenuation coefficient $\tilde{\mu}$.

transmission factor of activity in slice dx at depth x ; A is the activity distributed throughout the source; and t is the source thickness.

The effective activity of the source is obtained by integrating equation 2.18 over the thickness of the source. The result is

$$A_{\text{eff}} = A \frac{e^{-t/2}}{f_t} \quad (2.20)$$

where

$$f_t = \frac{\mu t / 2}{\sinh(\mu t / 2)} \quad 0 < f_t < 1. \quad (2.21)$$

In the above, A_{eff} is the effective activity of the source; f_t is a factor that modifies (increases) the $e^{-\mu t/2}$ transmission factor for a source in the middle of thickness t ; and A is the activity uniformly distributed throughout the source. The transmission factor $e^{-\mu t/2}$

does not adequately represent the attenuation of activity distributed uniformly throughout thickness t . The factor f_t is a correction factor that modifies $e^{-\mu t/2}$ by accounting for effects of attenuation throughout t . The magnitude of f_t is dependent on source thickness. The thicker the source, the larger the value of f_t . For thinner sources, the self-attenuation factor has less significance. For example, assuming that 140 keV photons are distributed uniformly in a 3 cm thick slab of water and using an effective attenuation coefficient of 0.140 cm^{-1} , the self-attenuation factor f_t of this source represents less than a 1% effect (Mintzer, Aarsvold et al. 2000).

If we replace the planar source in figure 2.8 with a volume source with a thickness t , we obtain the geometry of figure 2.12.




Figure 2.12. A source with thickness t embedded in an object with thickness T . The source has activity A . The depth of the source is defined as half the source $t/2$.

Conjugate imaging of a volumetric source is based on the same acquisitions of emission, transmission, and flood data as used for conjugate imaging of a planar source. However,

it is necessary to include the self-attenuation factor f_t in models and computations. The following discussion presents the application of classic conjugate imaging to the quantification of the activity in a volume source. The emission data of a thick source can be expressed as

$$E_L = C_L A_{eff} e^{-\mu(r-d-t/2)} \quad (2.22)$$

$$E_R = C_R A_{eff} e^{-\mu(d-t/2)}, \quad (2.23)$$

where

$$A_{eff} = A \frac{e^{-t/2}}{f_t}. \quad (2.24)$$

Here, E_L and E_R are the emission data acquired from the left and right cameras, respectively; C_L and C_R are the camera sensitivities of the left and right cameras, respectively; A_{eff} is the effective source activity; and f_t is the self-attenuation factor of the thick source. Substituting the value of A_{eff} (equation 2.24) into equations 2.22 and 2.23 yields the following emission data expressions:

$$E_L = C_L \frac{A}{f_t} e^{-\mu(r-d)} \quad (2.25)$$

$$E_R = C_R \frac{A}{f_t} e^{-\mu d}. \quad (2.26)$$

Taking the geometric mean of equations 2.25 and 2.26 and solving for A , the activity, we obtain:

$$A = \frac{\sqrt{E_L E_R}}{\sqrt{C_L C_R} e^{-\mu T}} f_t. \quad (2.27)$$

As mentioned previously, the self-attenuation factor of 140 keV photons distributed uniformly in a 3 cm thick slab of water with an attenuation coefficient of 0.140 cm^{-1} has less than a 1% affect on activity estimation. As a result, we assume the self-attenuation correction factor is one for sources less than or equal to 3 cm. Typically, the self-attenuation correction factor f_t is usually considered equal to one for source thickness less than half the object thickness, (Shung 1992).

The transmission data and flood data are measured using the same techniques discussed in the planar model section.

2.2.3 Sum of Two Sources

The sum of activity from two sources separated by a known distance, s , can be estimated by use of conjugate imaging (Sorenson 1974). We consider two sources of activity embedded in a uniformly attenuating object of thickness T . This is illustrated in figure 2.13. For this derivation, the sources are assumed to be thin enough relative to the attenuation coefficient that their thickness can be ignored, i.e. $f_t \approx 1$.

The emission data acquired from the left and right cameras are expressed as

$$E_L = C_L \left(A_1 e^{-\mu d_1} + A_2 e^{-\mu(d_1+s)} \right) \quad (2.28)$$

$$E_R = C_R \left(A_1 e^{-\mu(T-d_1)} + A_2 e^{-\mu(T-(d_1+s))} \right), \quad (2.29)$$




Figure 2.13. Two sources separated by a distance s in a uniformly attenuating object. The sources have activities A_1 and A_2 . The object has a thickness T . The source with activity A_1 has a depth d_1 . The source with activity A_2 has a depth $T - (d_1 + s)$.

where E_L and E_R are the emission data; C_L and C_R are the camera sensitivities of the left and right cameras, respectively; and A_1 and A_2 represent the activities of the two planar sources in Figure 2.13. Taking the geometric mean of equations 2.28 and 2.29 and solving for the sum of A_1 and A_2 yields

$$(A_1 + A_2) = \left(\frac{E_L E_R}{C_L C_R e^{-\mu T}} \right)^{1/2} \left(1 + \frac{2A_1 A_2}{(A_1 + A_2)^2} (\cosh \mu s - 1) \right)^{-1/2}. \quad (2.30)$$

By comparing equations 2.30 with equation 2.8, one finds that the second factor in equation 2.30 is the correction factor for the use of classic conjugate imaging in the estimations of the sum of the activities in two distinct sources. The correction factor is dependent on the source activity ratio and the distance between the sources. Assuming that the source activity ratio is $r = A_1/A_2$, the correction factor in equation 2.30 becomes

$$\left(1 + \frac{2r}{(r+1)^2} (\cosh \mu s - 1)\right)^{-1/2}. \quad (2.31)$$

This expression has a maximum value when $r = 1$. This implies the largest error in determining of the sum of two sources in this setting occurs when the source activities are approximately equal (i.e. $r = 1$). When $r = 1$, the correction term is

$$\left(1 + \frac{1}{2} (\cosh \mu s - 1)\right)^{-1/2}. \quad (2.32)$$

Equation 2.32 can be used to determine the upper bound for error in estimating the total activity in two sources separated by a distance s . For 140 keV photons in water, a case in which the attenuation coefficient is 0.140 cm^{-1} , the total activity of two sources separated by 2 cm is overestimated by 1%. However, for a separation of 6 cm, the total activity is overestimated by more than 10% (Mintzer, Aarsvold et al. 2000).

Conjugate imaging may be applied to estimate the total activity of more than two sources (Evans 1937; Sorenson 1974). Evans found that if one assumed the activity from multiple sources was placed at a single mean depth, only a small error was introduced into the estimated activity. He estimated the error for radium dial painter studies to be less than 4% and considered the error in small animal studies to be negligible (Evans 1937).

2.2.4 Applications

Robley Evans introduced conjugate counting in 1937 as a quantitative technique used in radium poisoning investigations. Evans used conjugate counting to estimate the amount of radium deposited in the bones of living people, specifically radium dial painters (Evans 1937). Since that time, classic conjugate imaging has been used to estimate the uptake of activity in tumors or normal organs for dosimetry calculations, therapy response estimations, and improved diagnostic assessment. However, with the advent of quantitative SPECT investigations, the employ of classic conjugate imaging has decreased such that its use is now limited.

Whole-body conjugate imaging is used to estimate biodistribution and dosimetry of radiopharmaceuticals (King and Farncombe 2003). Mozley used conjugate imaging to characterize the in vivo distribution of Tc-99m TRODAT (Mozley, Stubbs et al. 1998).

Several investigators have used conjugate imaging to estimate tumor uptake and uptake in normal surrounding tissues or organs. Thomas et. al. used conjugate imaging to estimate activity uptake in thyroid lesions to improve therapeutic treatment plans (Thomas, Maxon et al. 1976). Koral et al. applied conjugate imaging using pinhole collimators to the estimation of radiation dose to neck tumors (Koral, Adler et al. 1982). Hammond et al. presented conjugate imaging data on liver tumor uptake estimating error at less than 10% (Hammond, Moldofsky et al. 1984). Eary et al. used conjugate imaging to estimate the internal radiation dose to tumors and normal organs for radiolabeled antibody therapy (Eary, Appelbaum et al. 1989; Eary, Press et al. 1990). Eary et al. reported results with less than 6% error. Shulkin et al. used a modified conjugate imaging

technique to estimate the I-131 uptake in malignant pheochromocytomas; he reported less than 10% error (Shulkin, Sisson et al. 1988).

Sjogreen et al. used registered conjugate emission and transmission images to improve the quantitative measurements of organ uptake (Sjogreen, Ljungberg et al. 2001). Kojima et al. measured renal activity using conjugate imaging to evaluate the function of kidneys (Kojima, Ohyama et al. 2000). Yen et al. applied conjugate imaging to the determination of left ventricular volume, an important prognostic factor for patients with heart disease (Yen, Lim et al. 1994).

2.2.5 Extension of Classic Conjugate Imaging

Classic conjugate imaging, as described in this work, is the use of conjugate emission views and transmission measurements in the estimation of activity in a single source or the sum of activities in multiple sources embedded in a uniformly attenuating medium. The uniformly attenuating medium is assumed to contain no activity. The clinical application of standard conjugate imaging is useful for estimation of activity uptake in a single organ or tumor, or in multiple organs or tumors in dosimetric studies. These studies are usually performed under the assumptions that the body is a uniformly attenuating medium and uptake of activity in the body excluding the tumors or targeted organs is negligible. The classic conjugate imaging formulation is not sufficient for clinical cases in which organ or tumor uptake is surrounded by non-uniformly attenuating tissue or tissue with activity uptake or for clinical cases in which the estimation of the individual activities within multiple sources is desirable. Appropriate extensions of

classic conjugate imaging are necessary, if conjugate imaging is to be successfully applied in more complex settings.

2.2.5.1 Non-Uniformly Attenuating Medium

Sorenson and Thomas examined the estimation of activity in two overlapping sources embedded in a non-uniformly attenuating medium (Sorenson 1974; Thomas, Maxon et al. 1976). Figure 2.14 illustrates this geometry. Thomas designated one source to be the primary lesion and the other source to be the secondary lesion. In figure 2.14, the primary lesion is A_2 and the secondary lesion is A_4 . The activity in each lesion was assumed to be uniformly distributed throughout the lesion thickness. Thomas derived a formula for the estimation of the activity in the primary lesion (Thomas, Maxon et al. 1976).

If one assumes the two positions of acquisitions are left and right of an object, the activity in the primary lesion A_2 is expressed as

$$A_2 = \left(\frac{E_L E_R}{e^{-\mu_e t}} \right)^{1/2} \frac{f_2}{C} g(\alpha), \quad (2.33)$$

where

$$f_2 = \frac{\left(\frac{\mu_2 t_2}{2} \right)}{\sinh\left(\frac{\mu_2 t_2}{2} \right)},$$

$$g(\alpha) \equiv \left[1 + \left(\alpha \frac{f_2}{f_4} \right)^2 + 2\alpha \frac{f_2}{f_4} \cosh\left(\frac{\mu_2 t_2 + 2\mu_3 t_3 + \mu_4 t_4}{2} \right) \right]^{-1/2}, \quad (2.34)$$

and
$$\alpha = \frac{A_4}{A_2}, \quad (2.35)$$

where E_L and E_R are the left and right emission data; $e^{-\mu_e t}$ is the transmission factor through the object with thickness t ; f_2 is the self-attenuation correction factor for the primary lesion 1; and C is the system calibration factor, e.g. the camera sensitivities.

The term μ_e in equation 2.33 is the effective attenuation coefficient and is expressed as

$$\mu_e = \frac{1}{t} \sum_{i=1}^n \mu_i t_i, \quad (2.36)$$

where the index i represents the different regions of thickness t_i with attenuation coefficient, μ_i . The term $g(\alpha)$ is a correction factor that is dependent on the source activities ratio α and the thicknesses and attenuation coefficients of regions 2, 3, and 4.

Thomas compared his theoretical formalism with experimental data. The correction provided by the correction factor $g(\alpha)$ increases significantly with increasing source activities ratio. For a source activities ratio of 0.01, $g(\alpha)$ provides a correction of less than 5%. However, for a source activities ratio of 0.3, $g(\alpha)$ corrects the estimation activity by more than 20%. Thomas also reported the influence of source thickness on the source self-attenuation correction term f . For 140-keV photons the self-attenuation correction is less than 5% for source distributions less than 6-cm thick with attenuation coefficients within the experimental range reported. However, for larger source distributions the correction factor could exceed 10% (Thomas, Maxon et al. 1976).




Figure 2.14. Two activity distributions embedded in a non-uniformly attenuating medium. Activities A_2 and A_4 are uniformly distributed throughout regions 2 and 4. Activity A_2 is the primary lesion and activity A_4 is the secondary lesion.

2.2.5.2 Background Activity

Thomas et al. also examined estimation of activity of a single lesion surrounded by background activity in an attenuating medium (Thomas, Maxon et al. 1976). This is illustrated in figure 2.15. The primary lesion is region 3 and background activity is represented by regions 2 and 4.




Figure 2.15. A primary lesion surrounded by regions of background activity in a non-uniformly attenuating medium. A_3 is the primary lesion. A_2 and A_4 are the background activities surrounding the primary lesion.

Thomas expressed the activity of the primary lesion as

$$A_3 = \left(\frac{E_L E_R}{e^{-\mu_e t}} \right)^{1/2} \frac{f_3}{C} k(\gamma), \quad (2.37)$$

with

$$f_3 = \frac{\left(\frac{\mu_3 t_3}{2} \right)}{\sinh\left(\frac{\mu_3 t_3}{2} \right)}, \quad (2.38)$$

$$k(\gamma) = \left[1 + \left(\gamma_2 \frac{f_3}{f_2} \right)^2 + \left(\gamma_4 \frac{f_3}{f_4} \right)^2 + 2\gamma_2 \gamma_4 \frac{f_3^2}{f_2 f_4} \cosh\left(\frac{\mu_2 t_2 + 2\mu_3 t_3 + \mu_4 t_4}{2} \right) + 2\gamma_2 \frac{f_3}{f_4} \cosh\left(\frac{\mu_2 t_2 + \mu_3 t_3}{2} \right) + 2\gamma_4 \frac{f_3}{f_2} \cosh\left(\frac{\mu_3 t_3 + \mu_4 t_4}{2} \right) \right]^{-1/2} \quad (2.39)$$

and

$$\gamma_2 \equiv \frac{A_2}{A_3}, \quad (2.40)$$

$$\gamma_4 \equiv \frac{A_4}{A_3}. \quad (2.41)$$

In the equation 2.37, A_3 is the activity of the primary lesion; A_2 and A_4 are the activities in the volumes adjoining A_3 ; E_A and E_P are the anterior and posterior emission data, respectively; $e^{-\mu_e t}$ is the transmission factor through the object with thickness t ; f_3 is the self attenuation correction factor for primary source with thickness t_3 and attenuation coefficient μ_3 ; C is the system calibration factor, i.e. the camera sensitivity; and $k(\gamma)$ is the correction factor. The effective attenuation coefficient of the object is expressed as

$$\mu_e = \frac{1}{t} \sum_{i=1}^n \mu_i t_i. \quad (2.42)$$

In addition to the region thicknesses and attenuation coefficients, the source activity ratios of γ_2 and γ_4 must be known to yield the source activity.

Thomas compared his theoretical calculations with experimental data. The correction factor $k(\gamma)$ significantly increases with increasing source activities ratio. For a ratio of 0.15, $k(\gamma)$ is 20% or greater (Thomas, Maxon et al. 1976).

2.2.5.3 Multiple Sources

In 1981, Tsui formulated a general conjugate imaging equation for m regions of activity embedded in a medium with n attenuating regions where $m \leq n$. Each region has a uniform attenuation coefficient, μ_i , and a thickness, t_i , where $i = 1, 2, \dots, n$ (Tsui, Chen et al. 1981). The activity in the m regions is expressed as

$$A_{i(k)} = \rho_{i(k)} t_{i(k)}, \quad (2.43)$$

where $A_{i(k)}$ is the activity in source layer $i(k)$; $\rho_{i(k)}$ is the activity concentration in source layer $i(k)$; and $t_{i(k)}$ is the thickness of source layer $i(k)$; where $k = 1, 2, \dots, m$.

The emission data measured from the left and right cameras, E_L and E_R are expressed as

$$E_L = C_L \sum_{k=1}^m A_{i(k)} e^{\left(- \left[\frac{\mu_{i(k)} t_{i(k)}}{2} + \sum_{l=1}^{i(k)-1} \mu_l t_l \right] \right)} F_{i(k)} \quad (2.44)$$

$$E_R = C_R \sum_{k=1}^m A_{i(k)} e^{\left(- \left[\frac{\mu_{i(k)} t_{i(k)}}{2} + \sum_{l=i(k)+1}^{i(k)-1} \mu_l t_l \right] \right)} F_{i(k)}, \quad (2.45)$$

where C_L and C_R are the camera sensitivities of the left and right cameras, respectively, and $F_{i(k)}$ is the self-attenuation factor for source layer $i(k)$.

The geometric mean of the emission counts for m layers of activity embedded in n layers of attenuating medium is expressed as

$$\sqrt{E_L E_R} = \sqrt{C_L C_R} \left(- \sum_{l=1}^n \mu_l t_l / 2 \right) \left(\sum_{k=1}^m (A_{i(k)} F_{i(k)})^2 + 2 \sum_{\substack{k, k' \\ k \neq k'}}^m A_{i(k)} A_{i(k')} F_{i(k)} F_{i(k')} \times \cosh \left[\frac{1}{2} (\mu_{i(k)} t_{i(k)} + \mu_{i(k')} t_{i(k')}) + \sum_{\substack{k'' \\ k < k'' < k'}} \mu_{i(k'')} t_{i(k'')} \right] \right)^{1/2}, \quad (2.46)$$

where the indices k and k' denote different source layers and index k'' denotes a source layer between source layers k and k' .

For a single layer of radioactivity embedded in a medium with n attenuating regions, equation 2.47 simplifies to

$$\sqrt{E_1 E_2} = \sqrt{C_1 C_2} e^{\left(-\sum_{i=1}^n \mu_i t_i / 2\right)} \frac{A}{f}, \quad (2.47)$$

where

$$f = \frac{\sinh\left(\frac{\mu t}{2}\right)}{\frac{\mu t}{2}}. \quad (2.48)$$

Equation 2.47 is classic conjugate imaging of a thick source. The activity in the single region can be calculated if the detector/camera sensitivities, the attenuation coefficients, and the source thickness are known.

In the case of two or more aligning layers of activities embedded in a medium, the activity in each source can only be determined if the ratio of activity per unit area in each source layer and a reference source is known. The activity in the source layers can be determined through a system of simultaneous equations in terms of source activity per unit area, A_i .

The mathematical extensions provided by Thomas and Tsui require explicit knowledge of region thickness, attenuation coefficients, and activity ratios to determine individual activity distributions in an assembly of overlapping distributions. Clinically, knowledge of the activity ratios is usually not practical as such can generally be known only through use of invasive techniques.

2.3 Striatal Dopamine and Parkinsonism and Schizophrenia

Early single-photon emission imaging studies of the striatal dopaminergic system indicated that assessment of radiopharmaceutical binding to dopamine transporters and D2 dopamine receptors in the striata provide a reasonable measure of dopamine production and transmission. As a result, single-photon emission imaging performed during the last two decades has been an integral part of in vivo studies investigating abnormal dopamine production and transmission in the striatal dopaminergic system. Parkinsonism, a collection of motor disorders, and schizophrenia, a neuropsychiatric disorder, are associated with abnormal dopamine production in the striatal dopaminergic system.

Currently, diagnosis of Parkinsonism is primarily dependent on the expression of clinical symptoms and family history. Unfortunately, the various Parkinsonian disorders have very similar symptoms and the symptom presentations are usually in the later stages of the diseases. Parkinsonism imaging studies, especially Parkinson's disease imaging studies, have demonstrated success in aiding in diagnosing the disease, studying the progression of the disease (i.e. longitudinal studies), and assessing drug therapy. These studies have been qualitative and quantitative. Quantitative studies are more accurate, especially in longitudinal studies, as disease progression can be relatively slow. Consequently, clinical quantitative nuclear medicine techniques are being developed to provide a more accurate assessment of dopamine production as it pertains to Parkinsonism.

Single-photon emission studies of schizophrenia are performed to investigate the role of the striata in the pathology of the disorder as new hypotheses suggest that

excessive dopamine stimulation in the striata is not the primary reason for schizophrenia. Drug challenge imaging studies are also performed to evaluate the response of schizophrenic patients to neuroleptics. Many schizophrenia patients, who are clinically characterized to respond well to neuroleptics, do not do so. New neuroleptics are being produced to decrease side effects and increase symptom alleviation. Imaging studies are used to test and characterize these new drugs and their effects.

The motivation of the development of the modified conjugate imaging technique is the quantitative assessment of radiopharmaceutical uptake in the striata with direct clinical application to pathology studies and care management research for Parkinsonism and schizophrenia. In the following section, the anatomy and function of the striata are presented as well as Parkinsonism and schizophrenia. The role of striatal dopamine in the manifestation of these disorders is also discussed briefly.

2.3.1 Striata

The striata are parts of the basal ganglia, a collection of nuclei in the cerebral cortex. The basal ganglia are anatomically symmetric. One set of nuclei in the basal ganglia is located in the left hemisphere and a second is located in the right hemisphere. Collectively, the left striatum and right striatum are referred to as striata. In transaxial and coronal views of the brain, the striata appear as two aligned structures. The striata are the input portal of the basal ganglia. A striatum consists of a putamen and a caudate nucleus. However, for the scope of this dissertation, neither the putamen nor the caudate nucleus is referred to individually; they are only referred to collectively as a striatum. In addition to the striata, the basal ganglia consist of the globus pallidus, substantia nigra, and

subthalamic nucleus. The substantia nigra are each divided into two sections – substantia nigra pars compacta and substantia nigra pars reticulata. Figure 2.16 is a detailed illustration of the physiological structure of the basal ganglia (The HOPES Brain Tutorial 2003). The major tasks of the basal ganglia are to inhibit unwanted physical movement and to facilitate wanted physical movement. The basal ganglia receive input from the cerebral cortex and send information back to the cortex through the thalamus. There are two major pathways through the basal ganglia; one is direct, the other indirect. The direct and indirect pathways work together to legislate movement. The direct pathway facilitates movement and the indirect pathway inhibits movement. The amount of movement initiated is dependent on the activity in the direct and indirect pathways (Albin, Young et al. 1989; Mink 1996; Eidelberg 2003; Uhl 2003). Smooth movement is executed when the output from both pathways is in proper balance.

Peripheral, but crucial, to these pathways is the nigrostriatal pathway. The nigrostriatal pathway aids in maintaining the balance between the direct and indirect pathways through the release of the neurotransmitter dopamine. Neurons project from the substantia nigra pars compacta to the neurons in the striata. The neurons originating in the substantia nigra pars compacta synthesize dopamine and store it in vesicles. In response to stimulation, dopamine is released into the synaptic cleft (Booij, Tissingh et al. 1999). The dopamine crosses the synaptic cleft to bind with striatal dopamine receptors to stimulate the striatal neuron. The stimulated striatal neurons in turn send input to the direct and indirect pathways to modulate output from the basal ganglia (Albin, Young et al.; Mink). Figure 2.17 is a sketch of striatal neurons.

The release and production of dopamine are facilitated by dopamine transporters




Figure 2.16. Illustration of the basal ganglia. The basal ganglia consist of the caudate nucleus, putamen, globus pallidus, subthalamic nucleus and substantia nigra. Anatomically, the basal ganglia only consist of the caudate nucleus, putamen, and the globus pallidus. The subthalamic nucleus and the substantia nigra are included because of their role in the function of the basal ganglia. The caudate nucleus and the putamen are collectively called the neostriatum; it is often called the striatum. Throughout the remainder of this dissertation striatum will be used to refer to the caudate nucleus and the putamen. The substantia nigra is also divided into two components – the substantia nigra pars compacta (SNpc) and the substantia nigra reticulata (SNpr). (The HOPES Brain Tutorial 2003)




Figure 2.17. Striatal neurons. The pre-synaptic terminal is the end of a dopaminergic neuron that projects from the substantia nigra pars compacta. The dopaminergic neuron produces dopamine and releases dopamine into the synaptic cleft in response to an action potential. The dopamine released into the synaptic cleft migrates and eventually binds to the D₂ receptors located on the post-synaptic terminal of a striatal neuron. When enough dopamine has bound to the receptors, an action signal stimulates the striatal neuron. Dopamine left in the synaptic cleft is transported back to the pre-synaptic terminal of the dopaminergic neuron by the dopamine transporters for reuptake.

(DAT) and dopamine receptors. Dopamine transporters are located on the ends of the neurons projecting from the substantia nigra pars compacta. This location is called the pre-synaptic terminal. Dopamine transporters aid in the release of dopamine into the synaptic cleft and are responsible for transporting the unused dopamine in the synaptic cleft back to the pre-synaptic terminal (Bourguignon, Pauwel et al. 1997; Uhl 2003). Dopamine receptors are located at the post-synaptic terminal of the striatal neurons. These receptors are the input portals of the striatal neurons. Dopamine binding to the receptors serves as a stimulus to the striatal neurons. Dopamine receptors have been categorized into at least five subtypes, which are classified into two classes, D₁ and D₂. Receptors in the D₁ class are receptors D₁ and D₅. D₂ class receptors are D₂, D₃, and D₄. D₁ and D₂ receptors are the most abundant with their highest concentrations being in the striata (Bourguignon, Pauwel et al. 1997). D₁ dopamine is released to stimulate the striatal neurons to send excitatory input to the direct pathway. D₂ dopamine is released to stimulate striatal neurons to send inhibitory input to the indirect pathway (Cote and Crutcher 1991).

2.3.2 Parkinsonism

Parkinsonism is a group of motor disorders associated with an abnormal loss of dopaminergic neurons in the nigrostriatal pathway. The group of motor disorders under the Parkinsonism moniker is often divided into four categories: primary Parkinsonism, secondary Parkinsonism, Parkinsonism-Plus syndromes, and hereditary degenerative disorders. Primary Parkinsonism is the most common form of Parkinsonism and is commonly referred to as Parkinson's disease (PD). Until recently, the etiology of

Parkinson's disease was unknown. But research in recent years suggests gene mutations might be the cause of Parkinson's disease. Secondary Parkinsonism includes disorders with known causes. Examples of such disorders are vascular Parkinson's disease and drug-induced Parkinson's disease. Parkinsonism-Plus syndromes are disorders associated with other degenerative diseases, e.g., multiple system atrophy (MSA) and progressive supranuclear palsy (PSP). Hereditary degenerative disorders are hereditary, e.g., Huntington's disease.

It is estimated that approximately 500,000 to 1.5 million people suffer from some form of Parkinsonism. Parkinsonism generally affects elderly people, but has been known to affect young and middle-aged adults. In an incidence study of Parkinsonism performed in Olmstead County, Minnesota, the average annual incidence rate of Parkinsonism is reported to be 114.4 per 100,000 people (Bower, Maraganore et al. 1999). About 60-85% of all cases of Parkinsonism are Parkinson's disease. Secondary Parkinsonism accounts for about 8-10% of Parkinsonism; PSP and MSA are the next most common Parkinsonian syndromes (Albanese 2003; Mayo Clinic Staff). Pathological studies suggest that PSP represents 4.0-7.5% of all Parkinsonian disorders and MSA represents 4-22% (Stacy and Jankovic 1992; Testa, Monza et al. 2001). Based on the same epidemiology study in Olmstead, Minnesota, the average annual incidence of PSP and MSA for ages 50-99 years is 5.3 and 3.0 per 100,000 people, respectively (Bower, Maraganore et al. 1999).

Parkinsonism is the manifestation of the degeneration of dopaminergic neurons in the nigrostriatal pathway. Primary (Parkinson's disease) and secondary Parkinsonism are characterized by the degeneration of the dopamine-producing neurons projecting from the

substantia nigra pars compacta. This results in the loss of dopamine and dopamine transporters; the transporters are located on the pre-synaptic terminal of these neurons. Parkinsonism-Plus syndromes are associated with the loss of striatal neurons, dopamine-receiving neurons, as well as the dopamine-producing neurons in the substantia nigra pars compacta (Tatsch 2001; Seibyl 2003).

Currently, clinical diagnosis of the specific forms of Parkinsonism is primarily predicated on the presentation of a combination of motor-deficient symptoms. These motor features are tremors at rest, bradykinesia (difficulty initiating movement), rigidity, postural instability, and a shuffling gait (freezing phenomenon) (Fahn 2003). Two or more of these symptoms must be present for a diagnosis of Parkinsonism and tremors at rest or bradykinesia must be one of the symptoms. Although the etiologies of the various forms of Parkinsonism are different, the presentation of symptoms is similar, especially in the early symptomatic stages. In addition, psuedoparkinsonism, motor disorders not associated with the degeneration of the nigrostriatal dopaminergic system, present with Parkinsonism-like symptoms. This overlapping of symptom presentation results in a moderate number of misdiagnosed cases.

Pathological studies estimate that as many as 35 % of all PD cases clinically diagnosed have been erroneously characterized (Rajput, Rozdilsky et al. 1991; Seibyl 2003). Early stages of PSP and MSA are commonly mistaken for PD. It has been estimated that about 8% to 22% of clinically diagnosed PD cases are MSA cases (Schrag, Ben-Shlomo et al. 2002). Parkinson's disease and essential tremor are frequently mistaken for one another. Some studies suggest that up to 25% of idiopathic Parkinson's disease cases were actually essential tremor cases (Quinn 1994; Booi, Speelman et al.

2001). Conversely, patients in the early stages of idiopathic PD are often misdiagnosed with essential tremor, vascular PD, drug-induced PD, or Alzheimer's disease (Seibyl 2003). Misdiagnosis is an issue that has warranted concern due to the prognosis and treatment implications.

Drug therapy has been used to mitigate the error of misdiagnosis between the various forms of Parkinsonism and other motor disorders. The response to the therapy is dependent on the motor disorder. Patients with early stage Parkinson's disease respond well to the treatment. Patients with other forms of Parkinsonism or non-Parkinsonian motor disorders do not respond to the therapy. The most common drug therapy is levodopa, a dopamine-producing drug.

Due to the high incidence rate of Parkinson's disease most treatment options for Parkinsonism are directed at the management and care of Parkinson's disease. Patients with Parkinsonism-Plus syndromes have limited treatment options. The most common and widely clinically accepted Parkinson's disease treatment is dopamine replacement using levodopa. Levodopa therapy decreases the symptoms of Parkinson's disease. However, levodopa treatment is most beneficial in the early stages of Parkinson's disease and provides on average five years of symptom relief. About 60% of patients using levodopa therapy experience sudden episodes of abnormal involuntary movements (Fahn 2003). Dopamine agonists, drugs that readily bind to the dopamine receptors to simulate the neuron, have been used to manage Parkinson's disease. Unfortunately, some of the side effects of dopamine agonists are hallucinations, confusion, and psychosis. The severity of these side effects increases with age. Neuroprotective treatments, techniques that protect the remaining dopaminergic neurons from degeneration, have been

researched but no technique developed thus far has been able to markedly stop progression of the disease.

2.3.3 Schizophrenia

Schizophrenia is a chronic neuropsychiatric disorder widely associated with hyper-production of dopamine. Clinically, schizophrenia is divided into five main types based on symptom expression. These schizophrenic types are hebephrenic type, catatonic type, paranoid type, undifferentiated type, and residual type. Hebephrenic schizophrenia, also referred to as disorganized schizophrenia, is characterized by hallucinations, delusions, incoherence, inappropriate laughter, and social withdrawal. Catatonic schizophrenia is characterized by motor disturbances that include alternating between immobility and wild excitement, decreased sensitivity to painful stimulus, and a decline in personal hygiene (Mental Health Channel 2004). Paranoid schizophrenia is characterized by delusions of persecution, anxiety, and violence. Undifferentiated schizophrenia is a schizophrenic type that is characterized by the expression of schizophrenic symptoms from more than one type of schizophrenia. Residual schizophrenia is characterized by the diminished expression of prominent schizophrenic symptoms.

Schizophrenia affects approximately 1% of the United States population. Typically, the onset of schizophrenia occurs in the late adolescent stage or early twenties. However, schizophrenia has been known to affect people 50 years old and below; occurrences in people above age 50 are rare. Full recovery is extremely rare.

Schizophrenic symptoms worsen after the onset of schizophrenia in approximately 75% of the people afflicted with schizophrenia (Mental Health Channel 2004).

Clinically, schizophrenic symptoms are classified into two main groups - positive and negative. Positive symptoms are characterized as symptoms that generate distorted perceptions. Positive symptoms can be divided into the psychotic group and the disorganized group. Delusions and hallucinations make up the psychotic group. Disorganized thought and speech, inappropriate behavior and erratic behavior are examples of symptoms in the disorganized group (Andreasen 1995; Andreasen 1996). Negative symptoms are characterized by a reduction of normal function and behavioral deficits (Mental Health Channel 2004). Poor speech, lack of emotion, lack of energy, and a loss of interest are examples of negative symptoms. Schizophrenia is not defined by one symptom; schizophrenic symptoms are diverse and extensive. The World Health Organization and American Psychiatric Association compiled a symptomatic criterion-based system for diagnosing schizophrenia (Andreasen 1995). If a patient exhibits two symptoms negative or positive, the patient can be diagnosed with schizophrenia. In the event that a patient exhibits a typical schizophrenic symptom, schizophrenia can be diagnosed. In either case, the symptom or symptoms must not be from another existing medical condition.

Although the causes of schizophrenia are unknown, several hypotheses have been proposed. The predominant hypothesis is the dopamine hypothesis (Carlsson and Lindqvist 1963; Creese, Burt et al. 1976; Sedvall and Farde 1995; Mann 1996; Frederickson 1998). The dopamine hypothesis suggests that schizophrenia is the manifestation of hyperactivity of dopaminergic function in the brain, specifically the

striata. Initial support of the dopamine hypothesis was indirect and based on the pharmacological studies of neuroleptic drugs. Neuroleptic drugs used to treat schizophrenia alleviate symptoms by blocking D2 dopamine receptors (Davis, Kahn et al. 1991; Schwarz, Oertel et al. 1996; Harrison 2000; Depatie and Lal 2001; Strange 2001). Thus, it was deduced that people suffering from schizophrenia have excessive dopamine stimulation of dopamine D2 receptors. Post-mortem studies performed to investigate the hyperactivity of dopamine produced inconsistent results (Laruelle 2000). With the use of single-photon emission imaging, support the dopamine hypothesis now includes data that more directly reflects dopamine production and transmission in the brain of a schizophrenic patient. For the last couple decades, reports of imaging studies included data that support of the dopamine hypothesis. In 2000, Laruelle et al. reported increased dopamine transmission in schizophrenic patients in amphetamine-challenge imaging studies (Laruelle 2000). In 2000, Abi-Dargham et al. reported data that suggest increased D2 receptor stimulation in schizophrenia patients (Abi-Dargham, Rodenhiser et al. 2000). However, critics of the dopamine theory argue that the increase of dopamine activity is actually a response to long-term drug therapy (Toru, Kurumaji et al. 1994). Lomeña et al. performed an imaging study between neuroleptic-naïve patients (those that have never taken neuroleptics) and neuroleptic-free schizophrenia patients (those that have previously taken neuroleptics) and found no difference in D2 receptor density. This suggests that drug therapy does not increase dopamine activity (Lomena, Catafau et al. 2004).

Although the dopamine hypothesis is the most popular theory to date, dopamine hyperactivity in the striata is mainly associated with the psychotic symptoms of

schizophrenia. Laruelle et al. reported findings that demonstrated a possible relationship between dopamine hyperactivity in the striata and the negative symptoms of schizophrenia. However, it is widely accepted that negative symptoms are not influenced by the hyperactivity of striatal dopamine. Thus, alone, the dopamine hypothesis is not a sufficient hypothesis; it does not address the schizophrenic patient sub-group that has psychotic symptoms and that does not respond to neuroleptic drug therapy. It also does not offer a pathological explanation for the manifestation of negative symptoms. A hypothesis has been proposed by several investigators to explain the pathology of both negative symptoms and positive symptoms (Weinberger, Berman et al. 1988; Davis, Kahn et al. 1991; O'Donnell and Grace 1998). Their hypothesis suggests the underlying cause of negative symptoms is dopamine hypo-activity in the prefrontal cortex. Positive symptoms are associated with dopamine hyperactivity in the striata (Laruelle 2000).

Treatment of schizophrenia is continuous and requires steady administration of medication. Currently, antipsychotic drugs are used commonly in the management of schizophrenia. These drugs are divided into the traditional and the atypical category. Traditional antipsychotic drugs block dopamine D2 receptors. A common traditional antipsychotic drug is haloperidol, commercially known as Haldol. Traditional antipsychotic drugs are most effective in suppressing positive or psychotic symptoms. However, traditional antipsychotic drugs produce tardive dyskinesia, a parkinsonism-like disorder, due to the blocking of D2 receptors. Not all schizophrenia patients respond to traditional neuroleptics. Atypical neuroleptics help about 50% of schizophrenic patients that do not respond well to traditional neuroleptics. Patients that express positive symptoms are more likely to respond to neuroleptics than patients that express negative

symptoms. Atypical antipsychotic drugs block primarily serotonin, nonrepinephrine, histamine, and other dopamine receptors. These drugs suppress negative symptoms as well as positive symptoms. A common atypical neuroleptic is clozapine. A common side effect of atypical neuroleptics is agranulocytosis, the significant loss of white blood cells. New atypical drugs that work similar to atypical drugs but have milder side effects are available. These drugs work by blocking D2 receptors and serotonin receptors. Risperidone and Olanzapine are new atypical antipsychotic drugs.

2.4 Imaging the Striata

Critical to quantitative imaging studies is the selection and use of the appropriate radiopharmaceutical, especially for small tissue volumes, such as the striata. Selection of the right radiopharmaceutical is important because the properties of a radiopharmaceutical have a significant influence on the accuracy of measurements obtained using it. Imaging of the striata requires a radiopharmaceutical that participates in the physiological function of the nigrostriatal dopaminergic system without disruption. In addition, radiopharmaceuticals used to image the striatal dopaminergic system have a relatively high affinity for dopamine receptors or dopamine transporters, and have a low affinity for surrounding brain tissue (Halldin, Gulyas et al. 2001; Catafau 2003). These radiopharmaceuticals readily bind to the target area and unbind at a slow rate (low dissociation rate). This allows for adequate imaging time after equilibrium (steady-state condition) is reached to complete acceptable data acquisitions.

To date, most radiopharmaceuticals used for single-photon emission imaging studies of dopaminergic neurotransmission have been labeled with the radioisotope ^{123}I .

TRODAT, which is labeled with ^{99m}Tc , is the exception. Table 2.1 lists some radiopharmaceuticals used in dopaminergic investigations. The radiopharmaceuticals in bold print will be discussed in detail later in this chapter. Although ^{123}I is rather expensive, it readily binds to several pharmaceutical agents that demonstrate high affinity to the dopaminergic system in the striata. It is because of this chemical reactivity that ^{123}I is the principal radioisotope in striatal dopaminergic studies. ^{123}I has the moderately short half-life of 13.2 hours. It has a primary photon emission of 159 keV at an abundance of 83%. However, 3% of ^{123}I 's decays results in emissions of high-energy gammas with

Table 2.1. Radiopharmaceuticals used in single-photon emission imaging studies of dopaminergic neurotransmission in humans (Catafau 2003).

Target	Radioligand
D ₂ receptor (post-synaptic)	^{123}I - IBZM
	^{123}I -Epidepride
	^{123}I -BF
	^{123}I -Iodospiperone
	^{123}I -Iodolisuride
Dopamine transporter (pre-synaptic)	^{123}I - β-CIT
	^{123}I - FP-CIT
	^{123}I -IPT
	^{123}I -PE21
	^{99m}Tc - TRODAT-1

energies ranging from 248 to 783 keV. Figure 2.18 is a list of the energies and yields of the primary gamma emission and high-energy gamma emissions of ^{123}I . Due to the high-

energy emissions of ^{123}I , special considerations need to be made, if radiopharmaceutical uptake in the striata is to be accurately assessed. High-energy emissions further complicate the problem of obtaining quantitative measurements as they have the potential of contributing a significant number of counts to image data. These emissions can travel

Energy (keV)	% Yield
<i>Primary emission:</i>	
159.0	82.8
<i>High energy emissions:</i>	
248.0	0.0707
281.0	0.0786
346.3	0.1250
440.0	0.4250
505.3	0.3140
529.0	1.3800
538.5	0.3790
624.6	0.0828
687.9	0.0266
735.8	0.0612
783.6	0.0590
Total	3.0019

Figure 2.18. Percentage yield of ^{123}I decays. The distribution of high-energy source emissions from the lungs and GI tract region was simulated according to the percentage yield of ^{123}I decays as shown in the table.

through and/or scatter in the gamma camera's collimator and interact in the scintillation crystal depositing energy in the photopeak. Such occurrences are referred to as downscatter. Downscatter degrades image quality, and more importantly, produces an

error in the estimation of the activity in the striata. The influence of downscatter from shielding and collimator penetration can be significant, when imaging patients suffering from Parkinsonism, in which the activity uptakes in the striata are small. The extent of contribution is dependent on several factors – shielding, scatter correction techniques employed, camera parameters, and radiopharmaceutical biodistribution.

Approximately 85% of all imaging procedures in the United States are performed using a radiopharmaceutical labeled with ^{99m}Tc (Karesh 1996). As a result, nuclear medicine devices are designed to image ^{99m}Tc radiopharmaceuticals. Although ^{99m}Tc is the most common radioisotope used in nuclear medicine, it does not bind readily to presently used pharmaceutical agents exhibiting high affinity to the dopaminergic system. Consequently, most radiopharmaceuticals used in this context are not labeled with ^{99m}Tc . ^{99m}Tc TRODAT is the exception. Pharmacological studies are currently being performed to investigate the chemical reactivity of ^{99m}Tc with other dopaminergic pharmaceuticals, because integration of a ^{99m}Tc radiopharmaceutical into a clinical setting is easier than the integration of an ^{123}I labeled compound (Karesh 1996; Kung, Kung et al. 2003). ^{99m}Tc is readily available and relatively inexpensive. It has a half-life of 6 hours. Its principal photon emission is a 140 keV gamma with an abundance of 89%.

2.4.1 Non-specific Binding

In addition to striatal dopamine, non-neuronal dopamine can be found throughout the body as dopamine assists in regulation of gastrointestinal function, renal function, and respiratory function (Mezey, Eisenhofer et al. 1996; Pestana, Faria et al. 1997; Cabezas, Lezama et al. 2001; Cabezas, Israili et al. 2003; Tonini, Cipollina et al. 2004).

Consequently radiopharmaceutical uptake in the lungs, kidneys, gallbladder, and sections of the gastrointestinal tract is common. Radiopharmaceutical uptake in these organs is referred to as out-of- field of view (FOV) activity throughout the remainder of the dissertation. The degree of uptake depends on the specific radiopharmaceutical used. Out-of-FOV activity uptake of ^{123}I striatal radiopharmaceuticals poses as problem in quantitative studies because the higher energy photon emissions of ^{123}I are sufficiently energetic to escape the body, scatter, and deposit energy into the photopeak window. The inclusion of these counts is especially problematic, when small FOV cameras are being used, such as a dedicated imaging system. Since the out-of-FOV activity is typically in the lungs and GI tract, the higher energy photons will most likely penetrate the gamma camera from the bottom. For clinical gamma cameras in which the FOV is larger, counts due to interactions of higher energy photons will be toward the bottom of the FOV away from the striatal and background ROIs used in modified conjugate imaging. But due to the smaller dimensions of small FOV gamma cameras, counts from higher energy photons are more likely to be included in the striatal ROI and background ROI.

2.4.2 Pre- and Post-Synaptic Radiopharmaceuticals

The radiopharmaceuticals used in the imaging of the dopaminergic system in the striata are categorized as pre-synaptic or post-synaptic. Pre-synaptic radiopharmaceuticals, also known as dopamine transporter (DAT) ligands, bind to the dopamine receptors located at the pre-synaptic terminal of the dopaminergic neuron. (The pre-synaptic terminal is the end of a neuron that projects from the substantia nigra pars compacta to the striatum. See section 2.3.) These radiopharmaceuticals are commonly

used to confirm the presence of Parkinsonism when clinical symptoms are ambiguous. Often patients with essential tremor or drug-induced Parkinson's are diagnosed with some form of Parkinsonism because the symptoms are similar. Patients suffering from Parkinsonism, however, have fewer pre-synaptic dopaminergic nerve terminals due to degeneration of the dopaminergic neurons. Patients with essential tremor have normal amounts of pre-synaptic nerve terminals (Tatsch 2001). Thus, uptake of a dopaminergic radiopharmaceutical in a patient with Parkinsonism is considerably less than such uptake in a patient with essential tremor or drug-induced Parkinson's. DAT ligands are also used for early and pre-clinical diagnosis of Parkinsonism.

Post-synaptic radiopharmaceuticals bind to the D2 receptors located on the post-synaptic terminal of a striatal neuron. These radiopharmaceuticals are used to distinguish between Parkinson's disease and other Parkinsonism forms. Parkinson's disease affects the pre-synaptic dopaminergic nerve terminals but does not affect the post-synaptic D2 receptors. This is in contrast to what happens in other forms of Parkinsonism in which both the dopamine transporters and D2 receptors are affected (Tatsch 2001). Post-synaptic radiopharmaceuticals are also important in schizophrenia, as studies have found that some schizophrenic patients have excessive dopamine and D2 receptors.

There are several pre- and post-synaptic ligands available for dopaminergic imaging in humans. See table 2.1. The biodistributions of four radiopharmaceuticals are used in the characterization and simulation of the modified conjugate imaging technique using a novel dedicated conjugate imaging system. The ^{123}I radiopharmaceuticals, ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM, are chosen based on commercial availability and type of radiopharmaceutical – pre- or post-synaptic. Since ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I

IBZM are commercially available, European guidelines for nuclear medicine procedures exist (Tatsch, Asenbaum et al. 2002) . These guidelines are used in this study. ^{99m}Tc TRODAT is chosen due to its clinical promise and possible ease of integration into the clinic.

2.4.2.1 Biodistribution of ^{123}I β -CIT

Commercially known as DOPASCAN, ^{123}I β -CIT has a high specific binding in the striatum (Seibyl, Wallace et al. 1994; Bourguignon, Pauwel et al. 1997; Messa, Volonte et al. 1998; Seibyl, Marek et al. 1998; Catafau 2003; Kung, Kung et al. 2003). Seibyl reports the striatal-to-occipital ratios for the right and left striatum to be 9.3 ± 2.9 and 10.3 ± 3.5 , respectively (Seibyl, Wallace et al. 1994). The European Association of Nuclear Medicine Procedure Guidelines recommends that image data acquisition occur 18-24 hours post-injection. The recommended “standard” dose for an adult is 150-250 MBq – typically 185 MBq (5 mCi) (Tatsch, Asenbaum et al. 2002). The percent activity uptakes used in the simulation studies reported in this dissertation are taken from the biodistribution study performed by Seibyl et al (Seibyl, Wallace et al. 1994). The dose given in the study was 185 MBq within the recommended dose. The equilibrium time was taken to be approximately 21 hours, also within the recommended acquisition window. At the time of acquisition, the activity uptakes in the striata and brain are 2% and 10%, respectively of the initial activity administered. The percent uptake in the lungs, liver, and intestines is 21%, 25%, and 27%, respectively (Seibyl, Wallace et al. 1994). For the purposes of the simulation studies, the liver and small intestine uptakes are summed to get an overall activity uptake of 52%.

2.4.2.2 *Biodistribution of ^{123}I FP-CIT*

^{123}I FP-CIT, iodine-123-labelled N- ω -(fluoropropyl)-2 β -carbomethoxy-3 β -(4-iodophenyl) tropane is an N- ω -(fluoropropyl) analogue of β -CIT. It is commercially known as DATSCAN. ^{123}I FP-CIT has a relatively high specific binding in the striata (Booij, Sokole et al. 1998; Tissingh, Booij et al. 1998; Booij, Hemelaar et al. 1999; Benamer, Patterson et al. 2000; Booij, Speelman et al. 2001). Booij reports the ratios are 2.24 ± 0.32 ($P < 0.01$) and 0.89 ± 0.34 for healthy volunteers and Parkinson's patients, respectively, when the ratios are computed three hours post-injection (Booij, Hemelaar et al. 1999). Although the specific binding of ^{123}I FP-CIT is lower than that of ^{123}I β -CIT, ^{123}I FP-CIT has faster kinetics, making one-day imaging protocols possible, which is advantageous for clinical settings (Seibyl, Marek et al. 1998; Tissingh, Booij et al. 1998; Booij, Hemelaar et al. 1999; Tatsch and Frey 2002). However, the smaller specific binding ratio and fast kinetics of ^{123}I FP-CIT might decrease the accuracy in quantitative studies (Seibyl, Marek et al. 1998). The European Association of Nuclear Medicine Procedure Guidelines recommends that image data acquisition occurs 3-6 hours post-injection. The recommended standard dose for an adult is 150-250 MBq – typically 185 MBq (5 mCi) (Tatsch, Asenbaum et al. 2002). Figure 2.19 includes transaxial images of ^{123}I β -CIT and ^{123}I FP-CIT uptake in the striata of a healthy individual and a Parkinson's disease patient.

The percent activity uptakes used in the simulation studies reported in this dissertation were taken from the ^{123}I FP-CIT biodistribution study performed by Booij et al (Booij, Sokole et al. 1998). The simulated dose used is the recommended dose of 185 MBq. The equilibrium time is taken to be approximately 3 hours. At the time of

equilibrium the activity uptakes in the striata and brain are 0.8% and 2.2%, respectively, of the initial activity administered. The percent uptakes in the lungs, liver, and intestines are 9%, 18%, and 18%, respectively (Booij, Sokole et al. 1998). For the purposes of the simulation




Figure 2.19. ^{123}I β -CIT and ^{123}I FP-CIT striatal activity uptake. SPECT axial images in healthy control subject and Parkinson's disease patient. The same healthy subject and Parkinson's disease patient were used for both radiopharmaceutical studies. For ^{123}I β -CIT study, 6 mCi of activity are given to both subjects. For the ^{123}I FP-CIT study, 9 mCi of activity are given to both subjects. The striatal uptake for the Parkinson's disease patients is markedly less than that of the healthy individuals. Reprinted by permission of the Society of Nuclear Medicine from (Seibyl, Marek et al. 1998).

studies, the liver and intestine uptakes were summed to get an overall activity uptake of 36%.

2.4.2.3 Biodistribution of ^{99m}Tc TRODAT

Technetium-99m TRODAT, technetium-99m 2-[[2-[[[3-(4-chlorophenyl)-8-methyl-8-azabicyclo[3,2,1]oct-2-yl]-methyl](2-mercapto-ethyl)amino]ethyl]amino]ethanethiolato(3-)]-oxo-1R-(*exo-exo*)]-) is a relatively new pre-synaptic (DAT) ligand. Animal and human studies have shown that TRODAT readily binds to the dopamine transporters in the striata (Kung, Kim et al. 1996,2001; Meegalla, Plossl et al. 1997; Mozley, Stubbs et al. 1998; Fang, Wu et al. 2000; Tzen, Chin-Song et al. 2001). Huang performed a study on healthy controls and Stage I and II Parkinson's disease patients and reported specific uptake values in the striatum of 1.98 ± 0.24 , 1.62 ± 0.11 , and 1.22 ± 0.12 , for the healthy controls and Stage I and II Parkinson's disease patients, respectively. The transaxial slices in figure 2.20 represent specific uptakes in controls and Parkinson's disease patients (Huang, Lin et al. 2001). Huang chose to use the occipital region as the reference region. The occipital region has low uptake of TRODAT. Occipital data can be easily compared to striatum data from the same slice. Imaging took place 165 to 195 minutes post-injection (Huang, Lin et al. 2001). The administered dose was 740 MBq (20mCi) (Mozley, Schneider et al. 2000; Huang, Lin et al. 2001). The critical organs in healthy men and women controls were the liver, kidney, and the upper large intestine (Kung, Kim et al. 1996; Mozley, Stubbs et al. 1998).

Human biodistribution studies focus on reporting radiation dosimetry. Animal studies concentrate on absolute activity uptake in organs, as invasive procedures are required in obtaining uptake information. As a result, the activity uptake information used in this dissertation is from animal studies (Meegalla, Plossl et al. 1997; Fang, Wu et al. 2000). Fang et al. reported striata uptake as percent uptake per gram and the brain uptake

as percent uptake per organ. The striata percent uptake used in the reported simulation studies was determined by multiplying the striata uptake and mass of the simulated striata. Assuming equilibrium is reached at about 60 minutes post injection, the percent uptakes of initial activity used in the reported simulation studies in the striata and brain are 0.63% and 0.12%, respectively. The activity uptake in critical organs is not used in




Figure 2.20. ^{99m}Tc TRODAT-1 striatal activity uptake. SPECT axial images of Tc99m - TRODAT-1 uptake for healthy volunteer (A) and of patients with Hoehn and Yahr scale stage I (B) and stageII (C) Parkinson's disease. The Hoehn and Yahr scale (HYS) is a clinical classification system used to assess Parkinson's disease. The HYS divides Parkinson's disease into 6 stages. Stages I and II are classified as mild disease. Stage I is characterized by unilateral involvement. Stage II is characterized by bilateral disease involvement. A decrease in activity is apparent for the Parkinson's disease patients. Asymmetrical activity uptake is also apparent for the Parkinson's disease patients. Reprinted by permission of the Society of Nuclear Medicine from (Huang, Lin et al. 2001).

these studies because the photon emissions from these activity distributions are assumed to have negligible influence on the quantitative studies performed. The dose used in the studies was 740 MBq.

2.4.2.4 Biodistribution of ^{123}I IBZM

^{123}I IBZM, iodine-123 ((S)-(-)-2-hydroxy-3-iodo-6-methoxy-N [(1-ethyl-2-pyrrolidinyl)methyl] benzamide), is used to study dopamine systems, specifically those in Parkinsonism and schizophrenia patients (Tatsch and Frey 2002). The specific uptake ratio in the striata is relatively high (Kung, Alvai et al. 1990; Seibyl, Woods et al. 1992; Knable, Jones et al. 1995; Leslie, Abrams et al. 1996; Schwarz, Oertel et al. 1996; Leenders 2003; Lomena, Catafau et al. 2004). IBZM accumulates readily in the striatum, brain, large intestine, liver, spleen, and bladder. The critical organs are the wall of the large intestine, the spleen, and the wall of the urinary bladder (Kung, Alvai et al. 1990). The European Association of Nuclear Medicine reports the critical organs to be the thyroid and the bladder wall. The European Association of Nuclear Medicine Procedure Guidelines state the “standard” dose for an adult is 150-250 MBq – typically 185 MBq (5mCi) (Tatsch, Asenbaum et al. 2002). It is recommended that data acquisition begin 1.5-3 hours post-injection. Figure 2.20 is the image of ^{123}I IBZM activity uptake in the striata (Leslie, Abrams et al. 1996).




Figure 2.21 ^{123}I IBZM striatal activity uptake. Planar image of ^{123}I IBZM uptake in the striata of a healthy subject two-hour post-injection. Striatum-to-frontal cortex uptake ratio was 1.73. Reprinted by permission of the Society of Nuclear Medicine from (Leslie, Abrams et al. 1996).

The activity uptakes used in this dissertation are taken from the biodistribution study reported by Kung et al (Kung, Alvai et al. 1990). The percent uptakes reported are for two hours post-injection. The activity uptakes in the brain, striata, lungs, liver, small intestine, and gall bladder are 2%, 1.7%, 8.1%, 9.6%, 14.2%, and 3.8%, respectively. For the purposes of the simulation studies, the liver, small intestine, and gall bladder uptakes are summed to get an overall activity uptake of 27.5%. Kung reports a brain activity uptake of 3.72% two hours post injection (Kung, Alvai et al. 1990). This brain activity uptake is for the whole brain, including the striatum. The striatum uptake is determined by volume analysis and volume estimates provided in the text of the biodistribution study.

Chapter 3

CONJUGATE IMAGING OF TWO ALIGNED SOURCES

As discussed in chapter 2, classic conjugate imaging has been successfully applied to the estimation of activity in a single source and to the estimation of the sum of activities in multiple sources embedded in an attenuating medium. This estimation of activity is independent of source depth. Thomas and Tsui expanded classic conjugate imaging techniques to the estimation of activity in multiple sources aligned and embedded in an attenuating medium (Thomas, Maxon et al. 1976; Tsui, Chen et al. 1981). Their formulations, however, require *a priori* knowledge of activity ratios between the sources. Obtaining such activity ratios usually requires invasive procedures; a requirement that renders the technique clinically difficult at best.

Presented in this chapter is a mathematical formulation for estimation of the activities in each of two aligned sources embedded in a uniformly attenuating medium. The formulation does not require knowledge of the activity ratio between the sources. However, it does require that the depth of each source be known. The requirement of depth measurements introduces new variables as well as complexity as compared to analysis of conjugate imaging. Nonetheless, data acquisition is the same as that for classic conjugate imaging; two emission, two transmission, and two camera sensitivity data sets. The motivation for development of this extension is quantification of specific activity uptake in the striata in Parkinsonism and schizophrenia studies. Thus, the physical model used to formulate this extension is based on the anatomy of the brain, the

striata, and the specific and non-specific binding of radiopharmaceuticals used in imaging the striata.

3.1 Geometric Models

3.1.1 A Model of Striata and Head

Figure 3.1 is an illustration of a model representing activity uptake in the striata and in the surrounding brain tissue. The model consists of the skull, brain tissue, and the striata. The brain tissue has a total thickness T . The skull is represented by thicknesses t_0 and t_6 with attenuation coefficients μ_0 and μ_6 , respectively. The brain tissue is segmented into sections of thicknesses t_1 , t_2 , t_3 , t_4 , and t_5 with attenuation coefficients of μ_1 , μ_2 , μ_3 , μ_4 , and μ_5 , respectively. The striata are represented by thicknesses t_2 and t_4 . The nonspecific activities in the sections of brain tissue surrounding the striata are represented as A_1 , A_3 , and A_5 . The activities in the striatal regions are represented as A_2 and A_4 . Activity A_2 is the sum of the specific activity and the nonspecific activity in the relevant volume. Similarly, A_4 is a sum of specific and nonspecific activity

3.1.2 Simplified Model of Striata and Head

To obtain a more tractable model of striata and head, several assumptions are made about the characteristics of the brain and the striata of the model illustrated in figure 3.1. Assuming that the nonspecific activity concentration is constant throughout the brain tissue, the nonspecific activities in the sections of brain tissue are proportional to the brain thickness.




Figure 3.1. Seven-volume model of the head. Volumes 0 and 6 represent the skull. Regions 1, 3, and 5 represent the brain tissue excluding the striata. Regions 2 and 4 represent the striata.

The total nonspecific activity in the brain is the sum of A_1 , A_3 , A_5 , and the percentages of A_2 and A_4 representing the nonspecific activities in the striatal regions. The nonspecific activity in the brain tissue is referred to as background activity throughout the remainder of the dissertation.

The self-attenuation correction factor of a source with a thickness of 3 cm that emits 140 keV photons has less than 1% effect. Since the thickness of each striatum in humans is on average 3 cm, the activity distributions in the striatal regions are treated, in this model, as if the activity was concentrated in the middle of each striatal region (Hokama, Shenton et al. 1995; Mintzer, Aarsvold et al. 2000). The brain consists of soft tissue regions with similar attenuation coefficients. As a result, it is assumed the brain tissue is a single attenuating medium with a constant attenuation coefficient.

The application of these assumptions is illustrated in figure 3.2. The total brain thickness is denoted T . The skull is represented by thicknesses t_L and t_R ; the subscripts L and R denote the left and right sides, respectively. The attenuation coefficient for skull is

denoted μ_b . The attenuation coefficient for the brain tissue is denoted μ . The uniformly distributed background activity in the brain tissue is denoted A_B . The specifically bound activities in the striata are denoted A_1 and A_2 . The depth of A_1 is d_1 as measured from the left edge of brain tissue. The depth of A_2 is d_2 as measured from the right edge of brain tissue. The distance between the striatal activities is denoted s . The brain model depicted in figure 3.2 is the basis for the work presented in the dissertation. It is used to formulate the mathematics of conjugate imaging of two aligned sources.




Figure 3.2. Simplified seven-region model. The activity in sections 1, 3, and 5 in figure 3.1 represent non-specifically bound activity A_B . This activity is referred to as background activity and is uniformly distributed throughout the brain. Activities A_1 and A_2 represent the sum of specifically bound activity in the striata and background activity. The specific activity distributions in the striata are planar as the thickness of the striata is such that self-attenuation is negligible.

3.2 Parameters Required

Using conjugate imaging to estimate the activity in each of two aligned sources requires knowledge of the striatal depths d_1 and d_2 , the brain thickness T , the skull

thicknesses t_L and t_R , and the attenuation coefficient of the skull μ_b . Computed tomography or magnetic resonance imaging can be used to determine a head thickness, depths of the striata, and skull thicknesses. The skull attenuation factors can be determined alternatively from skull thickness measurements and an empirically determined bone attenuation coefficient. The techniques employed to obtain measurements of these parameters are not the focus of this discussion and will not be discussed further. For this discussion, it is assumed that the striatal depths d_1 and d_2 , the brain thickness T , the skull thicknesses t_L and t_R , and the attenuation coefficient of the skull μ_b have been determined.

3.3 Pixel-based Computations and Analyses

In conjugate-imaging data acquisition, each pixel in each emission, transmission, and flood projection represents the fraction of the imaged object that corresponds to the volume of the object from which the pixel records data. Thus, all activities, sensitivities, transmission factors and thicknesses are implicit functions of pixel location. The count data acquired from emission, transmission, and flood acquisitions are often referred to as maps. Conjugate imaging calculations are performed independently for each pixel within these maps. As a result, the conjugate imaging technique is implicitly pixel dependent. However, some desired quantities, e.g., total activity in each striatum, require that the values in a range of pixels delineated by regions of interest (ROIs) be summed or averaged. The ROIs used in the modified conjugate imaging technique are a striatal ROI and a background ROI. The striatal ROI is defined on an image or data set by delineating a region of inclusion of striatal pixels. The background ROI is a set of pixels defining brain tissue that does not include the striata.

All count maps are normalized and decay-corrected to yield count rate maps at a common time. The number of counts collected in a pixel between times t_1 and t_2 , relative to time t_0 , is

$$N = C \int_{t_1}^{t_2} A_0 e^{-t/\tau} dt, \quad (3.1)$$

where N is the number of counts collected; C is the pixel sensitivity in counts per unit time per unit activity; A_0 is the activity at time t_0 ; t_1 and t_2 are the start count acquisition time and end count acquisition time, respectively; and τ is the mean lifetime of the radioisotope. The normalized, decay corrected count rate at time t_0 represented by measured count N is

$$\dot{N}_o \equiv CA_0 = \frac{N}{\tau \left[e^{-t_1/\tau} - e^{-t_2/\tau} \right]}. \quad (3.2)$$

The denominator in equation 3.2 is referred to as the normalization factor. This factor is determined using a reference time t_0 . Each count map is divided by the appropriate normalization factor to obtain the corresponding decay-corrected count rate map.

3.4 Conjugate Imaging of Two Aligned Sources

As with classic conjugate imaging, two emission, two transmission, and two flood count map data sets are acquired in the conjugate imaging of two aligned sources.

However, the data analysis used in determining the activity above background activity (the striatal specific activity minus the non-specific activity) in each of the two aligned sources is quite different from that used in classic conjugate imaging. In this section, the data analysis in conjugate imaging of two aligned sources is discussed.

The emission data in conjugate imaging of two aligned sources, as in classic conjugate imaging, are measurements of the attenuated activity of the sources inside the object. The emission data are used to determine the activity in the striata above the background activity. The emission count maps of the model in figure 3.2 in the striatal ROI are expressed as

$$E_L = C_L k_L (A_1 e^{-\mu d_1} + A_2 e^{-\mu(d_1+s)} + B) \quad (3.3)$$

$$E_R = C_R k_R (A_1 e^{-\mu(d_2+s)} + A_2 e^{-\mu d_2} + B), \quad (3.4)$$

$$\text{where} \quad k_L = e^{-\mu_L} \quad (3.5)$$

$$k_L = e^{-\mu_L}, \quad (3.6)$$

E_L and E_R are the emission data from the left and right cameras, respectively; C_L and C_R are the camera sensitivities in the left and right cameras, respectively; k_L and k_R are the transmission factors through the left and right skull thicknesses, respectively; A_1 and A_2 are the specific-uptake activities in the left and right striata, respectively; μ is the effective attenuation coefficient of the brain tissue; d_1 and d_2 are the depths of the left and right striata, respectively; s is the distance between the two striata; and B is the effective background activity.

The effective background activity B is expressed as

$$B = \frac{A_B e^{-\mu T/2}}{f_T}, \quad (3.7)$$

where A_B is the true background activity, f_T is the self-attenuation correction factor of the brain, and $e^{-\mu T/2}$ is the transmission factor through one half the brain thickness T .

The measured activities from the left and right gamma cameras corrected for skull attenuation, respectively, are

$$A_L \equiv \frac{E_L}{C_L k_L} = (A_1 e^{-\mu d_1} + A_2 e^{-\mu(d_1+s)} + B) \quad (3.8)$$

$$A_R \equiv \frac{E_R}{C_R k_R} = (A_1 e^{-\mu(d_2+s)} + A_2 e^{-\mu d_2} + B). \quad (3.9)$$

The background activity in the striatal ROI is subtracted from the total activity in the striatal ROI to yield the sum of measured activities above background in each striatum. The specific activities as measured by the left and right gamma cameras within the striatal ROI, respectively, are thus defined as

$$\delta A_L \equiv (A_L - \langle B \rangle_{str}) = A_1 e^{-\mu d_1} + A_2 e^{-\mu(d_1+s)} \quad (3.10)$$

$$\delta A_R \equiv (A_R - \langle B \rangle_{str}) = A_1 e^{-\mu(d_2+s)} + A_2 e^{-\mu d_2}, \quad (3.11)$$

where $\langle B \rangle_{str}$ is the average background activity within the striatal ROI. Equations 3.10 and 3.11 are linear in A_1 and A_2 and can be solved for A_1 and A_2 yielding

$$A_1 = \frac{\delta A_L e^{\mu d_1} - \delta A_R e^{\mu(d_2-s)}}{1 - e^{-2\mu s}} \quad (3.12)$$

$$A_2 = \frac{\delta A_R e^{\mu d_2} - \delta A_L e^{\mu(d_1-s)}}{1 - e^{-2\mu s}} . \quad (3.13)$$

These expressions represent the activity in each pixel in the striatal ROI. These activities are summed over the striatal ROI to yield the total activity in each striatum. The total activity in each striatum is then

$$U_1 = \sum_{i,j \in str} A_1 \quad (3.14)$$

$$U_2 = \sum_{i,j \in str} A_2 , \quad (3.15)$$

where U_1 and U_2 are the total activity measured in each striatum; A_1 and A_2 are the activities of each striatum; and the summation indices i and j indicate pixel location.

The background activity in the striatal ROI is determined from the background activity in a ROI defined in a region on the acquired data that contains only background activity. In a background-only ROI where A_1 and A_2 equal zero, the emission count maps are

$$E_L = C_L k_L B \quad (3.16)$$

$$E_R = C_R k_R B , \quad (3.17)$$

where B is the effective background activity. Equation 3.7 describes effective background activity, taking into account self-attenuation. The measured background activities B_L and B_R in the background-only ROI as measured by the left and right gamma cameras, respectively, are

$$B_L \equiv \frac{E_L}{C_L k_L} \quad (3.18)$$

$$B_R \equiv \frac{E_R}{C_R k_R}. \quad (3.19)$$

The average background activity B for each pixel in the background-only ROI, is determined by taking the geometric mean of B_L and B_R

$$B = (B_L B_R)^{1/2}. \quad (3.20)$$

The background activity values determined for the pixels in the background-only ROI are then averaged to yield an average measured background activity $\langle B \rangle_{\text{bkgd}}$. From equation 3.7, the average activity in the background-only region can be obtained using

$$\langle A_B \rangle_{\text{bkgd}} = \left\langle \frac{B}{e^{-\mu \frac{T}{2}}} f_T \right\rangle_{\text{bkgd}}. \quad (3.21)$$

If the average brain tissue thickness imaged by the pixels in the striatal ROI is equal to that imaged by the pixels in the background-only ROI, the average measured background activity in the background-only ROI is equal to the averaged background activity in the striatal ROI. If the average brain tissue thickness imaged by the pixels in the striatal ROI $\langle T \rangle_{\text{striatal}}$ is different from that imaged by the pixels in the background ROI $\langle T \rangle_{\text{bkgd}}$, the average measured background activity in the striatal ROI is computed by scaling the average total background activity in the background ROI with the average thickness ratio. The approximate average background activity in the striatal ROI is expressed

$$\langle A_B \rangle_{\text{str}} = \langle A_B \rangle_{\text{bkgd}} \frac{\langle T \rangle_{\text{str}}}{\langle T \rangle_{\text{bkgd}}}, \quad (3.22)$$

where $\langle A_B \rangle_{\text{str}}$ is the average background activity in the striatal ROI; $\langle A_B \rangle_{\text{bkgd}}$ is the average background activity in the background-only ROI; $\langle T \rangle_{\text{str}}$ is the average brain tissue thickness within the striatal ROI; and $\langle T \rangle_{\text{bkgd}}$ is the average brain tissue thickness imaged by the background-only ROI. The measured average background activity in the striatal ROI is then

$$\langle B \rangle_{\text{str}} = \langle A_B \rangle_{\text{str}} \frac{e^{-\mu \langle T \rangle_{\text{str}} / 2}}{f_{\langle T \rangle_{\text{str}}}}. \quad (3.23)$$

The two flood count maps are acquired to calculate the camera sensitivities of the left and right gamma cameras C_L and C_R , respectively. The flood data are expressed as

$$F_L = C_L A_F \quad (3.24)$$

$$F_R = C_R A_F, \quad (3.25)$$

where F_L and F_R are the flood counts from the left and right gamma cameras; C_L and C_R are the camera sensitivities for the left and right cameras, respectively; A_F is the activity in the flood.

The transmission data are acquired for use in the calculations of the effective attenuation coefficient of the brain. This differs from the classic conjugate technique in that a map of effective attenuation coefficient μ_{eff} is desired instead of a map of transmission factor $e^{-\mu T}$ of the object. During acquisition of transmission data, the camera is detecting photons from the flood that are attenuated by the skull thicknesses t_L and t_R and the brain thickness T . The transmission data are acquired after the injection of the radiopharmaceutical. As a result, the transmission images will include counts representing transmission and emission gammas. The emission data are subtracted from the transmission/emission data to yield transmission data only. The data are

$$T_L = F_L k_L k_R e^{-\mu T} \quad (3.26)$$

$$T_R = F_R k_L k_R e^{-\mu T}, \quad (3.27)$$

where T_L and T_R are the transmission counts acquired by the left and right cameras, respectively. The transmission factors for the left and right skull thickness are

$$k_L = e^{-\mu_b t_L} \quad (3.28)$$

$$k_R = e^{-\mu_b t_R} . \quad (3.29)$$

Although not necessary, the following assumption is made to simplify some of the discussion. It is assumed

$$t_L = t_R = t_b , \quad (3.30)$$

that is, it is assumed that the left and right skull thicknesses are equal. Calculation of the geometric mean of the T_L and T_R of equations 3.27 and 3.28 yields μT (and thus μ) as

$$\mu T = -\frac{1}{2} \ln \left(\frac{T_L T_R}{F_L F_R} \right) - 2\mu_b t_b . \quad (3.31)$$

The skull transmission factors and the attenuation coefficient maps are then used in the emission data equations.

The sum of the two specific uptake activities A_1 and A_2 can also be determined using the same sets of acquired data – emission, transmission, and flood. This corresponds to use of classic conjugate imaging for determination of the sum of activity of multiple aligned sources.

As previously discussed in this section, the background activity in the striatal ROI $\langle B \rangle_{\text{str}}$ is subtracted from the striatal and background activity emission data in the striatal ROI. The resulting emission data in the striatal ROI as measured by the left and right gamma cameras, respectively, are

$$E'_L \equiv E_L - C_L k_L \langle B \rangle_{str} = C_L k_L (A_1 e^{-\mu d_1} + A_2 e^{-\mu(d_1+s)}) \quad (3.32)$$

$$E'_R \equiv E_R - C_R k_R B_{caud} = C_R k_R (A_1 e^{-\mu(d_2+s)} + A_2 e^{-\mu d_2}), \quad (3.33)$$

where E'_L and E'_R are the left and right specific activities above background activity in the striatal ROI.

The combination of Equations 3.32 and 3.33 leads to the expression of the sum of the specific uptake activities A_1 and A_2 , i.e.

$$(A_1 + A_2) = \left(\frac{E'_L E'_R}{C_L C_R k_L k_R e^{-\mu T}} \right)^{\frac{1}{2}} \left(1 + \frac{2A_1 A_2}{(A_1 + A_2)^2} (\cosh(\mu s) - 1) \right)^{-\frac{1}{2}} \quad (3.34)$$

with

$$k_L k_R e^{-\mu T} = \left(\frac{T_L T_R}{F_L F_R} \right)^{\frac{1}{2}}. \quad (3.35)$$

Equation 3.34 is essentially the classic conjugate imaging equation with a source separation correction factor (the second factor on the right-hand side of the equation). The expression is similar to equation 2.9 in Chapter 2. If the distance between two aligned sources, denoted as s , is reduced to zero, equation 3.34 reduces to the classic conjugate imaging equation for a single source. The magnitude of source separation influences the accuracy of the estimation of total activity. A large source separation results in an overestimation of the total activity.

The total striatal activity is obtained by summing the values of the pixels in the striatal ROI, which were determined using equation 3.34. The total activity in the striatal ROI is then expressed as

$$U_T = \sum_{i,j \in str} (A_1 + A_2), \quad (3.36)$$

where U_T is the total activity; A_1 and A_2 are the specific activity uptakes in the striata; and the summation indices i and j indicate the pixel location. The total activity calculation does not require precise knowledge of the brain tissue attenuation coefficient μ or of the depths of the striata.

3.5 Error Analysis

Regardless of technique, quantitative analysis of radiotracer uptake is associated with sources of error. These errors are due to counting statistics, depth dependent blurring, partial-volume effects, detection of scattered radiation, and possible error in measurement of the standard and/or administered activities. The conjugate imaging method presented in this chapter requires the estimation of several parameters. In addition to usual sources of error associated with estimating radiotracer uptake, errors in these estimated parameters may propagate error and further limit accuracy. The parameters of concern are the skull attenuation coefficient μ_b , the skull thickness t_b , the total brain-tissue thickness T , and the depths of the centers of striata, d_1 and d_2 .

Error in the estimation of $\mu_b t_b$ affects the determination of the average effective brain-tissue attenuation coefficient μ and calculation of the left and right specific activity

maps δA_L and δA_R . Error in T also effects the determination of μ . Errors in μ , δA_L , δA_R , d_1 , and d_2 result in error in the estimation of A_1 and A_2 and ultimately in the estimates of the total specific activity uptakes U_1 and U_2 .

The skull attenuation coefficient μ_b needs to be determined either by scaling the attenuation coefficient measured with tissue-equivalent material by a factor appropriate to the different compositions of tissue and skull matter or by actual experiments on samples of skull material. The thickness of the head, which is $T+2t_b$, can be measured externally. A computed tomography scan of the head is required to estimate accurately the dimension of the skull thickness t_b and the depths of the striata d_1 and d_2 . A database of these dimensions as a function of externally accessible measurements for primates or humans could be used. However, using dimensions from such a database might increase uncertainties in these quantities.

3.5.1 Uncertainty of Bone Attenuation Coefficient μ_b

The error in μ due to uncertainties in $\mu_b t_b$ can be estimated using,

$$\mu T = -\frac{1}{2} \ln \left(\frac{T_L T_R}{F_L F_R} \right) - 2\mu_b t_b. \quad (3.37)$$

This equation is used to determine the effective attenuation coefficient of brain tissue. (Refer to equation 3.31.) The last term $2\mu_b t_b$ represents and removes the attenuation of the skull. To estimate the error in μ due to uncertainties in $\mu_b t_b$ typical values for skull and brain thicknesses and an estimate of μ_b of primates and humans are used. The skull and

brain thicknesses are assumed to be 5 mm and 10 cm, respectively. Using attenuation measurements of 140 keV photons in water published by Thomas et al. the skull attenuation coefficient μ_b is determined to be approximately 1.3μ (Thomas, Maxon et al. 1976). The $2\mu_b t_b$ term results in about a 13% correction to μT using a 5 mm skull thickness, a 10 cm brain thickness, and a skull attenuation coefficient of 1.3μ . A 20% overestimate of $\mu_b t_b$ would cause about a 2.5% underestimate of μT . Error in the measurement of brain thickness T would also affect the accuracy of the value for μ . However, assuming that the determination of T to within 1 or 2% is possible, the total uncertainty in μ can be kept less than 5%.

Uncertainties in the specific-uptake activity maps δA_L and δA_R due to errors in $\mu_b t_b$ are estimated using the same values for the estimation of error in μ . A 20% overestimation of $\mu_b t_b$ causes less than 2% underestimates of δA_L and δA_R . Due to this uncertainty, we expect estimation of $\mu_b t_b$, T , and μ to result in an error of less than 2% in U_1 and U_2 .

3.5.2 Uncertainty of A_1 -to- A_2 Ratio

The uncertainty in the activity in each striatum A_1 and A_2 , due to errors in μ and errors in the depth measurements d_1 and d_2 , can be expressed as a function of the ratio of the two source activities

$$\frac{A_1}{A_2} = \frac{\delta A_L e^{\mu d_1} - \delta A_R e^{\mu(d_2-s)}}{\delta A_R e^{\mu d_2} - \delta A_L e^{\mu(d_1-s)}}. \quad (3.38)$$

where A_1 and A_2 are the activity uptakes in the striata and δA_L and δA_R are the left and right measured specific-uptake activities in the striatal ROI, respectively. Assuming d_1 and d_2 are equal, Equation 3.38 becomes

$$\frac{A_1}{A_2} = \frac{r - e^{-\mu s}}{1 - r e^{-\mu s}} \quad (3.39)$$

$$r = \frac{\delta A_L}{\delta A_R}, \quad (3.40)$$

where r is the ratio of measured specific activities in the striatal regions as detected by the left and right gamma cameras. Assuming the activities of the sources to be at least the background activity ($A_1 \geq 0$, $A_2 \geq 0$), the ratio of measured specific activities r is limited by the attenuation coefficient μ and distance between the two aligned source distributions s . The actual relationship is $e^{\mu s} \geq r \geq e^{-\mu s}$. For the values of $\mu \approx 0.14 \text{ cm}^{-1}$ and $s \approx 2 \text{ cm}$, r ranges from about 0.76 to 1.32. The ratio ranges from 0.8 to 1.2 for striatal uptake ratios of 1:5 to 5:1. For values of $\mu \approx 0.15 \text{ cm}^{-1}$ and $s \approx 2 \text{ cm}$, r ranges from about 0.7 to 1.2 for activity ratios of 1:5 to 5:1.

3.5.3 Uncertainty of Brain Tissue Attenuation Coefficient μ

The fractional uncertainty in A_1 due to the uncertainty in μ is

$$\frac{\Delta A_1}{A_1} = \left(\frac{d}{s} + \frac{e^{-\mu s}}{r - e^{-\mu s}} - \frac{2e^{-2\mu s}}{1 - e^{-2\mu s}} \right) s \Delta \mu, \quad (3.41)$$

where $\Delta A_1/A_1$ is the fractional uncertainty; d is the depth of each striatum; s is the separation distance between the striata; r is the ratio of the left and right measured specific-uptake activities; and $\Delta\mu$ is the uncertainty in μ . We have assumed that d_1 and d_2 are equal and have a value d . Assuming $T \approx 10$ cm, $t_b \approx 5$ mm, $d \approx 4$ cm, $s \approx 2$ cm, and $\mu \approx 0.142$ cm⁻¹, we find that a 2.5% underestimate in μ causes an underestimate of A_1 from less than 1% to about 7% for the values of r corresponding to activity ratios of A_1 to A_2 from 5:1 to 1:5.

Due to the symmetry of the problem, and thus of the expressions for A_1 and A_2 , the uncertainty in A_2 can be found by using the reciprocal value of r . The fractional uncertainty in A_2 due to the uncertainty in μ is

$$\frac{\Delta A_2}{A_2} = \left(\frac{d}{s} + \frac{re^{-\mu s}}{1 - re^{-\mu s}} - \frac{2e^{-2\mu s}}{1 - e^{-2\mu s}} \right) s \Delta\mu. \quad (3.42)$$

Assuming $T \approx 10$ cm, $t_b \approx 5$ mm, $d \approx 4$ cm, $s \approx 2$ cm, and $\mu \approx 0.142$ cm⁻¹, a 2.5% underestimate in μ causes an underestimate of A_2 from about 7% to about 2% for the values of r corresponding to activity ratios of A_1 and A_2 from 5:1 to 1:5. For activity ratios not equal to one, larger percentage errors occur in estimates of the smaller of the two activities.

3.5.4 Uncertainty of Position

This section discusses the uncertainty in A_1 and A_2 due to the errors in the measurements of s , the separation between the two striata, and the depths d_1 and d_2 of the striata. The total thickness of T is assumed to be constant and expressed as

$$T = d_1 + d_2 + s . \quad (3.43)$$

With T held constant, an error in one of the parameters will result in errors in one or both of the other parameters. The error analysis in this section exploits this relationship between s , d_1 , and d_2 .

The fractional uncertainty in A_1 due to the uncertainty in the separation distance s is

$$\frac{\nabla A_1}{A_1} = - \left(1 + \frac{4e^{-2\mu s}}{1 - e^{-2\mu s}} - \frac{2e^{-\mu s}}{r - e^{-\mu s}} \right) 2\mu \Delta s , \quad (3.44)$$

where $\Delta A_1/A_1$ is the fractional uncertainty; s is the separation distance between the striata; r is the ratio of the left and right measured specific-uptake activities; μ is the attenuation coefficient of the brain tissue; and Δs is the uncertainty in s . A measurement of s with an error of plus or minus one millimeter (± 1 mm) results in an error in A_1 of less than 0.5% to about 10% for the values of r corresponding to activity ratios from 5:1 to 1:5.

The uncertainty in A_2 due to error in the measurement of separation distance s can be found by using the reciprocal value of r . The fractional uncertainty in A_2 due to the uncertainty in μ can be expressed as

$$\frac{\Delta A_2}{A_2} = - \left(1 + \frac{4e^{-2\mu s}}{1 - e^{-2\mu s}} - \frac{2re^{-\mu s}}{1 - re^{-\mu s}} \right) 2\mu \Delta s. \quad (3.45)$$

A measurement of s with an error of plus or minus one millimeter (± 1 mm) results in an error in A_2 of less than 10% to about 0.5% for the values of r corresponding to activity ratios of A_1 and A_2 from 5:1 to 1:5.

The fractional uncertainty in A_1 due to an uncertainty in the d_1 measurement is

$$\frac{\Delta A_1}{A_1} = \left(1 + \frac{2e^{-2\mu s}}{1 - e^{-2\mu s}} \right) \mu \Delta d_1, \quad (3.46)$$

where $\Delta A_1/A_1$ is the fractional uncertainty; s is the separation distance between the striata; μ is the attenuation coefficient of the brain tissue; and Δd_1 is the uncertainty in d_1 . The fractional uncertainty of A_1 due to error in the estimation of the position of d_1 is independent of r . A measurement of d_1 with an error of plus or minus one millimeter (± 1 mm) translates to approximately a 5% error in A_1 .

The fractional uncertainty in A_1 due to an uncertainty in the d_2 measurement is

$$\frac{\Delta A_1}{A_1} = \left(\frac{1}{re^{\mu s} - 1} - \frac{e^{-2\mu s}}{1 - e^{-2\mu s}} \right) 2\mu \Delta d_2, \quad (3.47)$$

where $\Delta A_1/A_1$ is the fractional uncertainty; s is the separation distance between the striata; r is the ratio of the left and right measured specific-uptake activities; μ is the attenuation coefficient of the brain tissue; and Δd_2 is the uncertainty in d_2 . A measurement of d_2 with an error of plus or minus one millimeter (± 1 mm) results in an error in A_1 of about 1% to about 5% for the values of r corresponding to activity ratios from 5:1 to 1:5.

Due to the symmetry of the problem, the fractional uncertainty in A_2 due to measurement errors in d_2 is expressed as

$$\frac{\Delta A_2}{A_2} = \left(1 + \frac{2e^{-2\mu s}}{1 - e^{-2\mu s}} \right) \mu \Delta d_2. \quad (3.48)$$

The fractional uncertainty of A_2 due to error in the estimation of the position of d_2 is not dependent on the measured specific activities in the striatal regions, and thus, is independent of r . A measurement of d_2 with an error of plus or minus one millimeter (± 1 mm) results in about a 5% error in A_2 for activity ratios of A_1 and A_2 from 5:1 to 1:5.

Similarly, the fractional uncertainty in A_2 due to measurement errors in d_1 is expressed as

$$\frac{\Delta A_2}{A_2} = \left(\frac{1}{re^{\mu s} - 1} - \frac{e^{-2\mu s}}{1 - e^{-2\mu s}} \right) 2\mu \Delta d_1. \quad (3.49)$$

The fractional uncertainty in A_2 due to uncertainty in d_1 is dependent on r . A measurement of d_2 with an error of plus or minus one millimeter (± 1 mm) results in an

error in the A_1 of about 5% to about 25% for the values of r corresponding to activity ratios of A_1 and A_2 from 5:1 to 1:5.

Assuming A_L and A_R do not vary greatly over the striatal ROI, these fractional uncertainties are estimates of the fractional uncertainties in U_1 and U_2 . The fractional uncertainties in U_1 and U_2 due to positioning error are summarized for the values used throughout this error analysis section in table 3.1. For activity ratios not equal to one, the larger percentage error is associated with the estimation of the smaller of the two activities.

The above error analysis indicates that reasonable accuracy can be achieved for the case of approximately equal uptakes. Estimates of the smaller of the two activities become increasingly less certain as the ratio between the two activity uptakes deviates from one. Table 3.2 lists all the sources of error discussed and gives estimates of the worst-case uncertainty to be expected due to the assumptions used in this method.

Table 3.1 Positioning Errors

Type and magnitude of positioning error	% Error in U_1, U_2		
	$U_1/U_2 - 1$	$U_1/U_2 - 2$	$U_1/U_2 - 5$
$\pm 1\text{mm}$ error in separation s ; d_1, d_2 constant	$\pm 0.1\%, \pm 0.1\%$	$\pm 1.3\%, \pm 2.4\%$	$\pm 2.1\%, \pm 9.8\%$
$\pm 1\text{mm}$ error in depth d_1 ; d_2 constant	$\pm 5.1\%, \pm 4.9\%$	$\pm 5.1\%, \pm 9.9\%$	$\pm 5.1\%, \pm 25\%$
$\pm 1\text{mm}$ error in depth d_2 ; d_1 constant	$\pm 4.9\%, \pm 5.1\%$	$\pm 2.5\%, \pm 5.1\%$	$\pm 1.0\%, \pm 5.1\%$

Table 3.2 Summary of Errors

Source of error	Estimate of error	% Error in U_1, U_2		
		$U_1/U_2 - 1$	$U_1/U_2 - 2$	$U_1/U_2 - 5$
estimation of $\mu_b t_b$	$\pm 20\% \Rightarrow \pm 2.0\%$ in $\delta A_L, \delta A_L$	$\longleftrightarrow \pm 2.0\% \pm 2.0\% \longrightarrow$		
	$\Rightarrow \pm 2.5\%$ in μ	$\pm 1.7\% \pm 1.7\%$	$\pm 1.1\% \pm 3.0\%$	$\pm 0.7\% \pm 6.7\%$
measurement of T (± 2 mm)	$\pm 2.0\% \Rightarrow \pm 2.0\%$ in μ	$\pm 1.4\% \pm 1.4\%$	$\pm 0.9\% \pm 2.4\%$	$\pm 0.6\% \pm 5.3\%$
measurements of caudate positions, worst case error in depth $d_1, \pm 1$ mm; d_2 constant		$\pm 5.1\% \pm 4.9\%$	$\pm 5.1\% \pm 9.9\%$	$\pm 5.1\%, \pm 25\%$
Total uncertainty due to parameter measurement and estimation errors		$\pm 10\%, \pm 10\%$	$\pm 9.1\% \pm 17\%$	$\pm 8.4\% \pm 39\%$

Chapter 4

PHANTOM EXPERIMENTS

Phantom studies are conducted to characterize the accuracy of conjugate imaging applied to the estimation of two activities in two aligned sources, embedded within a uniformly attenuating medium. A quantitative SPECT study is also conducted to compare the accuracy of current clinical SPECT to that of the conjugate imaging technique. In addition, the sum of the two activities are calculated and compared to the sum of the activities calculated by standard conjugate imaging. The ratio of activity concentrations in each source and surrounding medium are varied to simulate the striatal activity uptake associated with Parkinson's disease patients. The ratios of the activity uptake in the striata are 1:1, 2:1, and 5:1. The 1:1 ratio represents the activity uptake of a normal individual. The 2:1 ratio represents the activity uptake of a Parkinson's disease patient presenting clinical symptoms. The 5:1 ratio represents the activity uptake of a severe Parkinson's disease patient. Since the background activity concentration at the time of imaging is relatively low for the relevant radiopharmaceuticals (as discussed in chapter 4), the ratio of striata activity concentration and background concentration is assumed and modeled to be 9:9:1 for the 1:1 striatal uptake ratio, 12:6:1 for the 2:1 striatal uptake ratio, and 15:3:1 for the 5:1 striatal uptake ratio. For these studies the source with the highest activity concentration is the left sphere and the source with the lower activity

concentration is the right sphere. A study is also performed with no background activity; the ratio of the two source activities is 5:1.

Although most radiopharmaceuticals used in imaging the dopaminergic system are tagged with ^{123}I , $^{99\text{m}}\text{Tc}$ pertechnetate is used in the studies. $^{99\text{m}}\text{Tc}$ pertechnetate is chosen because $^{99\text{m}}\text{Tc}$ photon emission spectrum is not as complex as that of ^{123}I . The high-energy photon emissions of ^{123}I complicate quantitative imaging calculations by introducing downscatter contributions to the photopeak of interest. Thus, foundational experimental assessment of the mathematical model presented in Chapter 3 is initiated with $^{99\text{m}}\text{Tc}$. Although part of the motivation for the use of $^{99\text{m}}\text{Tc}$ is its less complex photon emission spectrum, $^{99\text{m}}\text{Tc}$ is clinically relevant for dopaminergic studies as $^{99\text{m}}\text{Tc}$ TRODAT is used in the evaluation of striatal dopamine in Parkinson's disease studies.

4.1 Methods

4.1.1 Conjugate Imaging

4.1.1.1 *Equipment*

Imaging System. A Marconi IRIX triple-head gamma camera system is used for image data acquisition. Heads II and III are used to acquire conjugate image data. Heads II and III are positioned 34 cm apart at 180° separation and are equipped with low-energy high-resolution (LEHR) collimators. A 15% energy window is placed symmetrically about the 140 keV $^{99\text{m}}\text{Tc}$ photopeak. Images are acquired as 512 pixels x 512 pixels matrices with 1.18 mm x 1.18 mm pixels.

Phantom. A cylindrical tank with two spherical inserts, representative of the brain with striata, is used. When radioactivity is placed in the spherical inserts and the tank in

quantities relative to the activity concentration ratios, this simple phantom models the specific activity uptake in the striata and background uptake in the brain. The phantom wall thickness is 0.635 cm (0.25"). The phantom inner height and inner diameter are 30.48 cm (12") and 20.32 cm (8.0"), respectively. The outer volume of each sphere is 16.2 ml. The spheres were separated by 10.16 cm (4").

Flood/Transmission Source. A rectangular tank is used as the calibration flood and as the transmission source. Its inner dimensions are 58.75 cm x 43.0 cm x 2.54 cm. Self-attenuation of the flood source is accounted for by comparing the count rates of activity in a small syringe to the count rates of the same activity dispensed in the water-filled flood. The small syringe represents a source with negligible self-attenuation. The count rates are normalized to an experiment start time and the acquired data were decay corrected. The count rates from the syringe and the flood are determined from 256 pixels x 256 pixels images. Due to the fact that the physical dimensions of the flood image are larger than the 256 pixels x 256 pixels image, the activity imaged is a fraction of the total activity in the flood source. The pixel to flood size ratio is calculated and multiplied by the total number of pixels used (256 x 256) to determine the source fraction – the percentage of the flood activity in the 256 pixel x 256 pixel area. The count rate from the syringe is divided into the flood count rate scaled by the source fraction. The flood source self-attenuation factor is thus measured to be 85.3%.

4.1.1.2 Activity Concentrations

A weighing technique is used in the dispensing of the ^{99m}Tc pertechnetate solution into the spheres and phantom to attain the desired ratios. A pre-weighed tuberculin

syringe is filled with Tc-99m pertechnetate solution (~1.0 ml). The total activity (~4.5 mCi) is measured in a dose calibrator. Portions are dispensed by weight into the flood (~2.5 mCi), cylindrical tank (~1.0 mCi), and an additional volume of water (~25 ml). The further dilution of the solution is required to dispense more accurately much smaller amounts of activity (~5 to 30 μ Ci) into the spheres. After weighing the diluted activity, its activity concentration is calculated, and another syringe is weighed and filled. Portions of diluted solution are then dispensed by weight into the spheres. An analytic scale with better than 0.1 mg accuracy is used, resulting in better than 1% accuracy in this process of dividing the initial activity. Such accuracy is not achievable by separate dose-calibrator measurements of each quantity dispensed; basic accuracy and linearity of dose calibrators is typically no better than 2% within a single range, and is worse between ranges. The measurement of activities in different volumes can also result in inaccuracies due to different self-absorption and source-detector geometries.

4.1.1.3 Data Acquisition

The flood data are acquired for 30 minutes. The left and right transmission/emission data are acquired for 15 minutes each and the emission data for 15 minutes. Imaging times are adjusted to compensate approximately for decay; these adjustments ensured that all image data are approximately equivalent to that which would have resulted if all imaging had commenced at the same time.

Each data set requires a different imaging geometry. For the collection of flood data, the flood source is positioned midway between the cameras. For emission data acquisition, the phantom is placed midway between the cameras. When transmission data

are acquired by camera head II, the flood source is placed between the phantom and camera head III. For transmission data acquisition by camera head III, the flood source is placed between the phantom and camera head II. Figure 4.1 illustrates the experimental set-ups. Figure 4.2 shows transmission image acquisition on head III.

4.1.1.4 Data Processing

Software is developed in MATLAB to perform calculations for conjugate imaging of two aligned volumes and standard conjugate imaging. The code processes the image data and calculates the effective attenuation coefficient μ , the background activity concentration B , the activities above background in each of the two volumes A_1 and A_2 , and the sum of the activities above background A_T . The data are decay-corrected and normalized for acquisition times and durations. This converts the data into count rates at a common calibration time. All images acquired with head II are flipped horizontally for processing so that corresponding pixels in the opposing views coincided. The camera sensitivities are then calculated from the flood images. Emission and transmission rates are corrected for cylinder wall attenuation using a measured wall attenuation coefficient, μ_b ($\sim 0.118 \text{ cm}^{-1}$). The emission rates are subtracted from the transmission/emission rates to yield transmission rates. Due to the magnitude of the statistical fluctuations in the transmission data, the effective attenuation coefficient, μ , is determined using the fractional flood transmission averaged over a 37 pixel x 37 pixel source region of interest (ROI) centered on the source activities' projections. The effective background activity is determined by averaging emissions over a 37 pixel x 37 pixel background ROI. The source activities above background are then calculated and summed over the source ROI.




Figure 4.1. Experimental set-up. The top figure illustrates the acquisition of flood data, which are used to determine the cameras' sensitivities. The middle figure illustrates the acquisition of transmission data, which are used to determine the effective attenuation coefficient. The bottom figure illustrates acquisition of the emission data, which are used to determine the activity uptakes A_1 and A_2 .



Figure 4.2. The transmission data acquisition set-up. Shown are the cylinder-with-two-spheres phantom, the flood/transmission source, and the gamma camera system used in the experiments. The configuration shown is specifically that used to obtain transmission/emission image data from head III; the transmission/emission left data are acquired similarly, except the flood is on the right side of the phantom rather than the left. Flood data are acquired with the cylindrical phantom removed and the flood/transmission source centered between the two cameras. The emission data are acquired with the cylindrical phantom centered between heads II and III.

The background activity concentration is determined by correcting the effective background activity per pixel for self-attenuation and dividing by T and the pixel area.

Figure 4.3 shows flood, transmission/emission and transmission image maps.

The sum of the activities of the aligned volumes is estimated by summing the two estimates obtained from the above-outlined computations. The sum of the activities is also determined by applying a standard conjugate imaging analysis. In the case of two sources separated by a sufficiently large distance, the conjugate imaging calculations are modified to include a correction factor that accounts for the large distance. If the modifications are not made, the sum of the two activities would be overestimated. The correction factor depends on the product of the effective attenuation coefficient μ_{eff} and the distance between the volumes s and the ratio of the activities of the sources. Since the ratio is not known a priori, it is estimated by computing the ratio of the estimated source activities.

4.1.2 Quantitative SPECT

Two quantitative SPECT studies are performed. For the scope of this work, quantitative SPECT studies include attenuation correction and exclude scatter correction. This approach is used because most current clinical SPECT studies are not performed using scatter correction techniques. The two SPECT studies are performed with no background activity and a spheres activities ratio of 5:1. The results of these studies are compared to the conjugate imaging study performed with no background activity and a 5:1 spheres activities ratio.




Figure 4.3. Examples of flood, transmission/emission, emission, and transmission activity maps for heads II and III. The maps shown represent a region that covers the entire phantom. The white squares are representative of the source and background ROIs used in the computations. Flood activity maps (upper left). Transmission/emission activity maps (upper right). Emission activity maps (lower left). Transmission activity maps (lower right).

4.1.2.1 Equipment

Imaging System. The Marconi IRIX triple-head gamma camera system is used for SPECT data acquisition. Heads I, II, and III are used. They are separated by 120°. Each camera head is equipped with a LEHR collimator. A 15 % energy window is placed symmetrically about the 140 keV ^{99m}Tc - photopeak. Images were acquired as 128 x128 matrices with 4.72 mm x 4.72 mm pixels.

Phantom. The phantom used in the conjugate imaging study is used in the quantitative SPECT study. The reader is referred to the section 4.1.1.1 for a detailed description of the phantom.

4.1.2.2 Activity Concentrations

The weighing technique used in the conjugate imaging study to dispense the activity is used in the quantitative SPECT study as well. The reader is referred to section 4.1.1.2 for a detailed description of the weighing technique.

4.1.2.3 Data Acquisition

SPECT acquisition on the Marconi system involves a two-step process – emission and transmission acquisition. Transmission acquisition is performed for attenuation correction. The transmission scan is performed in two parts. The phantom is placed on the IRIX patient bed for the emission acquisition and for the first transmission acquisition. The phantom is removed from the field of view during a blank scan, the second transmission scan. The transmission and emission scans of the phantom also

includes the imaging of a point source of known activity, a calibration source for use in the quantitative analysis of the SPECT data.

The emission data are acquired with a continuous circular orbit with 120 projections per head over 360° at 1° per step and 8 seconds per stop. The Marconi IRIX system is equipped with BEACON attenuation-correction technology. This attenuation-correction technology uses two 10 mCi Ba-133 emission sources, with primary emission 356 keV. The first transmission scan is performed with camera heads I and II and the Ba-133 transmission sources. The data are acquired with a step-and-shoot circular orbit with 180 projections per head over 360° at 1° per step and 6 seconds per stop. The second scan is a static, blank, calibration scan of 5-minute acquisition. Camera heads I and II are used.

4.1.2.4 Data Processing

The collected SPECT data are processed into three data sets. These data sets are (1) filtered non-attenuation-corrected data (Filtered Non-AC), (2) filtered attenuation-corrected data (Filtered AC), and (3) non-filtered attenuation-corrected data (Non-Filtered AC). Filtered Non-AC data are reconstructed using ordered subsets maximum likelihood expectation maximization (OSEM) and post-filtered with the default smoothing filter. Attenuation correction is not applied. Filtered AC data are reconstructed using OSEM, corrected for attenuation, and post-filtered with the default smoothing filter. Non-Filtered AC data are reconstructed with OSEM and corrected for attenuation. No post-filtered is performed on Non-Filtered AC data. The activity in each sphere is determined by comparing the calibration source count information. ROIs are placed about the calibration source and the two spherical sources on image slices in which the sources

are visible. The size of the ROIs changes with each slice to account for the spherical volume of the sources. In other words, for each slice of the spherical volume, the circle representing a slice of the volume changes in size with the largest circle being defined halfway through the volume. The total counts for each sphere are determined by adding the counts in the ROIs from the slices. For the first SPECT study, a total of 15 projections are used for each sphere. For the second SPECT study, 16 projections are used for each sphere. The total counts from the calibration source are divided by the activity in the calibration sphere at the time of the study to estimate counts per unit activity. The total counts for each spherical source is divided by the counts per unit activity to determine the activity amount in each sphere.

4.2 Results

4.2.1 Conjugate Imaging

As previously mentioned, the left sphere contained more activity than the right sphere for activity uptake ratios not equal to one. Tables 4.1, 4.2, and 4.3 display the calculated estimates and corresponding percent errors for the uptake ratio studies that included background activity. Table 4.1 displays the effective attenuation coefficient and the estimated background activity concentration for each experiment. Table 4.2 shows the estimated activities in the left and right spheres. The estimates in these two tables are determined using the conjugate imaging algorithm for two aligned sources. Table 4.3 shows the sums of the estimated left and right activities listed in Table 4.2 and shows the estimates of the total activity as computed using the standard conjugate algorithm

Table 4.1. Effective Attenuation Coefficient and Background Activity Concentration

	Effective		%Error
Nominal Ratio	Attenuation Coefficient	Background Activity	Background Activity
	$\mu_{\text{eff}} (\text{cm}^{-1})$	Concentration ($\mu\text{Ci/cc}$)	Concentration
1:1	0.122	0.118	+ 13.5
2:1	0.117	0.119	+ 10.7
5:1	0.121	0.130	+ 13.0
1:5	0.128	0.135	+ 17.4

Table 4.2. Estimated Above Background Activity in Left and Right Spheres

	Left Sphere	% Error	Right Sphere	% Error
Nominal Ratio	Above-Background	Left Sphere	Above-Background	Right Sphere
	Activity (μCi)	Above -Background Activity	Activity (μCi)	Above -Background Activity
1:1	16.0	+ 11.1	14.0	+ 2.2
2:1	20.7	+ 2.5	8.14	- 12.8
5:1	31.8	+ 8.2	3.21	- 28.0
1:5	4.62	+ 3.6	33.0	+ 12.2

Table 4.3. Sum of Estimated Above Background Activities and Estimated Total Activity

	Left Sphere	% Error		
Nominal Ratio	Plus	Left Sphere	(Corrected) Classic Conjugate	% Error
	Right Sphere	Right Sphere	Total Activity (μCi)	(Corrected) Classic Conjugate
	Activity	Activity		Total Activity
1:1	30.0	+ 6.8	28.6	+ 1.8
2:1	28.8	+ 2.3	26.9	- 8.2
5:1	35.0	+ 3.4	30.5	- 9.9
1:5	37.6	+ 11.1	32.8	- 3.1

corrected for source separation. Standard conjugate imaging calculations, corrected for the separation of the two spheres, are used to estimate the sum of the two activities. The results of this computation are compared to the sum of the estimated source activities calculated from the modified conjugate technique. Corrected standard conjugate calculations yield total activity values within 10%. Total activities obtained by summing the estimates of source activities computed using modified conjugate calculations are within 11.1%.

Tables 4.4, 4.5, and 4.6 display the calculated estimates and percent errors for the uptake ratio studies that did not include background activity. Table 4.4 displays the calculated effective attenuation coefficient and calculated background activity concentration. Table 4.5 shows the estimated above background activity in the left and right spheres and associated percent errors. Table 4.6 compares the total activity estimates calculated from summing the individual activity estimates and calculated from the classic conjugate imaging technique with separation correction. Also listed in table 4.6 is the percent error associated with the classic imaging technique.

Table 4.4. No Background Activity Study: Estimates of Effective Attenuation Coefficient and Background Activity Concentration

Nominal Ratio	Effective	Background Activity
	Attenuation Coefficient	Concentration ($\mu\text{Ci/cc}$)
	$\mu_{\text{eff}} (\text{cm}^{-1})$	
5:1 *	0.127	0.003
5:1 **	0.123	0.005

* No ambient background subtraction

** Ambient background subtraction

Table 4.5. No Background Activity Study Estimated Above Background Activity in Left and Right Spheres

	Left Sphere	% Error	Right Sphere	% Error
Nominal Ratio	Above-Background	Left Sphere	Above-Background	Right Sphere
	Activity (μCi)	Above -Background Activity	Activity (μCi)	Above -Background Activity
5:1 [*]	34.0	+ 14.6	4.3	- 30.3
5:1 ^{**}	33.3	+ 12.2	3.6	- 42.1

* No ambient background subtraction

** Ambient background subtraction

Table 4.6. No Background Activity Study
Sum of Estimated Above Background Activities and Estimated Total Activity

Nominal Ratio	Left Sphere Plus Right Sphere Activity	% Error Left Sphere Plus Right Sphere Activity	(Corrected) Classic Conjugate Total Activity (μCi)	% Error (Corrected) Classic Conjugate Total Activity
5:1	38.3	+ 6.7	34.5	- 3.9
5:1	36.9	+ 2.8	33.1	- 7.8

* No ambient background subtraction

** Ambient background subtraction

4.2.2 Quantitative SPECT Study

As with the conjugate imaging study, the left sphere contains more activity than the right sphere for activity uptake ratios not equal to one. The Filtered Non-AC SPECT data underestimate the activity in the left sphere by about 70% in both SPECT studies. The activity in the right sphere is underestimated by about 73% and 74% using Filtered Non-AC SPECT data in SPECT studies 1 and 2, respectively. In SPECT study 1, the

Filtered AC SPECT data overestimate the activity in each sphere by about 12%. In SPECT study 2, the Filtered AC SPECT data overestimate the activity in the left sphere by about 11% and overestimates the activity in the right sphere by 10%. In SPECT study 1, the Non-Filtered AC SPECT data overestimate the activity in the left sphere by about 12% and overestimate the activity in the right sphere by about 16%. In SPECT study 2, the Non-Filtered AC SPECT data overestimate the activity in the left sphere by about 11% and overestimate the activity in the right sphere by about 14%. Tables 4.7 and 4.8 display the percent error for the three sets of data for SPECT studies 1 and 2.

Table 4.7. Percent Errors for quantitative SPECT Studies I

qSPECT 1				
Technique	Left Sphere Above-Background Activity (μCi)	% Error Left Sphere Above-Background Activity	Right Sphere Above-Background Activity (μCi)	% Error Right Sphere Above-Background Activity
Filtered Non-AC	6.3	-70.1	1.2	-72.7
Filtered AC	23.6	11.6	4.9	11.6
Non-Filtered AC	23.7	12.4	5.2	16.4

Table 4.8. Percent Errors for quantitative SPECT Studies II

qSPECT 2				
Technique	Left Sphere Above-Background Activity (μCi)	% Error Left Sphere Above-Background Activity	Right Sphere Above-Background Activity (μCi)	% Error Right Sphere Above-Background Activity
Filtered Non-AC	5.4	-70.4	1.0	-73.8
Filtered AC	20.2	10.8	4.2	9.7
Non-Filtered AC	20.2	11.2	4.3	14.2

4.3 Discussion

This phantom study is conducted to evaluate the application of a conjugate imaging technique to the estimation of activity in each of two aligned sources. Four conjugate imaging experiments are conducted using three different source-to-source ratios (nominal ratio values were 1:1, 2:1, and 5:1). Estimates of the source activities and the sums of the two activities are determined. Tables 4.1 – 4.6 summarize the results. The values in the tables associated with the nominal ratio 1:5 are the result of an extension of the experiment involving the 5:1 ratio.

The effective attenuation coefficient for each experiment is determined and used in estimation of the background concentration and the source activities. The estimates of the background concentrations are all overestimates (10.7% to 17.4%). Some statistical errors due to transmission maps with relatively low counts might be expected, but there may also be a systematic component to the error contributing to an overestimation of the effective attenuation coefficient.

Most of the estimates of the activities in the left and right spheres are within 13% of the known values. The 5:1 experiment resulted in a 28% underestimation of the smaller of the two sphere activities. An extension of this experiment is conducted to initiate investigation into the source of this error. The spheres used in the 5:1 experiment, together with their activities, are switched by rotating the base, upon which the spheres are mounted. The base is rotated until the right sphere is occupying the position of the left sphere and the left sphere is occupying the position of the right. Since the activity uptake ratios assume the sphere in the left position to be the reference source, the activity uptake ratio for this experiment is 1:5. The 1:5 experiment results in overestimates of the two activities that are less than 13% (3.6% and 12.2%). Throughout the computations, the source depths and separation distances are assumed known. A possible contributor to the large underestimation of the smaller activity in the 5:1 experiment may be a property of the phantom. The two spheres used as the sources are mounted on rods extending into the water in the cylinder. It is assumed the sources are separated by 10.1 cm, the separation of the mounting holes. It is possible that the rods, when attached to the cylinder cover, are not sufficiently perpendicular to the cover and thus the sphere positions with respect to the cylinder center are not the same as those of the mounting holes. If this is the case, the assumed values for the source depths and separation could have been incorrect for one or more experiments. Theoretical error analysis suggests that the conjugate imaging computations are sensitive to the accuracy of these values. A small mis-positioning of the larger-activity sphere could cause a significant fractional error in the estimation of the activity in the smaller-activity sphere. To corroborate the theoretical error analysis, the position of the left sphere used in the conjugate imaging calculations is moved in ± 1 mm

increments from its original position. The percent errors in the left and right activity estimations are calculated using each new source location. Since a mis-positioning error in the left sphere results in a bias between the percent errors for left and right activity estimates that increases with increasing activity uptake ratio, the left sphere is moved until the percent errors are relatively equal for all three ratios. It is estimated that the left sphere is approximately +3.6 mm from the putative position.

Most of the estimates computed using the data of these experiments are very good estimates given that some of the raw data images are count limited and no scatter correction is performed. The phantom used is a cylinder with a relatively large diameter and not representative of typical human or primate head sizes. Imaging of a smaller phantom with less attenuation would be more realistic and result in higher count transmission/emission images and thus higher count transmission maps. All the data are collected at rates designed to keep camera dead time errors below 1%. This may not have been appropriate, as the reduced counts in the data may have produced errors greater than the errors that would have resulted if more dead time had been tolerated. The strength of the flood/transmission source is a parameter that needs additional investigation.

The quantitative SPECT study is performed to compare the conjugate imaging calculations to clinically feasible techniques. As with the conjugate imaging studies, the SPECT study is performed with no scatter correction. Currently, scatter correction algorithms are not used widely in clinical settings. Filtered attenuation corrected SPECT provides the best estimates out of the three quantitative SPECT studies. The activity estimates in the left and right spheres are within 12% for both filtered attenuation corrected SPECT studies. The filtered attenuation corrected SPECT study and the 1:5

conjugate imaging study, where the mis-positioning error was minimized provided comparable results. This suggests that conjugate imaging of the two aligned sources is comparable in accuracy to current clinical capabilities and is a reasonable quantitative technique for estimating activity uptake in two aligned sources.

Chapter 5

SHIELDING CONSIDERATIONS FOR A NOVEL DEDICATED SYSTEM

This chapter presents the design and characterization of the shielding requirements for a novel dedicated mobile nuclear medicine system (Matthews, Mintzer et al. 1997 2000). The system is to perform modified conjugate imaging of ^{123}I and $^{99\text{m}}\text{Tc}$ radiopharmaceuticals in the striata, and comprises two small field-of-view (FOV) gamma cameras mounted on a mobile gantry. The FOV of each gamma camera, determined by the size of its scintillation crystal, is 10 cm x 10 cm. The conjugate imaging system is designed such that its two gamma cameras are positioned on opposite sides of a subject's head. Figure 5.1 is a schematic of the system. Figure 5.2 is a photograph of the system gantry. The design of the gamma cameras, the amount of shielding, and the payload of the gantry are interdependent in the system development. Shielding optimization is an important part of the design and construction of any gamma camera system. A practical balance between the thickness and weight of the shielding is necessary. The shielding must sufficiently reduce the background counts in the scintillation crystal, and the weight of the shielding must be within the limits of the gantry payload. In addition, the shielding design must take into account the mobility of the system and the emission spectra of the radiopharmaceuticals to be used in studies conducted with the novel system.

As discussed in chapter 2, most of the single-photon emitting radiopharmaceuticals used to image the dopaminergic system in the striata are labeled with ^{123}I . ^{123}I radiopharmaceuticals also pool readily in organs in the torso, specifically




Figure 5.1. Schematic of gantry of mobile dedicated conjugate imaging system.




Figure 5.2. Photograph of conjugate imaging gantry.

the lung, liver, and certain organs in the gastrointestinal tract. This generates activity distributions outside the FOV of the camera. This out-of-FOV activity could reduce the accuracy in estimations of activity in the striata in quantitative studies using small FOV gamma cameras due to the high-energy photon emissions of ^{123}I . Therefore, accurate quantification of focal uptake of ^{123}I radiopharmaceuticals in the striata requires that downscatter contribute minimally to image data used for uptake estimation. Reduction of downscatter contribution can occur with sufficient detector shielding and/or data processing that compensates for these events.

In this chapter, we describe two Monte Carlo simulation studies performed to determine and characterize an appropriate shielding configuration for the small FOV gamma cameras comprising the novel system. The first simulation study is performed to determine a shielding design that adequately shields high-energy gamma emissions originating from outside the FOV and is within the weight limits for the mobile conjugate imaging system. The second simulation study is performed to characterize the distribution of events across the scintillation crystal due to high-energy gamma emissions that penetrate the selected shielding design and deposit energy in the 159 keV photopeak window.

5.1 Monte Carlo Simulations

The simulations are performed using Monte Carlo N-Particle Transport code, version 4c (MCNP-4c) (Briesmeister). The simulations are performed in MCNP's photon transport mode. In this mode, only photons are transported and tracked. Electrons are generated but not transported. The thick target bremsstrahlung model (TTB) is used to

manage electron and bremsstrahlung generation and transport. In the TTB model, electrons are generated but are assumed to travel in the same direction as that of the incident photon and deposit their energy locally. Bremsstrahlung photons that are produced by the electrons were transported (Briesmeister 2000). The photon interactions that occurred during transport are modeled in MCNP's detailed physics mode. In detailed physics mode, photoelectric absorption, Compton scattering, coherent scattering and pair production are modeled. Fluorescence photons are transported and photoelectrons are treated with the TTB model for photoelectric absorption. The binding energies of atomic electrons are accounted for in Compton scattering (incoherent scattering) (Briesmeister 2000).

The simulation set-up consists of a model primate and one small FOV camera with the face of its collimator located 1 cm from the primate's head. Figure 5.3 is an illustration of the set-up used in the simulations. The activity placed in the striata represent in-FOV activity, and the activity in GI tract, liver, and lungs represent out-of-FOV activity. Simulations of activity in GI tract, liver and lungs omit the 159 keV gamma emission because its penetration of the shielding is negligible. The primary photon emission and the high-energy gammas are emitted from the appropriate sources in proportions corresponding to the relative yields of the decays of ^{123}I . To achieve acceptable statistical accuracy over most of the spectrum, with reasonable execution times, 1×10^8 source photons and 5×10^8 source photons are simulated and the results are recorded in tally bins of width 5-keV for the striatal sources and the lung and GIT sources, respectively.





Figure 5.3. Configuration of idealized primate and camera used in the simulations. The primate was modeled using water-filled shapes as shown, representing tissue-equivalent material. The activity was uniformly distributed throughout the spherical volumes representing the GI tract, liver and lungs. The camera was positioned with the collimator face located 1 cm from the primate's head.








5.1.1 Primate Model

The system is to be used to image primates as well as humans. Thus, a primate is modeled for the simulations. The primate is modeled with water-filled shapes to represent tissue equivalent material. The primate model is shown in Figure 5.3. The water is uniformly distributed in each shape. The head is a 10-cm diameter sphere; the neck a 5-cm diameter, 3-cm long cylinder; and the torso a 20-cm diameter, 30-cm long cylinder. The striatal spheres have a 2 cm diameter and are symmetric about the vertical bisector of the head sphere. The centers of the spheres are 1 cm superior the center of the head and 3.56 cm apart. The GI tract and the liver are represented by a 10-cm diameter sphere placed in the torso with the center of the sphere 20 cm from the top of the torso. For this study, the GI tract is assumed to include only the large and small intestines; the oral cavity, the pharynx, the esophagus, and the stomach were excluded. The lung is represented as a 10 cm diameter sphere located in the upper half of the torso with the center of the sphere 7 cm from the top of the torso. The modeled primate has no extremities.

5.1.2 Camera Model

The camera simulated consists of a lead collimator, a quartz light guide, photomultiplier tubes (PMTs), and a NaI scintillation crystal all encased in walls of tungsten, lead, and aluminum. The length of the camera is 26.27 cm. Figure 5.4 is a diagram of the simulated camera. The camera is modeled after the small gamma cameras in a novel conjugate - system. Figure 5.5 is a photograph of such a gamma camera. Figure 5.6 is a schematic of a small gamma camera. The gamma cameras have a high resolution hexagonal collimator



	Tungsten
	Lead
	22% Lead (collimator)
	Aluminum
	NaI(Tl)
	SiO₂ (lightguide)
	20% SiO₂ (PMTs)

Configurations Simulated			
	T_W	T_{Pb}	T_{Al}
1	1/4"	1/4"	0
2	3/8"	1/4"	1/8"
3	3/8"	3/8"	0
4	1/2"	3/8"	1/8"
5	3/4"	5/8"	1/8"

Figure 5.4. Detailed diagram of idealized camera components. Thicknesses of tungsten, lead, and aluminum wall components used in each simulation are listed. Lead distributed uniformly at 22% density within the collimator volume modeled a typical low-energy high-resolution collimator. The PMTs were modeled with SiO₂ uniformly distributed at 20% density within a 10 cm x 10 cm x 10 cm volume.

with hole size of 1.9 mm, septa thickness of 0.25 mm, hole length of 40 mm, and sensitivity of 250 cpm/ μ Ci. However, to ensure reasonable simulation execution times, a $4 \times 10 \times 10$ cm 22%-density uniform lead distribution is used as a small gamma camera. The gamma cameras have a high instead of a detailed collimator. The simulated NaI(Tl) scintillation crystal is 0.8 cm \times 10 cm \times 10 cm and positioned 0.2 cm behind the collimator to account for the presence of thin layers of material not modeled due to their low Z numbers. A 1.8-cm quartz light guide is placed 0.2 cm behind the NaI(Tl) crystal. The spacing is used to account for the presence of a thin layer of optical coupling compound between the light guide and crystal. The PMT glass envelopes are represented as a 10 cm \times 10 cm \times 10 cm volume of SiO₂ with a uniform density of 0.6 g/cm³.

Shielding for high-energy photons is composed usually of a highly dense material with a high atomic number (Z). The thickness of the shielding is inversely related to the material's Z number and density. The higher the Z number and density, the less shielding required. For the shielding design of the gamma cameras, tungsten and lead are used. Tungsten is used in the front half of the camera to provide shielding of the crystal. The back half of the camera is shielded with lead. Tungsten walls extended 10 cm along the length of the camera starting at the collimator face. These tungsten walls are intended to provide maximum practical shielding of the NaI(Tl) crystal. The back plate of the camera is tungsten. The remaining lateral shielding is lead with structural walls of aluminum.

Five different camera shielding configurations are simulated; each has different tungsten, lead, and aluminum wall thicknesses. Specifications of the configurations are



Figure 5.5. Small gamma camera for use in a novel dedicated conjugate imaging system. Photomultiplier tubes, scintillation crystal and collimator are shown next to the outer lead shielding.




Figure 5.6. Cross-sectional schematic of small gamma camera for use in the dedicated conjugate imaging system. The collimator is shown in front of the 0.6" scintillation crystal. The crystal is followed by the PMTs, which is followed by the positioning circuitry.

shown in figure 5.4. The thickness of tungsten is increased in increments of $\frac{1}{4}$ " from $\frac{1}{4}$ " to $\frac{3}{4}$ ". Due to lead's lack of hardness, the effect of including aluminum structural walls in the shielding design is investigated. Aluminum walls of $\frac{1}{8}$ " thickness are included for structural support in shielding designs two, four, and five. For shielding designs including aluminum walls, the thickness of lead is dictated by the thickness of the tungsten shielding and the aluminum walls.

5.2 Shielding Determination

To evaluate the count contribution from out-of-FOV activity as a function of shielding design, three simulation studies are performed. One simulation study models and investigates the detection of 159 keV photon emissions from the striatal sources. A second simulation study models and investigates the detection of downscatter from activity in the GI tract and the liver. The third simulation study models and investigates the detection of downscatter from activity in the lungs.

MCNP pulse-height tallies are collected to determine the event rates and the spectra of energy deposited in the NaI(Tl) crystal of the camera. The pulse-height tally results are fractions representing the number of energy-deposited events in each energy bin per source photon. Tallies are collected in 5 keV energy bins and a 20% energy window centered on 159 keV. These tallies are performed with and without MCNP's detector Gaussian energy broadening (GEB). Tallies performed without GEB record the energy as if deposited in an ideal detector. Tallies performed with GEB record the energy as if deposited in a detector with finite energy resolution. The camera simulated has 12%

energy resolution at full width at half maximum (FWHM) at 159 keV. The energy resolution of the simulated camera is a function of $E^{-1/2}$.

5.2.1 Data Analysis

The tallies collected in 5-keV bins are normalized by dividing by the bin size to yield plots of simulated pulse-height spectra. For ease of interpretation, the energy-window tallies of striata, GI tract, and lung activities are converted to tallies with units of cps/mCi. For the GI tract and lung results, the conversions include multiplication by 3.7×10^7 disintegrations/s/mCi and by 0.03 to reflect the fact that the emissions simulated represent only the 3% high-energy gamma emissions of ^{123}I . The tally results from the striata simulations are multiplied by 3.7×10^7 disintegrations/s/mCi and by 83% to represent the yield of the primary photon emission. The cps/mCi of striatal activity and out-of-FOV activity is then multiplied by the corresponding fraction of activity uptake for each radiopharmaceutical to yield the counts per second of activity.

Preliminary data analysis performed on the GI tract tally data assume the activity uptake in the GI tract and striata is 10 mCi and 10 μCi , respectively (Jangha, Mintzer et al. 2001). For each camera configuration, the counts per second recorded in the 159 keV, 20% energy window as a result of the detection of scattered high-energy emissions from a 10mCi GI source is determined through simulation. The counts per second due to high-energy emissions is compared to the counts per second due to striatal activity. The counts per second in the 159 keV 20% window from 10 μCi in the striata is not simulated, but is calculated to be 12.5 counts per second (Jangha, Mintzer et al. 2001) It is determined that

the out-of-FOV activity increase the striata count rate by 3% - 28% for shielding configurations 1-5, respectively.

However, biodistribution studies have shown at the time of imaging, i.e. at biokinetic equilibrium, the total activity uptake in the lungs on average is 10 times greater than that in striata. The combined activity uptake in the liver and GI tract on average is 25 times greater than that in the striata (Seibyl, Wallace et al. 1994). Subsequently, a more accurate simulation is performed using activity uptake values typical for ^{123}I radiopharmaceuticals used in striatal imaging studies and simulating striatal activity uptake and detection. Data analysis is performed on both GI and lung pulse-height tally data using activity uptake values typical of ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM. The activity uptake values are based on the administration of 5 mCi doses to healthy individuals (controls). Table 5.1 shows the percent uptake of the striata, GI tract, and lungs for an injection dose of 5 mCi of ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM, respectively. The last two columns in Table 5.1 are the percent uptakes of initial doses used to represent patients in the early (2:1) and late stages (5:1) of Parkinson's disease. These uptake ratios, which represent appropriate upper and lower bounds of ratios of activity uptake for Parkinson patients, are determined by subtracting the product of the percent dose for controls by the characteristic fractional loss of dopamine from the percent dose for controls. The reported fractional losses of dopamine in early- and late-stage Parkinson's disease patients are 50% and 75%, respectively (Marek, Seibyl et al. 1996; Seibyl 2003). It is assumed that there was a direct linear relationship between the percent loss and activity uptake.

For each camera configuration, the counts per second recorded in the 159 keV,

Table 5.1 Percent Uptake of Initial Dose in Critical Organs

Radiopharmaceutical	Controls			Parkinson's Patients	
				Early Stage	Late Stage
	striata	lungs	GI tract/liver	striata	striata
¹²³ I β-CIT	2.0%	21.0%	52.0%	1.0%	0.5%
¹²³ I FP-CIT	0.8%	9.0%	18.0%	0.4%	0.2%
¹²³ I IBZM	1.7%	8.12%	27.5%	0.8%	0.4%

20% energy window due to the detection of scattered high-energy emissions from GI activity and lung activity is determined. The count rates resultant from the GI activity and lung activity are added to determine the total count rate for out-of-FOV activity. The striatal count rate is determined in the 159 keV, 20% energy window. The percent increase in the striatal count rate is determined for each radiopharmaceutical.

It is assumed that the counts in the 159 keV, 20% energy window collected from striatal activity are concentrated in the physical location of the striatal projection. Since each striatal source has a diameter of 2 cm, the projected lateral area of the spheres is approximately 3.14 cm². Therefore, it is assumed that the projected striatal region can be encompassed by a square 4 cm² ROI drawn on the image data. A 4 cm² ROI is 4 % of the FOV of the gamma cameras. To determine the out-of-FOV contribution to the rate of emissions detected in a 4 cm² region, the total count rate across the face of the crystal due to activity uptakes in the GI tract and lungs are multiplied by this area fraction. The percent increase in striatal counts due to these contributions is determined by dividing the

out of FOV counts by the striatal counts. The percent increase in striatal counts due to only GI tract activity and to only lung activity are also determined.

5.2.2 Results

The normalized pulse-height spectra without GEB pictured in figure 5.7 represent ideal detector resolution. The spectra with GEB pictured in figure 5.8 show the spectra of out-of-FOV shielding penetration events as would be detected by the small FOV gamma camera. The spectra show characteristic tungsten and lead fluorescence x-ray peaks. A broad backscatter peak around 180 keV due to scattered 529 keV gamma emissions can be seen.

Figures 5.9 –5.11 show the percent increases due to out-of-FOV activity for ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM, respectively. Although the percent increase increases as the activity uptake ratio increased (total striatal activity decreased), the magnitude of increase is less than 2% for all configurations and all radiopharmaceuticals.




Figure 5.7. Pulse height spectrums as detected by ideal gamma camera. Pulse height spectrums from ^{123}I high-energy photons emitted from the GI tract and liver as seen by an ideal gamma camera (gamma camera with perfect energy resolution). K-shell x-ray peaks for tungsten (~ 70 keV) and lead (~ 88 keV) are shown. The 173 keV backscatter peak due to the 529 keV photon is also clearly shown.




Figure 5.8. Pulse height spectra as detected by a gamma camera with finite resolution. Pulse height spectra from ^{123}I high-energy photons emitted from the GI tract and liver as seen by a gamma camera with energy resolution of 12% at 159 keV. The iodine and tungsten fluorescence escape peaks that are clearly seen in figure 10.9 are not as clearly depicted. The backscatter peak for the 529 keV photon is also shown but not as clearly delineated as in figure 8.9. The soft peaks are due to finite energy resolution.




Figure 5.9. Comparison of percentage increase in striatal counts for ^{123}I β -CIT activity uptakes in the striata of controls and early, stage, and late-stage Parkinson's disease patients. The controls are considered to have normal striatal activity. It can be seen that as activity uptake in the striata decreases the significance of out-of-FOV count contribution increases. The percent increase is across the entire scintillation crystal.




Figure 5.10. Comparison of percentage increase in striatal counts for ^{123}I FP-CIT activity uptake in the striata of controls and early, stage, and late-stage Parkinson's disease patients. The controls are considered to have normal striatal activity. It can be seen that as activity uptake in the striata decreases the significance of out-of-FOV count contribution increases. The percent increase is across the entire scintillation crystal.




Figure 5.11. Comparison of percentage increase in striatal counts for ^{123}I IBZM activity uptake in the striata of controls and early, stage, and late-stage Parkinson's disease patients. The controls are considered to have normal striatal activity. It can be seen that as activity uptake in the striata decreases the significance of out-of-FOV count contribution increases. ^{123}I IBZM has the lowest percentage increase due to its low activity uptake in the lungs. The percent increase is across the entire scintillation crystal.

5.2.3 Discussion

In order to perform accurate conjugate imaging of two aligned sources, the activity outside the FOV must be appropriately shielded and/or appropriately accounted for. Five shielding configurations of thicknesses from $\frac{1}{4}$ to $\frac{3}{4}$ inch are investigated using MCNP-4c simulations. It is determined that configuration 4 with $\frac{1}{2}$ inch tungsten walls is a sufficient shielding design for the mobile conjugate imaging system. Although configuration 5 with $\frac{3}{4}$ inch tungsten walls shields better than configuration 4, configuration 4 provides comparable shielding to configuration 5 with about $\frac{3}{4}$ the weight.

The magnitude of the percent increase in striatal counts is dependent on the radiopharmaceutical used and the severity of the disease. The percent increase for patients in the late stages of Parkinson's disease is about 3-4 times the percent increase of controls. The differences in percent increases among the radiopharmaceuticals are a function of the ratio of uptake in the striata and out-of-FOV activity, specifically the lungs. Although the activity uptake in the GI tract is significantly larger than the activity uptake in the lungs for all three radiopharmaceuticals, the activity uptake in the lungs contributes more counts to the striatal ROI than the activity in the GI tract. Since the lungs are closer to the gamma cameras, photon emissions have a shorter distance to travel to reach the gamma cameras. For example, for ^{123}I IBZM, the activity uptake in the GI tract is much larger than the uptake in the lungs. As a result, the percent increase in the striatal counts due to ^{123}I IBZM is not as large as the percent increase due to ^{123}I β -CIT, which has a higher activity uptake in the lungs. However, the percent increase associated with ^{123}I FP-CIT is due to the relatively low activity uptake in the striata. Thus, an increase in striatal counts due to out-of-FOV activity results in a larger percent increase.

These studies assume that penetration events from the GI tract and lungs are uniformly distributed across the scintillation crystal. Under this assumption, background ROI subtraction used in the conjugate calculation of the striata activity uptake eliminate the penetration events due to GI tract and lung activity for all configurations simulated, with negligible additional statistical error. However, further analysis of the results reveals that the counts from the activity uptake in the GI tract and lungs are not uniformly distributed across the scintillation crystal. It is determined that there is approximately a 4:1 variation in the rates of penetration events across the crystal. Simulations are

performed to investigate the significance of non-uniform count distribution across the crystal on background subtraction and on the limits of quantitative accuracy associated with the ½ inch shielding configuration.

5.3 Shielding Characterization

Monte Carlo simulations are performed to determine the spatial distribution of shielding penetration events across the scintillation crystal for the shielding design with ½ inch tungsten walls. The code parameters, such as physics and mode, the simulation set-up, and the number and energy of photons as described in section 5.1 are used with the exception of collimator design. The collimator used in this study as well as in subsequent studies is hexagonal, with hole size 1.8 mm, septa thickness of 0.25 mm, hole length of 40 mm, and sensitivity of 250 cpm/ μ Ci. In addition, MCNP's ptrac feature is used to track source photons, scattered photons, and fluorescent photons that interacted in the scintillation crystal and to output the location of interaction events and energy deposited per event.

5.3.1 Data Analysis

Images of out-of-FOV activity and in-FOV striatal activity are generated from MCNP-4c ptrac outputs. The MCNP ptrac feature is used to track source photons, scattered photons, and fluorescent photons that interact in the scintillation crystal and to output the event location and energy deposited at each event. A code is written in Matlab to parse the ptrac files and to generate energy-windowed images. The Matlab code provides ROI selection capabilities and event filtering options based on interaction

location with respect to the gamma camera. The energy windows used are a 20% window centered on 159 keV and two 20% windows one centered on 129 keV, the other on 195 keV. The image size is 64 x 64 pixels, where the pixel size is 0.15625 cm.

The simulated images are generated based on energy-weighted centroiding, a concept used in the generation of images from an ideal Anger camera. Energy-weighted centroiding is analogous to calculation of the center of mass. Photons may interact and deposit energy at more than one location in the scintillation crystal. Energy-weighted centroiding calculates the location at which the total energy deposited by a photon can be considered to be concentrated. This location is called the energy-weighted centroid. For this application the energy-weighted centroid is expressed in x and y coordinates. The energy-weighted centroid is expressed as

$$x_E = \frac{\sum_i^N x_i E_i}{\sum_i^N E_i} \quad (5.1)$$

$$y_E = \frac{\sum_i^N y_i E_i}{\sum_i^N E_i} \quad (5.2)$$

where x_E and y_E are the x and y coordinates of the energy-weighted centroid, respectively; index i represents the i th energy-depositing interaction; N is the total number of energy-depositing interactions; x_i is the x coordinate for the i th interaction; y_i is the y coordinate for the i th interaction; and E_i is the energy deposited by the i th interaction. If the total energy deposited by a photon is within the 159 keV, 20% energy-

windowed image, a count is tallied in the pixel that corresponds to the location of the energy-weighted centroid. By using the energy-weighted centroid, cross-detection of an event by multiple photomultiplier tubes is ignored.

To determine a more accurate percent increase in the striatal counts due to out-of-FOV activity in the striatal region, a ROI with an area of 4.12 cm^2 (13 by 13 pixels) is placed around the pixels that correspond to the physical location of the striata on the striatal, GI tract, and lung images. The counts in the pixels within the ROI are summed and converted into cps/mCi. For the GI tract and lung images, in which only high-energy photons are simulated, the total ROI counts are multiplied by 3.7×10^7 disintegrations/s/mCi and by 0.03 to reflect the fact that the emissions simulated represent only 3% high-energy emissions of ^{123}I and are divided by the total number of particles simulated to yield cps/mCi. The cps/mCi in each pixel within the ROI in GI tract image are added to the cps/mCi in the corresponding pixel in the lungs image to yield a total out-of-FOV activity count rate per mCi. For the striatal images, the total ROI counts are multiplied by the 3.7×10^7 disintegrations/s/mCi; divided by the total number of simulated photons; and by 0.83 to reflect the fact that the emissions simulated represent 83% of ^{123}I gamma emissions. The cps/mCi in each pixel within the ROI in one striatal image are added to the cps/mCi in the corresponding pixel in the second striatal image to yield a total striatal count rate per mCi.

The calculated cps/mCi in each pixel in the striata, GI tract and lung images are multiplied by the percentage of 5 mCi of radiopharmaceutical taken up by the striata, GI tract and lung, respectively. The biodistributions of ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM are used. See table 5.1 for the percent uptakes in the striata, GI tract and lungs for

^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM. Then the count rates in each pixel in the striata, GI tract, and lung images are multiplied by 900 seconds to represent typical 15 minute planar images of activity uptake in the striata, GI tract and lungs. Assuming Poisson statistics, the error in the calculated counts was determined by taking the square root of the counts. The percent increase was determined for each radiopharmaceutical by dividing the total out-of-FOV activity counts by the total striatal counts.

It is noted that about 12% of counts due to striatal activity is displaced outside the 4 cm² (13 pixels x 13 pixels) ROI. As a result, the percent increase is determined for larger ROIs because the 4 cm² ROI does not account for partial volumes effects. A 9 cm² (19 pixels x 19 pixels) ROI and a 16 cm² (26 pixels x 26 pixels) ROI are also drawn around the striatal region to include displaced counts. To investigate the spatial distribution of counts due to out-of-FOV activity in non-striatal regions, background ROIs are drawn on the GI tract and lung count maps. However, the small FOV of the gamma camera limits the size of the background ROI. Thus, one background ROI is used in the investigations. The background ROI is about 3.125 cm by 2.03 cm (20 pixels x 13 pixels). The background ROI is placed below the striatal ROI; the distance between the ROIs depends on the size of the striatal ROI. For the 4 cm² striatal ROI, the background ROI is placed about 2.66 cm (17 pixels) from the bottom of the striatal ROI. For the 9 cm² striatal ROI, the background ROI is placed 2.34 cm (12 pixels) from the bottom of the striatal ROI. For the 16 cm² striatal ROI, the background ROI is placed about 1.25 cm (8 pixels) from the bottom of the striatal ROI. The counts in each ROI are assumed to be uniformly distributed. The counts within the striatal ROI and background ROI due to GI tract and lung activity are compared. The counts in the striatal ROI and background ROI

for the GI tract and lung count maps are compared to the counts in the striatal ROI and background ROI of the striatal count maps.

5.3.2 Results

Figure 5.12 is an image of the count distribution across the scintillation crystal from activity in the GI sphere for $\frac{1}{2}$ inch of shielding of tungsten, configuration 4. The white boxes represent the ROIs. The top ROI is the striatal ROI and the bottom ROI is the background ROI. There are more counts in the bottom half of the crystal. The counts in the background ROI due to activity in the GI tract are about 1.3 times larger than the counts in the striatal ROI due to activity in the GI tract. Figure 5.13 is an image of the count distribution across the scintillation crystal from activity in the lungs for $\frac{1}{2}$ inch of shielding of tungsten, configuration 4. There are more counts in the upper half of the scintillation crystal. The counts in the background ROI are about 40% of the counts in the striatal ROI.

Figure 5.14 shows the percent increases in counts in a 4 cm^2 striatal ROI due to the out-of-FOV activity uptake from 5 mCi (185 MBq) of ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM to control groups, patients in the early stages of Parkinson's disease, and patients in the late stages of Parkinson's disease. The percent increase due to out-of-FOV activity increases with the severity of Parkinson's disease for all three radiopharmaceuticals. The percent increases due to ^{123}I β -CIT are determined to be 0.44%, 0.87%, and 1.75% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. The percent increases due to ^{123}I FP-CIT are determined to be 0.45%, 0.90%, and 1.80% for controls, early-stage Parkinson's patients,




Figure 5.12. Count distribution across the face of the scintillation crystal from activity uptake in the GI tract and liver for a 159 keV energy-windowed image. The square is the striatal ROI. The lower rectangle is the background ROI. There are more counts in the bottom half of the image.




Figure 5.13. Count distribution across the face of the scintillation crystal from activity in the lungs for a 159 keV windowed image. The square is the striatal ROI. The lower rectangle is the background ROI. There are more counts in the upper half of the image. The striatal ROI contains approximately 2.5 times more counts than the background ROI.




Figure 5.14. Comparison of percentage increase in striatal counts in a 4 cm² ROI due to activity uptake of ¹²³I β-CIT, ¹²³I FP-CIT, and ¹²³I IBZM. Activity uptakes associated with control subjects, early and late stage Parkinson's disease patients are compared. It is shown that the percent increases in striatal counts due to out-of-FOV activity from ¹²³I radiopharmaceuticals increase as activity uptakes in the striata decrease.

and late-stage Parkinson's patients, respectively. The percent increases due to ¹²³I IBZM are determined to be 0.21%, 0.43%, and 0.86% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. However, due to partial volume effects and blurring, striatal counts are spread outside the 4 cm² area. As a result, larger ROIs are drawn around the striatal regions to include counts that are displaced. Figure 5.15 shows the percent increases in striatal counts in a 9 cm² ROI. The percent increase due to ¹²³I β-CIT are determined to be 0.77%, 1.54%, and 3.08% for controls, early-stage Parkinson's patients, and late-stage patients, respectively. The percent increases due to




Figure 5.15. Comparison of percent increases in striatal counts for a 9 cm² ROI due to activity uptakes of ¹²³I β-CIT, ¹²³I FP-CIT, and ¹²³I IBZM. Activity uptakes associated with control subjects, early and late stage Parkinson's disease patients are compared. The percent increase associated with the 9 cm² ROI is larger than that of the 4 cm² ROI.

¹²³I FP-CIT in the 9 cm² are determined to be 0.80%, 1.60%, and 3.19% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. The percent increases due to ¹²³I IBZM are determined to be 0.38%, 0.75%, and 1.50% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. Figure 5.16 shows the percent increases in striatal counts in a 16 cm² ROI. The percent increases due to ¹²³I β-CIT are determined to be 1.28%, 2.57%, and 5.13% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. The percent increases due to ¹²³I FP-CIT are determined to be 1.33%,




Figure 5.16. Comparison of percent increases in striatal counts for a 16 cm² ROI due to activity uptakes of ¹²³I β-CIT, ¹²³I FP-CIT, and ¹²³I IBZM. Activity uptakes associated with control subjects, early and late stage Parkinson's disease patients are compared. As seen with the smaller ROIs, the percent increase in striatal counts increases with decrease in activity uptake. The percent increases associated with the 16 cm² ROI is larger than that of 4 cm² and 9 cm² ROIs.

2.66%, and 5.32% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively. The percent increases due to ¹²³I IBZM are determined to be 0.63%, 1.25%, and 2.50% for controls, early-stage Parkinson's patients, and late-stage Parkinson's patients, respectively.

5.3.3 Discussion

In the context of conjugate imaging, subtraction of counts in a background region from counts in the striatal region is important. Thus, characterization of the spatial distribution of counts is important, specifically in the striatal and background ROIs. Preliminary simulations revealed that the spatial distribution of counts from out-of-FOV activity uptake in the GI tract was not uniformly distributed across the scintillation crystal. This prompted a second Monte Carlo simulation study. The second study is conducted to characterize the distribution of out-of-FOV counts across the scintillation crystal. This study is conducted using the shielding design with $\frac{1}{2}$ inch tungsten walls, as it was determined that this shielding design was the model of choice. In this second study, the percent increase in striatal counts is determined by using counts in a ROI placed around the projected striatal area. The percent increase determined using the counts encompassed in a ROI is larger than the percent increase determined under the assumption that out-of-FOV activity uniformly distributes counts across the scintillation crystal. The assumption of a uniform count distribution generates an underestimation of percent increase. The use of ROIs provides a more accurate estimation of percent increase in striatal counts due to out-of-FOV activity.

The simulation studies performed show that a large percentage of gamma photons leave the lungs at angles that allow them to travel through the collimator bypassing the $\frac{1}{2}$ inch thick tungsten shielding along the bottom of the camera and deposit energy in the upper portion of the scintillation crystal. The marked concentration of counts in the lower proportion of the scintillation crystal from activity uptake in the GI tract is due to the angle at which the photons escaped from the GI tract and impinge on the bottom of the

gamma camera. Although there is a significant difference in counts in the striatal ROI and the background ROI due to activity uptakes in the GI tract and lungs, the overall contributions of counts in the striatal ROI and background ROI due to out-of-FOV activity is negligible for ^{123}I β -CIT, ^{123}I FP-CIT, and ^{123}I IBZM. However, the influence of the high-energy penetration counts increases as the striatal activity uptake decreases. The magnitude of this influence is also dependent on the ROI size; the larger the ROI the larger the influence.

Chapter 6

EVALUATION AND CHARACTERIZATION OF THE MODIFIED CONJUGATE IMAGING TECHNIQUE

Monte Carlo simulations are performed to evaluate and characterize the modified conjugate imaging technique in the context of the novel dedicated imaging system. The biodistributions of radiopharmaceuticals are simulated and the photon emissions from the radiopharmaceuticals are detected using simulated gamma cameras. Estimates of the activity uptake in each of the striatum are determined using the modified conjugate imaging technique. To assess the range of the technique three activity uptake ratios, which represented three stages of Parkinsonism, are simulated.

The simulated conjugate imaging is similar to practical conjugate imaging in that both require data acquisition and image array generation before conjugate imaging calculations can be performed. However, data acquisition for simulated conjugate imaging is obtained through Monte Carlo simulations. The generation of image data is an explicit process, which included converting the results from the Monte Carlo simulations into image data arrays. These data arrays are then used in the conjugate imaging calculations.

6.1 Monte Carlo Simulations

6.1.1 Data Acquisition

Simulation studies are performed using photon emissions of ^{123}I and $^{99\text{m}}\text{Tc}$ using the primate/camera model discussed in chapter 5. To represent activity uptake of ^{123}I and $^{99\text{m}}\text{Tc}$ radiopharmaceuticals, activity is uniformly distributed throughout the striata, head, lungs, and GI tract. The left striatum, right striatum, head, lungs, and GI tract simulations are performed separately. For ^{123}I studies, two simulations are performed, one for 159 keV primary emissions and one for high-energy emissions ranging from 248 to 783 keV. The 159 keV and high-energy photons are emitted from the sources in proportion to their relative yields from ^{123}I . A simulation study using 140 keV photon emissions is performed for the $^{99\text{m}}\text{Tc}$ study. The 140 keV photons are emitted from the sources in proportion to their relative yields from $^{99\text{m}}\text{Tc}$. The parameters used in the simulations are described in section 5.2.1.

To represent acquired data from the second gamma camera Monte Carlo simulations are performed with the same geometry, but with a different random number sequence. This is done to ensure that the simulated photons used in generating data for the second camera does not undergo the same interactions as the photon used in generating the image data for the first gamma camera.

6.1.2 Image Generation

Images are generated from the ptrac output produced by the MCNP simulations. This process is described in section 5.3.1. In chapter 5, the images are only of scattered high-energy photon emissions from ^{123}I emitted from activity distributions in the lungs

and gastrointestinal tract. The images in this section are generated from the primary photon emissions and high-energy photon emissions, if applicable, from ^{123}I and $^{99\text{m}}\text{Tc}$ radiopharmaceuticals in the striata, brain, and lungs and gastrointestinal tract. The biodistributions of these radiopharmaceuticals used in chapter 5 are also used for these studies. Table 6.1 lists the biodistributions, which include the activity uptake in brain. The activity uptake in the striata in table 6.1 is based on the activity uptake of a healthy individual. The activity in the striata is varied to simulate the various stages of Parkinson's disease. As the severity of the disease increases, the activity uptake in the striata decreases. In this study, three striatal activity uptake ratios, 1:1, 2:1, and 5:1 are investigated. The absolute activity in the right sphere is decreased for the 2:1 and 5:1 ratios.

Table 6.1 Percent Uptake of Radiopharmaceuticals in Brain, Striata, and Critical Organs

Radiopharmaceutical	Striata	Brain	Lungs	Gastrointestinal Tract
^{123}I β -CIT	2	10	21	52
^{123}I FP-CIT	0.8	2.2	9	18
^{123}I IBZM	1.7	2	8.1	27.5
$^{99\text{m}}\text{Tc}$ TRODAT	0.63	0.12	-	-

- Percent uptake in that organ or system not used in simulation study.

^{123}I radiopharmaceutical biodistribution based on 5mCi of administered activity

$^{99\text{m}}\text{Tc}$ TRODAT biodistribution based on 20mCi of administered activity

Since the amount of activity that accumulates in the brain, the striata, and the critical organs differs, the pixel values in the image arrays for the head, the striata, and the out-of-FOV organs are scaled by the appropriate uptake fraction of administered activity. Once the primary and high-energy image arrays for each organ are scaled, the corresponding primary and high-energy image arrays are added together to generate data that reflect an emission image from an activity distribution in each organ. Then all the emission images from the organs are added together to simulate a complete emission image from all activity uptakes seen in actual patients. For the ^{123}I studies, the total injected dose is assumed to be 5 mCi. For the $^{99\text{m}}\text{Tc}$ study the total injected dose is assumed to be 20 mCi. The ^{123}I total dose is based on the recommendation of the European Association of Nuclear Medicine Procedure (Tatsch, Asenbaum et al. 2002). The $^{99\text{m}}\text{Tc}$ total dose is based on research studies (Meegalla, Plossl et al. 1997; Fang, Wu et al. 2000). The pixel values in the energy-windowed images are converted into counts per second per milliCurie (cps/mCi). A total imaging time of 15 minutes (900 seconds) is used to represent typical planar imaging protocols.

The photopeak window for the ^{123}I studies is a 20% window centered on 159 keV. The photopeak window for the $^{99\text{m}}\text{Tc}$ studies is a 20% window centered on 140 keV. Thus, the photopeak window used for the ^{123}I studies is 142 keV-175 keV and that used for the $^{99\text{m}}\text{Tc}$ study was 126 keV-154 keV. The image size is 64 pixels x 64 pixels. The image data arrays computed from the MCNP data simulated with a different random number sequence are flipped from left to right to represent the mirrored view as would be seen by the second camera. Figures 6.1-6.5 show the images for the emissions from a calibration source, activity in the brain, both striatum, and the lungs and GI tract as seen





Figure 6.1. Images of a calibration source in the 159 keV 20% energy window as seen by the left gamma camera. These images are from the ^{123}I β -CIT study. The left image is generated from the 159 keV emissions only. The right image is generated from the high-energy emissions only.





Figure 6.2. Images of brain activity in the 159 keV 20% energy window as seen as by the left gamma camera. These images are from the ^{123}I β -CIT study. The left image is generated from the 159 keV emissions only. The right image is generated from the high-energy emissions only.




Figure 6.3. Images of activity in the right striatum in the 159 keV 20% energy window as seen by the left gamma camera. These images are from the ^{123}I β -CIT study. The left image is generated from the 159 keV emissions only. The right image is generated from the high-energy emissions only.




Figure 6.4. Images of activity in the left striatum in the 159 keV 20% energy window as seen by the left gamma camera. These images are from the ^{123}I β -CIT study. The left image is generated from the 159 keV emissions only. The right image is generated from the high-energy emissions only.




Figure 6.5. Images of activity in the GI tract (left) and lungs (right) in the 159 keV 20% energy window as seen by the left gamma camera. These images are from the ^{123}I β -CIT study. These images are generated from the detection of high-energy emissions only.




Figure 6.6. ^{123}I β -CIT images 1:1 simulated data. These images are before conjugate imaging calculation. Image on the left is from the left simulated gamma camera. Image on the right is from the right simulated gamma camera.

by the left gamma camera for the case of ^{123}I β -CIT. Figure 6.6 is the sum of the calibration, brain activity uptake, striatal activity uptake, lungs activity uptake, and GI tract activity uptake image data. This represents a complete emission conjugate imaging images for the study.

6.1.3 Conjugate Imaging

The theoretical development of the conjugate imaging calculations for two aligned sources is discussed in chapter 3. In summary, conjugate imaging requires three data sets, emission, transmission, and flood data sets. The transmission data are used to determine the effective attenuation coefficient. The flood data are used to determine the camera sensitivities. The emission data are used to estimate the activity in the left and right striata above background using the estimated attenuation coefficient and the estimated camera sensitivities. A relative striatal uptake study is also performed to compare the accuracy of conjugate imaging estimations of striatal uptake ratios.

The conjugate imaging calculations performed in these studies require two regions of interest (ROIs), a striatal ROI and a background ROI. The conjugate imaging calculations are performed pixel by pixel within the striatal ROI. The background ROI is used to determine the background activity, the activity distributed in the brain. The activity in the background ROI is scaled by the ratio of average brain thickness in the striatal ROI and background ROI. The estimated background activity is subtracted from the estimated activity in the striatal ROI. The dimensions of these ROIs are determined in a ROI study. The ROI study is discussed in the following section.

Several deviations are made in order to simulate acceptable statistically accurate data within reasonable execution times. An effective attenuation coefficient, determined from transmission data, is not used in the conjugate imaging calculations. Instead, narrow-beam and broad-beam attenuation coefficients from published data are used. Use of a uniform attenuation coefficient is considered to be reasonable as Monte Carlo simulations have shown that transmission-based non-uniform attenuation correction is not crucial for brain studies (Van_Laere, Koole et al. 2000; Van_Laere, Koole et al. 2001; Kim, Varrone et al. 2003). A constant attenuation coefficient used in quantitative brain studies produce reasonably accurate results (Tsui 1996; Van_Laere, Koole et al. 2000; Van_Laere, Koole et al. 2001). The narrow-beam attenuation coefficients used are interpolated from the published data from the National Institute of Standards and Technology (NIST). The broad-beam attenuation coefficients used are from published striatal studies. For the ^{123}I studies the narrow-beam attenuation coefficient and the broad-beam attenuation coefficient are 0.1481 cm^{-1} and 12 cm^{-1} , respectively (NIST 2004 and Soret et al. 2003). The broad-beam attenuation coefficient used for ^{123}I is determined using a linear relationship between the narrow-beam and broad-beam attenuation coefficients of ^{123}I and $^{99\text{m}}\text{Tc}$ (Soret et al. 2003). The narrow beam attenuation and broad-beam attenuation coefficient used in the $^{99\text{m}}\text{Tc}$ studies are 0.1526 cm^{-1} and 0.12 cm^{-1} , respectively (NIST 2004)(Harris, Greer et al. 1984; Zaidi and Montandom 2002). The narrow-beam attenuation coefficient is used to represent a situation in which scatter photon detection is not taken into account. The broad-beam attenuation coefficient represented a situation in which scatter photon detection is taken into account.

A point source is used instead of a flood source to determine the sensitivity of the gamma cameras. The counts detected across the entire scintillation crystal are used to determine camera sensitivity. The simulated camera sensitivity data are generated into image data array using the same technique for organ-simulated data, excluding scaling the data for fractional activity uptake. Figure 6.1 shows an image of the calibration source emission. For the ^{123}I sensitivity studies, a total of 200 million photons are simulated from the point source, 100 million 159 keV photons and 100 million high-energy photons. A total of 100 million 140 keV photons are simulated for the $^{99\text{m}}\text{Tc}$ sensitivity studies. The values in the data arrays are converted into counts/decay and then summed over the entire face of the camera to yield the camera sensitivities.

Due to lack of counts, the background activity is averaged over the background ROI. After background subtraction, negative pixel values occur within the striatal ROI. To handle the negative pixel values a short iterative algorithm is used to sum all the negative pixel values. The sum of the negative pixel values is divided by the number of pixels with positive values. This value is subtracted from all positive value pixels and the negative values are replaced with zero. This process is repeated until no negative values were within the striatal ROI.

6.1.4 ROI Study

Although a crude ROI study is performed and described in chapter 5 for the shielding characterization, a MCNP study is performed to determine the best striatal and background ROI dimensions. Choosing an appropriate ROI size is of concern because counts from the striatal activity uptake are displaced due to blurring. An appropriate ROI

will include most of the counts from both striatum. The percent of counts contained in the ROIs depends on the proximity of the striatum to the gamma camera. Counts from the striatum nearest the gamma camera are consolidated to a smaller region than the counts from the far striatum.

The simulation set-up consists of a gamma camera and a sphere representing a striatum. The gamma camera has the same specifications described in the camera model section 5.1.2. The striatal sphere has a 2 cm diameter. To represent the left striatum, the near striatum, the striatal sphere is placed 4.22 cm from the gamma camera and 1 cm superior the gamma camera horizontal midline. To represent the right striatum, the far striatum, the striatal sphere is placed 6.0 cm from the gamma camera and 1 cm superior the gamma camera horizontal midline. For ^{123}I , 159 keV photons are emitted from the left and right striatum separately and detected by the both gamma cameras. For $^{99\text{m}}\text{Tc}$ 140 keV photons are similarly emitted and detected. The total counts in various ROIs are compared to the total image counts. Square ROIs are used. The ROI sizes used in this study are 13 pixels x 13 pixels, 15 pixels x 15 pixels, 17 pixels x 17 pixels, 19 pixels x 19 pixels, 26 pixels x 26 pixels, and 32 pixels x 32 pixels. The dimensions of a pixel are 0.15625 cm x 0.15625 cm. The 13 pixels x 13 pixels ROI is 2.03 cm x 2.03 cm; the side of the ROI is approximately the diameter of the striatal sphere. The ROI with 13 pixels x 13 pixels contain 89% of the total counts for the near striatum study and approximately 77% of the total counts for the far striatum study. Approximately 96% of the total counts from the near striatum study are within the ROI with 19 pixels x 19 pixels. Approximately 86% of the total counts from the far striatum study are within the 19 pixels x 19 pixels ROI. Although the percentage of total counts increases as the ROI increases, the 19 pixels

x19 pixels ROI is considered to be the best ROI size due to the small FOV of the gamma cameras. Larger ROIs increase cross talk between the striatum ROI and the background ROI, which is located below the striatal ROI. It should be noted that these count percents are determined in the absence of head attenuation and scatter.

The location of the background ROI is limited by the size of the small FOV gamma camera. As a result, background ROIs cannot be drawn too far from the striatal ROI. However, background ROIs too close to the striatal region might include counts from the activity uptake in the striata resulting in a possible overestimation of background activity. This overestimation of background activity reduces the amount of estimated activity in the striata. To determine the location of the background ROI, a background ROI is placed at different distances from the striatal ROI and the amount of background activity is estimated from the counts in the ROI. The counts in the ROI placed 4 pixels (0.625 cm) from the bottom edge of the gamma camera produce the background activity estimate closest to the known background activity. The background activity estimation using background ROIs ranging from 13 pixels x 13 pixels to 32 pixels x 32 pixels does not vary significantly. As result the size of the background ROI is kept at 13 pixels x 13 pixels ROI.

6.2 Results

6.2.1 Absolute Activity Estimation

6.2.1.1 13 pixels x 13 pixels ROI

Although it is determined in the ROI study that the 13 pixels x 13 pixels ROI does not contain adequate counts, conjugate imaging studies are performed using a 13 pixels x

13 pixels ROI for comparative purposes. Conjugate imaging studies for all ROI sizes are in Appendix A. Two conjugate imaging studies are performed for each radiopharmaceutical. The first study uses a narrow-beam attenuation coefficient and the second uses a broad-beam attenuation coefficient. Table 6.2 lists the percent error for the estimates of background activity per cm^2 for each activity uptake study. The percent error for the estimates of background activity, the brain activity uptake, for ^{123}I β -CIT study increases as the activity uptake increases. The percent errors for background activity estimates for the other radiopharmaceuticals decreases as the activity uptake increases. Table 6.3 lists the percent errors for the activity estimates of the left and right striata for the 1:1, 2:1, and 5:1 uptake ratios using the 13 pixels x 13 pixels ROI and the narrow-beam attenuation coefficients. For the ^{123}I radiopharmaceutical studies, the percent errors for right and left striata activity estimates for the 1:1 uptake ratio are relatively equal. However, as the activity uptake ratio increases, the difference between the percent errors increases. The percent error in the activity estimate for the left striatum remains relatively stable as the percent error associated with the activity estimate of the right striatum increases significantly. For the $^{99\text{m}}\text{Tc}$ TRODAT study, the difference between the percent errors also increases. While the percent error for the left striatum increases, the percent error for the right striatum decreases. The magnitude and direction of the percent errors for the left and right striatum in the $^{99\text{m}}\text{Tc}$ TRODAT study are considerably smaller than those for the ^{123}I radiopharmaceutical studies. The difference between the percent errors for the left and right striata is also smaller than the difference for the ^{123}I radiopharmaceutical studies. Table 6.4 lists the percent error in the activity estimation of the left and right striata for the 1:1, 2:1, and 5:1 uptake ratio using the 13 pixels x 13

Table 6.2. Percent Errors for Background Activity Estimates for Activity Uptakes

Background Activity Estimates			
ROI: 13 pixels x 13 pixels	Percent Error		
Radiopharmaceuticals	1:1	2:1	5:1
I-123 beta-CIT	-0.59%	-2.37%	-3.44%
I-123 FP-CIT	5.99%	2.82%	0.90%
I-123 IBZM	21.63%	14.62%	10.34%
Tc-99m TRODAT	24.41%	15.47%	8.73%

Table 6.3. Percent Errors for Activity Estimates for Narrow-beam Attenuation Study

Striatal Uptake Ratios						
ROI: 13 pixels x 13 pixels	Percent Error					
	1:1		2:1		5:1	
Radiopharmaceuticals	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
I-123 beta-CIT	-14.73%	-16.41%	-14.47%	-21.61%	-14.31%	-37.22%
I-123 FP-CIT	-13.32%	-14.07%	-13.06%	-16.94%	-12.91%	-25.53%
I-123 IBZM	-12.18%	-12.28%	-11.94%	-13.59%	-11.79%	-17.52%
Tc-99m TRODAT	6.23%	6.36%	6.61%	5.97%	6.85%	4.88%

*¹²³I Attenuation coefficient: 0.1481 cm⁻¹; ^{99m}Tc Attenuation coefficient: 0.1526 cm⁻¹

Table 6.4. Percent Errors for Activity Estimates for Broad-beam Attenuation Study

Striatal Uptake Ratios						
ROI: 13 pixels x 13 pixels	Percent Error					
	1:1		2:1		5:1	
Radiopharmaceuticals	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
I-123 beta-CIT *	-24.56%	-26.38%	-19.90%	-39.99%	-17.10%	-80.84%
I-123 FP-CIT*	-23.40%	-24.21%	-18.75%	-35.68%	-15.96%	-70.04%
I-123 IBZM*	-22.46%	-22.57%	-17.81%	-32.61%	-15.02%	-62.72%
Tc-99m TRODAT [^]	-8.07%	-7.92%	-1.69%	-20.33%	2.14%	-57.49%

*¹²³I Attenuation coefficient: 0.12cm⁻¹; ^{99m}Tc Attenuation coefficient: 0.12 cm⁻¹

pixels ROI and the broad-beam attenuation coefficients. The percent errors in this study are larger than the corresponding percent errors in the narrow-beam attenuation coefficient study. The activity estimates are underestimated for all radiopharmaceuticals. Similar to the narrow-beam attenuation coefficient study, the differences between the percent errors for the left and right striatal activity estimates increase as the activity uptakes increase. However, the magnitudes of the differences are considerably larger than the corresponding magnitudes in the broad-beam attenuation coefficient study. The percent errors in the activity estimates for the left striatum decrease as the percent errors associated with the activity estimates of the right striatum increase significantly for all radiopharmaceuticals.

6.2.1.2 19 pixels x 19 pixels ROI

The 19 pixels x 19 pixels ROI is chosen as the best ROI size because of the size of the conjugate imaging system FOV and adequate counts from both striata are contained in this ROI. The estimates of the background activity per cm^2 for the activity uptake studies are listed in Table 6.5. Similar to the 13 pixels x 13 pixels ROI study, the percent errors associated with the background activity estimates for the ^{123}I β -CIT study increase as the activity uptake ratios increase. The magnitudes of the percent errors are larger for the 19 pixels x 19 pixels ROI. The percent errors in the other radiopharmaceutical studies are smaller than the percent errors in the 13 pixels x 13 pixels ROI study and decrease as the activity uptake ratios increase. Table 6.6 lists the percent error for the estimates of activities in the left and right striata for the 1:1, 2:1, and 5:1 activity uptake ratios using the 19 pixels x 19 pixels ROI and the narrow-beam

Table 6.5. Percent Errors for Background Activity Estimates for Activity Uptakes

Background Activity Estimates			
ROI: 19 pixels x 19 pixels	Percent Error		
Radiopharmaceuticals	1:1	2:1	5:1
I-123 beta-CIT	-2.33%	-4.08%	-5.14%
I-123 FP-CIT	4.13%	1.02%	-0.88%
I-123 IBZM	19.49%	12.61%	8.40%
Tc-99m TRODAT	22.23%	13.44%	6.82%

Table 6.6. Percent Errors for Activity Estimates for Narrow-beam Attenuation Study

Striatal Uptake Ratios						
ROI: 19 pixels x 19 pixels	Percent Error					
	1:1		2:1		5:1	
Radiopharmaceuticals	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
I-123 beta-CIT	-9.23%	-3.43%	-10.14%	-1.35%	-10.69%	5.09%
I-123 FP-CIT	-6.52%	-3.06%	-7.44%	-0.64%	-7.99%	6.70%
I-123 IBZM	-4.43%	-2.65%	-5.38%	-0.10%	-5.94%	7.60%
Tc-99m TRODAT	15.77%	16.21%	13.88%	20.77%	12.77%	34.61%

* ^{123}I Attenuation coefficient: 0.1481 cm^{-1} ; $^{99\text{m}}\text{Tc}$ Attenuation coefficient: 0.1526 cm^{-1}

attenuation coefficient. The percent errors associated with the left and right striata for the ^{123}I radiopharmaceutical studies are within $\pm 11\%$ for all three uptake ratios. As the activity ratios increase, the activity estimates for the left striatum slightly decrease and the activity estimates for the right striata increase. Table 6.7 lists the percent errors for the estimates of activities in the left and right striata for the 1:1, 2:1, and 5:1 striatal uptake ratios using the 19 pixels by 19 pixels ROI and the broad-beam attenuation coefficient. In the ^{123}I radiopharmaceutical studies, the activity estimates for all activity uptake ratios are

Table 6.7 Percent Errors for Activity Estimates for Broad-beam Attenuation Study

Striatal Uptake Ratios						
ROI: 19 pixels x 19 pixels	Percent Error					
	1:1		2:1		5:1	
Radiopharmaceuticals	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
I-123 beta-CIT *	-20.45%	-14.11%	-16.66%	-20.83%	-14.39%	-40.81%
I-123 FP-CIT*	-17.81%	-14.03%	-14.03%	-20.70%	-11.77%	-40.65%
I-123 IBZM*	-15.80%	-13.85%	-12.03%	-20.61%	-9.77%	-40.84%
Tc-99m TRODAT [^]	0.18%	0.68%	4.60%	-7.42%	7.28%	-31.61%

* ¹²³I Attenuation coefficient: 0.12cm⁻¹; [^] ^{99m}Tc Attenuation coefficient: 0.12 cm⁻¹

underestimated. As the activity uptake ratios increase, the percent errors for the left striatum decrease and the percent errors for the right striatum increase. The percent errors for the right striatum are similar for all ¹²³I radiopharmaceutical studies. For the ^{99m}Tc TRODAT study, the percent errors for activity estimates for left striatum and right striatum increase as the activity uptake ratios increase. The percent errors for the ^{99m}Tc TRODAT are considerably smaller than the percent errors for the ¹²³I studies.

6.2.2 Relative Activity Estimation

Relative activity uptake ratios are also calculated using the activity estimates calculated from the 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient study and the 19 pixels x 19 pixels ROI and narrow-beam attenuation coefficient study described above. These studies are chosen because the activity estimates determined in these studies have relatively small percent errors. Figure 6.7 is a graph of the percent errors for the calculated relative uptake ratios using the 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient. The percent errors associated with the estimated




Figure 6.7. Relative uptake ratio for 13 pixels x 13 pixels ROI. Percent error in activity ratio between left and right striatum using the 13 pixels x 13 pixels ROI and a narrow-beam attenuation coefficient for the 1:1, 2:1, and 5:1 activity ratio.

activity ratios for the 1:1 and 2:1 ratios are less than 10% for all radiopharmaceutical studies. The percent error for the estimated activity ratio for the 5:1 ratio varies widely among the radiopharmaceuticals. The ^{123}I β -CIT study generates the worst percent error (36.5%) and $^{99\text{m}}\text{Tc}$ TRODAT generates the best percent error (1.9%).

Figure 6.8 is a graph of the percent errors for relative uptake ratios calculated using the 19 pixels x 19 pixels ROI and narrow-beam attenuation coefficient. The percent errors associated with the estimated activity ratios for the 1:1 and 2:1 activity uptake ratio are less than -9% for all radiopharmaceutical studies. Unlike the 13 pixels x 13 pixels




Figure 6.8. Relative uptake ratio for 19 pixels x 19 pixels ROI. Percent error in activity ratio between left and right striatum using the 19 pixels x 19 pixels ROI and a narrow-beam attenuation coefficient for 1:1, 2:1, and 5:1 activity ratios.

ROI study, the differences between percent errors among the radiopharmaceuticals are relatively small. The $^{99\text{m}}\text{Tc}$ TRODAT study generates the worst percent error of -16.2% and the ^{123}I IBZM study generates the best percent error of -12.6%

6.3 Discussion

Simulation studies are performed to assess the accuracy of activity estimates for each of the striata determined using the modified conjugate imaging technique and a proposed dedicated conjugate imaging system. The accuracies of the activity uptake

ratios between the left and right striata are also determined using the technique and imaging system. The biodistributions of four radiopharmaceuticals are used in the studies to evaluate the influence of biodistribution on the accuracy of estimation of activities in left and right striata and estimation of activity uptake ratios. The agents used are chosen because of their clinical viability. The activity uptake in the right sphere representing the right striata is varied to simulate unequal activity uptake in the left and right striata characteristic in Parkinson's disease.

The accuracy of the activity estimate for each striatum is influenced by ROI size, attenuation coefficient, and biodistribution of the radiopharmaceutical. ROI size and attenuation coefficient are correlated in their influence on accuracy. The activity estimates determined from the counts in the 13 pixels x 13 pixels ROI are on average smaller than activity estimates determined from counts in the 19 pixels x 19 pixels ROI due to the lack of counts. The extent of this underestimation of activity is dependent on the attenuation coefficient applied. The application of the broad-beam attenuation coefficient results in large underestimations of the activity uptakes for the 13 pixels x 13 pixels and 19 pixels x 19 pixels ROI. This indicates that the broad-beam attenuation coefficient chosen is too broad, overestimating the detection of scattered photons. The modified conjugate imaging technique uses a measured effective attenuation coefficient to account for attenuation and the inclusion of scattered photons. It is determined through error analysis of the modified technique that -2.5% error in the attenuation coefficient of the brain results in 1-7% underestimation errors in the activity estimates. Although the use of constant attenuation coefficients is warranted, the large percent errors for the activity estimates suggest that the broad-beam attenuation coefficient used is

underestimated by more than 2.5%. A measured effective attenuation coefficient might produce more accurate results, as the effects of geometry might thus be better modeled. The studies are performed without skull effects. Modeling of a skull will increase attenuation and scatter. An effective attenuation coefficient would include these effects. However, use of a transmission-based attenuation coefficient could introduce statistical error and might require an averaging technique as did the background activity estimations.

The differences in biodistribution influence the accuracies of the background and striatal activity estimates. When the percent uptake for the striata was larger or comparable to the percent uptake of brain, the percent errors for background activities are significantly overestimated. Scatter from the striatal region contribute counts to the background ROI, which result in an overestimation of background activity. Due to the large difference in volumes between the striata and the brain, comparable percent uptakes indicate that activity concentration in the brain is significantly smaller than that of the striata. As a result, a small amount of scatter contribution into the background ROI translates to a larger fractional increase. The magnitudes of the percent errors for the striatal estimates are also associated with the differences between the percent uptake of the striata and brain. Larger differences between the percent uptakes of striata and brain result in larger percent errors in activity estimates. The differences in magnitudes for the ^{123}I radiopharmaceuticals are less apparent for the broad-beam attenuation coefficient studies.

The $^{99\text{m}}\text{Tc}$ radiopharmaceutical studies usually results in an overestimation of activity, especially when the narrow-beam attenuation coefficient is applied. This is

contrary to the ^{123}I studies, which usually result in an underestimation of activity. This difference can be attributed to the difference in emission spectra between $^{99\text{m}}\text{Tc}$ and ^{123}I , which influences the detection and inclusion of scattered photons in calculations. The inclusion of counts from scattered photons decreases the estimated camera sensitivity by increasing the number of counts detected per known activity amount. Although the inclusion of counts from scattered photons also increases the estimates of activity, the influence of scattered photons in the calculation of camera sensitivity is greater. The relationship between the detection of scattered photons and camera sensitivity is dependent on the attenuation coefficient used. The relationship between the detection of scattered photons and activity estimates is also dependent on the attenuation coefficient used. This dependence is significant, when the attenuation coefficient used is not measured. Since the camera sensitivity for ^{123}I studies is determined using the primary and high-energy photon emissions of ^{123}I , the camera sensitivity is underestimated. This underestimation results in the underestimation of the activity estimates. This suggests that the underestimation of camera sensitivity has more influence than the inclusion of scattered photons in the estimation of activity estimates. This is supported by preliminary studies that show activity estimates determined using only the primary emission of ^{123}I are overestimated, when using narrow-beam attenuation coefficient. This is because narrow-beam attenuation coefficients do not account for the detection of scattered photons. As a result, the detected scattered photons are treated as primary photons. The influence of scatter is apparent in these simulation results. Accordingly, scatter correction is applied to this simulated data. The method and results of scatter correction are discussed in chapter 7.

6.4 Clinical Relevance

The clinical application of modified conjugate imaging using the dedicated conjugate imaging system is striatal imaging with specific application to the diagnosis and management of Parkinsonism and schizophrenia. Several simulation studies are performed where the ROI size and attenuation coefficient are changed. The 19 pixels x 19 pixels ROI is deemed more clinically applicable as the ROIs drawn in clinical settings are more likely to be drawn larger than the physical dimensions of the striata due to limited spatial resolution. Although the broad-beam attenuation coefficient would be considered a better representative of attenuation coefficients determined from transmission data, the broad-beam attenuation coefficient chosen was considered too broad. Thus, the 19 pixels x 19 pixels ROI study using the narrow-beam attenuation coefficient is used to determine the technique's efficacy. For absolute activity estimation, the modified conjugate imaging technique could be used to diagnosis early- and late-stage Parkinson's disease patients with an accuracy of $\pm 11\%$. This is limited to ^{123}I radiopharmaceuticals. $^{99\text{m}}\text{Tc}$ TRODAT could be used to diagnosis and differentiate between healthy and early-stage Parkinson's disease patients using the broad-beam attenuation coefficient with an accuracy of $\pm 8\%$. Relative uptake ratio estimations are useful in monitoring the progression of Parkinson's disease and schizophrenia. The modified conjugate imaging could monitor the onset of disease, which is often presented as an uptake ratio of 2:1 with an accuracy of $<10\%$ for all radiopharmaceuticals.

CHAPTER 7

SCATTER CORRECTIONS

Approximately 20%-40% of all counts detected are from scattered photon interactions, primarily Compton scattering events (Jaszczak, Greer et al. 1984; Rosenthal, Cullom et al. 1995). The detection of these scattered events degrades spatial resolution and inflates counts. For quantitative single-photon emission studies this count inflation leads to an overestimation of activity. Thus, an accurate quantitative emission study has to correct for scatter detection in addition to effects of attenuation. Conjugate imaging, classic or the new modified technique presented in this dissertation, accounts for attenuation by using transmission data. Classic conjugate imaging corrects for attenuation by using the transmission data to estimate $e^{-\mu t}$ where t is the thickness of the patient or some part of the patient. Modified conjugate imaging accounts for attenuation by using the transmission data to estimate an effective attenuation coefficient. However, conjugate imaging does not correct for the effects of scatter detection. In chapter 6, Monte Carlo simulations using the modified conjugate imaging and a model of a dedicated conjugate imaging system are performed to evaluate the accuracy of activity estimates in each striatum. In the MCNP model, the striata are modeled as left and right spheres. Following the initial simulation, a scatter correction technique was applied to selected data sets from the simulation studies presented in chapter 6. This is done to improve the quantitative accuracy of the activity estimates in the left and right striata. The triple energy window (TEW) technique is the scatter correction technique used in this simulation study. The

results from this study are evaluated and compared with the simulated results reported in the previous chapter. This chapter discusses the triple energy window technique and reports the results of application of the technique in this setting. This chapter concludes with a discussion section that evaluates the accuracy of activity estimations determined using scatter-corrected data and of the clinical value of the combination of the TEW technique and conjugate imaging of two aligned sources.

7.1 Triple Energy Window

The triple energy window (TEW) technique was introduced in 1991 by Ogawa et al. as a practical SPECT scatter technique (Ogawa, Harata et al. 1991). The TEW technique is an energy-distribution-based scatter correction technique. Energy-distribution-based scatter correction techniques use the number of counts in various energy windows to determine the amount of scatter in the photopeak window. TEW uses counts in two narrow energy windows adjacent to the photopeak window to estimate scatter in the photopeak window. Figure 7.1 illustrates example locations of the three energy windows used in TEW. TEW has become a popular scatter correction technique for SPECT studies using radiopharmaceuticals with secondary higher-energy photon emissions. TEW can account for downscatter in the photopeak from these higher-energy photons emissions, as well as scatter from primary photon emissions. The scatter estimate is defined as

$$C_{scat} = \left(\frac{C_{lower}}{w_s} + \frac{C_{upper}}{w_s} \right) \frac{w_m}{2} \quad (7.1)$$

where C_{lower} is the number of counts in the lower-energy window; C_{upper} is the number of




Figure 7.1. Locations of the energy windows used in standard TEW. C_{lower} are the counts in the lower energy window with a width of w_s . The counts in the lower window are from scattered photons with initial energy of that of the primary photon emissions or high-energy photon emissions. C_{upper} is the counts in the upper energy window with a width of w_s . The counts in the upper window are from downscattered high-energy photon emissions. C_{total} are the counts in the photopeak energy window with a width of w_m .

counts in the upper-energy window; w_s is the width of the sub-windows; w_m is the width of the photopeak window, and C_{scat} is the number of counts in the photopeak window that represent scattered photons. The scatter estimate is referred to as a trapezoidal approximation because the equation of the scatter estimate is that of the area of a trapezoid. The estimated scatter counts are then subtracted from the photopeak window to yield the primary photons

$$C_{prim} = C_{total} - C_{scat} \quad (7.2)$$

where C_{prim} is the counts of primary photons in the photopeak window (non-scattered photons) and C_{total} is the counts of primary and scattered photons in the photopeak window (primary and scattered photons). The estimation and subtraction of scatter is performed for each pixel in each SPECT projection.

For radiopharmaceutical studies with primary and secondary higher-energy photon emissions, such as ^{123}I imaging studies, the counts in the lower-energy window are from scattering of primary photons and scattered higher-energy photons. The upper-energy window contains counts from scattered high-energy photons. For $^{99\text{m}}\text{Tc}$ imaging studies, the higher-energy sub-window is set to zero, as $^{99\text{m}}\text{Tc}$ does not emit high-energy photons, i.e., photons with energies greater than 140 keV. Thus, the scatter estimate becomes

$$C_{scat} = \left(\frac{C_{lower}}{w_s} \right) \frac{w_m}{2} \quad (7.3)$$

This scatter estimate is known as a triangular approximation as the number of scattered counts in the photopeak window is estimated by the area of a triangle.

The TEW technique is not object or imaging system dependent. It does not require object and imaging system calibration to remove the scatter component adequately. The triple energy window method is emission spectrum dependent. The effectiveness of the TEW scatter compensation technique has been reported in quantitative imaging studies using ^{123}I and $^{99\text{m}}\text{Tc}$ activity distributions (King, Long et al. 1991; Ogawa, Harata et al. 1991; Buvat, Benali et al. 1994; Buvat, Rodriguea-Villafuerte et al. 1995; Haynor, Kaplan et al. 1995; Narita, Eberl et al. 1996; Pollard, Lewellen et al. 1996; Hashimoto, Sasaki et al. 1999; Leong, O' Connor et al. 1999; Buvat, Soret et al. 2001; Moore, Kijewski et al. 2001; Soret, Koulibaly et al. 2003; Vines, Ichise et al. 2003). Since TEW can be applied as a post-processing technique and can be applied directly to projections, it has been used for planar studies. (Ichihara, Ogawa et al. 1993; Okada, Imal et al. 1999; Deplon, Ferrer et al. 2003). TEW can be performed before an attenuation correction technique is applied (Ichihara, Ogawa et al. 1993).

7.2 Monte Carlo Simulations

Chapter 6 describes and presents results from simulated conjugate imaging studies. Estimated activity uptakes for the left and right striatal spheres are determined using several combinations of ROI sizes and attenuation coefficients. The best estimates are determined when the 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient and the 19 pixels x 19 pixels ROI and broad-beam attenuation coefficient are

used. Accordingly, the scatter correction studies presented in this chapter report results, which are determined using these parameters.

Similar to chapter 6 methods, simulation of data acquisition requires two steps. The first step is to perform MCNP simulations. The second acquisition step is processing the Monte Carlo output into image array data. Once the image data arrays are generated, the image data were scatter-corrected using the TEW technique. Then the activities in the left and right striata are estimated from the scatter-corrected image data using the modified conjugate imaging technique. Data acquisition, the simulated geometry, and the application of the modified conjugate imaging technique are discussed in the methods section in chapter 6. The application of TEW to the simulated data sets is described in the following section. Two simulation studies are performed using ^{123}I and $^{99\text{m}}\text{Tc}$ radiopharmaceuticals. The radiopharmaceuticals used are ^{123}I β -CIT, ^{123}I FP-CIT, ^{123}I IBZM, and $^{99\text{m}}\text{Tc}$ TRODAT. The percent uptakes of the administered radiopharmaceutical used in these studies are listed in Table 6.1. Narrow-beam and broad-beam attenuation coefficients are used in the ^{123}I and $^{99\text{m}}\text{Tc}$ studies. The selection of the attenuation coefficients is discussed in chapter 6. Two striatal ROI sizes are used in these studies. The first set of studies use a 13 pixels x 13 pixels ROI and a narrow-beam attenuation coefficient and the second use a 19 pixels x 19 pixels ROI and a broad-beam attenuation coefficient.

7.2.1 TEW Scatter Correction

Scatter correction is performed before conjugate imaging calculations are performed. Two scatter correction studies are performed. In the first study only the

emission data is scattered corrected. In the second study the camera sensitivity data and the emission data are scattered corrected. Typically, the lower and upper energy windows are narrow and abut the photopeak window to ensure that the scatter distribution in these side windows closely approximates the scatter distribution in the photopeak window. Unfortunately, narrow windows translate to low counts, which result in noisy scatter estimations (King and Farncombe 2003). In the ^{123}I studies, five sets of energy windows are investigated. Table 7.1 lists the upper and lower energy windows used in the ^{123}I studies. In study one, larger lower window of 12 keV lower energy window and a 5 keV upper energy window are used. Study 2 uses energy windows based on Monte Carlo spectral studies. Monte Carlo studies are performed to analyze the 159 keV primary component and the high-energy component of the ^{123}I spectrum. Study 3 uses a 5 keV and a 5 keV upper energy window. Study 4 uses overlapping energy windows. Study 5 uses a 12 keV lower energy window and no upper energy window.

Table 7.1. Scatter energy windows for ^{123}I studies.

Study	Lower Energy Window	Upper Energy Window
1	130-142	175-180
2	130-140	185-190
3	137-142	175-180
4	139-144	173-178
5	130-142	0

Figure 7.2 is the emission spectral data generated from the MCNP simulations of the ^{123}I emissions. The primary and high-energy components are graphed individually. The backscatter peak for the 529 keV photon emissions smears into the upper bound of photopeak window. (A backscattered 529 photon has a energy of 173 keV. The upper bound of the photopeak in this analysis is the 170-175 keV bin as the energy of the spectrum was binned into 5 keV bins.) The backscatter peak is relatively symmetric with a peak at 180-185 keV. As a result, the counts in the 170-175 keV window from the high energy emissions are equal to the total counts in the 185-190 keV window, which are from high energy emissions. Thus, the 185-190 keV window is considered to be a good representation of the high-energy scatter contribution in the 170-175 keV. To investigate the significance of the upper energy window, the last energy window set contains only an abutting lower energy window that is 13 keV wide.

For the $^{99\text{m}}\text{Tc}$ studies the upper energy window is omitted. The lower energy window is based on MCNP spectral data study of 140 keV photon emission of $^{99\text{m}}\text{Tc}$ in water. The width and placement of the lower energy window is varied. Figure 7.3 is the spectral data generated from the MCNP simulations of the $^{99\text{m}}\text{Tc}$ emissions. Table 7.2 lists the windows used in the $^{99\text{m}}\text{Tc}$ study. The lower energy window 121-126 keV used is based on previous studies (Zaidi 1996). An abutting energy window of 114-126 and an overlapping energy window of 123-128 keV are also used.




Figure 7.2. ^{123}I emission spectrum in water. ^{123}I emission spectrum in water as detected by a gamma camera from the conjugate imaging system simulated using MCNP 4c. The gamma camera was simulated with an energy resolution of 12% full width half maximum at 159 keV. The curve entitled Total represents the complete ^{123}I spectrum. The 159 keV photon emissions are represented by the Primary curve. The high-energy photon emissions are represented by the High curve. It is apparent from the graph that the high-energy photon emissions are responsible for the counts in the upper region of the Total spectrum. There is a symmetric curve about 180 keV in the high-energy photon spectrum. This curve is due to the backscattering of the 529 keV photon emissions. The selection of 185-190 keV scatter window was based on this backscatter curve. Due to the symmetric nature of the curve, the counts in the 185-190 window closely approximates the scatter counts in the corresponding window of 170-175 keV, which is inside the photopeak window.




Figure 7.3. ^{99m}Tc photon emissions in water. ^{99m}Tc photon emissions in water as detected by a gamma camera from the conjugate imaging system simulated using MCNP 4c. The gamma camera is simulated with an energy resolution of 12% full width half maximum at 140 keV.

Table 7.2. Scatter energy windows for ^{99m}Tc studies.

Study	Lower Energy Window	Upper Energy Window
1	121-126	none
2	114-126	none
3	123-128	none

7.3 Results

7.3.1 ^{123}I Radiopharmaceutical Studies

The results presented in the following sections are determined using TEW scatter energy windows listed in tables 7.2 and 7.3 for the ^{123}I and ^{99m}Tc studies, respectively.

The ROI and attenuation coefficient combination used are 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient and 19 pixels x 19 pixels ROI and broad-beam attenuation coefficients. Results from simulation studies performed with different combinations of ROI sizes and attenuation coefficients can be found in the Appendix B.

7.3.1.1 Scatter-Corrected Emission Only Study

Tables 7.3-7.12 lists the percent errors associated with the conjugate imaging calculations using TEW scatter corrected emission data from selected ^{123}I radiopharmaceuticals. The application of TEW to the emission data results in the underestimation of activity estimates for the left and right striatum. Tables 7.3-7.7 list percent errors from studies that used a 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient. The percent errors for the activity estimates in the left and right striata are approximately equal for each activity uptake ratio and each scatter

Table 7.3. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper:175-180

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-25.86%	-25.66%	-26.66%	-25.37%	-27.15%	-24.55%
¹²³ I FP-CIT	-24.74%	-25.38%	-25.54%	-24.79%	-26.02%	-23.09%
¹²³ I IBZM	-23.75%	-24.89%	-24.58%	-24.06%	-25.07%	-21.52%

Narrow-beam attenuation coefficient used: 0.1481 cm⁻¹

Table 7.4. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 185-190

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-21.70%	-21.76%	-22.11%	-23.20%	-22.33%	-27.38%
¹²³ I FP-CIT	-20.50%	-20.65%	-20.92%	-20.98%	-21.16%	-21.97%
¹²³ I IBZM	-19.44%	-19.65%	-19.91%	-19.25%	-20.16%	-17.96%

Narrow-beam attenuation coefficient used: 0.1481 cm⁻¹

Table 7.5. Percent Error in Striatal Activity Estimates; Lower:137-142, Upper:175-180

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-25.74%	-27.17%	-26.74%	-26.87%	-27.33%	-25.96%
¹²³ I FP-CIT	-25.00%	-26.95%	-26.00%	-26.43%	-26.57%	-24.77%
¹²³ I IBZM	-24.34%	-26.60%	-25.34%	-25.90%	-25.95%	-23.81%

Narrow-beam attenuation coefficient used: 0.1481 cm⁻¹

Table 7.6. Percent Error in Striatal Activity Estimates; Lower:139-144, Upper: 173-178

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-29.69%	-32.97%	-30.46%	-32.63%	-30.92%	-31.57%
¹²³ I FP-CIT	-29.86%	-32.51%	-30.67%	-31.80%	-31.14%	-29.52%
¹²³ I IBZM	-29.89%	-32.11%	-30.68%	-31.08%	-31.16%	-28.04%

Narrow-beam attenuation coefficient used: 0.1481 cm^{-1}

Table 7.7. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 000-000

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-21.81%	-15.79%	-21.49%	-20.27%	-21.30%	-33.72%
¹²³ I FP-CIT	-20.51%	-20.72%	-20.71%	-21.90%	-20.81%	-25.34%
¹²³ I IBZM	-19.33%	-19.55%	-19.57%	-19.82%	-19.68%	-20.49%

Narrow-beam attenuation coefficient used: 0.1481 cm^{-1}

correction study. Although the percent errors for the left and right activity estimates are close in magnitude, the percent errors do vary. The largest relative difference between percent errors for the left and right activity estimates is observed in the 5:1 uptake ratio. The percent errors for the left striatum increase with an increase in the activity uptake ratio (1:1 – 5:1). The changes in magnitude in the percent error for the right activity estimate are dependent on the radiopharmaceutical and scatter energy windows used. For

the studies using the 1) 130-142 keV and 175-180 keV, 2) 137-142 keV and 175-180 keV and 3) 139-144 keV and 173-178 keV scatter energy windows, the percent errors for the activity estimates in the right striatum decrease. The studies using the 1) 130-142 keV and 185-190 keV and 2) 130-142 keV and 000-000 keV scatter energy window produce similar percent errors, where the percent errors for the right activity estimates increase with the exception of ^{123}I IBZM and scatter energy windows of 130-142 keV and 185-190 keV. The most accurate activity estimates, defined by magnitude of percent error and relative difference between the percent error for the left and right activity estimates, are produced from the study using lower energy window 130-142 keV and high-energy window of 185-190 keV. The least accurate activity estimates are produced from the study using lower energy window of 139-144 keV and upper energy window of 173-178 keV.

Tables 7.8-7.12 list the percent errors from studies using a 19 pixels x 19 pixels ROI and a broad-beam attenuation coefficient 0.12 cm^{-1} . The percent errors in this study are larger than the 13 pixels x 13 pixels ROI study, with the exception of the percent errors associated with the ^{123}I β -CIT 5:1 uptake ratio. The differences in percent errors with respect to the radiopharmaceuticals in a given scatter energy window set are larger in this study. The percent errors for the left striatum decrease as the activity uptake ratios increase for all radiopharmaceuticals and all scatter energy window studies. The changes in direction (increase or decrease) and magnitudes of percent errors for the right striatum are dependent on the radiopharmaceutical and scatter energy window study. The magnitudes of changes in the percent errors associated with the activity estimate in right

Table 7.8. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper:175-180

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-28.31%	-24.02%	-26.86%	-22.59%	-25.96%	-18.56%
¹²³ I FP-CIT	-27.72%	-26.20%	-26.24%	-26.99%	-25.37%	-29.37%
¹²³ I IBZM	-26.94%	-27.40%	-25.53%	-29.67%	-24.66%	-36.37%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

Table 7.9. Percent Error in Striatal Activity Estimates; Lower:130-142, Upper: 185-190

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-25.10%	-20.11%	-23.00%	-20.57%	-21.73%	-21.47%
¹²³ I FP-CIT	-24.07%	-21.61%	-22.01%	-23.55%	-20.76%	-29.38%
¹²³ I IBZM	-23.08%	-22.39%	-21.10%	-25.46%	-19.86%	-34.47%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

Table 7.10. Percent Error in Striatal Activity Estimates; Lower:137-142, Upper:175-180

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-25.94%	-27.14%	-25.28%	-26.06%	-24.87%	-22.74%
¹²³ I FP-CIT	-26.41%	-29.01%	-25.77%	-29.76%	-25.35%	-31.78%
¹²³ I IBZM	-26.55%	-30.15%	-25.90%	-32.21%	-25.52%	-38.39%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

Table 7.11. Percent Error in Striatal Activity Estimates; Lower: 139-144, Upper: 173-178

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-32.21%	-30.66%	-31.70%	-26.94%	-31.34%	-15.75%
¹²³ I FP-CIT	-32.41%	-33.14%	-31.94%	-32.10%	-31.60%	-28.76%
¹²³ I IBZM	-32.52%	-34.86%	-31.99%	-35.59%	-31.67%	-37.83%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

Table 7.12. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 000-000

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-27.65%	-12.97%	-23.71%	-18.39%	-21.45%	-34.23%
¹²³ I FP-CIT	-24.82%	-21.89%	-22.15%	-25.68%	-20.51%	-36.96%
¹²³ I IBZM	-23.37%	-22.38%	-20.77%	-26.99%	-19.15%	-40.58%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

striatum are significantly larger than the changes in the percent errors for the activity estimate in left striatum. Similar to the 13 pixels x 13 pixels ROI study, the largest differences in percent errors between the left and right striata are for the 5:1 activity uptake ratio. The percent errors from studies using the scatter energy windows 139-144 keV and 173-178 keV are on average larger than the percent errors from the other energy window studies. Studies using the 130-142 keV and 185-190 keV energy windows

produce the most accurate activity estimates as a whole. It is important to note that the percent errors from both ROI studies are on the same magnitude.

7.3.1.2 Scatter-Corrected Emission and Camera Sensitivity Study

Tables 7.13-7.22 list percent errors associated with activity estimates determined using scattered corrected emission and camera sensitivity data. Tables 7.13-7.17 list the percent errors from studies that used the 13 pixels x13 pixels ROI and narrow-beam attenuation coefficient 0.1481 cm^{-1} . The percent errors for the uptake ratios are significantly smaller than the corresponding scatter correction study using the same energy windows but applied only to the emission data. Although the percent errors are smaller, the differences in percent errors among the radiopharmaceuticals and scatter energy window sets are larger. The percent errors for the left striatum increase for all scatter energy windows except the scatter energy window study using the lower energy window 130-142 keV only. The percent errors for the activity estimate for the right striatum do not follow a trend but instead are dependent on the radiopharmaceutical and scatter energy windows. The study using the lower energy window of 130-142 keV only produces the worst percent errors. The most accurate activity estimates, defined by the magnitude of percent error and relative difference between the percent error for the left and right activity estimates, are produced from the study using lower energy window 137-142 keV and high-energy window of 175-180 keV.

Table 7.13. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 175-180

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-2.52%	-2.26%	-3.57%	-1.88%	-4.22%	-0.81%
¹²³ I FP-CIT	-1.05%	-1.89%	-2.10%	-1.12%	-2.74%	1.12%
¹²³ I IBZM	0.26%	-1.24%	-0.84%	-0.16%	-1.48%	3.18%

Narrow-beam attenuation coefficient used: 0.1481 cm^{-1}

Table 7.14. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 185-190

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-3.75%	-3.83%	-4.26%	-5.59%	-4.52%	-10.73%
¹²³ I FP-CIT	-2.27%	-2.46%	-2.79%	-2.87%	-3.09%	-4.08%
¹²³ I IBZM	-0.98%	-1.23%	-1.55%	-0.75%	-1.86%	0.84%

Narrow-beam attenuation coefficient: 0.1481 cm^{-1}

Table 7.15. Percent Error in Striatal Activity Estimates Lower: 137-142, Upper: 175-180

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-1.55%	-3.45%	-2.88%	-3.06%	-3.67%	-1.84%
¹²³ I FP-CIT	-0.57%	-3.16%	-1.89%	-2.47%	-2.66%	-0.26%
¹²³ I IBZM	0.30%	-2.70%	-1.03%	-1.77%	-1.83%	1.00%

Narrow-beam attenuation coefficient: 0.1481 cm^{-1}

Table 7.16. Percent Error in Striatal Activity Estimates; Lower: 139-144, Upper: 173-178

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
^{123}I β -CIT	-0.30%	-4.95%	-1.40%	-4.47%	-2.04%	-2.97%
^{123}I FP-CIT	-0.54%	-4.31%	-1.70%	-3.29%	-2.36%	-0.07%
^{123}I IBZM	-0.59%	-3.73%	-1.71%	-2.28%	-2.39%	2.04%

Narrow-beam attenuation coefficient: 0.1481cm^{-1}

Table 7.17. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 000-000

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
^{123}I β -CIT	-14.56%	-7.99%	-14.21%	-12.88%	-14.00%	-27.57%
^{123}I FP-CIT	-11.88%	-12.12%	-12.10%	-13.42%	-12.22%	-17.24%
^{123}I IBZM	-10.58%	-10.82%	-10.84%	-11.12%	-10.97%	-11.86%

Narrow-beam attenuation coefficient: 0.1481cm^{-1}

Tables 7.18 -7.22 list the percent errors from studies that used the 19 pixels x19 pixels ROI and broad-beam attenuation coefficient 0.12 cm^{-1} . The percent errors in this study are also smaller than the corresponding study in which only the emission data are scatter

Table 7.18. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 175-180

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-5.74%	-0.11%	-3.83%	1.78%	-2.66%	7.08%
¹²³ I FP-CIT	-4.97%	-2.97%	-3.02%	-4.00%	-1.88%	-7.13%
¹²³ I IBZM	-3.95%	-4.55%	-2.09%	-7.53%	-0.95%	-16.34%

Broad-beam attenuation coefficient used: 0.12 cm⁻¹

Table 7.19. Percent Error of Striatal Uptakes; Lower: 130-142, Upper: 185-190

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-7.94%	-1.80%	-5.35%	-2.36%	-3.79%	-3.48%
¹²³ I FP-CIT	-6.67%	-3.64%	-4.13%	-6.03%	-2.59%	-13.20%
¹²³ I IBZM	-5.45%	-4.60%	-3.01%	-8.37%	-1.49%	-19.45%

Broad-beam attenuation coefficient: 0.12cm⁻¹

Table 7.20. Percent Error in Striatal Activity Estimates; Lower: 137-142, Upper: 175-180

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-1.82%	-3.42%	-0.95%	-1.98%	-0.40%	2.42%
¹²³ I FP-CIT	-2.45%	-5.89%	-1.59%	-6.88%	-1.03%	-9.56%
¹²³ I IBZM	-2.62%	-7.41%	-1.77%	-10.13%	-1.26%	-18.32%

Broad-beam attenuation coefficient: 0.12cm⁻¹

Table 7.21. Percent Error in Striatal Activity Estimates; Lower: 139-144, Upper: 173-178

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-3.88%	-1.68%	-3.15%	3.59%	-2.64%	19.47%
¹²³ I FP-CIT	-4.16%	-5.20%	-3.50%	-3.72%	-3.01%	1.02%
¹²³ I IBZM	-4.32%	-7.64%	-3.56%	-8.67%	-3.12%	-11.85%

Broad-beam attenuation coefficient: 0.12cm⁻¹

Table 7.22. Percent Error in Striatal Activity Estimates; Lower: 130-142, Upper: 000-000

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Radiopharmaceuticals						
¹²³ I β-CIT	-20.94%	-4.90%	-16.63%	-10.83%	-14.17%	-28.14%
¹²³ I FP-CIT	-16.66%	-13.41%	-13.71%	-17.61%	-11.88%	-30.12%
¹²³ I IBZM	-15.06%	-13.96%	-12.17%	-19.07%	-10.38%	-34.13%

Broad-beam attenuation coefficient: 0.12cm⁻¹

corrected. The percent errors for the left striatum decrease, while percent errors for the right striatum increase as the activity uptake ratios increase for all radiopharmaceuticals and scatter energy windows. The magnitudes of the percent errors in the 19 pixels x 19 pixels ROI study are higher than the percent errors in the 13 pixels x 13 pixels ROI study. The study using only a lower scatter energy window of 130-142 keV produces the least accurate activity estimates for all radiopharmaceuticals. There is no one scatter energy

window set that produces the best results for all three activity uptake ratios and for all radiopharmaceuticals. Reasonable estimates for the 1:1 and 2:1 activity uptake ratios are determined for all scatter energy window sets. For the 5:1 activity uptake ratio, the accuracy of the activity estimates is dependent on the scatter energy window set. The study with scatter energy windows 137-142 keV and 175-180 keV produces the most accurate activity estimates for ^{123}I β -CIT uptake. The study with scatter energy windows 139-144 keV and 173-178 keV produce the most accurate activity estimates for ^{123}I FP-CIT and ^{123}I IBZM uptake.

7.3.2 $^{99\text{m}}\text{Tc}$ Radiopharmaceutical Study

7.3.2.1 Scatter-Corrected Emission Only Study

Tables 7.23-7.24 list the percent errors from studies using scatter-corrected emission data only. Table 7.23 lists the percent errors associated with activity estimates determined using 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient 0.1526 cm^{-1} . Percent errors for the left striatum decrease, while the percent errors for the right striatum increase, as the activity uptake ratios increase. The study using the lower scatter energy window of 114-126 keV produces the most accurate activity estimates for both striata. The magnitude of increases in the percent errors for the right striatum is larger than the magnitude of decrease in the percent errors for the left striatum.

Table 7.24 lists the percent errors for studies using the 19 pixels x 19 pixels ROI and broad-beam attenuation coefficient of 0.12 cm^{-1} . Similar to the 13 pixels x 13 pixels ROI study, the percent errors for the activity estimates in the left striatum decrease, while the percent errors for the right striatum increase, as the activity ratio uptakes increase.

Table 7.23. ^{99m}Tc TRODAT Percent Error in Striatal Activity Estimates

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Scatter Energy Window						
121-126 keV	-4.27%	-7.10%	-4.16%	-8.51%	-4.08%	-12.63%
114-126 keV	-3.36%	-3.50%	-2.67%	-4.50%	-2.25%	-7.46%
123-128 keV	-11.23%	2.31%	-10.84%	-2.12%	-10.59%	-15.28%

Narrow-beam attenuation coefficient: 0.1526 cm^{-1} Table 7.24. ^{99m}Tc TRODAT Percent Error in Striatal Activity Estimates

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Scatter Energy Window						
121-126 keV	-11.71%	-12.66%	-8.09%	-19.99%	-5.88%	-41.78%
114-126 keV	-10.95%	-9.27%	-6.73%	-16.16%	-4.19%	-36.78%
123-128 keV	-20.04%	-1.10%	-15.62%	-10.96%	-12.91%	-40.35%

Broad-beam attenuation coefficient: 0.12 cm^{-1}

However, the magnitudes of the percent errors are significantly larger in the 19 pixels x 19 pixels ROI study. The most accurate activity estimates are determined using the 114-126 keV scatter energy window.

7.3.2.2 Scatter-Corrected Emission and Camera Sensitivity Study

Tables 7.25-7.26 list the percent errors from studies using scatter-corrected emission and camera sensitivity data. Table 7.25 lists the percent errors associated with activity estimates determined using 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient 0.1526 cm^{-1} . The activity estimates in this study are more accurate than the activity estimates determined using scatter-corrected emission data only. The exception is the study using the scatter energy window 123-128 keV for activity uptake ratios 1:1 and 2:1. For the studies using the scatter energy windows 121-126 keV, percent errors for the left and right striata increase as the activity uptake ratios increase. For the study using the scatter energy window 114-126 keV, the percent errors for the left striatum increase, while the percent error for the right striatum decrease for activity ratios 1:1 and 2:1. For the study using the scatter energy window 123-128 keV, the percent errors for the left striatum increase, while the percent errors for the right striatum decrease as the activity uptake ratios increase. The magnitudes of change in the percent errors for the right striatum are larger than those of the left striatum, especially for ^{123}I IBZM. The differences in the percent errors among the three studies are not as great as the differences in the studies using only scatter-corrected emission data. The study using the scatter energy window of 114-126 keV produces the most accurate activity estimates.

Table 7.26 lists the percent errors for studies using the 19 pixels x 19 pixels ROI and broad-beam attenuation coefficient of 0.12 cm^{-1} . The magnitudes of percent errors are less in this 19 pixels x 19 pixels ROI study than in the corresponding study using scatter corrected emission data only. The percent errors for the activity estimates in the left striatum decrease and the percent errors for the right striatum increase as the activity

Table 7.25. ^{99m}Tc TRODAT Percent Error in Striatal Activity Estimates

ROI: 13 pixels x 13 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Scatter Energy Window						
121-126 keV	3.09%	0.03%	3.20%	-1.49%	3.29%	-5.92%
114-126 keV	1.59%	1.44%	2.31%	0.39%	2.75%	-2.73%
123-128 keV	0.92%	16.31%	1.36%	11.28%	1.64%	-3.68%

Narrow-beam attenuation coefficient: 0.1526 cm^{-1} Table 7.26. ^{99m}Tc TRODAT Percent Error in Striatal Activity Estimates

ROI: 19 pixels x 19 pixels	Striatal Uptake Ratios					
Percent Error	1:1		2:1		5:1	
	Left Striatum	Right Striatum	Left Striatum	Right Striatum	Left Striatum	Right Striatum
Scatter Energy Window						
121-126 keV	-4.93%	-5.95%	-1.03%	-13.85%	1.35%	-37.31%
114-126 keV	-6.39%	-4.62%	-1.96%	-11.87%	0.71%	-33.55%
123-128 keV	-9.10%	12.43%	-4.07%	1.22%	-1.00%	-32.18%

Broad-beam attenuation coefficient: 0.12 cm^{-1}

uptake ratios increase. The most accurate activity estimates for the 1:1 and 2:1 uptake ratios are determined using the 114-126 keV scatter energy window. The most accurate activity estimates for the 5:1 uptake ratio is determined using the 123-128 keV scatter energy window.

7.4 Discussion

The triple energy window (TEW) technique is the scatter correction technique used to scatter correct data used in the modified conjugate imaging studies. It is chosen because of its applicability to planar images, its treatment of high-energy photon emission downscatter and easy implementation into the clinic. Conventional lower and upper scatter energy windows have been narrow and about the photopeak window (Ogawa, Harata et al. 1991; Ichihara, Ogawa et al. 1993; Rosenthal, Cullom et al. 1995; Narita, Eberl et al. 1996; Okada, Imal et al. 1999; Deplon, Ferrer et al. 2003; King and Farncombe 2003; Soret, Koulibaly et al. 2003). In addition to the standard scatter energy windows, a couple of unconventional scatter energy window sets are used in these studies. These are not equal in width and not adjacent to the photopeak window. A total of five energy window sets are used to scatter correct the image data in the ^{123}I studies. Three energy windows were used in the $^{99\text{m}}\text{Tc}$ study. Scatter correction techniques are usually applied to emission data only. In these studies, two scatter correction studies are compared; one in which scatter correction is applied only to the emission data and the second study in which scatter correction is applied to on the emission data and the camera sensitivity data. The effects of scatter correction on the conjugate imaging calculations are evaluated. Conjugate imaging calculations are performed with a 13 pixels x 13 pixels ROI and a narrow-beam attenuation coefficient and a 19 pixels x 19 pixels ROI and broad-beam attenuation coefficient were used for the ^{123}I and $^{99\text{m}}\text{Tc}$ studies.

It was anticipated that the application of scatter correction would further improve the accuracy of the estimations of activity in the left and right striatum. However, the magnitude of improvement is dependent on 1) which data sets scatter correction are

applied; 2) which scatter energy windows are used; 3) the size of the ROI; and 4) the radiopharmaceutical. Activity estimates determined using scatter-corrected emission and scatter-corrected camera sensitivity data are more accurate than activity estimates determined with scatter-corrected emission data only. Scatter correcting the emission data decreases the activity estimates, while scatter correcting the camera sensitivity counters this decrease by increasing activity estimates due to the increased camera sensitivity. Scattered corrected camera sensitivity increases activity estimates because a smaller number of counts now represent the known amount of activity, which translated into larger activity estimates.

Although the scatter-corrected emission and camera sensitivity data provide more accurate activity estimates, the scatter energy windows do influence the degree of accuracy. The scatter energy window selection is influenced by the ROI size and biodistribution of the radiopharmaceuticals. The accuracy of the activity estimates changes with ROI size because of the inclusion of more scatter. For radiopharmaceuticals, like ^{123}I β -CIT, with large brain to striatal activity concentration ratio, the activity estimates in the 19 pixels x 19 pixels ROI are, in general, not as underestimated as activity estimates determined using radiopharmaceuticals with small head to striatal activity concentration ratio, like ^{123}I FP-CIT, ^{123}I IBZM, and $^{99\text{m}}\text{Tc}$ TRODAT. For larger brain to striatal activity concentration ratios, a large fraction of counts due to scattered photons are from brain activity. Scatter correction removes these counts. However, for smaller brain to striatal activity concentration ratio, the fraction of counts from scattered photons from brain activity is small. As a result, scatter correction removes more primary counts. This is especially apparent in the 2:1 and 5:1 uptake ratios

where the overall activity amount is less. This is most likely compounded by the use of a too broad attenuation coefficient, which removes too large a fraction of counts.

Several studies are performed with various parameters. The best scatter correction results are produced when the 13 pixels x 13 pixels ROI, narrow-beam attenuation coefficient, and scatter corrected emission and camera sensitivity data are used. Although all scatter energy windows produce reasonable percent errors, the best scatter energy windows are 130-142 keV and 185-190 keV for uptake ratios 1:1 and 2:1. This is based on magnitude and relative difference of percent error between the left and right activity estimates. The larger lower scatter energy window of 130-142 keV provides a better scatter estimate because it introduces less noise due to the larger number of counts. The upper energy window of 185-190 keV is a better estimation of the high-energy scatter because unlike the 175-180 keV energy window it contains only scattered counts. The best scatter energy window for the 5:1 uptake ratio is the scatter energy windows 139-144 keV and 173-178 keV. For the ^{99m}Tc studies, the wider lower energy window 114-121 keV provides the most accurate activity estimates. In the wider window of 114-126 keV, the signal to noise ratio is smaller than 121-126 keV window providing for better estimation of scatter.

However, the use of the 13 pixels x 13 pixels ROI and narrow-beam attenuation coefficient does not model adequately the clinical setup. The ROI drawn in the clinic will most likely be drawn around the striatum highlighted from radiopharmaceutical. Due to blurring the striatum projected on the image will be larger than its physical size. This indicates that the ROI will not have the dimensions of the physical size of the striatum. The attenuation coefficient determined from the flood data will most likely be closer to

the value of the broad-beam attenuation coefficient even with scatter correction. Thus, the results determined from the 19 pixels x 19 pixels ROI and broad-beam attenuation coefficient most accurately model the clinical setting. The broad attenuation coefficient used in these studies overcompensates by removing too many counts. The activity estimates determined using these parameters are less accurate than the 13 pixels x 13 pixels ROI study, particularly for the striatum containing less activity (right striatum). Although the activities are less accurate, the percent errors are reasonable for the 1:1 and 2:1 uptake ratio for all radiopharmaceuticals and most scatter energy window sets. The percent errors for the activity estimates for both striata for the 5:1 uptake ratio are within 10% for the ^{123}I β -CIT and ^{123}I FP-CIT for most scatter energy windows.

7.5 Clinical Relevance

The application of TEW on emission data and camera sensitivity data improves the accuracy of the activity estimates for the left and right striata. The percent errors associated with the activity estimates for the left and right striata in the ^{123}I β -CIT and ^{123}I FP-CIT studies for the 1:1, 2:1, and 5:1 uptake ratios suggest that the modified conjugate imaging in conjunction with TEW is a feasible quantitative estimation technique for early to late Parkinson's disease striatal studies. ^{123}I IBZM and $^{99\text{m}}\text{Tc}$ TRODAT provide reasonable results for the healthy and early Parkinson's disease patients.

CHAPTER 8

CONCLUSIONS AND FUTURE WORK

8.1 Conclusions

In this dissertation, the development and characterization of conjugate imaging of two aligned activity distributions embedded in an attenuating medium, with specific application to the quantification of activity uptake in the striata, were reported. The actualization of this research involved the completion of six tasks. Tasks one two and three involved the model and general application of the technique. Tasks four, five and six were related to the specific application of striatal imaging and neurodegenerative diseases.

The development and assessment of conjugate imaging of the two aligned activity distributions with respect to the techniques general applicability were addressed by tasks one, two, and three. Task one was the derivation of a tractable mathematical model for conjugate imaging of two aligned activity distributions. Other extensions of conjugate imaging have been developed, but none have provided a practical means of estimating accurately two distinct aligned distributions. Task two was an error analysis of the mathematical model. The error analysis demonstrated the model was sensitive to the accuracy of the measured input parameters and the activity uptake ratio between the two activity distributions. The theoretic error was larger for a 5:1 activity uptake ratio than for a 1:1 ratio. Task three entailed the completion of a series of phantom studies to evaluate the practical application of the developed model. Software was developed to apply the

model to image data acquired using a phantom and a commercially available gamma camera system. The activity estimates obtained had accuracy consistent with that predicted by the mathematical model. The activity estimates were also comparable to activity estimates determined from quantitative analysis of attenuation-corrected SPECT. The results from the first three tasks indicated that the modified conjugate imaging technique presented can be used to estimate activity accurately in each of two aligned activity distributions with a diameter of 3cm or less separated by a distance s .

Tasks four, five, and six involved the specific application of conjugate imaging of activity uptake in the striata of a primate using a novel dedicated small field-of-view (FOV) imaging system and ^{123}I and $^{99\text{m}}\text{Tc}$ - labeled radiopharmaceuticals. These last three tasks utilized a Monte Carlo simulation model developed to simulate conjugate imaging of activity uptake in the striata of a primate using two small FOV gamma cameras. Task four was the designing and characterization of an appropriate shielding design for the small FOV gamma cameras comprising the dedicated conjugate imaging system. Task five was the use of the simulation model to assess the accuracy of the activity estimates in each striatum. To represent healthy, early-stage Parkinson's, and late-stage Parkinson's patients, three activity uptake ratios were investigated. In addition, four radiopharmaceuticals, two attenuation coefficients, and 2 regions of interest (ROIs) were investigated. The activity estimates obtained had accuracies consistent with the theoretic error analysis. However, the theoretic error analysis did not incorporate the influence of ROI size variability and the emission spectrum and biodistribution of radiopharmaceuticals. It was apparent from the differences in accuracies of activity estimates determined using different ROI sizes and different radiopharmaceuticals that

ROI size and choice of radiopharmaceuticals influence the accuracy of the activity estimates. Although the accuracies were consistent with the trends predicted by the theoretic error analysis, the influence of ROI size and properties of the radiopharmaceuticals introduced new correlations, which warrant further investigation. Task six was the investigation of the application of scatter correction to data acquired from conjugate imaging of two aligned activity distributions. The triple energy window (TEW) scatter correction technique was applied to the data obtained from the Monte Carlo simulations. It was determined that TEW scatter correction must be applied to emission and transmission data to improve accuracies of activity estimates. Even though, the scatter correction investigation was specific to the conjugate imaging of primate striata, it is probable key aspects of the investigation extrapolate to the general application of conjugate imaging of two aligned activity distributions. The results from tasks four through six indicated the viability of the modified conjugate imaging technique as a tool to evaluate striatal activity uptake.

8.2 Future Work

The modified conjugate imaging technique presented in this dissertation is a viable quantitative imaging technique to evaluate activity uptake in the striata. Its accuracy is comparable to that of diagnostic quantitative emission techniques currently used in clinical settings. A novel dedicated conjugate imaging system is capable of acquiring the data required for modified conjugate imaging of striatal uptake. A next step in the development of the imaging technique and imaging system is to refine their characterization.

The Monte Carlo simulations were performed with modifications. The attenuation coefficients and camera sensitivity were determined for one pixel and applied to all pixels instead of on a pixel-by-pixel basis. Future research could include performing Monte Carlo simulations using flood data to determine the effective attenuation coefficients and camera sensitivity for each pixel. An anthropomorphic phantom could be used in the Monte Carlo simulations to account for skull effects. Furthermore, the application of alternative scatter correction techniques that are applicable to planar data should be pursued.

Another research direction would include an investigation and clarification of the origins of the errors associated with the activity estimates. Error analysis of the modified conjugate imaging technique indicated that the accuracy of the activity estimates was sensitive to error in depth measurements. This sensitivity is exacerbated when the activity uptake ratio is not one to one. However, there was still a noticeable difference in the percent errors between the large and small activity striatum in the simulated studies, where the depth of each striatum was accurately known. A thorough investigation could be conducted to determine what other possible parameters might cause this percent error difference.

APPENDIX A

Listed below are the activity uptake estimates determined using the modified conjugate imaging technique for various ROI sizes and two attenuation coefficients. This data represent the complete set of data generated from the simulations.

Table A.1. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole			LStr (microCi)	RStr (microCi)	Sum (microCi)							
Non-SC				100	100	200							
Atten Coeff	0.1481												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15451	0.048653	6.3288	85.2662	83.5891	165.6842	171.4747	-14.73%	-16.41%	-17.16%	-14.26%	1.02	2.01%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15373	0.048653	6.2967	87.4493	91.456	173.5889	179.6554	-12.55%	-8.54%	-13.21%	-10.17%	0.96	-4.38%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15283	0.048653	6.2598	89.4319	94.6157	175.561	181.6964	-10.57%	-5.38%	-12.22%	-9.15%	0.95	-5.48%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1518	0.048653	6.2178	90.7658	96.5738	176.3087	182.4701	-9.23%	-3.43%	-11.85%	-8.76%	0.94	-6.01%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.14649	0.048653	6.0001	95.1331	101.4049	175.8076	181.9514	-4.87%	1.40%	-12.10%	-9.02%	0.94	-6.18%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.14099	0.048653	5.7749	99.3496	109.541	179.3776	185.6455	-0.65%	9.54%	-10.31%	-7.18%	0.91	-9.30%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15451	0.048653	6.3288	75.4421	73.6226	146.2603	149.6095	-24.56%	-26.38%	-26.87%	-25.20%	1.02	2.47%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15373	0.048653	6.2967	76.798	81.1624	153.2553	156.7646	-23.20%	-18.84%	-23.37%	-21.62%	0.95	-5.38%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.15283	0.048653	6.2598	78.4431	84.0752	155.0041	158.5534	-21.56%	-15.92%	-22.50%	-20.72%	0.93	-6.70%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1518	0.048653	6.2178	79.5545	85.8944	155.6747	159.2393	-20.45%	-14.11%	-22.16%	-20.38%	0.93	-7.38%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.14649	0.048653	6.0001	83.3828	90.3081	155.2638	158.819	-16.62%	-9.69%	-22.37%	-20.59%	0.92	-7.67%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.14099	0.048653	5.7749	86.7786	97.9962	158.4615	162.0895	-13.22%	-2.00%	-20.77%	-18.96%	0.89	-11.45%	

Table A.2. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy			LStr (microCi)	RStr (microCi)	Sum (microCi)							
Non-SC				100	100	200							
Atten Coeff	0.1481												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13637	0.066137	5.5859	104.2062	101.7587	202.3129	209.3835	4.21%	1.76%	1.16%	4.69%	1.02	2.41%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13568	0.066137	5.5576	107.3398	111.9397	212.8803	220.32	7.34%	11.94%	6.44%	10.16%	0.96	-4.11%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13489	0.066137	5.525	110.0653	116.7492	216.4444	224.0084	10.07%	16.75%	8.22%	12.00%	0.94	-5.73%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13398	0.066137	5.4879	112.2091	119.888	218.4146	226.0473	12.21%	19.89%	9.21%	13.02%	0.94	-6.41%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12929	0.066137	5.2958	120.4859	128.5718	222.0182	229.7768	20.49%	28.57%	11.01%	14.89%	0.94	-6.29%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12444	0.066137	5.097	127.5697	140.8909	229.0775	237.082	27.57%	40.89%	14.54%	18.54%	0.91	-9.45%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13637	0.066137	5.5859	92.2405	89.5869	178.5986	182.6884	-7.76%	-10.41%	-10.70%	-8.66%	1.03	2.96%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13568	0.066137	5.5576	94.2997	99.3133	187.95	192.2538	-5.70%	-0.69%	-6.03%	-3.87%	0.95	-5.05%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13489	0.066137	5.525	96.5148	103.7759	191.1084	195.4844	-3.49%	3.78%	-4.45%	-2.26%	0.93	-7.00%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13398	0.066137	5.4879	98.3026	106.6877	192.8646	197.2807	-1.70%	6.69%	-3.57%	-1.36%	0.92	-7.86%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12929	0.066137	5.2958	105.6174	114.5459	196.1178	200.6084	5.62%	14.55%	-1.94%	0.30%	0.92	-7.79%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12444	0.066137	5.097	111.4544	126.1396	202.4468	207.0819	11.45%	26.14%	1.22%	3.54%	0.88	-11.64%	

Table A.3. ¹²³I β-CIT 2:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole			LStr (microCi)	RStr (microCi)	Sum (microCi)								
Non-SC				100	50	150								
Atten Coeff	0.1481													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15175	0.048813	6.2155	85.5289	39.1927	120.6074	124.7845	-14.47%	-21.61%	-19.60%	-16.81%	2.18	9.11%		
Atten Coeff	0.1481													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15098	0.048813	6.184	86.8469	45.5504	126.0203	130.3969	-13.15%	-8.90%	-15.99%	-13.07%	1.91	-4.67%		
Atten Coeff	0.1481													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15009	0.048813	6.1477	88.4166	48.1324	127.137	131.5553	-11.58%	-3.74%	-15.24%	-12.30%	1.84	-8.15%		
Atten Coeff	0.1481													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.14908	0.048813	6.1065	89.8561	49.3254	127.4433	131.8731	-10.14%	-1.35%	-15.04%	-12.08%	1.82	-8.91%		
Atten Coeff	0.1481													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.14386	0.048813	5.8927	94.348	52.8252	126.3031	130.6953	-5.65%	5.65%	-15.80%	-12.87%	1.79	-10.70%		
Atten Coeff	0.1481													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13846	0.048813	5.6715	98.372	60.532	129.5579	134.0702	-1.63%	21.06%	-13.63%	-10.62%	1.63	-18.74%		
Atten Coeff	0.12													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15175	0.048813	6.2155	80.1	30.003	106.4551	108.8708	-19.90%	-39.99%	-29.03%	-27.42%	2.67	33.49%		
Atten Coeff	0.12													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15098	0.048813	6.184	80.7552	36.1433	111.2495	113.7809	-19.24%	-27.71%	-25.83%	-24.15%	2.23	11.72%		
Atten Coeff	0.12													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.15009	0.048813	6.1477	82.05	38.5289	112.2425	114.7983	-17.95%	-22.94%	-25.17%	-23.47%	2.13	6.48%		
Atten Coeff	0.12													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.14908	0.048813	6.1065	83.342	39.5832	112.5226	115.0852	-16.66%	-20.83%	-24.98%	-23.28%	2.11	5.27%		
Atten Coeff	0.12													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.14386	0.048813	5.8927	87.4238	42.6673	111.542	114.0835	-12.68%	-14.67%	-25.64%	-23.94%	2.05	2.45%		
Atten Coeff	0.12													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13846	0.048813	5.6715	90.6714	49.9445	114.458	117.07	-9.33%	-0.11%	-23.69%	-21.95%	1.82	-9.23%		

Table A.4. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy			LStr (microCi)	RStr (microCi)	Sum (microCi)								
Non-SC				100	50	150								
Atten Coeff	0.1481													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.1356	0.066243	5.5541	103.8559	49.309	148.3487	153.4904	3.86%	-1.38%	-1.10%	2.33%	2.11	5.31%		
Atten Coeff	0.1481													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13491	0.066243	5.526	105.7571	57.9126	155.9671	161.3879	5.76%	15.83%	3.98%	7.59%	1.83	-8.69%		
Atten Coeff	0.1481													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13412	0.066243	5.4935	107.7927	62.2475	158.5941	164.1112	7.79%	24.50%	5.73%	9.41%	1.73	-13.42%		
Atten Coeff	0.1481													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13322	0.066243	5.4567	109.8481	64.7632	159.9531	165.5195	9.85%	29.53%	6.64%	10.35%	1.70	-15.19%		
Atten Coeff	0.1481													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12856	0.066243	5.2656	117.0373	73.0175	162.7306	168.3992	17.04%	46.04%	8.49%	12.27%	1.60	-19.86%		
Atten Coeff	0.1481													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12373	0.066243	5.068	123.159	85.7935	169.6449	175.5628	23.16%	71.59%	13.10%	17.04%	1.44	-28.22%		
Atten Coeff	0.12													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.1356	0.066243	5.5541	97.0942	38.1207	130.9442	133.9177	-2.91%	-23.76%	-12.70%	-10.72%	2.55	27.35%		
Atten Coeff	0.12													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13491	0.066243	5.526	96.0969	46.4178	137.6906	140.6261	-1.90%	-7.16%	-8.21%	-6.12%	2.11	5.67%		
Atten Coeff	0.12													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13412	0.066243	5.4935	99.6781	50.4827	140.0217	143.2132	-0.32%	0.97%	-6.65%	-4.52%	1.97	-1.27%		
Atten Coeff	0.12													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13322	0.066243	5.4567	101.4409	52.7893	141.2386	144.459	1.44%	5.58%	-5.84%	-3.69%	1.92	-3.92%		
Atten Coeff	0.12													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12856	0.066243	5.2656	107.719	60.332	143.7582	147.0393	7.72%	20.66%	-4.16%	-1.97%	1.79	-10.73%		
Atten Coeff	0.12													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12373	0.066243	5.068	112.5486	72.4754	149.957	153.3846	12.55%	44.95%	-0.03%	2.26%	1.55	-22.35%		

Table A.5. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole			LStr (microCi)	RStr (microCi)	Sum (microCi)							
Non-SC				100	20	120							
Atten Coeff	0.1481												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.15007	0.048936	6.1469	85.6892	12.5564	92.5433	95.6605	-14.31%	-37.22%	-22.88%	-20.28%	6.82	36.49%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14931	0.048936	6.1157	86.491	18.0099	96.4345	99.7109	-13.51%	-9.95%	-19.64%	-16.91%	4.80	-3.95%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14843	0.048936	6.0799	87.8202	20.2495	97.0538	100.3698	-12.18%	1.25%	-19.12%	-16.37%	4.34	-13.26%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14744	0.048936	6.0391	89.3109	21.0184	97.0659	100.3746	-10.69%	5.09%	-19.11%	-16.35%	4.25	-15.02%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14228	0.048936	5.8276	93.8563	23.7841	95.3168	98.5735	-6.14%	18.92%	-20.57%	-17.86%	3.95	-21.08%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.13694	0.048936	5.6089	97.7888	31.1779	98.3914	101.7735	-2.21%	55.89%	-18.01%	-15.19%	3.14	-37.27%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.15007	0.048936	6.1469	82.8973	3.8328	81.6578	83.4597	-17.10%	-80.84%	-31.95%	-30.45%	21.63	332.57%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14931	0.048936	6.1157	83.1346	9.1343	85.1099	87.0043	-16.87%	-54.33%	-29.08%	-27.50%	9.10	82.03%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14843	0.048936	6.0799	84.226	11.207	85.664	87.5757	-15.77%	-43.97%	-28.61%	-27.02%	7.52	50.31%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14744	0.048936	6.0391	85.6109	11.838	85.6839	87.5974	-14.39%	-40.81%	-28.60%	-27.00%	7.23	44.64%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.14228	0.048936	5.8276	89.8171	14.1908	84.1665	86.0505	-10.18%	-29.05%	-29.86%	-28.29%	6.33	26.58%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio	
0.13694	0.048936	5.6089	93.0064	21.1626	86.924	88.8815	-6.99%	5.81%	-27.56%	-25.93%	4.39	-12.10%	

Table A.6. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy			LStr (microCi)	RStr (microCi)	Sum (microCi)								
Non-SC				100	20	120								
Atten Coeff	0.1481													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13506	0.066405	5.5319	103.6582	17.8496	114.8941	118.7787	3.66%	-10.75%	-4.25%	-1.02%	5.81	16.15%		
Atten Coeff	0.1481													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13437	0.066405	5.5039	104.8183	25.5156	120.7397	124.8592	4.82%	27.58%	0.62%	4.05%	4.11	-17.84%		
Atten Coeff	0.1481													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13358	0.066405	5.4716	106.441	29.5745	122.7208	126.9212	6.44%	47.87%	2.27%	5.77%	3.60	-28.02%		
Atten Coeff	0.1481													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13269	0.066405	5.4349	108.4209	31.7219	123.7374	127.9789	8.42%	58.61%	3.11%	6.65%	3.42	-31.64%		
Atten Coeff	0.1481													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12804	0.066405	5.2446	115.131	39.7152	126.0584	130.399	15.13%	98.58%	5.05%	8.67%	2.90	-42.02%		
Atten Coeff	0.1481													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12324	0.066405	5.0477	120.7654	52.6776	132.8607	137.4614	20.77%	163.39%	10.72%	14.55%	2.29	-54.15%		
Atten Coeff	0.12													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13506	0.066405	5.5319	100.0177	7.2503	101.3861	103.6318	0.02%	-63.75%	-15.51%	-13.64%	13.79	175.90%		
Atten Coeff	0.12													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13437	0.066405	5.5039	100.3849	14.698	106.5695	108.951	0.38%	-26.51%	-11.19%	-9.21%	6.83	36.60%		
Atten Coeff	0.12													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13358	0.066405	5.4716	101.5857	18.5328	108.331	110.7602	1.59%	-7.34%	-9.72%	-7.70%	5.48	9.63%		
Atten Coeff	0.12													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.13269	0.066405	5.4349	103.311	20.4838	109.2459	111.6993	3.31%	2.42%	-8.96%	-6.92%	5.04	0.87%		
Atten Coeff	0.12													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12804	0.066405	5.2446	109.1381	27.8163	111.3622	113.8745	9.14%	39.08%	-7.20%	-5.10%	3.92	-21.53%		
Atten Coeff	0.12													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	%Err Ratio		
0.12324	0.066405	5.0477	113.4617	40.193	117.463	120.128	13.46%	100.97%	-2.11%	0.11%	2.82	-43.54%		

Table A.7. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			40	40	80										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036241	0.010658	1.4844	34.6736	34.3725	68.5012	70.8953	-13.32%	-14.07%	-14.37%	-11.38%	1.01	0.88%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036058	0.010658	1.4769	36.151	36.928	72.0297	74.547	-9.62%	-7.68%	-9.96%	-6.82%	0.98	-2.10%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035846	0.010658	1.4683	36.9593	38.0219	73.1168	75.6721	-7.60%	-4.95%	-8.60%	-5.41%	0.97	-2.79%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035606	0.010658	1.4584	37.3938	38.7755	73.6726	76.2474	-6.52%	-3.06%	-7.91%	-4.69%	0.96	-3.56%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034359	0.010658	1.4073	38.7877	40.4335	74.2187	76.8124	-3.03%	1.08%	-7.23%	-3.98%	0.96	-4.07%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.033069	0.010658	1.3545	39.9703	42.5285	75.3107	77.9426	-0.07%	6.32%	-5.86%	-2.57%	0.94	-6.02%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036241	0.010658	1.4844	30.5402	30.3141	60.473	61.8578	-23.40%	-24.21%	-24.41%	-22.68%	1.01	1.08%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036058	0.010658	1.4769	31.8373	32.686	63.595	65.0513	-20.41%	-18.29%	-20.51%	-18.69%	0.97	-2.60%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035846	0.010658	1.4683	32.527	33.6817	64.5578	66.0361	-18.68%	-15.80%	-19.30%	-17.45%	0.97	-3.43%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035606	0.010658	1.4584	32.8775	34.3869	65.0519	66.5415	-17.81%	-14.03%	-18.69%	-16.62%	0.96	-4.39%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034359	0.010658	1.4073	34.0893	35.9039	65.5447	67.0456	-14.78%	-10.24%	-18.07%	-16.19%	0.95	-5.05%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.033069	0.010658	1.3545	35.0605	37.8735	66.523	68.0462	-12.35%	-5.32%	-16.85%	-14.94%	0.93	-7.43%			

Table A.8. ¹²³I FP -CIT 1:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			40	40	80										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030445	0.014555	1.247	41.9472	41.4942	82.8219	85.7164	4.87%	3.74%	3.53%	7.15%	1.01	1.09%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030291	0.014555	1.2407	43.8116	44.7141	87.2975	90.3484	9.53%	11.79%	9.12%	12.94%	0.98	-2.02%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030113	0.014555	1.2334	44.8447	46.2407	88.8816	91.9879	12.11%	15.60%	11.10%	14.98%	0.97	-3.02%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029911	0.014555	1.2252	45.4758	47.3035	89.7924	92.9305	13.69%	18.26%	12.24%	16.16%	0.96	-3.86%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028864	0.014555	1.1823	47.7272	49.9105	91.3593	94.5521	19.32%	24.78%	14.20%	18.19%	0.96	-4.37%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02778	0.014555	1.1379	49.6075	52.9052	93.2049	96.462	24.02%	32.26%	16.51%	20.58%	0.94	-6.23%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030445	0.014555	1.247	37.0766	36.5866	73.1158	74.7902	-7.31%	-8.53%	-8.61%	-6.51%	1.01	1.34%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030291	0.014555	1.2407	38.5881	39.5746	77.0763	78.8413	-3.53%	-1.06%	-3.65%	-1.45%	0.98	-2.49%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030113	0.014555	1.2334	39.4571	40.9735	78.4794	80.2765	-1.36%	2.43%	-1.90%	0.35%	0.96	-3.70%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029911	0.014555	1.2252	39.9698	41.966	79.2887	81.1044	-0.08%	4.92%	-0.89%	1.38%	0.95	-4.76%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028864	0.014555	1.1823	41.9364	44.3413	80.6928	82.5406	4.84%	10.65%	0.87%	3.18%	0.95	-5.42%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02778	0.014555	1.1379	43.5124	47.1444	82.3475	84.233	8.78%	17.86%	2.93%	5.29%	0.92	-7.70%			

Table A.9. ¹²³I FP -CIT 2:1 Activity Uptake Ratio Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)												
Non-SC			40	20	60												
Atten Coeff	0.1481																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.035157	0.010689	1.44	34.7742	16.6124	50.6044	52.3685		-13.06%	-16.94%	-15.66%	-12.74%	2.09	4.66%				
Atten Coeff	0.1481																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034979	0.010689	1.4328	35.9077	18.5607	53.1984	55.0454		-10.23%	-7.20%	-11.34%	-8.26%	1.93	-3.27%				
Atten Coeff	0.1481																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034774	0.010689	1.4243	36.5502	19.4209	53.9498	55.8238		-8.62%	-2.90%	-10.08%	-6.96%	1.88	-5.90%				
Atten Coeff	0.1481																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034541	0.010689	1.4148	37.0245	19.8725	54.2613	56.1466		-7.44%	-0.64%	-9.56%	-6.42%	1.86	-6.84%				
Atten Coeff	0.1481																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.033332	0.010689	1.3653	38.4472	20.9854	54.4487	56.3412		-3.88%	4.93%	-9.25%	-6.10%	1.83	-8.40%				
Atten Coeff	0.1481																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.03208	0.010689	1.314	39.5894	22.8257	55.3782	57.3048		-1.03%	14.13%	-7.70%	-4.49%	1.73	-13.28%				
Atten Coeff	0.12																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.035157	0.010689	1.44	32.4993	12.8649	44.669	45.6835		-18.75%	-35.68%	-25.55%	-23.86%	2.53	26.31%				
Atten Coeff	0.12																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034979	0.010689	1.4328	33.4179	14.6736	46.965	48.0333		-16.46%	-26.63%	-21.73%	-19.94%	2.28	13.87%				
Atten Coeff	0.12																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034774	0.010689	1.4243	33.9672	15.4559	47.6308	48.7148		-15.08%	-22.72%	-20.62%	-18.81%	2.20	9.88%				
Atten Coeff	0.12																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.034541	0.010689	1.4148	34.387	15.8594	47.9086	48.9992		-14.03%	-20.70%	-20.15%	-18.33%	2.17	8.41%				
Atten Coeff	0.12																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.033332	0.010689	1.3653	35.6813	16.8343	48.0631	49.1781		-10.80%	-15.63%	-19.86%	-18.04%	2.12	5.96%				
Atten Coeff	0.12																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.03208	0.010689	1.314	36.6376	18.5542	48.9164	50.0315		-8.41%	-7.23%	-18.47%	-16.61%	1.97	-1.27%				

Table A.10. ¹²³I FP -CIT 2:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			40	20	60										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030164	0.014556	1.2355	41.8015	20.5113	61.4189	63.547	4.50%	2.56%	2.36%	5.91%	2.04	1.90%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030011	0.014556	1.2292	43.1736	23.0972	64.7592	67.009	7.93%	15.49%	7.93%	11.68%	1.87	-6.54%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029635	0.014556	1.222	43.9278	24.4321	65.9453	68.2379	9.82%	22.16%	9.91%	13.73%	1.80	-10.10%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029635	0.014556	1.2138	44.5099	25.243	66.5775	68.8929	11.27%	26.22%	10.96%	14.82%	1.76	-11.84%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028597	0.014556	1.1713	46.4115	27.6227	67.7339	70.0916	16.03%	38.11%	12.89%	16.82%	1.68	-15.99%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027524	0.014556	1.1274	47.8528	30.741	69.4314	71.851	19.63%	53.71%	15.72%	19.75%	1.56	-22.17%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030164	0.014556	1.2355	39.0131	15.9973	54.2144	55.4463	-2.47%	-20.01%	-9.64%	-7.59%	2.44	21.94%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.030011	0.014556	1.2292	40.1024	18.4109	57.1728	58.4742	0.26%	-7.95%	-4.71%	-2.54%	2.18	8.91%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029835	0.014556	1.222	40.7151	19.6492	58.2239	59.5501	1.79%	-1.75%	-2.96%	-0.75%	2.07	3.60%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029635	0.014556	1.2138	41.2046	20.3982	58.7865	60.1261	3.01%	1.99%	-2.02%	0.21%	2.02	1.00%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028597	0.014556	1.1713	42.8466	22.5842	59.8261	61.1905	7.12%	12.92%	-0.29%	1.98%	1.90	-5.14%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027524	0.014556	1.1274	43.9728	25.5542	61.349	62.7497	9.93%	27.77%	2.25%	4.58%	1.72	-13.96%			

Table A.11. ¹²³I FP -CIT 5:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)												
Non-SC			40	8	48												
Atten Coeff	0.1481																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034499	0.010723	1.4131	34.8357	5.958	39.4517	40.7854	-12.91%	-25.53%	-17.81%	-15.03%	5.85	16.94%					
Atten Coeff	0.1481																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034324	0.010723	1.4059	35.7633	7.5425	41.4438	42.8522	-10.59%	-5.72%	-13.66%	-10.72%	4.74	-5.17%					
Atten Coeff	0.1481																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034122	0.010723	1.3977	36.3061	8.2636	42.0021	43.4322	-9.23%	3.30%	-12.50%	-9.52%	4.39	-12.13%					
Atten Coeff	0.1481																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.033894	0.010723	1.3883	36.8044	8.5358	42.16	43.5964	-7.99%	6.70%	-12.17%	-9.17%	4.31	-13.76%					
Atten Coeff	0.1481																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.032707	0.010723	1.3397	38.2464	9.3412	42.1035	43.5403	-4.38%	16.77%	-12.28%	-9.29%	4.09	-18.11%					
Atten Coeff	0.1481																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.031479	0.010723	1.2894	39.3774	11.0428	42.9503	44.4215	-1.56%	38.04%	-10.52%	-7.46%	3.57	-28.68%					
Atten Coeff	0.12																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034499	0.010723	1.4131	33.6158	2.3968	34.8143	35.5853	-15.96%	-70.04%	-27.47%	-25.86%	14.03	180.51%					
Atten Coeff	0.12																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034324	0.010723	1.4059	34.3677	3.8681	36.5785	37.3928	-14.08%	-51.65%	-23.79%	-22.10%	8.88	77.70%					
Atten Coeff	0.12																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034122	0.010723	1.3977	34.8323	4.5234	37.0738	37.9009	-12.92%	-43.46%	-22.76%	-21.04%	7.70	54.01%					
Atten Coeff	0.12																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.033894	0.010723	1.3883	35.2938	4.7477	37.2159	38.0465	-11.77%	-40.65%	-22.47%	-20.74%	7.43	48.68%					
Atten Coeff	0.12																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.032707	0.010723	1.3397	36.6372	5.4167	37.1743	38.0054	-8.41%	-32.29%	-22.55%	-20.82%	6.76	35.27%					
Atten Coeff	0.12																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.031479	0.010723	1.2894	37.5962	6.9995	37.9341	38.7854	-6.01%	-12.51%	-20.97%	-19.20%	5.37	7.43%					

Table A.12. ¹²³I FP -CIT 5:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			40	8	48										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02997	0.014589	1.2276	41.7178	7.926	48.1021	49.7323	4.29%	-0.92%	0.21%	3.61%	5.26	5.27%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029819	0.014589	1.2214	42.7956	10.133	50.7317	52.4609	6.99%	26.66%	5.69%	9.29%	4.22	-15.53%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029644	0.014589	1.2142	43.3839	11.3544	51.6766	53.4426	8.46%	41.93%	7.66%	11.34%	3.82	-23.58%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029445	0.014589	1.2061	43.9408	12.0155	52.1512	53.9356	9.85%	50.19%	8.65%	12.37%	3.66	-26.86%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028414	0.014589	1.1638	45.6491	14.2687	53.1233	54.9477	14.12%	78.36%	10.67%	14.47%	3.20	-36.02%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027347	0.014589	1.1201	46.6396	17.461	54.7213	56.6092	17.10%	118.26%	14.00%	17.94%	2.88	-46.35%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02997	0.014589	1.2276	40.1782	3.6477	42.4497	43.3923	0.45%	-54.40%	-11.56%	-9.60%	11.01	120.29%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029819	0.014589	1.2214	41.0153	5.7179	44.7787	45.7786	2.54%	-28.53%	-6.71%	-4.63%	7.17	43.46%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029644	0.014589	1.2142	41.4755	6.8611	45.6169	46.6382	3.69%	-14.24%	-4.96%	-2.84%	6.05	20.90%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029445	0.014589	1.2061	41.955	7.4651	46.0404	47.0725	4.89%	-6.69%	-4.08%	-1.93%	5.62	62.02%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028414	0.014589	1.1638	43.4175	9.5456	46.9173	47.9729	8.54%	19.32%	-2.26%	-0.06%	4.55	-9.03%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027347	0.014589	1.1201	44.2862	12.6144	48.3524	49.4453	10.72%	57.68%	0.73%	3.01%	3.51	-29.78%			

Table A.13. ¹²³I IBZM 1:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			84	84	168										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038468	0.010021	1.5757	73.7714	73.6867	146.803	151.9336	-12.18%	-12.26%	-12.62%	-9.56%	1.00	0.11%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038274	0.010021	1.5677	77.7528	78.2086	154.8244	160.2354	-7.44%	-6.89%	-7.84%	-4.62%	0.99	-0.58%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038049	0.010021	1.5585	79.4953	80.2366	157.5951	163.1028	-5.36%	-4.48%	-6.19%	-2.92%	0.99	-0.92%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.037794	0.010021	1.548	80.2753	81.7702	159.1376	164.6992	-4.43%	-2.65%	-5.28%	-1.96%	0.98	-1.83%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036471	0.010021	1.4938	82.6633	84.891	161.6658	167.3158	-1.59%	1.06%	-3.77%	-0.41%	0.97	-2.62%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035102	0.010021	1.4378	84.5832	87.4921	163.594	169.3113	0.69%	4.16%	-2.62%	0.78%	0.97	-3.32%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038468	0.010021	1.5757	65.1337	65.0446	129.5998	132.5676	-22.46%	-22.57%	-22.86%	-21.09%	1.00	0.14%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038274	0.010021	1.5677	68.5995	69.103	136.6977	139.828	-18.33%	-17.73%	-18.63%	-16.77%	0.99	-0.73%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.038049	0.010021	1.5585	70.1193	70.9238	139.1516	142.3381	-16.52%	-15.57%	-17.17%	-15.27%	0.99	-1.13%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.037794	0.010021	1.548	70.7311	72.3651	140.52	143.7378	-15.80%	-13.85%	-16.36%	-14.44%	0.98	-2.26%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.036471	0.010021	1.4938	72.7821	75.2304	142.7752	146.0446	-13.35%	-10.44%	-15.01%	-13.07%	0.97	-3.25%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.035102	0.010021	1.4378	74.4365	77.6295	144.5003	147.8092	-11.39%	-7.58%	-13.99%	-12.02%	0.96	-4.11%			

Table A.14. ¹²³I IBZM 1:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	pmv		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			84	84	168										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02917	0.013699	1.1948	88.706	88.5104	176.4347	182.6009	5.60%	5.37%	5.02%	8.69%	1.00	0.22%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029023	0.013699	1.1888	93.5122	94.0716	186.2798	192.79	11.32%	11.99%	10.88%	14.76%	0.99	-0.59%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028852	0.013699	1.1818	95.6119	96.7175	189.9188	196.5562	13.82%	15.14%	13.05%	17.00%	0.99	-1.14%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028659	0.013699	1.1739	96.6221	98.6828	191.9792	198.6886	15.03%	17.48%	14.27%	18.27%	0.98	-2.09%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027655	0.013699	1.1328	99.9662	102.8303	195.6668	202.505	19.01%	22.42%	16.47%	20.54%	0.97	-2.79%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026617	0.013699	1.0902	102.5003	106.3384	198.3641	205.2965	22.02%	26.59%	18.07%	22.20%	0.96	-3.61%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02917	0.013699	1.1948	78.3288	78.121	155.7595	159.3263	-6.75%	-7.00%	-7.29%	-5.16%	1.00	0.27%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.029023	0.013699	1.1888	82.5032	83.1205	164.4715	168.2378	-1.78%	-1.05%	-2.10%	0.14%	0.99	-0.74%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028852	0.013699	1.1818	84.3149	85.5136	167.695	171.5351	0.37%	1.80%	-0.18%	2.10%	0.99	-1.40%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028659	0.013699	1.1739	85.11	87.3595	169.5229	173.4048	1.32%	4.00%	0.91%	3.22%	0.97	-2.57%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027655	0.013699	1.1328	88.0064	91.1511	172.816	176.7733	4.77%	8.51%	2.87%	5.22%	0.97	-3.45%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026617	0.013699	1.0902	90.1826	94.3998	175.2338	179.2464	7.36%	12.38%	4.31%	6.69%	0.96	-4.47%			

Table A.15. ¹²³I IBZM 2:1 Activity Uptake Ratio, Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)												
Non-SC			84	42	126												
Atten Coeff	0.1481																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.036254	0.009956	1.4849	73.9724	36.2906	109.2666	113.077	-11.94%	-13.59%	-13.26%	-10.26%	2.04	1.92%					
Atten Coeff	0.1481																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.03607	0.009956	1.4774	77.2293	39.5292	115.2432	119.2436	-8.06%	-5.88%	-8.54%	-5.36%	1.95	-2.31%					
Atten Coeff	0.1481																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.035858	0.009956	1.4688	78.6175	41.0638	117.2906	121.3635	-6.41%	-2.23%	-6.91%	-3.68%	1.91	-4.27%					
Atten Coeff	0.1481																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.035618	0.009956	1.4589	79.4786	41.9566	118.3479	122.4585	-5.38%	-0.10%	-6.07%	-2.81%	1.89	-5.28%					
Atten Coeff	0.1481																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034371	0.009956	1.4078	81.9196	43.8974	120.0319	124.2022	-2.48%	4.52%	-4.74%	-1.43%	1.87	-6.69%					
Atten Coeff	0.1481																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.033081	0.009956	1.355	83.7275	45.9473	121.5307	125.755	-0.32%	9.40%	-3.55%	-0.19%	1.82	-8.89%					
Atten Coeff	0.12																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.036254	0.009956	1.4849	69.0376	28.304	96.4714	98.6636	-17.81%	-32.61%	-23.44%	-21.70%	2.44	21.96%					
Atten Coeff	0.12																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.03607	0.009956	1.4774	71.9165	31.1725	101.7425	104.0565	-14.39%	-25.78%	-19.25%	-17.42%	2.31	15.35%					
Atten Coeff	0.12																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.035858	0.009956	1.4688	73.1356	32.5428	103.5559	105.912	-12.93%	-22.52%	-17.81%	-15.94%	2.25	12.37%					
Atten Coeff	0.12																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.035618	0.009956	1.4589	73.891	33.3444	104.4943	106.8722	-12.03%	-20.61%	-17.07%	-15.18%	2.22	10.80%					
Atten Coeff	0.12																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.034371	0.009956	1.4078	76.1114	35.0361	105.999	108.412	-9.39%	-16.58%	-15.67%	-13.96%	2.17	8.62%					
Atten Coeff	0.12																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err					
0.033081	0.009956	1.355	77.7089	36.8988	107.3413	109.7859	-7.49%	-12.15%	-14.81%	-12.87%	2.11	5.30%					

Table A.16. ¹²³I IBZM 2:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			84	42	126										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02865	0.013585	1.1735	88.3875	44.3293	131.576	136.1414	5.22%	5.55%	4.43%	8.05%	1.99	-0.31%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028505	0.013585	1.1676	92.1562	48.551	138.9253	143.7503	9.71%	15.60%	10.26%	14.09%	1.90	-5.09%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028338	0.013585	1.1607	93.6638	50.7911	141.6452	146.5672	11.50%	20.93%	12.42%	16.32%	1.84	-7.80%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028148	0.013585	1.1529	94.5687	52.2221	143.1627	148.1392	12.58%	24.34%	13.62%	17.57%	1.81	-9.46%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027162	0.013585	1.1126	97.1623	55.8903	146.0272	151.107	15.67%	33.07%	15.89%	19.93%	1.74	-13.08%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026143	0.013585	1.0708	98.8151	59.6014	148.2597	153.4213	17.64%	41.91%	17.67%	21.76%	1.66	-17.10%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.02865	0.013585	1.1735	82.3946	34.7688	116.148	118.7885	-1.91%	-17.21%	-7.82%	-5.72%	2.37	18.49%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028505	0.013585	1.1676	85.6792	38.5557	122.6523	125.4434	2.00%	-8.20%	-2.66%	-0.44%	2.22	11.11%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028338	0.013585	1.1607	86.9468	40.6085	125.062	127.9094	3.51%	-3.31%	-0.74%	1.52%	2.14	7.05%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028148	0.013585	1.1529	87.6906	41.9386	126.409	129.288	4.39%	-0.15%	0.32%	2.61%	2.09	4.55%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027162	0.013585	1.1126	89.889	45.332	128.9691	131.9086	7.01%	7.93%	2.36%	4.69%	1.98	-0.85%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026143	0.013585	1.0708	91.1651	48.8732	130.9719	133.9595	8.53%	16.36%	3.95%	6.32%	1.87	-6.73%			

Table A.17. ¹²³I IBZM 5:1 Activity Uptake Ratio Primary and High-Energy Emissions

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			84	16.8	100.8										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0349	0.009972	1.4295	74.0972	13.6571	86.006	88.9199	-11.79%	-17.52%	-14.68%	-11.79%	5.35	6.94%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034723	0.009972	1.4222	76.9210	16.3267	90.7325	93.8164	-8.43%	-2.82%	-9.99%	-6.93%	4.71	-5.77%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034519	0.009972	1.4139	78.0987	17.5662	92.3549	95.4966	-7.03%	4.56%	-8.38%	-5.26%	4.45	-11.08%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034288	0.009972	1.4044	79.0095	18.0764	93.0617	96.231	-5.94%	7.60%	-7.68%	-4.53%	4.37	-12.58%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.033087	0.009972	1.3552	81.4713	19.3347	94.2651	97.4788	-3.01%	15.09%	-6.48%	-3.29%	4.21	-15.73%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.031845	0.009972	1.3044	83.1912	21.0889	95.5084	98.7704	-0.96%	25.53%	-5.25%	-2.01%	3.94	-21.10%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0349	0.009972	1.4295	71.3836	6.2633	75.901	77.5856	-15.02%	-62.72%	-24.70%	-23.03%	11.40	127.94%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034723	0.009972	1.4222	73.9119	8.4187	80.0844	81.8676	-12.01%	-49.89%	-20.55%	-18.78%	8.78	75.59%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034519	0.009972	1.4139	74.9524	9.5194	81.5219	83.3399	-10.77%	-43.34%	-19.13%	-17.32%	7.87	57.47%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.034288	0.009972	1.4044	75.7949	9.939	82.1498	83.9827	-9.77%	-40.84%	-18.50%	-16.68%	7.63	52.52%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.033087	0.009972	1.3552	78.1034	10.9534	83.2276	85.0865	-7.02%	-34.80%	-17.43%	-15.59%	7.13	42.61%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.031845	0.009972	1.3044	79.6431	12.5305	84.3424	86.2296	-5.19%	-25.41%	-16.33%	-14.45%	6.36	27.12%			

Table A.18. ¹²³I IBZM 5:1 Activity Uptake Ratio, Primary 159 keV Emissions

Data Type	prmy		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			84	16.8	100.8										
Atten Coeff	0.1481														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028295	0.013578	1.1589	88.2037	17.8278	103.8001	107.3236	5.00%	6.12%	2.96%	6.47%	4.95	-1.05%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028151	0.013578	1.1531	91.3529	21.2474	109.7024	113.4398	8.75%	26.47%	8.83%	12.54%	4.30	-14.01%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027986	0.013578	1.1463	92.5082	23.2462	111.9157	115.736	10.13%	38.37%	11.03%	14.82%	3.98	-20.41%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027798	0.013578	1.1386	93.3519	24.3592	113.0501	116.9132	11.13%	45.00%	12.15%	15.99%	3.83	-23.35%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026825	0.013578	1.0988	95.4726	27.7555	115.4039	119.3592	13.66%	65.21%	14.49%	18.41%	3.44	-31.20%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.025818	0.013578	1.0575	96.5843	31.6064	117.3849	121.4209	14.98%	88.13%	16.45%	20.46%	3.06	-38.88%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028295	0.013578	1.1589	84.8406	8.7655	91.6069	93.644	1.00%	-47.82%	-9.12%	-7.10%	9.68	93.58%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.028151	0.013578	1.1531	87.5939	11.8246	96.8318	98.993	4.28%	-29.62%	-3.94%	-1.79%	7.41	48.16%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027986	0.013578	1.1463	88.5376	13.6751	98.7938	101.0032	5.40%	-18.60%	-1.99%	0.20%	6.47	29.49%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.027798	0.013578	1.1386	89.2527	14.6981	99.8017	102.036	6.25%	-12.51%	-0.99%	1.23%	6.07	21.45%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.026825	0.013578	1.0988	91.0088	17.8702	101.909	104.1973	8.34%	6.37%	1.10%	3.37%	5.09	1.86%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.025818	0.013578	1.0575	91.7317	21.6051	103.6889	106.0247	9.20%	28.60%	2.87%	5.18%	4.25	-15.08%			

Table A.19. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			126.4	126.4	252.8										
Atten Coeff	0.1526														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0092817	0.005544	0.36018	134.2712	134.4403	267.5392	277.4698	6.23%	6.36%	5.83%	9.76%	0.9967	-0.13%			
Atten Coeff	0.1526														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0092347	0.005544	0.37825	142.7976	140.5125	281.7836	292.2429	12.97%	11.16%	11.47%	15.60%	1.0163	1.63%			
Atten Coeff	0.1526														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0091805	0.005544	0.37603	144.9782	144.4505	287.2131	297.8739	14.70%	14.28%	13.61%	17.83%	1.0037	0.37%			
Atten Coeff	0.1526														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0091189	0.005544	0.37351	146.3315	146.8877	290.4317	301.212	15.77%	16.21%	14.89%	19.15%	0.9962	-0.38%			
Atten Coeff	0.1526														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0087996	0.005544	0.36043	150.3009	151.2339	296.3029	307.3011	18.91%	19.65%	17.21%	21.56%	0.9938	-0.62%			
Atten Coeff	0.1526														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0084693	0.005544	0.3469	153.1292	152.903	298.444	309.5217	21.15%	20.97%	18.06%	22.44%	1.0015	0.15%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0092817	0.005544	0.36018	116.2045	116.3931	231.5815	236.8845	-8.07%	-7.92%	-8.39%	-6.30%	0.9984	-0.16%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0092347	0.005544	0.37825	123.882	121.3855	243.9431	249.5292	-1.99%	-3.97%	-3.50%	-1.29%	1.0206	2.06%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0091805	0.005544	0.37603	125.5768	125.0115	248.6653	254.3596	-0.65%	-1.10%	-1.64%	0.62%	1.0045	0.45%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0091189	0.005544	0.37351	126.6328	127.2562	251.467	257.2254	0.18%	0.68%	-0.53%	1.75%	0.9951	-0.49%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0087996	0.005544	0.36043	130.0647	131.1101	256.6061	262.4822	2.90%	3.73%	1.51%	3.83%	0.9920	-0.80%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0084693	0.005544	0.3469	132.6866	132.4578	258.4916	264.4109	4.97%	4.79%	2.25%	4.59%	1.0017	0.17%			

Table A.20. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)												
Non-SC			126.4	63.2	189.6												
Atten Coeff	0.1526																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0086145	0.004774	0.35285	134.7488	66.975	199.9628	207.3274		6.61%	5.97%	5.47%	9.35%	2.012	0.60%				
Atten Coeff	0.1526																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0085708	0.004774	0.35106	141.6489	71.0809	210.6389	218.3983		12.06%	12.47%	11.10%	15.19%	1.993	-0.36%				
Atten Coeff	0.1526																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0085205	0.004774	0.349	143.1434	74.1871	214.7109	222.6255		13.25%	17.38%	13.24%	17.42%	1.929	-3.53%				
Atten Coeff	0.1526																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0084634	0.004774	0.34666	143.9434	76.3289	217.1691	225.178		13.88%	20.77%	14.54%	18.76%	1.886	-5.71%				
Atten Coeff	0.1526																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0081671	0.004774	0.33452	146.3141	80.2091	221.3619	229.531		15.75%	26.91%	16.75%	21.06%	1.824	-8.79%				
Atten Coeff	0.1526																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0078605	0.004774	0.32197	147.8364	82.1916	222.9692	231.2		16.96%	30.05%	17.60%	21.94%	1.799	-10.07%				
Atten Coeff	0.12																
Striata ROI Size	13 13																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0086145	0.004774	0.35285	124.2577	50.3542	173.0693	177.0014		-1.69%	-20.33%	-8.72%	-6.64%	2.468	23.38%				
Atten Coeff	0.12																
Striata ROI Size	15 15																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0085708	0.004774	0.35106	130.5595	53.6044	182.334	186.4774		3.29%	-15.18%	-3.83%	-1.65%	2.436	21.78%				
Atten Coeff	0.12																
Striata ROI Size	17 17																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0085205	0.004774	0.349	131.674	56.4906	185.8762	190.103		4.17%	-10.62%	-1.96%	0.27%	2.331	16.55%				
Atten Coeff	0.12																
Striata ROI Size	19 19																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0084634	0.004774	0.34666	132.2186	58.5083	188.0181	192.2956		4.60%	-7.42%	-0.83%	1.42%	2.260	12.99%				
Atten Coeff	0.12																
Striata ROI Size	26 26																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0081671	0.004774	0.33452	134.1143	62.0895	191.6893	196.0533		6.10%	-1.76%	1.10%	3.40%	2.160	8.00%				
Atten Coeff	0.12																
Striata ROI Size	32 32																
Bkgd ROI Size	13 13																
Units:	microCi																
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor		LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err				
0.0078605	0.004774	0.32197	135.393	63.9053	193.1059	197.5035		7.11%	1.12%	1.85%	4.17%	2.119	5.93%				

Table A.21. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)										
Non-SC			126.40	25.28	151.68										
Atten Coeff	0.1526														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0081116	0.004338	0.33225	135.0529	26.5133	158.0251	163.7083	6.85%	4.86%	4.18%	7.93%	5.094	1.86%			
Atten Coeff	0.1526														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0080705	0.004338	0.33057	140.9830	29.4446	166.4949	172.4927	11.54%	16.47%	9.77%	13.72%	4.788	-4.24%			
Atten Coeff	0.1526														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0080232	0.004338	0.32863	142.0722	32.0572	169.699	175.8252	12.40%	26.81%	11.88%	15.92%	4.432	-11.36%			
Atten Coeff	0.1526														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0079694	0.004338	0.32642	142.5470	34.0285	171.6355	177.8413	12.77%	34.61%	13.16%	17.25%	4.189	-16.22%			
Atten Coeff	0.1526														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0076903	0.004338	0.315	143.9812	37.6585	174.5309	180.8582	13.91%	48.97%	15.07%	19.24%	3.823	-23.53%			
Atten Coeff	0.1526														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0074017	0.004338	0.30317	144.7553	39.8597	175.7526	182.1342	14.52%	57.67%	15.87%	20.08%	3.632	-27.37%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0081116	0.004338	0.33225	129.1048	10.7460	136.7296	139.7626	2.14%	-57.49%	-9.86%	-7.86%	12.014	140.28%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0080705	0.004338	0.33057	134.5862	12.9555	144.0804	147.2818	6.48%	-48.75%	-5.01%	-2.90%	10.388	107.77%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0080232	0.004338	0.32863	135.3580	15.4025	146.8693	150.1397	7.09%	-39.07%	-3.17%	-1.02%	8.788	75.76%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0079694	0.004338	0.32642	135.6017	17.2900	148.5593	151.8726	7.28%	-31.61%	-2.06%	0.13%	7.843	56.86%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0076903	0.004338	0.315	136.5961	20.7337	151.0983	154.4772	8.07%	-17.96%	-0.38%	1.84%	6.588	31.76%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.0074017	0.004338	0.30317	137.1026	22.8558	152.1781	155.5865	8.47%	-9.59%	0.33%	2.58%	5.999	19.97%			

APPENDIX B

Listed below are the activity uptake estimates determined using the modified conjugate imaging technique and triple energy window scatter correction for various ROI sizes and two attenuation coefficients. This data represents the complete set of data generated from the scatter-corrected emission only and scattered-corrected emission and camera sensitivity simulations. Tables B.1 - B.55 list the data from the study using scatter-corrected emission data only. Tables B.56 - B.109 list the data from the study using scatter-corrected emission and camera sensitivity data. The activity estimates listed for the ^{123}I radiopharmaceuticals were determined using the primary 159 keV and high-energy emissions.

Table B.1. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13_h16			
SC			100	100	200				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								175-180			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.083855	0.055389	3.4347	74.1404	74.3441	143.9079	148.9374	-25.86%	-25.66%	-28.05%	-25.53%	0.997	-0.27%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.08343	0.055389	3.4173	76.8105	81.1271	150.3839	155.6394	-23.19%	-18.67%	-24.81%	-22.18%	0.947	-5.32%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.082941	0.055389	3.3973	79.6584	84.1088	153.3662	158.7259	-20.34%	-15.89%	-23.32%	-20.64%	0.947	-5.29%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.082384	0.055389	3.3745	81.6334	85.5625	154.2168	159.6063	-18.37%	-14.44%	-22.89%	-20.20%	0.954	-4.59%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0795	0.055389	3.2563	89.0941	91.3286	155.9446	161.3946	-10.91%	-8.67%	-22.03%	-19.30%	0.976	-2.45%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076516	0.055389	3.1341	92.8595	103.8389	161.1257	166.7555	-7.14%	3.84%	-19.44%	-16.62%	0.894	-10.57%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.083855	0.055389	3.4347	65.4388	65.6434	127.037	129.9461	-34.56%	-34.36%	-36.48%	-35.03%	0.997	-0.31%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.08343	0.055389	3.4173	67.3856	72.064	132.7685	135.8087	-32.61%	-27.94%	-33.62%	-32.10%	0.935	-6.49%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.082941	0.055389	3.3973	69.8995	74.7177	135.4111	138.5118	-30.10%	-25.28%	-32.29%	-30.74%	0.936	-6.45%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.082384	0.055389	3.3745	71.6938	75.9796	136.1761	139.2943	-28.31%	-24.02%	-31.91%	-30.35%	0.944	-5.64%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0795	0.055389	3.2563	78.5222	80.998	137.7632	140.9178	-21.48%	-19.00%	-31.12%	-29.54%	0.969	-3.06%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076516	0.055389	3.1341	81.0314	93.1266	142.4094	145.6697	-18.97%	-6.87%	-28.80%	-27.17%	0.870	-12.99%

Table B.2. ¹²³I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13_hi6			
SC			100	50	150					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083841	0.055674	3.4341	73.3427	37.3158	105.2874	108.9415	-26.66%	-25.37%	-29.81%	-27.37%	1.965	-1.73%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083416	0.055674	3.4167	75.0037	43.3171	110.0469	113.875	-25.00%	-13.37%	-26.64%	-24.08%	1.732	-13.42%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082927	0.055674	3.3967	77.3796	46.1569	112.5789	116.4972	-22.62%	-7.69%	-24.95%	-22.34%	1.676	-16.18%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082371	0.055674	3.3739	79.2624	47.3619	113.1254	117.063	-20.74%	-5.28%	-24.58%	-21.96%	1.674	-16.32%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079487	0.055674	3.2558	86.3	53.3679	115.0158	119.0216	-13.70%	6.74%	-23.32%	-20.65%	1.617	-19.15%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076503	0.055674	3.1336	89.6957	66.1626	120.296	124.4949	-10.30%	32.33%	-19.80%	-17.00%	1.356	-32.22%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083841	0.055674	3.4341	68.3244	29.3643	92.9348	95.0481	-31.68%	-41.27%	-38.04%	-36.63%	2.327	16.34%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083416	0.055674	3.4167	69.3553	35.1171	97.1503	99.3646	-30.64%	-29.77%	-35.23%	-33.76%	1.975	-1.25%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082927	0.055674	3.3967	71.4172	37.6787	99.3951	101.6617	-28.58%	-24.64%	-33.74%	-32.23%	1.895	-5.23%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082371	0.055674	3.3739	73.1442	38.7061	99.8912	102.1693	-26.86%	-22.59%	-33.41%	-31.89%	1.890	-5.51%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079487	0.055674	3.2558	79.5294	44.0023	101.6224	103.9414	-20.47%	-12.00%	-32.25%	-30.71%	1.807	-9.63%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076503	0.055674	3.1336	81.64	56.4503	106.358	108.7904	-18.36%	12.90%	-29.09%	-27.47%	1.446	-27.69%	

Table B.3. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13_hi6			
SC			100	20	120					130-142			
Atten Coeff	0.1481									142-175			
Sriata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083888	0.055795	3.436	72.8538	15.0893	80.9689	83.7188	-27.15%	-24.55%	-32.53%	-30.23%	4.828	-3.44%	
Atten Coeff	0.1481												
Sriata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083463	0.055795	3.4186	73.9078	20.6158	84.6647	87.5627	-26.09%	3.08%	-29.45%	-27.03%	3.585	-28.30%	
Atten Coeff	0.1481												
Sriata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082973	0.055795	3.3986	76.0095	23.3575	86.8765	89.8576	-23.99%	16.79%	-27.60%	-25.12%	3.254	-34.92%	
Atten Coeff	0.1481												
Sriata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082417	0.055795	3.3758	77.8586	24.3888	87.1702	90.1633	-22.14%	21.94%	-27.36%	-24.86%	3.192	-36.15%	
Atten Coeff	0.1481												
Sriata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079531	0.055795	3.2576	84.6818	30.486	88.771	91.8304	-15.32%	52.43%	-26.02%	-23.47%	2.778	-44.45%	
Atten Coeff	0.1481												
Sriata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076546	0.055795	3.1353	87.8394	43.4443	94.1371	97.4051	-12.16%	117.22%	-21.55%	-18.83%	2.022	-59.56%	
Atten Coeff	0.12												
Sriata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083888	0.055795	3.436	70.0469	7.5885	71.4497	73.0394	-29.95%	-62.06%	-40.46%	-39.13%	9.231	84.61%	
Atten Coeff	0.12												
Sriata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.083463	0.055795	3.4186	70.5275	12.9352	74.7286	76.4043	-29.47%	-35.32%	-37.73%	-36.33%	5.452	9.05%	
Atten Coeff	0.12												
Sriata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082973	0.055795	3.3986	72.3284	15.4279	76.6916	78.4157	-27.67%	-22.86%	-36.09%	-34.65%	4.688	-6.24%	
Atten Coeff	0.12												
Sriata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.082417	0.055795	3.3758	74.0385	16.2884	76.9648	78.6961	-25.96%	-18.56%	-35.86%	-34.42%	4.545	-9.09%	
Atten Coeff	0.12												
Sriata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079531	0.055795	3.2576	80.202	21.6944	78.4421	80.2132	-19.80%	8.47%	-34.63%	-33.16%	3.697	-26.06%	
Atten Coeff	0.12												
Sriata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076546	0.055795	3.1353	82.0596	34.3267	83.2554	85.1491	-17.94%	71.63%	-30.62%	-29.04%	2.391	-52.19%	

Table B.4. ^{123}I FP -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018554	0.012061	0.75999	30.1034	29.8497	59.0006	61.0626	-24.74%	-25.38%	-26.25%	-23.67%	1.008	0.85%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01846	0.012061	0.75614	31.4976	31.9687	61.8299	63.9898	-21.26%	-20.08%	-22.71%	-20.01%	0.985	-1.47%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018352	0.012061	0.7517	32.3678	32.8031	62.8247	65.0203	-19.08%	-17.99%	-21.47%	-18.72%	0.987	-1.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018229	0.012061	0.74666	32.804	33.3602	63.1126	65.3183	-17.99%	-16.60%	-21.11%	-18.35%	0.983	-1.67%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017591	0.012061	0.72052	34.5348	34.7787	63.5662	65.7878	-13.66%	-13.05%	-20.54%	-17.77%	0.993	-0.70%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016931	0.012061	0.69347	35.2913	37.6621	64.6479	66.9071	-11.77%	-5.84%	-19.19%	-16.37%	0.937	-6.29%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018554	0.012061	0.75999	26.6026	26.3243	52.0852	53.2779	-33.49%	-34.19%	-34.89%	-33.40%	1.011	1.06%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01846	0.012061	0.75614	27.7623	28.2734	54.5875	55.8375	-30.59%	-29.32%	-31.77%	-30.20%	0.982	-1.81%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018352	0.012061	0.7517	28.5385	29.0082	55.4697	56.7399	-28.65%	-27.48%	-30.66%	-29.08%	0.984	-1.62%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018229	0.012061	0.74666	28.9116	29.5192	55.7271	57.0032	-27.72%	-26.20%	-30.34%	-28.75%	0.979	-2.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017591	0.012061	0.72052	30.4873	30.7643	56.1421	57.4277	-23.78%	-23.09%	-29.82%	-28.22%	0.991	-0.90%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016931	0.012061	0.69347	30.9559	33.5704	57.1135	58.4213	-22.61%	-16.07%	-28.61%	-26.97%	0.922	-7.79%	

Table B.5. ¹²³I FP -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	20	60					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018543	0.012179	0.75953	29.7839	15.0413	43.6984	45.2148	-25.54%	-24.79%	-27.17%	-24.64%	1.980	-0.99%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018449	0.012179	0.75569	30.7775	16.846	45.8318	47.4246	-23.06%	-15.77%	-23.61%	-20.96%	1.827	-8.65%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018341	0.012179	0.75125	31.4642	17.6212	46.5998	48.22	-21.34%	-11.89%	-22.33%	-19.63%	1.786	-10.72%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018218	0.012179	0.74622	31.8697	18.0746	46.8213	48.4496	-20.33%	-9.63%	-21.96%	-19.25%	1.763	-11.84%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01758	0.012179	0.72009	33.4539	19.5448	47.3366	48.9637	-16.37%	-2.26%	-21.11%	-18.36%	1.712	-14.42%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01692	0.012179	0.69306	34.0604	22.5497	48.4749	50.1648	-14.85%	12.75%	-19.21%	-16.39%	1.510	-24.48%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018543	0.012179	0.75953	27.7565	11.8152	38.5727	39.4497	-30.61%	-40.92%	-35.71%	-34.25%	2.349	17.46%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018449	0.012179	0.75569	28.5526	13.4954	40.4609	41.3822	-28.62%	-32.52%	-32.57%	-31.03%	2.116	5.79%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018341	0.012179	0.75125	29.1534	14.1905	41.1416	42.0788	-27.12%	-29.05%	-31.43%	-29.87%	2.054	2.72%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018218	0.012179	0.74622	29.5059	14.603	41.3406	42.2626	-26.24%	-26.99%	-31.10%	-29.53%	2.021	1.03%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01758	0.012179	0.72009	30.9313	15.913	41.8104	42.7636	-22.67%	-20.44%	-30.32%	-28.73%	1.944	-2.81%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01692	0.012179	0.69306	31.2378	18.855	42.8328	43.8112	-21.91%	-5.73%	-28.61%	-26.98%	1.657	-17.16%	

Table B.6. ¹²³I FP -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	8	48					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018551	0.012238	0.75985	29.5903	6.153	34.1873	35.3495	-26.02%	-23.09%	-28.78%	-26.36%	4.8091	-3.82%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018457	0.012238	0.756	30.342	7.7689	35.8746	37.0998	-24.15%	-2.89%	-25.26%	-22.71%	3.9056	-21.89%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018349	0.012238	0.75157	30.9156	8.5086	36.5279	37.7777	-22.71%	6.36%	-23.90%	-21.30%	3.6335	-27.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018226	0.012238	0.74653	31.2977	8.8995	36.6563	37.9117	-21.76%	11.24%	-23.63%	-21.02%	3.5168	-29.66%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017588	0.012238	0.72039	32.7605	10.4316	37.1237	38.3991	-18.10%	30.40%	-22.66%	-20.00%	3.1405	-37.19%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016927	0.012238	0.69335	33.2605	13.5051	38.2639	39.5886	-16.85%	68.81%	-20.28%	-17.53%	2.4628	-50.74%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018551	0.012238	0.75985	28.4471	3.1069	30.1698	30.8411	-28.88%	-61.16%	-37.15%	-35.75%	9.1561	83.12%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018457	0.012238	0.756	29.0236	4.6257	31.6641	32.3725	-27.44%	-42.18%	-34.03%	-32.56%	6.2744	25.49%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018349	0.012238	0.75157	29.5164	5.2972	32.2437	32.9665	-26.21%	-33.79%	-32.83%	-31.32%	5.5721	11.44%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018226	0.012238	0.74653	29.8518	5.6507	32.3607	33.0868	-25.37%	-29.37%	-32.58%	-31.07%	5.2828	5.66%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017588	0.012238	0.72039	31.1511	7.0335	32.789	33.527	-22.12%	-12.08%	-31.69%	-30.15%	4.429	-11.42%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016927	0.012238	0.69335	31.3444	10.0543	33.8141	34.58	-21.64%	25.68%	-29.55%	-27.96%	3.1175	-37.65%	

Table B.7. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC			84	84	168					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016821	0.011411	0.68901	64.0534	63.095	125.8882	130.2678	-23.75%	-24.69%	-25.07%	-22.45%	1.0152	1.52%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016736	0.011411	0.68552	67.4984	66.7136	132.1326	136.7505	-19.64%	-20.58%	-21.35%	-18.60%	1.0118	1.18%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016638	0.011411	0.6815	68.9822	67.9791	134.0298	138.714	-17.88%	-19.07%	-20.22%	-17.43%	1.0148	1.48%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016526	0.011411	0.67692	69.457	69.0986	134.7396	139.4486	-17.31%	-17.74%	-19.80%	-16.99%	1.0052	0.52%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015948	0.011411	0.65322	71.4023	70.9995	135.6249	140.3648	-15.00%	-15.48%	-19.27%	-16.45%	1.0057	0.57%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015349	0.011411	0.6287	72.0849	74.1379	136.7005	141.478	-14.18%	-11.74%	-18.63%	-15.79%	0.9723	-2.77%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016821	0.011411	0.68901	56.6445	55.6035	111.1346	113.6795	-32.57%	-33.81%	-33.85%	-32.33%	1.0187	1.87%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016736	0.011411	0.68552	59.6739	58.8234	116.6579	119.3293	-28.96%	-29.97%	-30.56%	-28.97%	1.0145	1.45%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016638	0.011411	0.6815	61.0143	59.9205	118.3393	121.0492	-27.36%	-28.67%	-29.56%	-27.95%	1.0183	1.83%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016526	0.011411	0.67692	61.3696	60.9813	118.97	121.6943	-26.94%	-27.40%	-29.18%	-27.56%	1.0064	0.64%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015948	0.011411	0.65322	63.1082	62.6863	119.7705	122.5131	-24.87%	-25.37%	-28.71%	-27.08%	1.0067	0.67%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015349	0.011411	0.6287	63.4789	65.7538	120.737	123.5017	-24.43%	-21.72%	-28.13%	-26.49%	0.9654	-3.46%	

Table B.8. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I ₀₁₃ hi6			
SC			84	42	126					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016966	0.011368	0.69493	63.3536	31.6946	93.8127	97.0679	-24.58%	-24.06%	-25.55%	-22.96%	1.986	-0.68%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01688	0.011368	0.69141	65.9479	34.8415	98.6549	102.0814	-21.49%	-17.04%	-21.70%	-18.98%	1.893	-5.36%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016781	0.011368	0.68735	67.0297	35.9759	100.1213	103.5997	-20.20%	-14.34%	-20.54%	-17.78%	1.863	-6.84%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016669	0.011368	0.68274	67.4174	36.87	100.6216	104.1186	-19.74%	-12.21%	-20.14%	-17.37%	1.829	-8.57%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015481	0.011368	0.63411	69.2871	42.1595	102.4668	106.0337	-17.52%	0.38%	-18.68%	-15.85%	1.643	-17.83%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016966	0.011368	0.69493	59.0492	25.0364	82.8113	84.694	-29.70%	-40.39%	-34.28%	-32.78%	2.359	17.93%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01688	0.011368	0.69141	61.3065	27.6816	87.095	89.0769	-27.02%	-34.09%	-30.86%	-29.30%	2.215	10.74%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016781	0.011368	0.68735	62.2631	28.6898	88.395	90.4071	-25.88%	-31.69%	-29.85%	-28.25%	2.170	8.51%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016669	0.011368	0.68274	62.554	29.5385	88.8403	90.8633	-25.53%	-29.67%	-29.49%	-27.89%	2.118	5.89%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016095	0.011368	0.65884	63.9119	31.3297	89.5714	91.6122	-23.91%	-25.41%	-28.91%	-27.29%	2.040	2.00%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015481	0.011368	0.63411	63.8822	34.6356	90.505	92.5698	-23.95%	-17.53%	-28.17%	-26.53%	1.844	-7.78%	

Table B.9. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I _{o13} hi6			
SC			84	16.8	100.8					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016996	0.011466	0.69614	62.9423	13.1847	73.9054	76.4163	-25.07%	-21.52%	-26.68%	-24.19%	4.7739	-4.52%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016909	0.011466	0.69261	65.0286	15.7322	77.0708	80.5263	-22.59%	-6.36%	-22.75%	-20.11%	4.1335	-17.33%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01681	0.011466	0.68855	65.8733	16.7908	79.0397	81.7386	-21.58%	-0.05%	-21.59%	-18.91%	3.9232	-21.54%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016696	0.011466	0.68393	66.215	17.5529	79.4079	82.1224	-21.17%	4.48%	-21.22%	-18.53%	3.7723	-24.55%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016113	0.011466	0.65996	67.5731	19.5705	80.1476	82.8941	-19.56%	16.49%	-20.49%	-17.76%	3.4528	-30.94%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015508	0.011466	0.63521	67.6876	23.0272	81.242	84.0381	-19.42%	37.07%	-19.40%	-16.63%	2.9395	-41.21%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016996	0.011466	0.69614	60.4997	6.7052	65.2229	66.6745	-27.98%	-60.09%	-35.29%	-33.85%	9.0228	80.46%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016909	0.011466	0.69261	62.2968	9.0088	68.7321	70.2676	-25.84%	-46.38%	-31.81%	-30.29%	6.915	38.30%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01681	0.011466	0.68855	63.026	9.9662	69.7693	71.3301	-24.97%	-40.68%	-30.78%	-29.24%	6.324	26.48%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016696	0.011466	0.68393	63.2842	10.6899	70.0955	71.6684	-24.66%	-36.37%	-30.46%	-28.90%	5.92	18.40%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016113	0.011466	0.65996	64.4288	12.567	70.7706	72.3594	-23.30%	-25.20%	-29.79%	-28.21%	5.1268	2.54%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015508	0.011466	0.63521	64.1988	16.0096	71.757	73.3751	-23.57%	-4.70%	-28.81%	-27.21%	4.01	-19.80%	

Table B.10. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13_hi6			
SC			100	100	200				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								185-190			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.093488	0.055585	3.8293	78.301	78.241	152.8461	158.1879	-21.70%	-21.76%	-23.58%	-20.91%	1.001	0.08%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.093014	0.055585	3.8099	81.3178	84.7305	159.873	165.4602	-18.68%	-15.27%	-20.06%	-17.27%	0.960	-4.03%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.092469	0.055585	3.7875	84.2377	87.059	162.2464	167.9166	-15.76%	-12.94%	-18.88%	-16.04%	0.968	-3.24%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.091848	0.055585	3.7621	85.3354	89.9252	163.7306	169.4525	-14.66%	-10.07%	-18.13%	-15.27%	0.949	-5.10%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.088632	0.055585	3.6304	92.0973	94.9595	164.2006	169.9391	-7.90%	-5.04%	-17.90%	-15.03%	0.970	-3.01%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.085306	0.055585	3.4941	97.314	105.3731	169.1096	175.0191	-2.69%	5.37%	-15.45%	-12.49%	0.924	-7.65%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.093488	0.055585	3.8293	69.1332	69.0571	134.9222	138.0119	-30.87%	-30.94%	-32.54%	-30.99%	1.001	0.11%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.093014	0.055585	3.8099	71.4514	75.1528	141.1394	144.3713	-28.55%	-24.85%	-29.43%	-27.81%	0.951	-4.93%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.092469	0.055585	3.7875	74.1026	77.1538	143.2424	146.5225	-25.90%	-22.85%	-28.38%	-26.74%	0.960	-3.95%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.091848	0.055585	3.7621	74.8964	79.8889	144.5675	147.8778	-25.10%	-20.11%	-27.72%	-26.06%	0.938	-6.25%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.088632	0.055585	3.6304	81.0963	84.2661	145.044	148.3654	-18.90%	-15.73%	-27.48%	-25.82%	0.962	-3.76%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.085306	0.055585	3.4941	85.2586	94.1516	149.4421	152.8638	-14.74%	-5.85%	-25.28%	-23.57%	0.906	-9.45%

Table B.11. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I ¹³¹ h ¹⁶			
SC			100	50	150					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.093381	0.055876	3.8249	77.8912	38.4021	111.6741	115.5475	-22.11%	-23.20%	-25.55%	-22.97%	2.0283	1.42%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092908	0.055876	3.8055	79.9107	44.011	116.9005	120.9638	-20.09%	-11.98%	-22.07%	-19.36%	1.8157	-9.22%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092363	0.055876	3.7832	82.3241	46.1967	118.8787	123.0121	-17.68%	-7.61%	-20.75%	-17.99%	1.782	-10.90%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.091743	0.055876	3.7578	83.3339	48.8061	120.1026	124.2817	-16.67%	-2.39%	-19.93%	-17.15%	1.7074	-14.63%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.088531	0.055876	3.6262	89.765	53.9384	120.6646	124.8653	-10.24%	7.88%	-19.56%	-16.76%	1.6642	-16.79%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.085208	0.055876	3.4901	94.6035	64.7687	125.8035	130.1912	-5.40%	29.54%	-16.13%	-13.21%	1.4606	-26.97%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.093381	0.055876	3.8249	72.6783	29.9806	98.5671	100.8071	-27.32%	-40.04%	-34.29%	-32.80%	2.4242	21.21%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092908	0.055876	3.8055	74.1014	35.3105	103.1947	105.5449	-25.90%	-29.38%	-31.20%	-29.64%	2.0986	4.93%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092363	0.055876	3.7832	76.2669	37.2221	104.9487	107.3395	-23.73%	-25.56%	-30.03%	-28.44%	2.049	2.45%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.091743	0.055876	3.7578	76.9978	39.7146	106.0444	108.4621	-23.00%	-20.57%	-29.30%	-27.69%	1.9388	-3.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.088531	0.055876	3.6262	82.8594	44.2201	106.6015	109.0331	-17.14%	-11.56%	-28.93%	-27.31%	1.8738	-6.31%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.085208	0.055876	3.4901	86.6026	54.556	111.2048	113.7461	-13.40%	9.11%	-25.86%	-24.17%	1.5874	-20.63%	

Table B.12. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13_hi6			
SC			100	20	120					130-142			
Atten Coeff	0.1481									142-175			
Sriata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.093145	0.056239	3.8152	77.6741	14.5244	86.0261	88.9404	-22.33%	-27.38%	-28.31%	-25.88%	5.35	6.96%	
Atten Coeff	0.1481												
Sriata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092674	0.056239	3.7959	79.0989	19.6169	90.2117	93.2902	-20.90%	-1.92%	-24.82%	-22.26%	4.03	-19.36%	
Atten Coeff	0.1481												
Sriata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.09213	0.056239	3.7736	81.211	21.7375	91.8099	94.9491	-18.79%	8.69%	-23.49%	-20.88%	3.74	-25.28%	
Atten Coeff	0.1481												
Sriata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.091512	0.056239	3.7483	82.1509	24.2323	92.821	96.0032	-17.85%	21.16%	-22.65%	-20.00%	3.39	-32.20%	
Atten Coeff	0.1481												
Sriata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.088308	0.056239	3.6171	88.3663	29.558	93.3387	96.5499	-11.63%	47.79%	-22.22%	-19.54%	2.99	-40.21%	
Atten Coeff	0.1481												
Sriata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.084993	0.056239	3.4813	93.171	40.4673	98.6211	102.0357	-6.83%	102.34%	-17.82%	-14.97%	2.30	-53.95%	
Atten Coeff	0.12												
Sriata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.093145	0.056239	3.8152	74.831	6.5576	75.9078	77.5924	-25.17%	-67.21%	-36.74%	-35.34%	11.41	128.23%	
Atten Coeff	0.12												
Sriata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.092674	0.056239	3.7959	75.7197	11.4397	79.6188	81.3988	-24.28%	-42.80%	-33.65%	-32.17%	6.62	32.38%	
Atten Coeff	0.12												
Sriata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.09213	0.056239	3.7736	77.5942	13.3174	81.0381	82.8534	-22.41%	-33.41%	-32.47%	-30.96%	5.83	16.53%	
Atten Coeff	0.12												
Sriata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.091512	0.056239	3.7483	78.2666	15.7052	81.947	83.7876	-21.73%	-21.47%	-31.71%	-30.18%	4.98	-0.33%	
Atten Coeff	0.12												
Sriata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.088308	0.056239	3.6171	83.8955	20.4211	82.467	84.3258	-16.10%	2.11%	-31.28%	-29.73%	4.11	-17.83%	
Atten Coeff	0.12												
Sriata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.084993	0.056239	3.4813	87.5964	30.8374	87.1998	89.1779	-12.40%	54.19%	-27.33%	-25.69%	2.84	-43.19%	

Table B.13. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Iot3 hi6			
SC			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02062	0.012193	0.8446	31.8019	31.7418	62.7771	64.9711	-20.50%	-20.65%	-21.53%	-18.79%	1.002	0.19%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020516	0.012193	0.84033	33.287	33.8596	65.7848	68.0839	-16.78%	-15.35%	-17.77%	-14.90%	0.983	-1.69%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020395	0.012193	0.8354	34.2262	34.5482	66.6953	69.0262	-14.43%	-13.63%	-16.63%	-13.72%	0.991	-0.93%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020259	0.012193	0.82979	34.4961	35.4052	67.1971	69.5455	-13.76%	-11.49%	-16.00%	-13.07%	0.974	-2.57%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.019549	0.012193	0.80074	36.0957	36.7237	67.4493	69.8065	-9.76%	-8.19%	-15.69%	-12.74%	0.983	-1.71%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018816	0.012193	0.77068	37.1447	39.2215	68.5265	70.9213	-7.14%	-1.95%	-14.34%	-11.35%	0.947	-5.30%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02062	0.012193	0.8446	28.0816	28.0136	55.4174	56.6864	-29.80%	-29.97%	-30.73%	-29.14%	1.002	0.24%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020516	0.012193	0.84033	29.3312	29.9522	58.078	59.4079	-26.67%	-25.12%	-27.40%	-25.74%	0.979	-2.07%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020395	0.012193	0.8354	30.19	30.5361	58.8846	60.233	-24.53%	-23.66%	-26.39%	-24.71%	0.989	-1.13%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020259	0.012193	0.82979	30.3705	31.3579	59.331	60.6897	-24.07%	-21.61%	-25.84%	-24.14%	0.969	-3.15%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.019549	0.012193	0.80074	31.8245	32.5198	59.5687	60.9328	-20.44%	-18.70%	-25.54%	-23.83%	0.979	-2.14%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018816	0.012193	0.77068	32.6229	34.9106	60.5357	61.9218	-18.44%	-12.72%	-24.33%	-22.60%	0.934	-6.55%	

Table B.14. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	20	60					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020599	0.01228	0.84372	31.6338	15.8038	46.5336	48.1481	-20.92%	-20.98%	-22.44%	-19.75%	2.002	0.08%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020494	0.01228	0.83945	32.7194	17.5684	48.7951	50.4903	-18.20%	-12.16%	-18.67%	-15.85%	1.862	-6.88%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020374	0.01228	0.83453	33.4538	18.1999	49.4994	51.2196	-16.37%	-9.00%	-17.50%	-14.63%	1.838	-8.09%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020238	0.01228	0.82693	33.6805	18.9578	49.8891	51.6239	-15.80%	-5.21%	-16.85%	-13.96%	1.777	-11.17%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.019529	0.01228	0.79991	35.1243	20.311	50.1523	51.8971	-12.19%	1.56%	-16.41%	-13.50%	1.729	-13.53%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018796	0.01228	0.76988	36.0769	22.8659	51.2374	53.0226	-9.81%	14.33%	-14.60%	-11.63%	1.578	-21.11%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020599	0.01228	0.84372	29.496	12.3808	41.0744	42.0081	-26.26%	-38.10%	-31.54%	-29.99%	2.382	19.12%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020494	0.01228	0.83945	30.3867	14.0123	43.0754	44.056	-24.03%	-29.94%	-28.21%	-26.57%	2.169	8.43%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020374	0.01228	0.83453	31.0491	14.5606	43.6996	44.6946	-22.38%	-27.20%	-27.17%	-25.51%	2.132	6.62%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020238	0.01228	0.82693	31.1962	15.2896	44.0472	45.0508	-22.01%	-23.55%	-26.59%	-24.92%	2.040	2.02%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.019529	0.01228	0.79991	32.4925	16.5006	44.2948	45.3045	-18.77%	-17.50%	-26.18%	-24.49%	1.969	-1.54%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018796	0.01228	0.76988	33.188	18.9581	45.2685	46.3019	-17.03%	-5.21%	-24.55%	-22.83%	1.751	-12.47%	

Table B.15. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW $\ln 13_{hi6}$			
SC			40	8	48				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								185-190			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020575	0.012348	0.84277	31.5348	6.2427	36.3909	37.6255	-21.16%	-21.97%	-24.19%	-21.61%	5.051	1.03%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020471	0.012348	0.8385	32.3805	7.7964	38.191	39.4933	-19.05%	-2.55%	-20.44%	-17.72%	4.153	-16.93%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020351	0.012348	0.83358	32.9937	8.3939	38.7717	40.0956	-17.52%	4.92%	-19.23%	-16.47%	3.931	-21.39%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020215	0.012348	0.82799	33.1971	9.0909	39.0926	40.4301	-17.01%	13.64%	-18.56%	-15.77%	3.652	-26.97%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.019507	0.012348	0.799	34.5613	10.4581	39.3497	40.6998	-13.60%	30.73%	-18.02%	-15.21%	3.305	-33.91%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018775	0.012348	0.76901	35.453	13.0632	40.4649	41.8605	-11.37%	63.29%	-15.70%	-12.79%	2.714	-45.72%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020575	0.012348	0.84277	30.3462	3.0026	32.1132	32.8269	-24.13%	-62.47%	-33.10%	-31.61%	10.107	102.13%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020471	0.012348	0.8385	31.0215	4.4507	33.707	34.46	-22.45%	-44.37%	-29.78%	-28.21%	6.970	39.40%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020351	0.012348	0.83358	31.5678	4.9776	34.222	34.9876	-21.08%	-37.78%	-28.70%	-27.11%	6.342	26.84%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020215	0.012348	0.82799	31.6976	5.6495	34.5093	35.2829	-20.76%	-29.38%	-28.11%	-26.49%	5.611	12.21%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.019507	0.012348	0.799	32.9135	6.8819	34.7523	35.5335	-17.72%	-13.98%	-27.60%	-25.97%	4.783	-4.35%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018775	0.012348	0.76901	33.543	9.3953	35.7535	36.5614	-16.14%	17.44%	-25.51%	-23.83%	3.570	-28.60%

Table B.16. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW I ¹²³ I _{hi6}				
SC			84	84	168				130-142				
Atten Coeff	0.1481								142-175				
Striata ROI Size	13 13								185-190				
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018781	0.011653	0.76925	67.67	67.4978	134.3319	139.0267	-19.44%	-19.65%	-20.04%	-17.25%	1.003	0.26%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018685	0.011653	0.76536	71.2243	71.2652	140.9261	145.8513	-15.21%	-15.16%	-16.12%	-13.18%	0.999	-0.06%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018576	0.011653	0.76087	72.883	72.3779	142.7893	147.7796	-13.23%	-13.84%	-15.01%	-12.04%	1.007	0.70%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018451	0.011653	0.75576	73.2336	73.7715	143.6699	148.691	-12.82%	-12.18%	-14.48%	-11.49%	0.993	-0.73%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017805	0.011653	0.7293	75.0901	75.664	144.4694	149.5185	-10.61%	-9.92%	-14.01%	-11.00%	0.992	-0.76%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017137	0.011653	0.70193	76.033	78.4868	145.5446	150.6311	-9.48%	-6.56%	-13.37%	-10.34%	0.969	-3.13%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018781	0.011653	0.76925	59.7573	59.5685	118.5872	121.3028	-28.86%	-29.09%	-29.41%	-27.80%	1.003	0.32%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018685	0.011653	0.76536	62.8793	62.9243	124.4205	127.2696	-25.14%	-25.09%	-25.94%	-24.24%	0.999	-0.07%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018576	0.011653	0.76087	64.4069	63.8535	126.0717	128.9586	-23.33%	-23.98%	-24.96%	-23.24%	1.009	0.87%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018451	0.011653	0.75576	64.6136	65.1961	126.8537	129.7586	-23.08%	-22.39%	-24.49%	-22.76%	0.991	-0.89%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017805	0.011653	0.7293	66.2641	66.9037	127.5794	130.5009	-21.11%	-20.35%	-24.06%	-22.32%	0.990	-0.96%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017137	0.011653	0.70193	66.926	69.6292	128.5462	131.4898	-20.33%	-17.11%	-23.48%	-21.73%	0.961	-3.88%	

Table B.17. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I ¹²³ I _{h16}			
SC			84	42	126					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018962	0.011508	0.77668	67.2775	33.913	100.0967	103.5701	-19.91%	-19.25%	-20.56%	-17.80%	1.984	-0.81%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018866	0.011508	0.77274	69.9901	36.9239	105.1599	108.8123	-16.69%	-12.09%	-16.54%	-13.64%	1.895	-5.24%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018755	0.011508	0.76821	71.1932	37.904	106.5691	110.271	-15.25%	-9.75%	-15.42%	-12.48%	1.878	-6.09%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018629	0.011508	0.76306	71.434	39.0762	107.2332	110.96	-14.96%	-6.96%	-14.89%	-11.94%	1.828	-8.60%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017977	0.011508	0.73634	72.8707	40.9742	107.9349	111.6882	-13.25%	-2.44%	-14.34%	-11.36%	1.778	-11.08%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017302	0.011508	0.7087	73.4504	43.9434	108.9381	112.7302	-12.56%	4.63%	-13.54%	-10.53%	1.671	-16.43%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018962	0.011508	0.77668	62.7001	26.63	88.3571	90.366	-25.36%	-36.60%	-29.88%	-28.28%	2.354	17.72%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018866	0.011508	0.77274	65.0574	29.3278	92.8368	94.9494	-22.55%	-30.17%	-26.32%	-24.64%	2.218	10.91%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018755	0.011508	0.76821	66.1588	30.1707	94.0858	96.2272	-21.24%	-28.17%	-25.33%	-23.63%	2.193	9.64%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.018629	0.011508	0.76306	66.2783	31.307	94.6784	96.8323	-21.10%	-25.46%	-24.86%	-23.15%	2.117	5.85%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017977	0.011508	0.73634	67.5113	33.0612	95.3151	97.4867	-19.63%	-21.28%	-24.35%	-22.63%	2.042	2.10%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017302	0.011508	0.7087	67.7947	35.9691	96.2178	98.4123	-19.29%	-14.36%	-23.64%	-21.90%	1.885	-5.76%	

Table B.18. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW I ¹²³ I _{hi6}			
SC			84	16.8	100.8				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								185-190			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018948	0.011601	0.77612	67.0622	13.7826	78.9032	81.5826	-20.16%	-17.96%	-21.72%	-19.06%	4.866	-2.69%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018852	0.011601	0.77219	69.2597	16.3465	83.0196	85.8489	-17.55%	-2.70%	-17.64%	-14.83%	4.237	-15.26%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018742	0.011601	0.76766	70.2146	17.2527	84.1094	86.9788	-16.41%	2.69%	-16.56%	-13.71%	4.070	-18.60%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018616	0.011601	0.76251	70.4003	18.2969	84.6379	87.5297	-16.19%	8.92%	-16.03%	-13.16%	3.847	-23.06%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.017964	0.011601	0.73581	71.6286	20.2414	85.2709	88.1911	-14.73%	20.48%	-15.41%	-12.51%	3.539	-29.23%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01729	0.011601	0.70819	72.0657	23.3165	86.3176	89.2846	-14.21%	38.79%	-14.37%	-11.42%	3.091	-38.18%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018948	0.011601	0.77612	64.4838	6.8851	69.6322	71.1813	-23.23%	-59.02%	-30.92%	-29.38%	9.366	87.31%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018852	0.011601	0.77219	66.3873	9.1942	73.2754	74.9114	-20.97%	-45.27%	-27.31%	-25.68%	7.221	44.41%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018742	0.011601	0.76766	67.2413	9.9902	74.2423	75.9016	-19.95%	-40.53%	-26.35%	-24.70%	6.731	34.61%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.018616	0.011601	0.76251	67.3183	11.0085	74.7136	76.366	-19.86%	-34.47%	-25.88%	-24.22%	6.115	22.30%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.017964	0.011601	0.73581	68.3404	12.8268	75.2923	76.9816	-18.64%	-23.65%	-25.31%	-23.63%	5.328	6.56%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01729	0.011601	0.70819	68.4699	15.8549	76.2347	77.9515	-18.49%	-5.63%	-24.37%	-22.67%	4.319	-13.63%

Table B.19. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC			100	100	200				137-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								175-180			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080967	0.056577	3.3164	74.2607	72.8325	142.0695	147.0347	-25.74%	-27.17%	-28.97%	-26.48%	1.020	1.96%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080557	0.056577	3.2996	76.8137	78.9971	147.3112	152.4595	-23.19%	-21.00%	-26.34%	-23.77%	0.972	-2.76%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080084	0.056577	3.2803	80.6494	81.952	150.8184	156.0893	-19.35%	-18.05%	-24.59%	-21.96%	0.984	-1.59%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.079547	0.056577	3.2583	83.728	82.5991	151.87	157.1776	-16.27%	-17.40%	-24.07%	-21.41%	1.014	1.37%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076762	0.056577	3.1442	87.5684	95.0654	156.0819	161.5361	-12.43%	-4.93%	-21.96%	-19.23%	0.921	-7.89%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.073881	0.056577	3.0262	92.4776	107.8255	161.4137	167.0525	-7.52%	7.83%	-19.29%	-16.47%	0.858	-14.23%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080967	0.056577	3.3164	65.7111	64.1456	125.4175	128.2895	-34.29%	-35.85%	-37.29%	-35.86%	1.024	2.44%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080557	0.056577	3.2996	67.6073	69.968	130.058	133.0362	-32.39%	-30.03%	-34.97%	-33.48%	0.966	-3.37%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.080084	0.056577	3.2803	71.1037	72.4924	133.1667	136.2161	-28.90%	-27.51%	-33.42%	-31.89%	0.981	-1.92%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.079547	0.056577	3.2583	74.0625	72.8554	134.1127	137.1837	-25.94%	-27.14%	-32.94%	-31.41%	1.017	1.66%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076762	0.056577	3.1442	76.6469	84.8609	137.9022	141.0597	-23.35%	-15.14%	-31.05%	-29.47%	0.903	-9.68%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.073881	0.056577	3.0262	80.2627	97.1505	142.6918	145.9579	-19.74%	-2.85%	-28.65%	-27.02%	0.826	-17.38%

Table B.20. ¹²³I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC			100	50	150					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080812	0.056576	3.3101	73.2582	36.5637	103.8198	107.4217	-26.74%	-26.87%	-30.79%	-28.39%	2.004	0.18%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080403	0.056576	3.2933	74.6719	42.1656	107.3484	111.0812	-25.33%	-15.67%	-28.43%	-25.95%	1.771	-11.45%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079931	0.056576	3.274	77.7915	45.2579	110.5136	114.3584	-22.21%	-9.48%	-26.32%	-23.76%	1.719	-14.06%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079395	0.056576	3.252	80.6887	45.7402	111.3007	115.1714	-19.31%	-8.52%	-25.80%	-23.22%	1.764	-11.80%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076615	0.056576	3.1382	84.4915	57.9204	115.5734	119.604	-15.51%	15.84%	-22.95%	-20.26%	1.459	-27.06%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.073739	0.056576	3.0204	89.3098	70.7458	121.0455	125.273	-10.69%	41.49%	-19.30%	-16.48%	1.262	-36.88%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080812	0.056576	3.3101	68.3183	28.635	91.6426	93.7257	-31.68%	-42.73%	-38.90%	-37.52%	2.386	19.29%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080403	0.056576	3.2933	69.1479	34.0185	94.7702	96.9295	-30.85%	-31.96%	-36.82%	-35.38%	2.033	1.63%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079931	0.056576	3.274	71.9231	36.7512	97.5773	99.8016	-28.08%	-26.50%	-34.95%	-33.47%	1.957	-2.15%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079395	0.056576	3.252	74.7182	36.9707	98.288	100.5276	-25.28%	-26.06%	-34.47%	-32.98%	2.021	1.05%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076615	0.056576	3.1382	77.2917	48.6959	102.1336	104.4674	-22.71%	-2.61%	-31.91%	-30.36%	1.587	-20.64%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.073739	0.056576	3.0204	80.808	61.0583	107.0486	109.4982	-19.19%	22.12%	-28.63%	-27.00%	1.323	-33.83%	

Table B.21. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC			100	20	120					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080666	0.056633	3.3041	72.6669	14.8082	79.4694	82.1673	-27.33%	-25.96%	-33.78%	-31.53%	4.907	-1.86%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080257	0.056633	3.2873	73.4021	20.0716	82.1053	84.9142	-26.60%	0.36%	-31.58%	-29.24%	3.657	-26.86%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079786	0.056633	3.268	76.0934	23.2519	84.999	87.9152	-23.91%	16.26%	-29.17%	-26.74%	3.273	-34.55%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079251	0.056633	3.2461	78.8834	23.6385	85.5087	88.4414	-21.12%	18.19%	-28.74%	-26.30%	3.337	-33.26%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076476	0.056633	3.1325	82.7574	35.5915	89.7141	92.8188	-17.24%	77.96%	-25.24%	-22.65%	2.325	-53.50%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.073606	0.056633	3.0149	87.4998	48.5743	95.3436	98.6603	-12.50%	142.87%	-20.55%	-17.78%	1.801	-63.97%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080666	0.056633	3.3041	69.8919	7.3339	70.1308	71.6906	-30.11%	-63.33%	-41.56%	-40.26%	9.530	90.60%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.080257	0.056633	3.2873	70.087	12.4527	72.4732	74.0976	-29.91%	-37.74%	-39.61%	-38.25%	5.628	12.57%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079786	0.056633	3.268	72.4302	15.3158	75.0417	76.7285	-27.57%	-23.42%	-37.47%	-36.06%	4.729	-5.42%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.079251	0.056633	3.2461	75.1288	15.4516	75.5067	77.2034	-24.87%	-22.74%	-37.08%	-35.66%	4.862	-2.76%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.076476	0.056633	3.1325	77.7945	26.9429	79.2968	81.0948	-22.21%	34.71%	-33.92%	-32.42%	2.887	-42.25%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.073606	0.056633	3.0149	81.2196	39.4661	84.3533	86.276	-18.78%	97.33%	-29.71%	-28.10%	2.058	-58.84%	

Table B.22. ^{123}I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	40	80					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01803	0.01251	0.73849	30.002	29.2188	58.1945	60.2283	-25.00%	-26.95%	-27.26%	-24.71%	1.027	2.68%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017938	0.01251	0.73475	31.3582	31.161	60.7314	62.8539	-21.60%	-22.10%	-24.09%	-21.43%	1.006	0.63%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017833	0.01251	0.73044	32.5238	31.9418	61.8416	64.0029	-18.69%	-20.15%	-22.70%	-20.00%	1.018	1.82%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017713	0.01251	0.72554	33.2208	32.2598	62.1561	64.3284	-16.95%	-19.35%	-22.30%	-19.59%	1.030	2.98%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017093	0.01251	0.70014	34.0103	35.3407	63.1657	65.3732	-14.97%	-11.65%	-21.04%	-18.28%	0.962	-3.76%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016452	0.01251	0.67386	34.9745	38.3456	64.2646	66.5102	-12.56%	-4.14%	-19.67%	-16.86%	0.912	-8.79%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01803	0.01251	0.73849	26.5665	25.7145	51.3744	52.5508	-33.58%	-35.71%	-35.78%	-34.31%	1.033	3.31%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017938	0.01251	0.73475	27.7072	27.4929	53.6185	54.8464	-30.73%	-31.27%	-32.98%	-31.44%	1.008	0.78%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017833	0.01251	0.73044	28.7809	28.1447	54.6019	55.8523	-28.05%	-29.64%	-31.75%	-30.18%	1.023	2.26%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017713	0.01251	0.72554	29.4342	28.3955	54.8836	56.1404	-26.41%	-29.01%	-31.40%	-29.82%	1.037	3.66%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017093	0.01251	0.70014	29.9163	31.3767	55.7923	57.0698	-25.21%	-21.56%	-30.26%	-28.66%	0.953	-4.65%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016452	0.01251	0.67386	30.5787	34.2887	56.7817	58.0817	-23.55%	-14.28%	-29.02%	-27.40%	0.892	-10.82%	

Table B.23. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Jo13 hi6			
SC			40	20	60					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01796	0.012474	0.73563	29.6018	14.7136	43.109	44.6045	-26.00%	-26.43%	-28.15%	-25.66%	2.012	0.59%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017869	0.012474	0.7319	30.503	16.4303	44.9659	46.5283	-23.74%	-17.85%	-25.06%	-22.45%	1.857	-7.17%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017764	0.012474	0.72761	31.3794	17.2695	45.911	47.5068	-21.55%	-13.65%	-23.48%	-20.82%	1.817	-9.15%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017645	0.012474	0.72273	32.0017	17.5231	46.0956	47.6977	-20.00%	-12.38%	-23.17%	-20.50%	1.826	-8.69%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017027	0.012474	0.69742	32.816	20.45	47.0873	48.7274	-17.96%	2.25%	-21.52%	-18.79%	1.605	-19.77%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016388	0.012474	0.67124	33.7279	23.5245	48.2649	49.9486	-15.68%	17.62%	-19.56%	-16.75%	1.434	-28.31%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01796	0.012474	0.73563	27.61	11.5122	38.0536	38.9185	-30.98%	-42.44%	-36.58%	-35.14%	2.398	19.92%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017869	0.012474	0.7319	28.3246	13.1148	39.697	40.6007	-29.19%	-34.43%	-33.84%	-32.33%	2.160	7.99%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017764	0.012474	0.72761	29.1071	13.8534	40.5348	41.4579	-27.23%	-30.73%	-32.44%	-30.90%	2.101	5.05%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017645	0.012474	0.72273	29.6927	14.0483	40.7013	41.6282	-25.77%	-29.76%	-32.16%	-30.62%	2.114	5.68%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017027	0.012474	0.69742	30.213	16.8753	41.5943	42.5436	-24.47%	-15.62%	-30.68%	-29.09%	1.790	-10.48%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016388	0.012474	0.67124	30.8148	19.8623	42.6539	43.6288	-22.96%	-0.69%	-28.91%	-27.29%	1.551	-22.43%	

Table B.24. ^{123}I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC			40	8	48					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017866	0.012521	0.7318	29.3704	6.0188	33.6437	34.7861	-26.57%	-24.77%	-29.91%	-27.53%	4.880	-2.40%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017776	0.012521	0.72809	30.0011	7.6032	35.0594	36.2564	-25.00%	-4.96%	-26.96%	-24.47%	3.946	-21.08%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017671	0.012521	0.72382	30.7075	8.48	35.9062	37.1349	-23.23%	6.00%	-25.20%	-22.64%	3.621	-27.58%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017553	0.012521	0.71897	31.2856	8.7007	36.0448	37.2786	-21.79%	8.76%	-24.91%	-22.34%	3.596	-28.08%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016938	0.012521	0.69379	32.1046	11.5761	37.0576	38.3347	-19.74%	44.70%	-22.80%	-20.14%	2.773	-44.53%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016303	0.012521	0.66775	33.0375	14.6541	38.3062	39.633	-17.41%	83.18%	-20.20%	-17.43%	2.254	-54.91%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017866	0.012521	0.7318	28.2438	2.9983	29.6911	30.3516	-29.39%	-62.52%	-38.14%	-36.77%	9.420	88.40%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017776	0.012521	0.72809	28.7049	4.498	30.9454	31.6377	-28.24%	-43.78%	-35.53%	-34.09%	6.382	27.63%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017671	0.012521	0.72382	29.3159	5.2907	31.6967	32.4073	-26.71%	-33.87%	-33.97%	-32.48%	5.541	10.82%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.017553	0.012521	0.71897	29.861	5.4579	31.8227	32.5364	-25.35%	-31.78%	-33.70%	-32.22%	5.471	9.42%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016938	0.012521	0.69379	30.3905	8.2335	32.7356	33.4748	-24.02%	2.92%	-31.80%	-30.26%	3.691	-26.18%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016303	0.012521	0.66775	31.0119	11.222	33.859	34.6274	-22.47%	40.28%	-29.46%	-27.86%	2.763	-44.73%	

Table B.25. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13_hi6			
SC			84	84	168				137-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								175-180			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.017076	0.01189	0.69945	63.5519	61.6535	123.8875	128.2171	-24.34%	-26.60%	-26.26%	-23.68%	1.031	3.08%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01699	0.01189	0.69591	66.9004	64.98	129.5384	134.0656	-20.36%	-22.64%	-22.89%	-20.20%	1.030	2.96%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01689	0.01189	0.69182	68.8184	66.0569	131.4545	136.0485	-18.07%	-21.36%	-21.75%	-19.02%	1.042	4.18%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016777	0.01189	0.68718	69.5545	66.7572	131.9962	136.6092	-17.20%	-20.53%	-21.43%	-18.69%	1.042	4.19%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01619	0.01189	0.66312	70.1462	70.5174	133.2463	137.9031	-16.49%	-16.05%	-20.69%	-17.91%	0.995	-0.53%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015582	0.01189	0.63823	70.9651	73.7489	134.3219	139.0161	-15.52%	-12.20%	-20.05%	-17.25%	0.962	-3.77%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.017076	0.01189	0.69945	56.2961	54.2366	109.3688	111.8733	-32.98%	-35.43%	-34.90%	-33.41%	1.038	3.80%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01699	0.01189	0.69591	59.2595	57.1791	114.3673	116.9862	-29.45%	-31.93%	-31.92%	-30.37%	1.036	3.64%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01689	0.01189	0.69182	61.0468	58.0476	116.0645	118.7223	-27.33%	-30.90%	-30.91%	-29.33%	1.052	5.17%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016777	0.01189	0.68718	61.7019	58.6699	116.5472	119.216	-26.55%	-30.15%	-30.63%	-29.04%	1.052	5.17%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01619	0.01189	0.66312	61.9224	62.3444	117.672	120.3666	-26.28%	-25.78%	-29.96%	-28.35%	0.993	-0.68%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015582	0.01189	0.63823	62.4185	65.4958	118.6411	121.3578	-25.69%	-22.03%	-29.38%	-27.76%	0.953	-4.70%

Table B.26. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW, lo13, hi6			
SC			84	42	126				137-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								175-180			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016926	0.011669	0.6933	62.7104	31.1203	92.2861	95.4876	-25.34%	-25.90%	-26.76%	-24.22%	2.015	0.75%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016841	0.011669	0.68979	65.1024	33.9719	96.6788	100.0359	-22.50%	-19.11%	-23.27%	-20.61%	1.916	-4.18%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016742	0.011669	0.68574	66.4106	35.1746	98.2904	101.7044	-20.94%	-16.25%	-21.99%	-19.28%	1.888	-5.60%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01663	0.011669	0.68115	66.9856	35.7387	98.6951	102.1236	-20.26%	-14.91%	-21.67%	-18.95%	1.874	-6.28%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016047	0.011669	0.6573	67.601	39.1749	99.8729	103.3476	-19.52%	-6.73%	-20.74%	-17.98%	1.726	-13.72%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015445	0.011669	0.63262	68.2643	42.51	101.0119	104.5304	-18.73%	1.21%	-19.83%	-17.04%	1.606	-19.71%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016926	0.011669	0.6933	58.4934	24.3412	81.4642	83.3158	-30.37%	-42.04%	-35.35%	-33.88%	2.403	20.15%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016841	0.011669	0.68979	60.5611	26.9133	85.3504	87.2923	-27.90%	-35.92%	-32.26%	-30.72%	2.250	12.51%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016742	0.011669	0.68574	61.7343	27.9665	86.779	88.7539	-26.51%	-33.41%	-31.13%	-29.56%	2.207	10.37%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01663	0.011669	0.68115	62.2433	28.4711	87.1405	89.1239	-25.90%	-32.21%	-30.84%	-29.27%	2.186	9.31%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016047	0.011669	0.6573	62.5162	31.8242	88.202	90.2126	-25.58%	-24.23%	-30.00%	-28.40%	1.964	-1.78%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015445	0.011669	0.63262	62.8434	35.0965	89.2278	91.2642	-25.19%	-16.44%	-29.18%	-27.57%	1.791	-10.47%

Table B.27. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW_Io13_hi6			
SC			84	16.8	100.8					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016843	0.011598	0.68991	62.2037	12.7996	72.6252	75.0916	-25.95%	-23.81%	-27.95%	-25.50%	4.860	-2.80%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016758	0.011598	0.68641	64.0211	15.366	76.1992	78.7972	-23.78%	-8.54%	-24.41%	-21.83%	4.166	-16.67%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01666	0.011598	0.68239	64.9639	16.6429	77.579	80.2285	-22.66%	-0.94%	-23.04%	-20.41%	3.903	-21.93%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016548	0.011598	0.67781	65.441	17.1264	77.8606	80.5214	-22.09%	1.94%	-22.76%	-20.12%	3.821	-23.58%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015969	0.011598	0.65408	66.0407	20.3925	78.9224	81.6315	-21.38%	21.38%	-21.70%	-19.02%	3.238	-35.23%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015369	0.011598	0.62953	66.6805	23.7519	80.1551	82.917	-20.62%	41.38%	-20.48%	-17.74%	2.807	-43.85%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016843	0.011598	0.68991	59.8103	6.4033	64.0939	65.5198	-28.80%	-61.89%	-36.41%	-35.00%	9.341	86.81%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016758	0.011598	0.68641	61.3397	8.753	67.2574	68.7597	-26.98%	-47.90%	-33.28%	-31.79%	7.008	40.16%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01666	0.011598	0.68239	62.1452	9.9157	68.4819	70.0142	-26.02%	-40.98%	-32.06%	-30.54%	6.267	25.35%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016548	0.011598	0.67781	62.5651	10.3512	68.7346	70.2735	-25.52%	-38.39%	-31.81%	-30.28%	6.044	20.88%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015969	0.011598	0.65408	62.8381	13.5376	69.6938	71.2612	-25.19%	-19.42%	-30.86%	-29.30%	4.642	-7.17%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015369	0.011598	0.62953	63.1376	16.8372	70.8041	72.4026	-24.84%	0.22%	-29.76%	-28.17%	3.750	-25.00%	

Table B.28. 123I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo6 hi6			
SC			100	100	200				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076438	0.056568	3.1309	70.315	67.0328	133.1668	137.8206	-29.69%	-32.97%	-33.42%	-31.09%	1.049	4.90%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076051	0.056568	3.115	72.0317	74.5487	139.2743	144.1417	-27.97%	-25.45%	-30.36%	-27.93%	0.966	-3.38%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075605	0.056568	3.0968	74.0063	78.5345	141.5847	146.5327	-25.99%	-21.47%	-29.21%	-26.73%	0.942	-5.77%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075098	0.056568	3.076	76.9422	78.2963	141.469	146.4131	-23.06%	-21.70%	-29.27%	-26.79%	0.983	-1.73%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072468	0.056568	2.9683	77.342	91.1672	143.588	148.6037	-22.66%	-8.83%	-28.21%	-25.70%	0.848	-15.16%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069748	0.056568	2.8569	83.5691	99.4709	147.6695	152.8273	-16.43%	-0.53%	-26.17%	-23.59%	0.840	-15.99%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076438	0.056568	3.1309	62.4056	58.8493	117.5604	120.2524	-37.59%	-41.15%	-41.22%	-39.87%	1.060	6.04%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076051	0.056568	3.115	63.3341	66.0947	122.9699	125.7857	-36.67%	-33.91%	-38.52%	-37.11%	0.958	-4.18%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075605	0.056568	3.0968	64.8914	69.8262	125.0238	127.8866	-35.11%	-30.17%	-37.49%	-36.06%	0.929	-7.07%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075098	0.056568	3.076	67.7858	69.3391	124.9328	127.7937	-32.21%	-30.66%	-37.53%	-36.10%	0.978	-2.24%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072468	0.056568	2.9683	66.9104	82.1085	126.8591	129.7626	-33.09%	-17.89%	-36.57%	-35.12%	0.815	-18.51%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069748	0.056568	2.8569	72.2838	89.8277	130.538	133.5255	-27.72%	-10.17%	-34.73%	-33.24%	0.805	-19.53%

Table B.29. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW, lo13 hi6			
SC			100	50	150				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076338	0.056583	3.1268	69.5364	33.6856	97.9054	101.3001	-30.46%	-32.63%	-34.73%	-32.47%	2.064	3.21%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075952	0.056583	3.111	70.0774	40.691	102.413	105.9761	-29.92%	-18.62%	-31.72%	-29.35%	1.722	-13.89%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075506	0.056583	3.0927	71.4338	44.7751	104.2846	107.917	-28.57%	-10.45%	-30.48%	-28.06%	1.595	-20.23%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075	0.056583	3.072	74.0735	44.5923	103.9774	107.5973	-25.93%	-10.82%	-30.68%	-28.27%	1.661	-16.94%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072374	0.056583	2.9644	74.0189	57.6958	106.3963	110.1117	-25.98%	15.39%	-29.07%	-26.59%	1.283	-35.85%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069657	0.056583	2.8532	80.0422	66.1128	110.4397	114.2979	-19.96%	32.23%	-26.37%	-23.80%	1.211	-39.47%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076338	0.056583	3.1268	64.9457	26.1818	86.4228	88.3862	-35.05%	-47.64%	-42.38%	-41.08%	2.481	24.03%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075952	0.056583	3.111	64.7672	33.0431	90.4199	92.4811	-35.23%	-33.91%	-39.72%	-38.35%	1.960	-2.00%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075506	0.056583	3.0927	65.7051	36.9331	92.0862	94.1878	-34.29%	-26.13%	-38.61%	-37.21%	1.779	-11.05%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075	0.056583	3.072	68.3021	36.5298	91.825	93.9195	-31.70%	-26.94%	-38.78%	-37.39%	1.870	-6.51%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072374	0.056583	2.9644	66.953	49.5735	94.0219	96.1731	-33.05%	-0.85%	-37.32%	-35.88%	1.351	-32.47%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069657	0.056583	2.8532	72.1098	57.4254	97.6661	99.9017	-27.89%	14.85%	-34.89%	-33.40%	1.256	-37.21%

Table B.30. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo6 hi6			
SC			100	20	120				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076201	0.056657	3.1212	69.0843	13.6861	75.9607	78.5377	-30.92%	-31.57%	-36.70%	-34.55%	5.048	0.96%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075815	0.056657	3.1054	68.9271	20.387	79.4175	82.1403	-31.07%	1.94%	-33.82%	-31.55%	3.381	-32.38%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075371	0.056657	3.0872	69.916	24.5356	81.0126	83.802	-30.08%	22.68%	-32.49%	-30.17%	2.850	-43.01%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.074865	0.056657	3.0665	72.4045	24.3823	80.7621	83.5407	-27.60%	21.91%	-32.70%	-30.38%	2.970	-40.61%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072244	0.056657	2.9591	72.1179	37.5975	82.9301	85.8117	-27.88%	87.99%	-30.89%	-28.49%	1.918	-61.64%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069532	0.056657	2.848	78.0532	46.1566	87.1278	90.1619	-21.95%	130.78%	-27.39%	-24.87%	1.691	-66.18%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.076201	0.056657	3.1212	66.4832	6.5893	67.0346	68.5245	-33.52%	-67.05%	-44.14%	-42.90%	10.090	101.79%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075815	0.056657	3.1054	65.6476	13.2211	70.1071	71.6818	-34.35%	-33.89%	-41.58%	-40.27%	4.965	-0.69%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.075371	0.056657	3.0872	66.2168	17.2112	71.5307	73.1443	-33.78%	-13.94%	-40.39%	-39.05%	3.847	-23.05%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.074865	0.056657	3.0665	68.6627	16.8507	71.3203	72.9279	-31.34%	-15.75%	-40.57%	-39.23%	4.075	-18.50%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.072244	0.056657	2.9591	67.0731	30.0257	73.298	74.9667	-32.93%	50.13%	-38.92%	-37.53%	2.234	-55.32%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.069532	0.056657	2.848	72.1283	38.0247	77.0788	78.8374	-27.87%	90.12%	-35.77%	-34.30%	1.897	-62.06%

Table B.31. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC			40	40	80					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0169	0.012608	0.6922	28.0577	26.9943	54.1057	55.9965	-29.86%	-32.51%	-32.37%	-30.00%	1.039	3.94%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016814	0.012608	0.6887	29.2077	29.2017	56.7585	58.7422	-26.98%	-27.00%	-29.05%	-26.57%	1.000	0.02%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016715	0.012608	0.68466	29.8534	30.2562	57.5734	59.5855	-25.37%	-24.36%	-28.03%	-25.52%	0.987	-1.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016603	0.012608	0.68007	30.592	30.3007	57.6582	59.6733	-23.52%	-24.25%	-27.93%	-25.41%	1.010	0.96%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016022	0.012608	0.65625	30.7571	33.3368	58.3212	60.3592	-23.11%	-16.66%	-27.10%	-24.55%	0.923	-7.74%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01542	0.012608	0.63162	32.0703	35.3649	59.1633	61.2306	-19.82%	-11.59%	-26.05%	-23.46%	0.907	-9.32%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0169	0.012608	0.6922	24.8762	23.7252	47.7655	48.8593	-37.81%	-40.69%	-40.29%	-38.93%	1.049	4.85%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016814	0.012608	0.6887	25.7849	25.7877	50.113	51.2606	-35.54%	-35.53%	-37.36%	-35.92%	1.000	-0.01%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016715	0.012608	0.68466	26.3181	26.7625	50.8362	52.0003	-34.20%	-33.09%	-36.45%	-35.00%	0.983	-1.66%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016603	0.012608	0.68007	27.0363	26.7421	50.914	52.0799	-32.41%	-33.14%	-36.36%	-34.90%	1.011	1.10%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016022	0.012608	0.65625	26.8996	29.7467	51.5131	52.6926	-32.75%	-25.63%	-35.61%	-34.13%	0.904	-9.57%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01542	0.012608	0.63162	28.0061	31.6512	52.2738	53.4706	-29.98%	-20.87%	-34.66%	-33.16%	0.885	-11.52%	

Table B.32. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, I _o 13, h16			
SC			40	20	60					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01695	0.012514	0.69427	27.7314	13.6406	40.2525	41.6486	-30.67%	-31.80%	-32.91%	-30.59%	2.033	1.65%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016864	0.012514	0.69075	28.4066	15.6387	42.2239	43.6916	-28.98%	-21.81%	-29.63%	-27.18%	1.816	-9.18%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016765	0.012514	0.6867	28.7983	16.7268	42.859	44.3501	-28.00%	-16.37%	-28.57%	-26.08%	1.722	-13.92%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016653	0.012514	0.6821	29.4214	16.7788	42.8722	44.3633	-26.45%	-16.11%	-28.55%	-26.06%	1.753	-12.33%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01607	0.012514	0.65821	29.4099	19.8281	43.3836	44.8963	-26.48%	-0.86%	-27.69%	-25.17%	1.483	-25.84%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015466	0.012514	0.63351	30.6545	21.8727	44.26	45.8044	-23.36%	9.36%	-26.23%	-23.66%	1.401	-29.93%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01695	0.012514	0.69427	25.8791	10.645	35.5325	36.3399	-35.30%	-46.78%	-40.78%	-39.43%	2.431	21.56%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016864	0.012514	0.69075	26.3411	12.5493	37.2781	38.127	-34.15%	-37.25%	-37.87%	-36.46%	2.099	4.95%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016765	0.012514	0.6867	26.6205	13.5827	37.8426	38.7053	-33.45%	-32.09%	-36.93%	-35.49%	1.960	-2.01%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016653	0.012514	0.6821	27.224	13.5809	37.857	38.7197	-31.94%	-32.10%	-36.91%	-35.47%	2.005	0.23%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01607	0.012514	0.65821	26.9107	16.6164	38.3228	39.1984	-32.72%	-16.92%	-36.13%	-34.67%	1.620	-19.02%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015466	0.012514	0.63351	27.946	18.5446	39.1134	40.0077	-30.14%	-7.28%	-34.81%	-33.32%	1.507	-24.65%	

Table B.33. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Iot3_h16			
SC			40	8	48					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016921	0.012517	0.69309	27.5454	5.6384	31.5428	32.6138	-31.14%	-29.52%	-34.29%	-32.05%	4.885	-2.29%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016835	0.012517	0.68958	27.9387	7.5144	33.0874	34.2189	-30.15%	-6.07%	-31.07%	-28.71%	3.718	-25.64%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016737	0.012517	0.68553	28.1837	8.625	33.5943	34.747	-29.54%	7.81%	-30.01%	-27.61%	3.268	-34.65%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016624	0.012517	0.68093	28.7447	8.684	33.566	34.7175	-28.14%	8.55%	-30.07%	-27.67%	3.310	-33.80%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016042	0.012517	0.65709	28.6393	11.7711	34.0622	35.2399	-28.40%	47.14%	-29.04%	-26.58%	2.433	-51.34%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01544	0.012517	0.63243	29.8345	13.8481	34.9373	36.1486	-25.41%	73.10%	-27.21%	-24.69%	2.154	-56.91%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016921	0.012517	0.69309	26.4895	2.8058	27.8374	28.4566	-33.78%	-64.93%	-42.01%	-40.72%	9.441	88.82%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016835	0.012517	0.68958	26.6861	4.6182	29.2068	29.8611	-33.28%	-42.27%	-39.15%	-37.79%	5.778	15.57%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016737	0.012517	0.68553	26.8187	5.6882	29.6585	30.3252	-32.95%	-28.90%	-38.21%	-36.82%	4.715	-5.70%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016624	0.012517	0.68093	27.3603	5.6995	29.636	30.3021	-31.60%	-28.76%	-38.26%	-36.87%	4.800	-3.99%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016042	0.012517	0.65709	26.9503	8.7814	30.0901	30.7717	-32.62%	9.77%	-37.31%	-35.89%	3.069	-38.62%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.01544	0.012517	0.63243	27.9329	10.7465	30.8796	31.5811	-30.17%	34.33%	-35.67%	-34.21%	2.599	-48.01%	

Table B.34. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW I ₀₆ hi6			
SC			84	84	168				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016084	0.011772	0.6588	58.8902	57.03	114.4599	118.4601	-29.89%	-32.11%	-31.87%	-29.49%	1.033	3.26%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016003	0.011772	0.65547	61.9544	60.5052	120.0836	124.2804	-26.24%	-27.97%	-28.52%	-26.02%	1.024	2.40%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015909	0.011772	0.65162	63.0778	61.9185	121.5539	125.802	-24.91%	-26.29%	-27.65%	-25.12%	1.019	1.87%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015802	0.011772	0.64725	63.9955	62.1574	121.8129	126.07	-23.81%	-26.00%	-27.49%	-24.96%	1.030	2.96%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015249	0.011772	0.62459	64.2209	65.3951	122.3346	126.6101	-23.55%	-22.15%	-27.18%	-24.64%	0.982	-1.80%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.014676	0.011772	0.60114	65.2599	67.7285	123.193	127.4983	-22.31%	-19.37%	-26.67%	-24.11%	0.964	-3.64%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016084	0.011772	0.6588	52.1748	50.1626	101.0476	103.3614	-37.89%	-40.28%	-39.85%	-38.48%	1.040	4.01%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016003	0.011772	0.65547	54.84	53.2832	106.0225	108.4504	-34.71%	-36.57%	-36.89%	-35.45%	1.029	2.92%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015909	0.011772	0.65162	55.8097	54.5628	107.3256	109.7833	-33.56%	-35.04%	-36.12%	-34.65%	1.023	2.29%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015802	0.011772	0.64725	56.6844	54.7167	107.5574	110.0203	-32.52%	-34.86%	-35.98%	-34.51%	1.036	3.60%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015249	0.011772	0.62459	56.5888	57.9142	108.032	110.5059	-32.63%	-31.05%	-35.70%	-34.22%	0.977	-2.29%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.014676	0.011772	0.60114	57.3931	60.1477	108.8064	111.2979	-31.67%	-28.40%	-35.23%	-33.75%	0.954	-4.58%

Table B.35. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC			84	42	126					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016037	0.011676	0.65686	58.2293	28.9462	85.5186	88.4853	-30.68%	-31.08%	-32.13%	-29.77%	2.012	0.58%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015955	0.011676	0.65354	60.303	31.9893	89.7461	92.8634	-28.21%	-23.84%	-28.77%	-26.30%	1.885	-5.75%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015862	0.011676	0.6497	60.9015	33.4812	90.9061	94.0658	-27.50%	-20.28%	-27.85%	-25.34%	1.819	-9.05%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015755	0.011676	0.64534	61.5896	33.741	91.0747	94.2401	-26.68%	-19.66%	-27.72%	-25.21%	1.825	-8.73%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015204	0.011676	0.62275	61.4731	37.0461	91.5322	94.7188	-26.82%	-11.80%	-27.36%	-24.83%	1.659	-17.03%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.014633	0.011676	0.59937	62.363	39.4421	92.349	95.5665	-25.76%	-6.09%	-26.71%	-24.15%	1.581	-20.94%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.016037	0.011676	0.65686	54.3088	22.6517	75.4909	77.2067	-35.35%	-46.07%	-40.09%	-38.72%	2.398	19.88%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015955	0.011676	0.65354	56.0414	25.4466	79.2318	81.0351	-33.28%	-39.41%	-37.12%	-35.69%	2.202	10.12%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015862	0.011676	0.6497	56.4858	26.856	80.2615	82.0893	-32.76%	-36.06%	-36.30%	-34.85%	2.103	5.16%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015204	0.011676	0.62275	57.1318	27.0534	80.4134	82.2446	-31.99%	-35.59%	-36.18%	-34.73%	2.112	5.59%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.015204	0.011676	0.62275	56.6871	30.3541	80.8318	82.6756	-32.52%	-27.73%	-35.85%	-34.38%	1.868	-6.62%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.014633	0.011676	0.59937	57.3339	32.667	81.5691	83.4312	-31.75%	-22.22%	-35.26%	-33.78%	1.755	-12.24%	

Table B.36. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo6 hi6			
SC			84	16.8	100.8				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016042	0.011602	0.65707	57.8276	12.0899	67.391	69.6804	-31.16%	-28.04%	-33.14%	-30.87%	4.783	-4.34%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01596	0.011602	0.65374	59.3039	14.8727	70.7697	73.1854	-29.40%	-11.47%	-29.79%	-27.40%	3.987	-20.25%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015867	0.011602	0.64991	59.5863	16.4093	71.6624	74.1145	-29.06%	-2.33%	-28.91%	-26.47%	3.631	-27.37%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01576	0.011602	0.64555	60.1385	16.6795	71.7533	74.2092	-28.41%	-0.72%	-28.82%	-26.38%	3.606	-27.89%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015209	0.011602	0.62294	59.8319	20.01	72.1234	74.6044	-28.77%	19.11%	-28.45%	-25.99%	2.990	-40.20%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.014638	0.011602	0.59956	60.6373	22.4472	72.9585	75.4749	-27.81%	33.61%	-27.62%	-25.12%	2.701	-45.97%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.016042	0.011602	0.65707	55.5845	6.1397	59.4751	60.7988	-33.83%	-63.45%	-41.00%	-39.68%	9.053	81.07%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01596	0.011602	0.65374	56.755	8.7384	62.467	63.8639	-32.43%	-47.99%	-38.03%	-36.64%	6.495	29.90%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015867	0.011602	0.64991	56.8834	10.2235	63.2613	64.6795	-32.28%	-39.15%	-37.24%	-35.83%	5.564	11.28%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.01576	0.011602	0.64555	57.3944	10.4444	63.3447	64.7651	-31.67%	-37.83%	-37.16%	-35.75%	5.495	9.90%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.015209	0.011602	0.62294	56.7582	13.7895	63.6869	65.1222	-32.43%	-17.92%	-36.82%	-35.39%	4.116	-17.68%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.014638	0.011602	0.59956	57.3153	16.1524	64.4409	65.8971	-31.77%	-3.85%	-36.07%	-34.63%	3.548	-29.03%

Table B.37. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW	lo13	hi10			
SC			100	100	200					130-142					
Atten Coeff	0.1481									142-175					
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12229	0.054287	5.0091	78.1917	84.2076	159.0344	164.592	-21.81%	-15.79%	-20.48%	-17.70%	0.929	-7.14%			
Atten Coeff	0.1481														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12167	0.054287	4.9837	80.5542	91.4762	166.3394	172.1512	-19.45%	-8.52%	-16.83%	-13.92%	0.881	-11.94%			
Atten Coeff	0.1481														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12096	0.054287	4.9545	82.9754	94.0088	168.6161	174.5076	-17.02%	-5.99%	-15.69%	-12.75%	0.883	-11.74%			
Atten Coeff	0.1481														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12015	0.054287	4.9213	83.4723	97.005	169.5579	175.4816	-16.53%	-3.00%	-15.22%	-12.26%	0.860	-13.95%			
Atten Coeff	0.1481														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.11594	0.054287	4.749	88.1457	102.6442	169.187	175.0975	-11.85%	2.64%	-15.41%	-12.45%	0.859	-14.13%			
Atten Coeff	0.1481														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.11159	0.054287	4.5707	92.6631	111.984	173.5471	179.6085	-7.34%	11.98%	-13.23%	-10.20%	0.827	-17.25%			
Atten Coeff	0.12														
Striata ROI Size	13 13														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12229	0.054287	5.0091	68.4332	74.9282	140.3867	143.6012	-31.57%	-25.07%	-29.81%	-28.20%	0.913	-8.67%			
Atten Coeff	0.12														
Striata ROI Size	15 15														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12167	0.054287	4.9837	70.0291	81.8565	146.8497	150.2116	-29.97%	-18.14%	-26.58%	-24.89%	0.856	-14.45%			
Atten Coeff	0.12														
Striata ROI Size	17 17														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12096	0.054287	4.9545	72.167	84.1099	148.8677	152.2758	-27.83%	-15.89%	-25.57%	-23.86%	0.858	-14.20%			
Atten Coeff	0.12														
Striata ROI Size	19 19														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.12015	0.054287	4.9213	72.353	87.0344	149.7113	153.1383	-27.65%	-12.97%	-25.14%	-23.43%	0.831	-16.87%			
Atten Coeff	0.12														
Striata ROI Size	26 26														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			
0.11594	0.054287	4.749	76.4149	92.2183	149.4324	152.8529	-23.59%	-7.78%	-25.28%	-23.57%	0.829	-17.14%			
Atten Coeff	0.12														
Striata ROI Size	32 32														
Bkgd ROI Size	13 13														
Units:	microCi														
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err			

Table B.38. ¹²³I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13_h0			
SC			100	50	150				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11921	0.05459	4.8827	78.5089	39.8647	113.925	117.8788	-21.49%	-20.27%	-24.05%	-21.41%	1.969	-1.53%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1186	0.05459	4.8579	80.0302	45.634	118.8821	123.0167	-19.97%	-8.73%	-20.75%	-17.99%	1.754	-12.31%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11791	0.05459	4.8294	82.0671	47.6066	120.4433	124.6334	-17.93%	-4.79%	-19.70%	-16.91%	1.724	-13.81%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11712	0.05459	4.797	82.729	49.8546	120.9572	125.1679	-17.27%	-0.29%	-19.36%	-16.55%	1.659	-17.03%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11301	0.05459	4.6291	87.6732	54.2897	120.1824	124.3682	-12.33%	8.58%	-19.88%	-17.09%	1.615	-19.25%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10877	0.05459	4.4553	92.1576	63.3017	124.5007	128.8431	-7.84%	26.60%	-17.00%	-14.10%	1.456	-27.21%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11921	0.05459	4.8827	73.1393	31.3558	100.5548	102.8412	-26.86%	-37.29%	-32.96%	-31.44%	2.333	16.63%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1186	0.05459	4.8579	74.0562	36.893	104.9446	107.3361	-25.94%	-26.21%	-30.04%	-28.44%	2.007	0.37%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11791	0.05459	4.8294	75.8704	38.6336	106.3296	108.7532	-24.13%	-22.73%	-29.11%	-27.50%	1.964	-1.81%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11712	0.05459	4.797	76.2936	40.8042	106.7954	109.2313	-23.71%	-18.39%	-28.80%	-27.18%	1.870	-6.51%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11301	0.05459	4.6291	80.7407	44.7696	106.1565	108.579	-19.26%	-10.46%	-29.23%	-27.61%	1.803	-9.83%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10877	0.05459	4.4553	84.3095	53.3212	110.023	112.5374	-15.69%	6.64%	-26.65%	-24.98%	1.581	-20.94%

Table B.39. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13_h0			
SC			100	20	120				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11735	0.05476	4.8067	78.6999	13.2566	85.6572	88.5517	-21.30%	-33.72%	-28.62%	-26.21%	5.937	18.73%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11676	0.05476	4.7824	79.714	18.1274	89.231	92.2693	-20.29%	-9.36%	-25.64%	-23.11%	4.397	-12.05%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11607	0.05476	4.7543	81.4959	19.7867	90.1821	93.258	-18.50%	-1.07%	-24.85%	-22.29%	4.119	-17.63%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11529	0.05476	4.7224	82.1813	21.6415	90.3685	93.4577	-17.82%	8.21%	-24.69%	-22.12%	3.797	-24.05%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11126	0.05476	4.5571	87.1522	25.508	89.1781	92.2362	-12.85%	27.54%	-25.68%	-23.14%	3.417	-31.67%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10708	0.05476	4.386	91.6223	34.3167	93.5302	96.7579	-8.38%	71.58%	-22.06%	-19.37%	2.670	-46.60%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11735	0.05476	4.8067	75.9637	5.2106	75.5796	77.2527	-24.04%	-73.95%	-37.02%	-35.62%	14.579	191.57%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11676	0.05476	4.7824	76.4707	9.9149	78.7501	80.5067	-23.53%	-50.43%	-34.37%	-32.91%	7.713	54.25%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11607	0.05476	4.7543	78.0644	11.3725	79.596	81.3746	-21.94%	-43.14%	-33.67%	-32.19%	6.864	37.29%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11529	0.05476	4.7224	78.55	13.1532	79.7738	81.5604	-21.45%	-34.23%	-33.52%	-32.03%	5.972	19.44%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11126	0.05476	4.5571	83.0791	16.5513	78.7657	80.5353	-16.92%	-17.24%	-34.36%	-32.89%	5.019	0.39%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10708	0.05476	4.386	86.654	24.8987	82.6633	84.5323	-13.35%	24.49%	-31.11%	-29.56%	3.480	-30.39%

Table B.40. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Jo13_hi0			
SC			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025969	0.012362	1.0637	31.7964	31.7102	62.7441	64.937	-20.51%	-20.72%	-21.57%	-18.83%	1.003	0.27%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025837	0.012362	1.0583	33.1671	33.8171	65.6618	67.9566	-17.08%	-15.46%	-17.92%	-15.05%	0.981	-1.92%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025686	0.012362	1.0521	34.015	34.4553	66.4425	68.7646	-14.96%	-13.86%	-16.95%	-14.04%	0.987	-1.28%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025514	0.012362	1.045	34.1786	35.2579	66.8	69.1345	-14.55%	-11.86%	-16.50%	-13.58%	0.969	-3.06%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02462	0.012362	1.0084	35.4654	36.4602	66.7302	69.0623	-11.34%	-8.85%	-16.59%	-13.67%	0.973	-2.73%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023696	0.012362	0.97059	36.4432	38.6802	67.6592	70.0237	-8.89%	-3.30%	-15.43%	-12.47%	0.942	-5.78%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025969	0.012362	1.0637	28.0789	27.9836	55.3884	56.6567	-29.80%	-30.04%	-30.76%	-29.18%	1.003	0.34%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025837	0.012362	1.0583	29.2169	29.9228	57.9694	59.2969	-26.96%	-25.19%	-27.54%	-25.88%	0.976	-2.36%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025686	0.012362	1.0521	29.9906	30.4663	58.661	60.0042	-25.02%	-23.83%	-26.67%	-24.99%	0.984	-1.56%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025514	0.012362	1.045	30.0717	31.2448	58.9794	60.3299	-24.82%	-21.89%	-26.28%	-24.59%	0.962	-3.75%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02462	0.012362	1.0084	31.2257	32.3225	58.929	60.2784	-21.94%	-19.19%	-26.34%	-24.65%	0.966	-3.39%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023696	0.012362	0.97059	31.9795	34.4427	59.7623	61.1308	-20.05%	-13.89%	-25.30%	-23.59%	0.928	-7.15%	

Table B.41. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Jo13_hi0			
SC			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025969	0.012362	1.0637	31.7964	31.7102	62.7441	64.937	-20.51%	-20.72%	-21.57%	-18.83%	1.003	0.27%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025837	0.012362	1.0583	33.1671	33.8171	65.6618	67.9566	-17.08%	-15.46%	-17.92%	-15.05%	0.981	-1.92%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025686	0.012362	1.0521	34.015	34.4553	66.4425	68.7646	-14.96%	-13.86%	-16.95%	-14.04%	0.987	-1.28%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025514	0.012362	1.045	34.1786	35.2579	66.8	69.1345	-14.55%	-11.86%	-16.50%	-13.58%	0.969	-3.06%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02462	0.012362	1.0084	35.4654	36.4602	66.7302	69.0623	-11.34%	-8.85%	-16.59%	-13.67%	0.973	-2.73%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023696	0.012362	0.97059	36.4432	38.6802	67.6592	70.0237	-8.89%	-3.30%	-15.43%	-12.47%	0.942	-5.78%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025969	0.012362	1.0637	28.0789	27.9836	55.3884	56.6567	-29.80%	-30.04%	-30.76%	-29.18%	1.003	0.34%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025837	0.012362	1.0583	29.2169	29.9228	57.9694	59.2969	-26.96%	-25.19%	-27.54%	-25.88%	0.976	-2.36%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025686	0.012362	1.0521	29.9906	30.4663	58.661	60.0042	-25.02%	-23.83%	-26.67%	-24.99%	0.984	-1.56%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025514	0.012362	1.045	30.0717	31.2448	58.9794	60.3299	-24.82%	-21.89%	-26.28%	-24.59%	0.962	-3.75%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02462	0.012362	1.0084	31.2257	32.3225	58.929	60.2784	-21.94%	-19.19%	-26.34%	-24.65%	0.966	-3.39%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023696	0.012362	0.97059	31.9795	34.4427	59.7623	61.1308	-20.05%	-13.89%	-25.30%	-23.59%	0.928	-7.15%	

Table B.42. ^{123}I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Jo13_h10			
SC			40	20	60					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025536	0.012396	1.0459	31.7171	15.6206	46.4387	48.0494	-20.71%	-21.90%	-22.60%	-19.92%	2.030	1.52%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025406	0.012396	1.0407	32.7318	17.3249	48.5846	50.2722	-18.17%	-13.38%	-19.03%	-16.21%	1.889	-5.54%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025257	0.012396	1.0345	33.4142	17.8591	49.1569	50.8646	-16.46%	-10.70%	-18.07%	-15.23%	1.871	-6.45%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025088	0.012396	1.0276	33.5779	18.5158	49.398	51.1151	-16.06%	-7.42%	-17.67%	-14.81%	1.813	-9.33%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02421	0.012396	0.99163	34.8486	19.5876	49.3297	51.0451	-12.88%	-2.06%	-17.78%	-14.92%	1.779	-11.04%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023301	0.012396	0.9544	35.8415	21.7258	50.2739	52.0244	-10.40%	8.63%	-16.21%	-13.29%	1.650	-17.51%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025536	0.012396	1.0459	29.5958	12.1928	40.9906	41.9221	-26.01%	-39.04%	-31.68%	-30.13%	2.427	21.37%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025406	0.012396	1.0407	30.4227	13.7718	42.8894	43.8654	-23.94%	-31.14%	-28.52%	-26.89%	2.209	10.45%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025257	0.012396	1.0345	31.0433	14.2295	43.3964	44.3842	-22.39%	-28.85%	-27.67%	-26.03%	2.182	9.08%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025088	0.012396	1.0276	31.1385	14.8649	43.6123	44.6056	-22.15%	-25.68%	-27.31%	-25.66%	2.095	4.74%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02421	0.012396	0.99163	32.2893	15.8145	43.5631	44.5557	-19.28%	-20.93%	-27.39%	-25.74%	2.042	2.09%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023301	0.012396	0.9544	33.0657	17.8517	44.4099	45.423	-17.34%	-10.74%	-25.98%	-24.30%	1.852	-7.39%	

Table B.43. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW I ₀₁₃ hi6			
SC			40	8	48				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025239	0.012468	1.0338	31.6759	5.9729	36.2399	37.4679	-20.81%	-25.34%	-24.50%	-21.94%	5.303	6.07%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025111	0.012468	1.0286	32.4781	7.4382	37.9277	39.2193	-18.80%	-7.02%	-20.98%	-18.29%	4.366	-12.67%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024964	0.012468	1.0225	33.0654	7.9102	38.3698	39.6779	-17.34%	-1.12%	-20.06%	-17.34%	4.180	-16.40%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024796	0.012468	1.0157	33.2367	8.4788	38.5455	39.862	-16.91%	5.99%	-19.70%	-16.95%	3.920	-21.60%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023928	0.012468	0.9801	34.5228	9.4726	38.4579	39.7741	-13.69%	18.41%	-19.88%	-17.14%	3.644	-27.11%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02303	0.012468	0.94331	35.5459	11.5703	39.4131	40.7683	-11.14%	44.63%	-17.89%	-15.07%	3.072	-38.56%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025239	0.012468	1.0338	30.5114	2.7236	31.9793	32.6892	-23.72%	-65.96%	-33.38%	-31.90%	11.203	124.05%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025111	0.012468	1.0286	31.1529	4.0888	33.4737	34.2205	-22.12%	-48.89%	-30.26%	-28.71%	7.619	52.38%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024964	0.012468	1.0225	31.6857	4.4951	33.866	34.6224	-20.79%	-43.81%	-29.45%	-27.87%	7.049	40.98%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024796	0.012468	1.0157	31.7967	5.0432	34.0246	34.7859	-20.51%	-36.96%	-29.12%	-27.53%	6.305	26.10%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023928	0.012468	0.9801	32.9706	5.9137	33.9592	34.7207	-17.57%	-26.08%	-29.25%	-27.67%	5.575	11.51%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02303	0.012468	0.94331	33.781	7.9071	34.8165	35.6008	-15.55%	-1.16%	-27.47%	-25.83%	4.272	-14.56%

Table B.44. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW, lo13, hi6			
SC			84	84	168				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026029	0.011882	1.0662	67.7616	67.5753	134.474	139.1738	-19.33%	-19.55%	-19.96%	-17.16%	1.003	0.28%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025997	0.011882	1.0608	71.1875	71.3448	140.9399	145.8656	-15.25%	-15.07%	-16.11%	-13.18%	0.998	-0.22%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025745	0.011882	1.0545	72.7302	72.4403	142.6562	147.6419	-13.42%	-13.76%	-15.09%	-12.12%	1.004	0.40%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025573	0.011882	1.0475	72.9822	73.7521	143.3888	148.4001	-13.12%	-12.20%	-14.65%	-11.67%	0.990	-1.04%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024677	0.011882	1.0108	74.5456	75.4737	143.905	148.9343	-11.26%	-10.15%	-14.34%	-11.35%	0.988	-1.23%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023751	0.011882	0.97284	75.3863	78.081	144.8661	149.9289	-10.25%	-7.05%	-13.77%	-10.76%	0.965	-3.45%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026029	0.011882	1.0662	59.8391	59.6359	118.7126	121.431	-28.76%	-29.00%	-29.34%	-27.72%	1.003	0.34%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025897	0.011882	1.0608	62.8346	63.0065	124.4323	127.2817	-25.20%	-24.99%	-25.93%	-24.24%	0.997	-0.27%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025745	0.011882	1.0545	64.2493	63.9303	125.9528	128.8371	-23.51%	-23.89%	-25.03%	-23.31%	1.005	0.50%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025573	0.011882	1.0475	64.3669	65.2018	126.6035	129.5026	-23.37%	-22.38%	-24.64%	-22.92%	0.987	-1.28%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.024677	0.011882	1.0108	65.7425	66.7686	127.0747	129.9846	-21.74%	-20.51%	-24.36%	-22.63%	0.985	-1.54%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023751	0.011882	0.97284	66.3223	69.2873	127.938	130.8676	-21.04%	-17.52%	-23.85%	-22.10%	0.957	-4.28%

Table B.45. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I ¹²³ I hi6			
SC			84	42	126					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025307	0.011662	1.0366	67.5633	33.6752	100.1128	103.586	-19.57%	-19.82%	-20.55%	-17.79%	2.006	0.32%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025179	0.011662	1.0313	70.2315	36.5859	105.0458	108.6933	-16.39%	-12.89%	-16.63%	-13.74%	1.920	-4.02%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025031	0.011662	1.0253	71.4137	37.4521	106.3065	109.9982	-14.98%	-10.83%	-15.63%	-12.70%	1.907	-4.66%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024864	0.011662	1.0184	71.6527	38.4461	106.8294	110.5409	-14.70%	-8.46%	-15.21%	-12.27%	1.864	-6.81%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023993	0.011662	0.98275	73.1509	39.8106	107.2096	110.9355	-12.92%	-5.21%	-14.91%	-11.96%	1.837	-8.13%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023092	0.011662	0.94586	73.9931	42.1404	108.0685	111.8274	-11.91%	0.33%	-14.23%	-11.25%	1.756	-12.21%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025307	0.011662	1.0366	63.0041	26.3684	88.371	90.3798	-25.00%	-37.22%	-29.86%	-28.27%	2.389	19.47%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025179	0.011662	1.0313	65.3376	28.9708	92.7354	94.8452	-22.22%	-31.02%	-26.40%	-24.73%	2.255	12.76%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025031	0.011662	1.0253	66.4196	29.7046	93.8527	95.9881	-20.93%	-29.27%	-25.51%	-23.82%	2.236	11.80%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024864	0.011662	1.0184	66.555	30.665	94.3179	96.4649	-20.77%	-26.99%	-25.14%	-23.44%	2.170	8.52%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023993	0.011662	0.98275	67.9009	31.8836	94.6681	96.8238	-19.17%	-24.09%	-24.87%	-23.16%	2.130	6.48%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023092	0.011662	0.94586	68.5085	34.1264	95.4402	97.6153	-18.44%	-18.75%	-24.25%	-22.53%	2.007	0.37%	

Table B.46. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only;
Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW I ¹²³ I h16			
SC			84	16.8	100.8					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024742	0.011734	1.0134	67.4662	13.3574	78.8298	81.5042	-19.68%	-20.49%	-21.80%	-19.14%	5.051	1.02%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024617	0.011734	1.0083	69.6861	15.7603	82.81	85.6291	-17.04%	-6.19%	-17.85%	-15.05%	4.422	-11.57%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024472	0.011734	1.0024	70.6606	16.4961	83.759	86.6126	-15.88%	-1.81%	-16.91%	-14.07%	4.283	-14.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024308	0.011734	0.99565	70.9025	17.307	84.1544	87.025	-15.59%	3.02%	-16.51%	-13.67%	4.097	-18.06%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023457	0.011734	0.96079	72.3913	18.5079	84.4683	87.3535	-13.82%	10.17%	-16.20%	-13.34%	3.911	-21.77%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022576	0.011734	0.92473	73.2666	20.7018	85.3166	88.2389	-12.78%	23.23%	-15.36%	-12.46%	3.539	-29.22%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024742	0.011734	1.0134	64.9225	6.4274	69.5666	71.1127	-22.71%	-61.74%	-30.99%	-29.45%	10.101	102.02%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024617	0.011734	1.0083	66.8643	8.5757	73.0891	74.7192	-20.40%	-48.95%	-27.49%	-25.87%	7.797	55.94%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024472	0.011734	1.0024	67.7543	9.2019	73.931	75.5812	-19.34%	-45.23%	-26.66%	-25.02%	7.363	47.26%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024308	0.011734	0.99565	67.9102	9.9823	74.284	75.944	-19.15%	-40.58%	-26.31%	-24.66%	6.803	36.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023457	0.011734	0.96079	69.2628	11.0391	74.5757	76.2446	-17.54%	-34.29%	-26.02%	-24.36%	6.274	25.49%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022576	0.011734	0.92473	69.9166	13.1431	75.3388	77.0294	-16.77%	-21.77%	-25.26%	-23.58%	5.320	6.39%	

Table B.47. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	126.4	252.8				114-126			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024309	0.004264	0.099571	122.1541	121.9801	242.5674	251.5711	-3.36%	-3.50%	-4.05%	-0.49%	1.001	0.14%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024186	0.004264	0.099067	129.0229	126.97	253.6408	263.0555	2.08%	0.45%	0.33%	4.06%	1.016	1.62%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024044	0.004264	0.098486	129.9295	130.3973	256.9172	266.4536	2.79%	3.16%	1.63%	5.40%	0.996	-0.36%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0023883	0.004264	0.097825	130.2773	132.199	258.1007	267.6809	3.07%	4.59%	2.10%	5.89%	0.985	-1.45%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0023047	0.004264	0.0944	131.4717	134.292	258.559	268.1562	4.01%	6.24%	2.28%	6.07%	0.979	-2.10%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0022182	0.004264	0.090857	132.8004	134.8142	258.4471	268.0401	5.06%	6.66%	2.23%	6.03%	0.985	-1.49%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024309	0.004264	0.099571	105.751	105.5683	209.9616	214.7696	-16.34%	-16.48%	-16.95%	-15.04%	1.002	0.17%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024186	0.004264	0.099067	111.9197	109.6901	219.5693	224.5973	-11.46%	-13.22%	-13.15%	-11.16%	1.020	2.03%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024044	0.004264	0.098486	112.4226	112.9561	222.4192	227.5125	-11.06%	-10.64%	-12.02%	-10.00%	0.995	-0.47%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0023883	0.004264	0.097825	112.563	114.6884	223.4514	228.5683	-10.95%	-9.27%	-11.61%	-9.59%	0.981	-1.85%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0023047	0.004264	0.0944	113.5098	116.635	223.8621	228.9883	-10.20%	-7.73%	-11.45%	-9.42%	0.973	-2.68%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0022182	0.004264	0.090857	114.7807	117.012	223.7734	228.8977	-9.19%	-7.43%	-11.48%	-9.46%	0.981	-1.91%

Table B.48. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW				
SC			126.4	63.2	189.6				114-126				
Atten Coeff	0.1526								126-154				
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025919	0.003843	0.10617	123.0246	60.3577	181.258	187.9319	-2.67%	-4.50%	-4.40%	-0.88%	2.038	1.91%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025788	0.003843	0.10563	128.4413	63.8925	189.5985	196.5816	1.61%	1.10%	0.00%	3.68%	2.010	0.51%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025637	0.003843	0.10501	128.791	66.8838	192.1715	199.2556	1.89%	5.83%	1.36%	5.09%	1.926	-3.72%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025465	0.003843	0.1043	128.4585	68.9017	192.9749	200.0933	1.63%	9.02%	1.78%	5.53%	1.864	-6.78%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024573	0.003843	0.10065	128.3028	71.8116	193.2599	200.3948	1.51%	13.63%	1.93%	5.69%	1.787	-10.67%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0023651	0.003843	0.096874	128.9437	72.9093	193.2381	200.3736	2.01%	15.36%	1.92%	5.68%	1.769	-11.57%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025919	0.003843	0.10617	113.5313	45.2017	156.8771	160.4404	-10.18%	-28.48%	-17.26%	-15.38%	2.512	25.58%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025788	0.003843	0.10563	118.4389	48.062	164.1139	167.8427	-6.30%	-23.95%	-13.44%	-11.48%	2.464	23.21%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025637	0.003843	0.10501	118.4434	50.9631	166.3534	170.1365	-6.29%	-19.36%	-12.26%	-10.27%	2.324	16.21%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025465	0.003843	0.1043	117.8892	52.9865	167.0567	170.8582	-6.73%	-16.16%	-11.89%	-9.88%	2.225	11.24%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024573	0.003843	0.10065	117.4057	55.8947	167.316	171.1266	-7.12%	-11.56%	-11.75%	-9.74%	2.100	5.02%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0023651	0.003843	0.096874	117.9158	56.9304	167.3044	171.1156	-6.71%	-9.92%	-11.76%	-9.75%	2.071	3.56%	

Table B.49. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	25.28	151.68				114-126			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026355	0.003576	0.10795	123.5559	23.3933	142.9879	148.1255	-2.25%	-7.46%	-5.73%	-2.34%	5.282	5.63%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026222	0.003576	0.1074	128.1041	26.0578	149.5934	154.9787	1.35%	3.08%	-1.38%	2.17%	4.916	-1.68%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026068	0.003576	0.10677	128.1207	28.7919	151.6351	157.1088	1.36%	13.89%	-0.03%	3.58%	4.450	-11.00%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0025893	0.003576	0.10606	127.3838	30.942	152.1716	157.6766	0.78%	22.40%	0.32%	3.95%	4.117	-17.66%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024986	0.003576	0.10234	126.4165	34.3706	152.2803	157.807	0.01%	35.96%	0.40%	4.04%	3.678	-26.44%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024048	0.003576	0.098502	126.7029	35.8124	152.3755	157.9125	0.24%	41.66%	0.46%	4.11%	3.538	-29.24%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026355	0.003576	0.10795	118.2072	8.9896	123.7161	126.4579	-6.48%	-64.44%	-18.44%	-16.63%	13.149	162.99%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026222	0.003576	0.1074	122.3606	11.0954	129.4488	132.3232	-3.20%	-56.11%	-14.66%	-12.76%	11.028	120.56%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0026068	0.003576	0.10677	122.0673	13.7813	131.2287	134.1506	-3.43%	-45.49%	-13.48%	-11.56%	8.857	77.15%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0025893	0.003576	0.10606	121.0997	15.9817	131.7023	134.6412	-4.19%	-36.78%	-13.17%	-11.23%	7.577	51.55%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024986	0.003576	0.10234	119.7535	19.4957	131.8113	134.7622	-5.26%	-22.88%	-13.10%	-11.15%	6.143	22.85%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0024048	0.003576	0.098502	119.8643	20.9203	131.9031	134.8597	-5.17%	-17.25%	-13.04%	-11.09%	5.730	14.59%

Table B.50. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	126.4	252.8				121-126			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0036038	0.004847	0.14761	121.009	117.4222	236.3864	245.1606	-4.27%	-7.10%	-6.49%	-3.02%	1.031	3.05%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035856	0.004847	0.14686	127.1651	122.9043	246.9302	256.0956	0.61%	-2.77%	-2.32%	1.30%	1.035	3.47%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035645	0.004847	0.146	128.3139	126.0127	250.0377	259.3186	1.51%	-0.31%	-1.09%	2.58%	1.018	1.83%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035406	0.004847	0.14502	128.7564	127.6435	251.1871	260.5108	1.86%	0.98%	-0.64%	3.05%	1.009	0.87%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0034166	0.004847	0.13995	131.1375	128.3484	251.0173	260.3345	3.75%	1.54%	-0.71%	2.98%	1.022	2.17%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0032884	0.004847	0.13469	132.5675	128.9327	251.3074	260.6354	4.88%	2.00%	-0.59%	3.10%	1.028	2.82%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0036038	0.004847	0.14761	105.1495	101.2334	204.6106	209.296	-16.81%	-19.91%	-19.06%	-17.21%	1.039	3.87%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035856	0.004847	0.14686	110.5663	105.9156	213.7574	218.6522	-12.53%	-16.21%	-15.44%	-13.51%	1.044	4.39%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035645	0.004847	0.146	111.3396	108.8437	216.4594	221.4161	-11.91%	-13.89%	-14.38%	-12.41%	1.023	2.29%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035406	0.004847	0.14502	111.5933	110.3971	217.4615	222.4412	-11.71%	-12.66%	-13.98%	-12.01%	1.011	1.08%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0034166	0.004847	0.13995	113.8699	110.8404	217.3282	222.3049	-9.91%	-12.31%	-14.03%	-12.06%	1.027	2.73%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0032884	0.004847	0.13469	115.231	111.2657	217.591	222.5736	-8.84%	-11.97%	-13.93%	-11.96%	1.036	3.56%

Table B.51. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW				
SC			126.4	63.2	189.6				121-126				
Atten Coeff	0.1526								126-154				
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037983	0.004339	0.15558	121.1368	57.8199	176.4159	182.9077	-4.16%	-8.51%	-6.95%	-3.53%	2.095	4.75%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037791	0.004339	0.15479	126.078	61.6881	184.3318	191.1186	-0.25%	-2.39%	-2.78%	0.80%	2.044	2.19%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037569	0.004339	0.15388	126.617	64.3986	186.6488	193.5265	0.17%	1.90%	-1.56%	2.07%	1.966	-1.69%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037317	0.004339	0.15285	126.3883	66.1928	187.4416	194.3527	-0.01%	4.74%	-1.14%	2.51%	1.909	-4.53%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036011	0.004339	0.1475	127.2687	67.7016	187.2641	194.171	0.69%	7.12%	-1.23%	2.41%	1.880	-6.01%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0034659	0.004339	0.14196	128.0684	68.5877	187.4767	194.3925	1.32%	8.52%	-1.12%	2.53%	1.867	-6.64%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037983	0.004339	0.15558	111.9736	42.9284	152.684	156.15	-11.41%	-32.08%	-19.47%	-17.64%	2.608	30.42%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037791	0.004339	0.15479	116.3784	46.1678	159.5514	163.1752	-7.93%	-26.95%	-15.85%	-13.94%	2.521	26.04%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037569	0.004339	0.15388	116.5984	48.7733	161.5675	165.2402	-7.75%	-22.83%	-14.79%	-12.85%	2.391	19.53%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037317	0.004339	0.15285	116.1726	50.5651	162.262	165.9527	-8.09%	-19.99%	-14.42%	-12.47%	2.297	14.87%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036011	0.004339	0.1475	116.8687	51.9745	162.1189	165.8076	-7.54%	-17.76%	-14.49%	-12.55%	2.249	12.43%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0034659	0.004339	0.14196	117.5585	52.7795	162.3117	166.0054	-6.99%	-16.49%	-14.39%	-12.44%	2.227	11.37%	

Table B.52. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	25.28	151.68				121-126			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0037448	0.004045	0.15339	121.2426	22.0873	138.9688	143.957	-4.08%	-12.63%	-8.38%	-5.09%	5.489	9.78%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0037258	0.004045	0.15261	125.4622	24.9974	145.2992	150.5271	-0.74%	-1.12%	-4.21%	-0.76%	5.019	0.38%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.003704	0.004045	0.15172	125.6449	27.479	147.1754	152.4842	-0.60%	8.70%	-2.97%	0.53%	4.572	-8.55%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0036791	0.004045	0.1507	125.0281	29.3821	147.7494	153.0895	-1.09%	16.23%	-2.59%	0.93%	4.255	-14.90%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035503	0.004045	0.14542	125.0803	31.419	147.5087	152.8507	-1.04%	24.28%	-2.75%	0.77%	3.981	-20.38%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0034171	0.004045	0.13996	125.5566	32.5575	147.7424	153.0983	-0.67%	28.79%	-2.60%	0.94%	3.856	-22.87%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0037448	0.004045	0.15339	116.0932	7.9704	120.2359	122.8978	-8.15%	-68.47%	-20.73%	-18.98%	14.566	191.31%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0037258	0.004045	0.15261	119.8974	10.353	125.7297	128.5199	-5.14%	-59.05%	-17.11%	-15.27%	11.581	131.62%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.003704	0.004045	0.15172	119.7939	12.7737	127.3649	130.1987	-5.23%	-49.47%	-16.03%	-14.16%	9.378	87.56%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0036791	0.004045	0.1507	118.9734	14.7179	127.8714	130.7221	-5.88%	-41.78%	-15.70%	-13.82%	8.084	61.67%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0035503	0.004045	0.14542	118.7873	16.7448	127.6729	130.5248	-6.02%	-33.76%	-15.83%	-13.95%	7.094	41.88%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0034171	0.004045	0.13996	119.1217	17.842	127.8845	130.7441	-5.76%	-29.42%	-15.69%	-13.80%	6.676	33.53%

Table B.53. ^{99m}Tc TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	126.4	252.8				123-128			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.004918	0.004983	0.20144	112.2099	129.3195	239.8634	248.7636	-11.23%	2.31%	-5.12%	-1.60%	0.8677	-13.23%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048931	0.004983	0.20042	117.9607	135.8845	251.3698	260.697	-6.68%	7.50%	-0.57%	3.12%	0.8681	-13.19%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048644	0.004983	0.19925	119.374	139.2292	255.065	264.5287	-5.56%	10.15%	0.90%	4.64%	0.8574	-14.26%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048318	0.004983	0.19791	119.595	141.5103	256.6621	266.1843	-5.38%	11.95%	1.53%	5.29%	0.8451	-15.49%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0046626	0.004983	0.19098	122.1236	144.4359	259.2899	268.9096	-3.38%	14.27%	2.57%	6.37%	0.8455	-15.45%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0044876	0.004983	0.18381	124.9645	144.465	260.1057	269.7569	-1.14%	14.29%	2.89%	6.71%	0.865	-13.50%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.004918	0.004983	0.20144	95.1927	113.8771	207.6252	212.378	-24.69%	-9.91%	-17.87%	-15.99%	0.8359	-16.41%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048931	0.004983	0.20042	100.0941	119.6635	217.61	222.5914	-20.81%	-5.33%	-13.92%	-11.95%	0.8365	-16.35%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048644	0.004983	0.19925	101.1047	122.7903	220.8234	225.878	-20.01%	-2.86%	-12.65%	-10.65%	0.8234	-17.66%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0048318	0.004983	0.19791	101.065	125.0096	222.2143	227.3003	-20.04%	-1.10%	-12.10%	-10.09%	0.8085	-19.15%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0046626	0.004983	0.19098	103.2335	127.6267	224.5192	229.658	-18.33%	0.97%	-11.19%	-9.15%	0.8089	-19.11%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0044876	0.004983	0.18381	106.0533	127.3478	225.2428	230.3988	-16.10%	0.75%	-10.90%	-8.86%	0.8328	-16.72%

Table B.54. ^{99m}Tc TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW			
SC			126.4	63.2	189.6				123-128			
Atten Coeff	0.1526								126-154			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0042104	0.004278	0.17246	112.6936	61.8616	172.3912	178.7529	-10.84%	-2.12%	-9.08%	-5.72%	1.822	-8.91%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041891	0.004278	0.17159	116.8185	66.4623	180.4283	187.0909	-7.58%	5.16%	-4.84%	-1.32%	1.758	-12.12%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041645	0.004278	0.17058	117.5502	68.9727	182.8445	189.6	-7.00%	9.13%	-3.56%	0.00%	1.704	-14.78%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041366	0.004278	0.16943	117.2186	70.9621	183.7536	190.5462	-7.26%	12.28%	-3.08%	0.50%	1.652	-17.41%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0039918	0.004278	0.1635	118.1725	73.4321	184.9008	191.7389	-6.51%	16.19%	-2.48%	1.13%	1.609	-19.54%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0038419	0.004278	0.15737	119.7367	73.8036	185.3014	192.1535	-5.27%	16.78%	-2.27%	1.35%	1.622	-18.88%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0042104	0.004278	0.17246	103.2513	47.8447	149.2107	152.6075	-18.31%	-24.30%	-21.30%	-19.51%	2.158	7.90%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041891	0.004278	0.17159	106.7774	51.8908	156.1857	159.7436	-15.52%	-17.89%	-17.62%	-15.75%	2.058	2.89%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041645	0.004278	0.17058	107.2118	54.2758	158.288	161.8958	-15.18%	-14.12%	-16.51%	-14.61%	1.975	-1.23%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0041366	0.004278	0.16943	106.6615	56.2714	159.082	162.7098	-15.62%	-10.96%	-16.10%	-14.18%	1.895	-5.23%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0039918	0.004278	0.1635	107.3163	58.6233	160.0922	163.7447	-15.10%	-7.24%	-15.56%	-13.64%	1.831	-8.47%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0038419	0.004278	0.15737	108.8188	58.8368	160.4503	164.1106	-13.91%	-6.90%	-15.37%	-13.44%	1.850	-7.52%

Table B.55. ^{99m}Tc TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission Data Only; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW				
SC			126.4	25.28	151.68				123-128				
Atten Coeff	0.1526								126-154				
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036083	0.003891	0.1478	113.0128	21.418	130.1735	134.851	-10.59%	-15.28%	-14.18%	-11.10%	5.277	5.53%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.00359	0.003891	0.14705	116.1713	24.85	136.0223	140.9258	-8.09%	-1.70%	-10.32%	-7.09%	4.675	-6.50%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003569	0.003891	0.14619	116.5041	26.8705	137.5274	142.4954	-7.83%	6.29%	-9.33%	-6.06%	4.336	-13.28%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003545	0.003891	0.14521	115.8585	28.6934	137.9892	142.9838	-8.34%	13.50%	-9.03%	-5.73%	4.038	-19.24%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0034209	0.003891	0.14012	115.9431	30.9426	138.2309	143.2453	-8.27%	22.40%	-8.87%	-5.56%	3.747	-25.06%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0032925	0.003891	0.13486	116.8147	31.5656	138.4035	143.4265	-7.58%	24.86%	-8.75%	-5.44%	3.701	-25.99%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036083	0.003891	0.1478	108.1115	8.2522	112.632	115.1283	-14.47%	-67.36%	-25.74%	-24.10%	13.101	162.02%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.00359	0.003891	0.14705	110.8206	11.2628	117.7099	120.3273	-12.33%	-55.45%	-22.40%	-20.67%	9.840	96.79%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003569	0.003891	0.14619	110.9181	13.2121	119.0217	121.6738	-12.25%	-47.74%	-21.53%	-19.78%	8.395	67.90%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003545	0.003891	0.14521	110.0776	15.0801	119.4288	122.0952	-12.91%	-40.35%	-21.26%	-19.50%	7.300	45.99%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0034209	0.003891	0.14012	109.8959	17.3142	119.6503	122.3276	-13.06%	-31.51%	-21.12%	-19.35%	6.347	26.94%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0032925	0.003891	0.13486	110.6786	17.8582	119.8085	122.4906	-12.44%	-29.36%	-21.01%	-19.24%	6.198	23.95%	

Table B.56. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6		
SC-Total			100	100	200					130-142		
Atten Coeff	0.1481									142-175		
Striata ROI Size	13 13									175-180		
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11025	0.072623	4.5158	97.4762	97.7439	189.2031	195.8155	-2.52%	-2.26%	-5.40%	-2.09%	0.997	-0.27%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10969	0.072623	4.4929	100.9866	106.6619	197.7173	204.627	0.99%	6.66%	-1.14%	2.31%	0.947	-5.32%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10905	0.072623	4.4665	104.731	110.5821	201.6383	208.685	4.73%	10.58%	0.82%	4.34%	0.947	-5.29%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10831	0.072623	4.4366	107.3276	112.4934	202.7567	209.8425	7.33%	12.49%	1.38%	4.92%	0.954	-4.59%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10452	0.072623	4.2812	117.1366	120.0744	205.0283	212.1937	17.14%	20.07%	2.51%	6.10%	0.976	-2.45%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1006	0.072623	4.1205	122.0872	136.5223	211.8402	219.242	22.09%	36.52%	5.92%	9.62%	0.894	-10.57%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11025	0.072623	4.5158	86.0357	86.3047	167.022	170.8467	-13.96%	-13.70%	-16.49%	-14.58%	0.997	-0.31%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10969	0.072623	4.4929	88.5953	94.7462	174.5575	178.5546	-11.40%	-5.25%	-12.72%	-10.72%	0.935	-6.49%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10905	0.072623	4.4665	91.9005	98.2352	178.0319	182.1085	-8.10%	-1.76%	-10.98%	-8.95%	0.936	-6.45%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10831	0.072623	4.4366	94.2695	99.8942	179.0377	183.1373	-5.74%	-0.11%	-10.48%	-8.43%	0.944	-5.64%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10452	0.072623	4.2812	103.2372	106.4922	181.1243	185.2719	3.24%	6.49%	-9.44%	-7.36%	0.969	-3.06%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1006	0.072623	4.1205	106.5361	122.4382	187.2329	191.5194	6.54%	22.44%	-6.38%	-4.24%	0.870	-12.99%

Table B.57. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			100	50	150					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11023	0.073198	4.515	96.4273	49.061	138.4267	143.231	-3.57%	-1.88%	-7.72%	-4.51%	1.965	-1.73%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10967	0.073198	4.4922	98.6111	56.9512	144.6842	149.7173	-1.39%	13.90%	-3.54%	-0.19%	1.732	-13.42%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10903	0.073198	4.4658	101.7349	60.6848	148.0132	153.1647	1.73%	21.37%	-1.32%	2.11%	1.676	-16.18%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1083	0.073198	4.4359	104.2104	62.2691	148.7317	153.9066	4.21%	24.54%	-0.85%	2.61%	1.674	-16.32%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10451	0.073198	4.2805	113.463	70.1655	151.2172	156.4837	13.46%	40.33%	0.81%	4.32%	1.617	-19.15%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10058	0.073198	4.1199	117.9275	86.9873	158.1593	163.6797	17.93%	73.97%	5.44%	9.12%	1.356	-32.22%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11023	0.073198	4.515	89.8296	38.6067	122.1861	124.9645	-10.17%	-22.79%	-18.54%	-16.69%	2.327	16.34%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10967	0.073198	4.4922	91.1849	46.1702	127.7284	130.6397	-8.82%	-7.66%	-14.85%	-12.91%	1.975	-1.25%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10903	0.073198	4.4658	93.8959	49.538	130.6797	133.6598	-6.10%	-0.92%	-12.88%	-10.89%	1.895	-5.23%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1083	0.073198	4.4359	96.1664	50.8889	131.3321	134.3272	-3.83%	1.78%	-12.45%	-10.45%	1.890	-5.51%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10451	0.073198	4.2805	104.5614	57.852	133.6082	136.657	4.56%	15.70%	-10.93%	-8.90%	1.807	-9.63%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10058	0.073198	4.1199	107.3363	74.2181	139.8342	143.0322	7.34%	48.44%	-6.78%	-4.65%	1.446	-27.69%	

Table B.58. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW to13 hi6			
SC-Total			100	20	120					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11029	0.073357	4.5175	95.7846	19.8386	106.4539	110.0694	-4.22%	-0.81%	-11.29%	-8.28%	4.828	-3.44%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10973	0.073357	4.4947	97.1704	27.1046	111.3131	115.1231	-2.83%	35.52%	-7.24%	-4.06%	3.585	-28.30%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10909	0.073357	4.4683	99.9336	30.7093	114.221	118.1404	-0.07%	53.55%	-4.82%	-1.55%	3.254	-34.92%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10836	0.073357	4.4383	102.3647	32.0652	114.6071	118.5423	2.36%	60.33%	-4.49%	-1.21%	3.192	-36.15%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10456	0.073357	4.2829	111.3355	40.0815	116.7118	120.7341	11.34%	100.41%	-2.74%	0.61%	2.778	-44.45%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10064	0.073357	4.1222	115.4869	57.1184	123.7668	128.0634	15.49%	185.59%	3.14%	6.72%	2.022	-59.56%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11029	0.073357	4.5175	92.0942	9.977	93.9386	96.0266	-7.91%	-50.12%	-21.72%	-19.98%	9.231	84.61%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10973	0.073357	4.4947	92.7261	17.0065	98.2495	100.4527	-7.27%	-14.97%	-18.13%	-16.29%	5.452	9.05%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10909	0.073357	4.4683	95.0939	20.2839	100.8304	103.0971	-4.91%	1.42%	-15.97%	-14.09%	4.688	-6.24%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10836	0.073357	4.4383	97.3422	21.4152	101.1895	103.4658	-2.66%	7.08%	-15.68%	-13.78%	4.545	-9.09%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10456	0.073357	4.2829	105.4457	28.5228	103.1318	105.4604	5.45%	42.61%	-14.06%	-12.12%	3.697	-26.06%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10064	0.073357	4.1222	107.888	45.131	109.4602	111.9499	7.89%	125.66%	-8.78%	-6.71%	2.391	-52.19%	

Table B.59. ¹²³I FP-CIT, 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC-Total			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024394	0.015857	0.9992	39.5765	39.2449	77.5711	80.2821	-1.05%	-1.89%	-3.04%	0.35%	1.009	0.85%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024271	0.015857	0.99414	41.4114	42.0309	81.2897	84.1306	3.53%	5.08%	1.61%	5.16%	0.985	-1.47%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024129	0.015857	0.9883	42.5556	43.1279	82.5988	85.4855	6.39%	7.82%	3.25%	6.86%	0.987	-1.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023967	0.015857	0.98167	43.1291	43.8603	82.9773	85.8773	7.82%	9.65%	3.72%	7.35%	0.983	-1.67%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023128	0.015857	0.9473	45.4047	45.7254	83.5737	86.4946	13.51%	14.31%	4.47%	8.12%	0.993	-0.70%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022259	0.015857	0.91175	46.3993	49.5162	84.9959	87.9661	16.00%	23.79%	6.24%	9.96%	0.937	-6.29%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024394	0.015857	0.9992	34.9758	34.61	68.479	70.0472	-12.56%	-13.48%	-14.40%	-12.44%	1.011	1.06%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024271	0.015857	0.99414	36.5005	37.1725	71.7689	73.4124	-8.75%	-7.07%	-10.29%	-8.23%	0.982	-1.81%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024129	0.015857	0.9883	37.5211	38.1385	72.9289	74.5989	-6.20%	-4.65%	-8.84%	-6.75%	0.984	-1.62%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023967	0.015857	0.98167	38.0115	38.8104	73.2672	74.945	-4.97%	-2.97%	-8.42%	-6.32%	0.979	-2.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023128	0.015857	0.9473	40.0832	40.4474	73.8129	75.5032	0.21%	1.12%	-7.73%	-5.62%	0.991	-0.90%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022259	0.015857	0.91175	40.6993	44.1367	75.0901	76.8094	1.75%	10.34%	-6.14%	-3.99%	0.922	-7.79%	

Table B.60. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	20	60					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02438	0.016012	0.9986	39.1585	19.7756	57.4525	59.4462	-2.10%	-1.12%	-4.25%	-0.92%	1.980	-0.99%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024256	0.016012	0.99354	40.4648	22.1483	60.2574	62.3516	1.16%	10.74%	0.43%	3.92%	1.827	-8.65%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024114	0.016012	0.98771	41.3676	23.1675	61.2671	63.3972	3.42%	15.84%	2.11%	5.66%	1.786	-10.72%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023952	0.016012	0.98109	41.9007	23.7636	61.5584	63.6992	4.75%	18.82%	2.60%	6.17%	1.763	-11.84%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023114	0.016012	0.94674	43.9835	25.6965	62.2358	64.4013	9.96%	28.48%	3.73%	7.34%	1.712	-14.42%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022246	0.016012	0.9112	44.781	29.6472	63.7324	65.9542	11.95%	48.24%	6.22%	9.92%	1.510	-24.48%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02438	0.016012	0.9986	36.4928	15.5341	50.7135	51.8665	-8.77%	-22.33%	-15.48%	-13.56%	2.349	17.46%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024256	0.016012	0.99354	37.5395	17.7431	53.1999	54.4073	-6.15%	-11.28%	-11.34%	-9.32%	2.116	5.79%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024114	0.016012	0.98771	38.3294	18.657	54.091	55.3232	-4.18%	-6.72%	-9.85%	-7.79%	2.054	2.72%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023952	0.016012	0.98109	38.7929	19.1994	54.3526	55.591	-3.02%	-4.00%	-9.41%	-7.35%	2.021	1.03%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023114	0.016012	0.94674	40.667	20.9216	54.9703	56.2235	1.67%	4.61%	-8.38%	-6.29%	1.944	-2.81%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022246	0.016012	0.9112	41.0699	24.7897	56.3144	57.6008	2.67%	23.95%	-6.14%	-4.00%	1.657	-17.16%	

Table B.61. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	8	48					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02439	0.01609	0.99901	38.9039	8.0897	44.9478	46.4745	-2.74%	1.12%	-6.36%	-3.18%	4.8091	-3.82%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024266	0.01609	0.99395	39.8922	10.2141	47.1661	48.7769	-0.27%	27.68%	-1.74%	1.62%	3.9056	-21.89%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024124	0.01609	0.98812	40.6463	11.1867	48.0251	49.6683	1.62%	39.83%	0.05%	3.48%	3.6334	-27.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023962	0.01609	0.98149	41.1486	11.7006	48.1939	49.8445	2.87%	46.26%	0.40%	3.84%	3.5168	-29.66%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023123	0.01609	0.94713	43.0718	13.715	48.8085	50.4853	7.68%	71.44%	1.68%	5.18%	3.1405	-37.19%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022255	0.01609	0.91158	43.7293	17.7558	50.3075	52.0465	9.32%	121.95%	4.81%	8.43%	2.4628	-50.74%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02439	0.01609	0.99901	37.4009	4.0848	39.6658	40.5484	-6.50%	-48.94%	-17.36%	-15.52%	9.1561	83.12%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024266	0.01609	0.99395	38.1588	6.0817	41.6303	42.5618	-4.60%	-23.98%	-13.27%	-11.33%	6.2744	127.44%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024124	0.01609	0.98812	38.8067	6.9644	42.3925	43.3428	-2.98%	-12.95%	-11.68%	-9.70%	5.5722	57.22%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023962	0.01609	0.98149	39.2477	7.4293	42.5462	43.5008	-1.88%	-7.13%	-11.36%	-9.37%	5.2628	28.28%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023123	0.01609	0.94713	40.9559	9.2473	43.1094	44.0797	2.39%	15.59%	-10.19%	-8.17%	4.429	-11.42%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022255	0.01609	0.91158	41.2101	13.2188	44.4571	45.464	3.03%	65.24%	-7.38%	-5.28%	3.1175	-37.65%	

Table B.62. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			84	84	168					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022116	0.015003	0.90587	84.2143	82.9543	165.5116	171.296	0.26%	-1.24%	-1.48%	1.96%	1.0152	1.52%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022004	0.015003	0.90128	88.7436	87.7118	173.7215	179.7928	5.65%	4.42%	3.41%	7.02%	1.0118	1.18%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021875	0.015003	0.896	90.6944	89.3756	176.2158	182.3744	7.97%	6.40%	4.89%	8.56%	1.0148	1.48%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021728	0.015003	0.88999	91.3186	90.8474	177.1491	183.3402	8.71%	8.15%	5.45%	9.13%	1.0052	0.52%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020967	0.015003	0.85883	93.8762	93.3467	178.3129	184.5448	11.76%	11.13%	6.14%	9.85%	1.0057	0.57%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02018	0.015003	0.82659	94.7736	97.4728	179.7271	186.0083	12.83%	16.04%	6.98%	10.72%	0.9723	-2.77%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022116	0.015003	0.90587	74.4734	73.1047	146.1143	149.4602	-11.34%	-12.97%	-13.03%	-11.04%	1.0187	1.87%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022004	0.015003	0.90128	78.4563	77.3381	153.3761	156.8883	-6.60%	-7.93%	-8.70%	-6.61%	1.0145	1.45%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021875	0.015003	0.896	80.2186	78.7805	155.5867	159.1495	-4.50%	-6.21%	-7.39%	-5.27%	1.0183	1.83%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021728	0.015003	0.88999	80.6858	80.1752	156.4159	159.9977	-3.95%	-4.55%	-6.90%	-4.76%	1.0064	0.64%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020967	0.015003	0.85883	82.9716	82.4169	157.4683	161.0742	-1.22%	-1.88%	-6.27%	-4.12%	1.0067	0.67%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02018	0.015003	0.82659	83.459	86.4498	158.739	162.374	-0.64%	2.92%	-5.51%	-3.35%	0.9654	-3.46%	

Table B.63. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC-Total			84	42	126				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								175-180			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022306	0.014946	0.91366	83.2942	41.9334	123.3403	127.6202	-0.84%	-0.16%	-2.11%	1.29%	1.986	-0.68%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022193	0.014946	0.90903	86.7051	45.8079	129.7066	134.2116	3.22%	9.07%	2.94%	6.52%	1.893	-5.36%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022063	0.014946	0.9037	88.1274	47.2993	131.6346	136.2079	4.91%	12.62%	4.47%	8.10%	1.863	-6.84%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021915	0.014946	0.89764	88.6372	48.4749	132.2923	136.8901	5.52%	15.42%	4.99%	8.64%	1.829	-8.57%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021148	0.014946	0.86621	90.6971	51.0396	133.3536	137.9908	7.97%	21.52%	5.84%	9.52%	1.777	-11.15%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020354	0.014946	0.83369	91.0952	55.4292	134.7171	139.4079	8.45%	31.97%	6.92%	10.64%	1.643	-17.83%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022306	0.014946	0.91366	77.635	32.9166	108.8762	111.3515	-7.58%	-21.63%	-13.59%	-11.63%	2.359	17.93%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022193	0.014946	0.90903	80.6028	36.3944	114.5082	117.114	-4.04%	-13.35%	-9.12%	-7.05%	2.215	10.74%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022063	0.014946	0.9037	81.8605	37.7199	116.2174	118.8628	-2.55%	-10.19%	-7.76%	-5.66%	2.170	8.51%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021915	0.014946	0.89764	82.2429	38.8357	116.8029	119.4626	-2.09%	-7.53%	-7.30%	-5.19%	2.118	5.89%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021148	0.014946	0.86621	84.0282	41.1907	117.7641	120.4472	0.03%	-1.93%	-6.54%	-4.41%	2.040	2.00%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.020354	0.014946	0.83369	83.9891	45.5372	118.9915	121.7062	-0.01%	8.42%	-5.56%	-3.41%	1.844	-7.78%

Table B.64. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			84	16.8	100.8					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022345	0.015075	0.91525	82.7534	17.3346	97.1672	100.4684	-1.48%	3.18%	-3.60%	-0.33%	4.7739	-4.52%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022232	0.015075	0.91061	85.4964	20.6839	102.3807	105.872	1.78%	23.12%	1.57%	5.03%	4.1335	-17.33%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022101	0.015075	0.90527	86.607	22.0758	103.9175	107.4659	3.10%	31.40%	3.09%	6.61%	3.9232	-21.54%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021953	0.015075	0.8992	87.0563	23.0777	104.4016	107.9705	3.64%	37.37%	3.57%	7.11%	3.7723	-24.55%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021184	0.015075	0.86771	88.8418	25.7303	105.3742	108.9851	5.76%	53.16%	4.54%	8.12%	3.4528	-30.94%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020389	0.015075	0.83514	88.9924	30.275	106.813	110.4892	5.94%	80.21%	5.97%	9.61%	2.9395	-41.21%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022345	0.015075	0.91525	79.542	8.8157	85.7519	87.6604	-5.31%	-47.53%	-14.93%	-13.04%	9.0228	80.46%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022232	0.015075	0.91061	81.9034	11.8444	90.3656	92.3645	-2.50%	-29.50%	-10.35%	-8.35%	6.9149	38.30%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022101	0.015075	0.90527	82.8635	13.1031	91.7293	93.7813	-1.35%	-22.01%	-9.00%	-6.96%	6.324	26.48%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021953	0.015075	0.8992	83.2029	14.0545	92.162	94.226	-0.95%	-16.34%	-8.57%	-6.52%	5.92	18.40%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021184	0.015075	0.86771	84.7078	16.5225	93.0457	95.1346	0.84%	-1.65%	-7.69%	-5.62%	5.1268	2.54%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020389	0.015075	0.83514	84.4055	21.0487	94.3426	96.47	0.48%	25.29%	-6.41%	-4.30%	4.01	-19.80%	

Table B.65. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			100	100	200					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11492	0.068326	4.707	96.2485	96.1746	187.8801	194.4463	-3.75%	-3.83%	-6.06%	-2.78%	1.001	0.08%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11433	0.068326	4.6831	99.9567	104.1517	196.5176	203.3855	-0.04%	4.15%	-1.74%	1.69%	0.960	-4.03%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11366	0.068326	4.6556	103.5459	107.0138	199.435	206.4049	3.55%	7.01%	-0.28%	3.20%	0.968	-3.24%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1129	0.068326	4.6244	104.8952	110.537	201.2594	208.2928	4.90%	10.54%	0.63%	4.15%	0.949	-5.10%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10895	0.068326	4.4625	113.207	116.7252	201.8372	208.891	13.21%	16.73%	0.92%	4.45%	0.970	-3.01%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10486	0.068326	4.295	119.6194	129.5257	207.8714	215.1353	19.62%	29.53%	3.94%	7.57%	0.924	-7.65%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11492	0.068326	4.707	84.9793	84.8857	165.8479	169.6457	-15.02%	-15.11%	-17.08%	-15.18%	1.001	0.11%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11433	0.068326	4.6831	87.8288	92.3787	173.49	177.4627	-12.17%	-7.62%	-13.26%	-11.27%	0.951	-4.93%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11366	0.068326	4.6556	91.0877	94.8383	176.0751	180.107	-8.91%	-5.16%	-11.96%	-9.95%	0.960	-3.95%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1129	0.068326	4.6244	92.0635	98.2002	177.7039	181.773	-7.94%	-1.80%	-11.15%	-9.11%	0.938	-6.25%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10895	0.068326	4.4625	99.6844	103.5808	178.2897	182.3723	-0.32%	3.58%	-10.86%	-8.81%	0.962	-3.76%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10486	0.068326	4.295	104.8008	115.7322	183.6958	187.9018	4.80%	15.73%	-8.15%	-6.05%	0.906	-9.45%	

Table B.66. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13 hi6			
SC-Total			100	50	150					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11478	0.068684	4.7016	95.7447	47.2043	137.271	142.0323	-4.26%	-5.59%	-8.49%	-5.31%	2.0283	1.42%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1142	0.068684	4.6778	98.2271	54.0987	143.6953	148.69	-1.77%	6.20%	-4.20%	-0.67%	1.8157	-9.21%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11353	0.068684	4.6503	101.1937	56.7855	146.127	151.2078	1.19%	13.57%	-2.58%	0.81%	1.782	-10.90%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11277	0.068684	4.6191	102.435	59.9929	147.6314	152.7684	2.44%	19.99%	-1.58%	1.85%	1.7075	-14.63%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10882	0.068684	4.4574	110.34	66.3016	148.3222	153.4858	10.34%	32.60%	-1.12%	2.32%	1.6642	-16.79%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10474	0.068684	4.2901	116.2677	79.6144	154.639	160.0324	16.29%	59.23%	3.09%	6.69%	1.4606	-26.97%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11478	0.068684	4.7016	89.337	36.8525	121.1597	123.9131	-10.66%	-26.30%	-19.23%	-17.39%	2.4242	21.21%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1142	0.068684	4.6778	91.0863	43.4041	126.8481	129.7369	-8.91%	-13.19%	-15.43%	-13.51%	2.0986	4.93%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11353	0.068684	4.6503	93.7481	45.7538	129.004	131.9429	-6.25%	-8.49%	-14.00%	-12.04%	2.049	2.45%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11277	0.068684	4.6191	94.6466	48.8177	130.3509	133.3227	-5.35%	-2.36%	-13.10%	-11.12%	1.9388	-3.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10882	0.068684	4.4574	101.8516	54.3559	131.0358	134.0247	1.85%	8.71%	-12.64%	-10.65%	1.8738	-6.31%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10474	0.068684	4.2901	106.4529	67.0608	136.6941	139.8179	6.45%	34.12%	-8.87%	-6.79%	1.5874	-20.63%	

Table B.67. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			100	20	120					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1145	0.06913	4.6897	95.4778	17.8536	105.7442	109.3265	-4.52%	-10.73%	-11.88%	-8.89%	5.35	6.96%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11392	0.06913	4.666	97.2293	24.1133	110.8892	114.6733	-2.77%	20.57%	-7.59%	-4.44%	4.03	-19.36%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11325	0.06913	4.6386	99.8254	26.72	112.8537	116.7125	-0.17%	33.60%	-5.96%	-2.74%	3.74	-25.28%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11249	0.06913	4.6075	100.9808	29.7867	114.0965	118.0082	0.98%	48.93%	-4.92%	-1.66%	3.39	-32.20%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10855	0.06913	4.4462	108.6208	36.333	114.7329	118.6802	8.62%	81.67%	-4.39%	-1.10%	2.99	-40.21%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10447	0.06913	4.2793	114.5268	49.7428	121.2261	125.4234	14.53%	148.71%	1.02%	4.52%	2.30	-53.95%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1145	0.06913	4.6897	91.9831	8.0606	93.3067	95.3775	-8.02%	-59.70%	-22.24%	-20.52%	11.41	128.23%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11392	0.06913	4.666	93.0754	14.0618	97.8683	100.0563	-6.92%	-29.69%	-18.44%	-16.62%	6.62	32.38%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11325	0.06913	4.6386	95.3796	16.3699	99.6129	101.8443	-4.62%	-18.15%	-16.99%	-15.13%	5.83	16.53%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11249	0.06913	4.6075	96.2061	19.305	100.7302	102.9927	-3.79%	-3.48%	-16.06%	-14.17%	4.98	-0.33%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10855	0.06913	4.4462	103.1252	25.1019	101.3693	103.6542	3.13%	25.51%	-15.53%	-13.62%	4.11	-17.83%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10447	0.06913	4.2793	107.6744	37.9057	107.1869	109.6185	7.67%	89.53%	-10.68%	-8.65%	2.84	-43.19%	

Table B.68. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025347	0.014987	1.0382	39.0912	39.0173	77.1663	79.8632	-2.27%	-2.46%	-3.54%	-0.17%	1.002	0.19%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025218	0.014987	1.0329	40.9168	41.6206	80.8634	83.6895	2.29%	4.05%	1.08%	4.61%	0.983	-1.69%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02507	0.014987	1.0269	42.0712	42.467	81.9826	84.8478	5.18%	6.17%	2.48%	6.06%	0.991	-0.93%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024902	0.014987	1.02	42.403	43.5204	82.5993	85.4861	6.01%	8.80%	3.25%	6.86%	0.974	-2.57%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02403	0.014987	0.98428	44.3692	45.1411	82.9094	85.8069	10.92%	12.65%	3.64%	7.26%	0.983	-1.71%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023128	0.014987	0.94733	45.6587	48.2115	84.2336	87.1773	14.15%	20.53%	5.29%	8.97%	0.947	-5.30%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025347	0.014987	1.0382	34.5183	34.4347	68.1196	69.6795	-13.70%	-13.91%	-14.85%	-12.90%	1.002	0.24%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025218	0.014987	1.0329	36.0542	36.8176	71.3901	73.0249	-9.86%	-7.96%	-10.76%	-8.72%	0.979	-2.07%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02507	0.014987	1.0269	37.1098	37.5354	72.3816	74.0391	-7.23%	-6.16%	-9.52%	-7.45%	0.989	-1.13%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024902	0.014987	1.02	37.3317	38.5455	72.9303	74.6004	-6.67%	-3.64%	-8.84%	-6.75%	0.969	-3.15%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02403	0.014987	0.98428	39.119	39.9737	73.2225	74.8993	-2.20%	-0.07%	-6.47%	-6.38%	0.979	-2.14%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023128	0.014987	0.94733	40.1004	42.9124	74.4111	76.115	0.25%	7.28%	-6.99%	-4.86%	0.934	-6.55%	

Table B.69. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	20	60					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02532	0.015095	1.0371	38.8847	19.4262	57.1996	59.1841	-2.79%	-2.87%	-4.67%	-1.36%	2.002	0.08%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025192	0.015095	1.0319	40.219	21.5953	59.9794	62.0633	0.55%	7.98%	-0.03%	3.44%	1.862	-6.88%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025044	0.015095	1.0258	41.1218	22.3715	60.8452	62.9596	2.80%	11.86%	1.41%	4.93%	1.838	-8.09%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024876	0.015095	1.0189	41.4004	23.3032	61.3242	63.4566	3.50%	16.52%	2.21%	5.76%	1.777	-11.17%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024005	0.015095	0.98325	43.1752	24.9665	61.6477	63.7925	7.94%	24.83%	2.75%	6.32%	1.729	-13.53%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023104	0.015095	0.94635	44.3461	28.107	62.9816	65.176	10.87%	40.54%	4.97%	8.63%	1.578	-21.11%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02532	0.015095	1.0371	36.2568	15.2186	50.4891	51.6368	-9.36%	-23.91%	-15.85%	-13.94%	2.382	19.12%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025192	0.015095	1.0319	37.3517	17.224	52.9488	54.1541	-6.62%	-13.88%	-11.75%	-9.74%	2.169	8.43%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025044	0.015095	1.0258	38.1659	17.898	53.716	54.939	-4.59%	-10.51%	-10.47%	-8.44%	2.132	6.62%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024876	0.015095	1.0189	38.3468	18.7942	54.1433	55.3769	-4.13%	-6.03%	-9.76%	-7.71%	2.040	2.02%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024005	0.015095	0.98325	39.9402	20.2828	54.4476	55.6888	-0.15%	1.41%	-9.25%	-7.19%	1.969	-1.54%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023104	0.015095	0.94635	40.7951	23.3034	55.6445	56.9148	1.99%	16.52%	-7.26%	-5.14%	1.751	-12.47%	

Table B.70. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	8	48					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025292	0.015179	1.0359	38.7629	7.6736	44.7321	46.2496	-3.09%	-4.08%	-6.81%	-3.65%	5.051	1.03%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025163	0.015179	1.0307	39.8024	9.5834	46.9448	48.5456	-0.49%	19.79%	-2.20%	1.14%	4.153	-16.93%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025016	0.015179	1.0246	40.5562	10.3178	47.6586	49.2859	1.39%	28.97%	-0.71%	2.68%	3.931	-21.39%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024848	0.015179	1.0178	40.8062	11.1746	48.053	49.6971	2.02%	39.68%	0.11%	3.54%	3.652	-26.97%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023978	0.015179	0.96214	42.4831	12.8552	48.3691	50.0267	6.21%	60.69%	0.77%	4.23%	3.305	-33.91%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023078	0.015179	0.94528	43.5793	16.0575	49.7398	51.4553	8.95%	100.72%	3.62%	7.20%	2.714	-45.72%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025292	0.015179	1.0359	37.3019	3.6908	39.4738	40.3511	-6.75%	-53.87%	-17.76%	-15.94%	10.107	102.13%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025163	0.015179	1.0307	38.132	5.4708	41.433	42.3586	-4.67%	-31.62%	-13.68%	-11.75%	6.970	197.01%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025016	0.015179	1.0246	38.8034	6.1185	42.066	43.0071	-2.99%	-23.52%	-12.36%	-10.40%	6.342	134.20%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024848	0.015179	1.0178	38.9631	6.9444	42.4192	43.3701	-2.59%	-13.20%	-11.63%	-9.65%	5.611	61.07%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023978	0.015179	0.96214	40.4577	8.4593	42.7179	43.6781	1.14%	5.74%	-11.00%	-9.00%	4.783	-4.35%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023078	0.015179	0.94528	41.2314	11.5488	43.9486	44.9417	3.08%	44.36%	-8.44%	-6.37%	3.570	-28.60%	

Table B.71. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC-Total			84	84	168				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								185-190			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023085	0.014324	0.94557	83.1807	82.969	165.1223	170.8931	-0.98%	-1.23%	-1.71%	1.72%	1.003	0.26%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022968	0.014324	0.94079	87.5497	87.6	173.2278	179.282	4.23%	4.29%	3.11%	6.72%	0.999	-0.06%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022834	0.014324	0.93527	89.5886	88.9677	175.5181	181.6523	6.65%	5.91%	4.48%	8.13%	1.007	0.70%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02268	0.014324	0.92899	90.0195	90.6807	176.6006	182.7726	7.17%	7.95%	5.12%	8.79%	0.993	-0.73%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021886	0.014324	0.89647	92.3016	93.007	177.5834	183.7897	9.88%	10.72%	5.70%	9.40%	0.992	-0.76%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021065	0.014324	0.86282	93.4606	96.4768	178.905	185.1574	11.26%	14.85%	6.49%	10.21%	0.969	-3.13%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023085	0.014324	0.94557	73.4543	73.2222	145.7687	149.1067	-12.55%	-12.83%	-13.23%	-11.25%	1.003	0.32%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022968	0.014324	0.94079	77.2919	77.3473	152.939	156.4412	-7.99%	-7.92%	-8.96%	-6.88%	0.999	-0.07%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022834	0.014324	0.93527	79.1697	78.4894	154.9686	158.5173	-5.75%	-6.56%	-7.76%	-5.64%	1.009	0.87%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02268	0.014324	0.92899	79.4238	80.1397	155.9299	159.5006	-5.45%	-4.60%	-7.18%	-5.06%	0.991	-0.89%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021886	0.014324	0.89647	81.4525	82.2388	156.822	160.4131	-3.03%	-2.10%	-6.65%	-4.52%	0.990	-0.96%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021065	0.014324	0.86282	82.2661	85.5889	158.0104	161.6287	-2.06%	1.89%	-5.95%	-3.79%	0.961	-3.88%

Table B.72. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Jo13 hi6			
SC-Total			84	42	126					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									185-190			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023308	0.014146	0.9547	82.6982	41.6862	123.0399	127.3095	-1.55%	-0.75%	-2.35%	1.04%	1.984	-0.81%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02319	0.014146	0.94986	86.0204	45.3872	129.2637	133.7532	2.41%	8.06%	2.59%	6.15%	1.895	-5.24%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023054	0.014146	0.94429	87.5114	46.5921	130.9959	135.5463	4.18%	10.93%	3.97%	7.58%	1.878	-6.09%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022899	0.014146	0.93796	87.8075	48.0329	131.8122	136.3932	4.53%	14.36%	4.61%	8.25%	1.828	-8.60%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022096	0.014146	0.90512	89.5735	50.366	132.6748	137.2883	6.64%	19.92%	5.30%	8.96%	1.778	-11.08%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021268	0.014146	0.87114	90.286	54.0157	133.9078	138.5692	7.48%	28.61%	6.28%	9.98%	1.671	-16.43%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023308	0.014146	0.9547	77.0717	32.7339	108.6096	111.0789	-8.25%	-22.06%	-13.80%	-11.84%	2.354	17.72%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02319	0.014146	0.94986	79.9692	36.0501	114.115	116.7128	-4.80%	-14.17%	-9.43%	-7.37%	2.218	10.91%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023054	0.014146	0.94429	81.3231	37.0862	115.6513	118.2835	-3.19%	-11.70%	-8.21%	-6.12%	2.193	9.64%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022899	0.014146	0.93796	81.47	38.4828	116.3773	119.0273	-3.01%	-8.37%	-7.64%	-5.53%	2.117	5.85%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022096	0.014146	0.90512	82.9857	40.6392	117.1624	119.8318	-1.21%	-3.24%	-7.01%	-4.90%	2.042	2.10%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021268	0.014146	0.87114	83.334	44.2137	118.2719	120.9695	-0.79%	5.27%	-6.13%	-3.99%	1.885	-5.76%	

Table B.73. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 185-190 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW I ¹²³ I hi6			
SC-Total			84	16.8	100.8				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								185-190			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023291	0.01426	0.95401	82.4336	16.9417	96.9687	100.2822	-1.86%	0.84%	-3.78%	-0.51%	4.866	-2.69%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023173	0.01426	0.94918	85.1348	20.0933	102.0485	105.5263	1.35%	19.60%	1.24%	4.69%	4.237	-15.26%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023037	0.01426	0.94361	86.3095	21.2072	103.3882	106.9153	2.75%	26.23%	2.57%	6.07%	4.070	-18.60%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022883	0.01426	0.93728	86.5368	22.4933	104.0379	107.5925	3.02%	33.89%	3.21%	6.74%	3.847	-23.06%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022082	0.01426	0.90447	88.0466	24.881	104.8159	108.4055	4.82%	48.10%	3.96%	7.55%	3.539	-29.23%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021253	0.01426	0.87052	88.5839	28.6609	106.1025	109.7496	5.46%	70.60%	5.26%	8.88%	3.091	-38.18%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023291	0.01426	0.95401	79.2641	8.4632	85.5927	87.4968	-5.64%	-49.62%	-15.09%	-13.20%	9.366	87.31%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023173	0.01426	0.94918	81.604	11.3016	90.0709	92.0819	-2.85%	-32.73%	-10.64%	-8.65%	7.221	44.41%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023037	0.01426	0.94361	82.6537	12.2801	91.2595	93.2991	-1.60%	-26.90%	-9.46%	-7.44%	6.731	34.61%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022883	0.01426	0.93728	82.7484	13.5317	91.8368	93.8945	-1.49%	-19.45%	-8.89%	-6.65%	6.115	22.30%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022082	0.01426	0.90447	84.0048	15.7668	92.5501	94.6266	0.01%	-6.15%	-8.18%	-6.12%	5.328	6.56%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021253	0.01426	0.87052	84.1639	19.489	93.7085	95.8188	0.20%	16.01%	-7.04%	-4.94%	4.319	-13.63%

Table B.74. ^{123}I β -CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole									TEW	lo13	hi6		
SC-Total										137-142				
Atten Coeff	0.1481									142-175				
Striata ROI Size	13 13									175-180				
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10734	0.075004	4.3966	98.4473	96.554	188.3414	194.9237	-1.55%	-3.45%	-5.83%	-2.54%	1.020	1.96%		
Atten Coeff	0.1481													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10679	0.075004	4.3743	101.8318	104.7264	195.2903	202.1154	1.83%	4.73%	-2.35%	1.06%	0.972	-2.76%		
Atten Coeff	0.1481													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10617	0.075004	4.3486	106.9168	108.6437	199.9397	206.9274	6.92%	8.64%	-0.03%	3.46%	0.984	-1.59%		
Atten Coeff	0.1481													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10546	0.075004	4.3195	110.9981	109.5016	201.3338	208.3702	11.00%	9.50%	0.67%	4.19%	1.014	1.37%		
Atten Coeff	0.1481													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10176	0.075004	4.1682	116.0893	126.0281	206.9176	214.1482	16.09%	26.03%	3.46%	7.07%	0.921	-7.89%		
Atten Coeff	0.1481													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.097943	0.075004	4.0118	122.5974	142.9441	213.986	221.4612	22.60%	42.94%	6.99%	10.73%	0.858	-14.23%		
Atten Coeff	0.12													
Striata ROI Size	13 13													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10734	0.075004	4.3966	87.1131	85.0378	166.2659	170.0732	-12.89%	-14.96%	-16.87%	-14.96%	1.024	2.44%		
Atten Coeff	0.12													
Striata ROI Size	15 15													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10679	0.075004	4.3743	89.6269	92.7564	172.4178	176.366	-10.37%	-7.24%	-13.79%	-11.82%	0.966	-3.37%		
Atten Coeff	0.12													
Striata ROI Size	17 17													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10617	0.075004	4.3486	94.2621	96.103	176.5389	180.5815	-5.74%	-3.90%	-11.73%	-9.71%	0.981	-1.92%		
Atten Coeff	0.12													
Striata ROI Size	19 19													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10546	0.075004	4.3195	98.1846	96.5843	177.793	181.8643	-1.82%	-3.42%	-11.10%	-9.07%	1.017	1.66%		
Atten Coeff	0.12													
Striata ROI Size	26 26													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.10176	0.075004	4.1682	101.6107	112.5	182.8168	187.0027	1.61%	12.50%	-8.59%	-6.50%	0.903	-9.68%		
Atten Coeff	0.12													
Striata ROI Size	32 32													
Bkgd ROI Size	13 13													
Units:	microCi													
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err		
0.097943	0.075004	4.0118	106.4041	128.7924	189.1664	193.4962	6.40%	28.79%	-5.42%	-3.25%	0.826	-17.38%		

Table B.75. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			100	50	150					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10713	0.075002	4.3682	97.1183	48.4724	137.6338	142.4088	-2.88%	-3.06%	-8.24%	-5.06%	2.004	0.18%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10659	0.075002	4.3659	98.9925	55.8989	142.3116	147.2603	-1.01%	11.80%	-5.13%	-1.83%	1.771	-11.45%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10596	0.075002	4.3403	103.128	59.9983	146.5078	151.6048	3.13%	20.00%	-2.33%	1.07%	1.719	-14.06%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10525	0.075002	4.3112	106.9689	60.6377	147.5511	152.6826	6.97%	21.28%	-1.63%	1.79%	1.764	-11.80%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10157	0.075002	4.1603	112.0103	76.785	153.2155	158.5589	12.01%	53.57%	2.14%	5.71%	1.459	-27.06%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.097756	0.075002	4.0041	118.3979	93.7876	160.4699	166.0743	18.40%	87.58%	6.98%	10.72%	1.262	-36.88%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10713	0.075002	4.3682	90.5695	37.9614	121.4904	124.252	-9.43%	-24.08%	-19.01%	-17.17%	2.366	19.29%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10659	0.075002	4.3659	91.6693	45.0982	125.6368	128.4993	-8.33%	-9.80%	-16.24%	-14.33%	2.033	1.63%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10596	0.075002	4.3403	95.3484	48.721	129.3581	132.3068	-4.65%	-2.56%	-13.76%	-11.80%	1.957	-2.15%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10525	0.075002	4.3112	99.0539	49.012	130.3002	133.2693	-0.95%	-1.98%	-13.13%	-11.15%	2.021	1.05%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10157	0.075002	4.1603	102.4655	64.5561	135.3984	138.4923	2.47%	29.11%	-9.73%	-7.67%	1.587	-20.64%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.097756	0.075002	4.0041	107.1271	80.9449	141.9142	145.1617	7.13%	61.89%	-5.39%	-3.23%	1.323	-33.83%	

Table B.76. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			100	20	120					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10694	0.075079	4.3802	96.3345	19.6312	105.3525	108.9291	-3.67%	-1.84%	-12.21%	-9.23%	4.907	-1.86%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1064	0.075079	4.358	97.3091	26.609	108.8499	112.5706	-2.69%	33.05%	-9.29%	-6.19%	3.657	-26.86%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10577	0.075079	4.3324	100.8769	30.825	112.683	116.549	0.88%	54.13%	-6.10%	-2.88%	3.273	-34.55%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10506	0.075079	4.3034	104.5757	31.3376	113.3587	117.2467	4.58%	56.69%	-5.53%	-2.29%	3.337	-33.26%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10138	0.075079	4.1527	109.7113	47.1835	118.9339	123.0498	9.71%	135.92%	-0.89%	2.54%	2.325	-53.50%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.097579	0.075079	3.9968	115.9983	64.3948	126.3968	130.7938	16.00%	221.97%	5.33%	8.99%	1.801	-63.97%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10694	0.075079	4.3802	92.6555	9.7226	92.9723	95.0401	-7.34%	-51.39%	-22.52%	-20.80%	9.530	90.60%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1064	0.075079	4.358	92.9142	16.5086	96.0777	98.2311	-7.09%	-17.46%	-19.94%	-18.14%	5.628	12.56%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10577	0.075079	4.3324	96.0206	20.3042	99.4827	101.7188	-3.98%	1.52%	-17.10%	-15.23%	4.729	-5.42%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10506	0.075079	4.3034	99.5981	20.4841	100.0991	102.3484	-0.40%	2.42%	-16.58%	-14.71%	4.862	-2.76%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10138	0.075079	4.1527	103.132	35.7181	105.1237	107.5072	3.13%	78.59%	-12.40%	-10.41%	2.887	-42.25%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.097579	0.075079	3.9968	107.6727	52.3202	111.8271	114.376	7.67%	161.60%	-6.81%	-4.69%	2.058	-58.84%	

Table B.77. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	40	80					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023902	0.016585	0.97902	39.7736	38.7354	77.1484	79.8446	-0.57%	-3.16%	-3.56%	-0.19%	1.027	2.68%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023781	0.016585	0.97406	41.5715	41.3101	80.5116	83.3254	3.93%	3.28%	0.64%	4.16%	1.006	0.63%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023641	0.016585	0.96834	43.1168	42.3451	81.9833	84.8485	7.79%	5.86%	2.48%	6.06%	1.018	1.82%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023483	0.016585	0.96185	44.0408	42.7668	82.4003	85.28	10.10%	6.92%	3.00%	6.60%	1.030	2.98%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02266	0.016585	0.92817	45.0874	46.8511	83.7387	86.6651	12.72%	17.13%	4.67%	8.33%	0.962	-3.76%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02181	0.016585	0.89333	46.3656	50.8347	85.1955	88.1725	15.91%	27.09%	6.49%	10.22%	0.912	-8.79%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023902	0.016585	0.97902	35.2192	34.0897	68.107	69.6665	-11.95%	-14.78%	-14.87%	-12.92%	1.033	3.31%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023781	0.016585	0.97406	36.7313	36.4473	71.082	72.7098	-8.17%	-8.88%	-11.15%	-9.11%	1.008	0.78%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023641	0.016585	0.96834	38.1548	37.3114	72.3657	74.0433	-4.61%	-6.72%	-9.52%	-7.45%	1.023	2.26%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023483	0.016585	0.96185	39.0209	37.6438	72.7591	74.4252	-2.45%	-5.89%	-9.05%	-6.97%	1.037	3.66%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02266	0.016585	0.92817	39.66	41.596	73.9637	75.6574	-0.85%	3.99%	-7.55%	-5.43%	0.953	-4.65%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02181	0.016585	0.89333	40.5381	45.4564	75.2754	76.9988	1.35%	13.64%	-5.91%	-3.75%	0.892	-10.82%	

Table B.78. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data, Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	20	60					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023809	0.016537	0.97522	39.2431	19.5058	57.1495	59.1321	-1.89%	-2.47%	-4.75%	-1.45%	2.012	0.59%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023688	0.016537	0.97028	40.4378	21.7816	59.6113	61.6825	1.09%	8.91%	-0.65%	2.80%	1.857	-7.17%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023549	0.016537	0.96459	41.5996	22.8942	60.8642	62.9797	4.00%	14.47%	1.44%	4.97%	1.817	-9.15%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023392	0.016537	0.95812	42.4246	23.2303	61.1089	63.2328	6.06%	16.15%	1.85%	5.39%	1.826	-8.69%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022572	0.016537	0.92457	43.5041	27.1105	62.4235	64.5978	8.76%	35.55%	4.04%	7.66%	1.605	-19.77%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021725	0.016537	0.88987	44.7131	31.1864	63.9847	66.2168	11.78%	55.93%	6.64%	10.36%	1.434	-28.31%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023809	0.016537	0.97522	36.6026	15.2618	50.4476	51.5942	-8.49%	-23.69%	-15.92%	-14.01%	2.398	19.92%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023688	0.016537	0.97028	37.5498	17.3863	52.6263	53.8243	-6.13%	-13.07%	-12.29%	-10.29%	2.160	7.99%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023549	0.016537	0.96459	38.5872	18.3654	53.7369	54.9606	-3.53%	-8.17%	-10.44%	-8.40%	2.101	5.05%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023392	0.016537	0.95812	39.3636	18.6238	53.9577	55.1865	-1.59%	-6.88%	-10.07%	-8.02%	2.114	5.68%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022572	0.016537	0.92457	40.0533	22.3716	55.1415	56.3999	0.13%	11.86%	-8.10%	-6.00%	1.790	-10.48%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021725	0.016537	0.88987	40.8512	26.3315	56.5462	57.8387	2.13%	31.66%	-5.76%	-3.60%	1.551	-22.43%	

Table B.79. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW, lo13, hi6			
SC-Total			40	8	48					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023685	0.016599	0.97015	38.9363	7.9791	44.6014	46.1159	-2.66%	-0.26%	-7.06%	-3.93%	4.880	-2.40%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023565	0.016599	0.96523	39.7724	10.0795	46.4782	48.0651	-0.57%	25.99%	-3.17%	0.14%	3.946	-21.08%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023427	0.016599	0.95957	40.7088	11.2419	47.6008	49.2297	1.77%	40.52%	-0.83%	2.56%	3.621	-27.58%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02327	0.016599	0.95313	41.4753	11.5345	47.7846	49.4203	3.69%	44.18%	-0.45%	2.96%	3.596	-28.08%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022455	0.016599	0.91976	42.561	15.3464	49.1272	50.8202	6.40%	91.83%	2.35%	5.88%	2.773	-44.53%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021612	0.016599	0.88524	43.7977	19.4269	50.7824	52.5414	9.49%	142.84%	5.80%	9.46%	2.254	-54.91%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023685	0.016599	0.97015	37.4428	3.9749	39.3615	40.2371	-6.39%	-50.31%	-18.00%	-16.17%	9.420	88.40%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023565	0.016599	0.96523	38.0541	5.963	41.0243	41.942	-4.86%	-25.46%	-14.53%	-12.62%	6.382	27.63%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023427	0.016599	0.95957	38.864	7.0139	42.0203	42.9624	-2.84%	-12.33%	-12.46%	-10.50%	5.541	10.82%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02327	0.016599	0.95313	39.5867	7.2366	42.1874	43.1334	-1.03%	-9.56%	-12.11%	-10.14%	5.471	9.42%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022455	0.016599	0.91976	40.2887	10.9151	43.3976	44.3775	0.72%	36.44%	-9.59%	-7.55%	3.691	-26.18%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021612	0.016599	0.88524	41.1125	14.877	44.8869	45.9055	2.78%	85.96%	-6.49%	-4.36%	2.763	-44.73%	

Table B.80. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13_hi6			
SC-Total			84	84	168					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022638	0.015763	0.92726	84.2507	81.7339	164.2374	169.9773	0.30%	-2.70%	-2.24%	1.18%	1.031	3.08%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022524	0.015763	0.92256	88.6898	86.1439	171.7289	177.7306	5.58%	2.55%	2.22%	5.79%	1.030	2.96%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022391	0.015763	0.91715	91.2325	87.5716	174.269	180.3594	8.61%	4.25%	3.73%	7.36%	1.042	4.18%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022241	0.015763	0.911	92.2083	88.5	174.9872	181.1026	9.77%	5.36%	4.16%	7.80%	1.042	4.19%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021462	0.015763	0.8791	92.9928	93.4848	176.6444	182.8179	10.71%	11.29%	5.15%	8.82%	0.995	-0.53%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020657	0.015763	0.8461	94.0784	97.7689	178.0703	184.2935	12.00%	16.39%	5.99%	9.70%	0.962	-3.77%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022638	0.015763	0.92726	74.6316	71.9014	144.9901	148.3102	-11.15%	-14.40%	-13.70%	-11.72%	1.038	3.80%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022524	0.015763	0.92256	78.5603	75.8023	151.6166	155.0885	-6.48%	-9.76%	-9.75%	-7.69%	1.036	3.64%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022391	0.015763	0.91715	80.9297	76.9536	153.8666	157.39	-3.66%	-8.39%	-8.41%	-6.32%	1.052	5.17%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022241	0.015763	0.911	81.7981	77.7786	154.5064	158.0444	-2.62%	-7.41%	-8.03%	-5.93%	1.052	5.17%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021462	0.015763	0.8791	82.0905	82.6499	155.9977	159.5689	-2.27%	-1.61%	-7.14%	-5.02%	0.993	-0.68%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020657	0.015763	0.8461	82.7481	86.8276	157.2824	160.8839	-1.49%	3.37%	-6.38%	-4.24%	0.953	-4.70%	

Table B.81. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data, Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			84	42	126					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022439	0.015469	0.91911	83.135	41.2562	122.3436	126.5878	-1.03%	-1.77%	-2.90%	0.47%	2.015	0.75%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022326	0.015469	0.91446	86.3062	45.0365	128.1669	132.6175	2.75%	7.23%	1.72%	5.25%	1.916	-4.18%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022195	0.015469	0.90909	88.0405	46.6309	130.3034	134.8293	4.81%	11.03%	3.42%	7.01%	1.888	-5.60%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022046	0.015469	0.90299	88.8027	47.3787	130.8399	135.3852	5.72%	12.81%	3.84%	7.45%	1.874	-6.28%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021274	0.015469	0.87138	89.6186	51.9341	132.4014	137.0078	6.69%	23.65%	5.08%	8.74%	1.726	-13.72%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020475	0.015469	0.83867	90.4979	56.3555	133.9113	138.5758	7.74%	34.18%	6.28%	9.98%	1.606	-19.71%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022439	0.015469	0.91911	77.5447	32.2691	107.9969	110.4516	-7.68%	-23.17%	-14.29%	-12.34%	2.403	20.15%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022326	0.015469	0.91446	80.2857	35.679	113.1489	115.7232	-4.42%	-15.05%	-10.20%	-8.16%	2.250	12.51%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022195	0.015469	0.90909	81.841	37.0752	115.0428	117.6609	-2.57%	-11.73%	-8.70%	-6.62%	2.207	10.37%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022046	0.015469	0.90299	82.5159	37.7441	115.522	118.1514	-1.77%	-10.13%	-8.32%	-6.23%	2.186	9.31%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021274	0.015469	0.87138	82.8777	42.1893	116.9293	119.5947	-1.34%	0.45%	-7.20%	-5.08%	1.964	-1.78%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020475	0.015469	0.83867	83.3115	46.5274	118.2892	120.9888	-0.82%	10.78%	-6.12%	-3.98%	1.791	-10.47%	

Table B.82. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data, Energy Windows Used: 137-142 keV, 142-175 keV, 175-180 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			84	16.8	100.8					137-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									175-180			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022329	0.015376	0.91461	82.4633	16.9684	96.2792	99.5488	-1.83%	1.00%	-4.48%	-1.24%	4.860	-2.80%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022216	0.015376	0.90998	84.8727	20.3706	101.0172	104.4614	1.04%	21.25%	0.22%	3.63%	4.166	-16.67%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022086	0.015376	0.90464	86.1226	22.0635	102.8464	106.3589	2.53%	31.33%	2.03%	5.51%	3.903	-21.93%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021938	0.015376	0.89657	86.755	22.7045	103.2197	106.7471	3.28%	35.15%	2.40%	5.90%	3.821	-23.58%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02117	0.015376	0.86711	87.55	27.0343	104.6274	108.2188	4.23%	60.92%	3.80%	7.36%	3.238	-35.23%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020375	0.015376	0.83456	88.3982	31.4879	106.2615	109.923	5.24%	87.43%	5.42%	9.05%	2.807	-43.85%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022329	0.015376	0.91461	79.2904	8.4888	84.9692	86.8596	-5.61%	-49.47%	-15.71%	-13.83%	9.341	86.81%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022216	0.015376	0.90998	81.318	11.6038	89.1631	91.1547	-3.19%	-30.93%	-11.54%	-9.57%	7.008	40.16%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022086	0.015376	0.90464	82.3858	13.1453	90.7864	92.8177	-1.92%	-21.75%	-9.93%	-7.92%	6.267	25.35%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021938	0.015376	0.89657	82.9424	13.7226	91.1214	93.1614	-1.26%	-18.32%	-9.60%	-7.58%	6.044	20.88%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02117	0.015376	0.86711	83.3044	17.9468	92.393	94.4708	-0.83%	6.83%	-8.34%	-6.28%	4.642	-7.17%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020375	0.015376	0.83456	83.7014	22.321	93.8649	95.984	-0.36%	32.86%	-6.88%	-4.78%	3.750	-25.00%	

Table B.83. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC-Total			100	100	200					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10836	0.080209	4.4394	99.7016	95.0476	188.8207	195.4196	-0.30%	-4.95%	-5.59%	-2.29%	1.049	4.90%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10783	0.080209	4.4169	102.1366	105.7046	197.4808	204.3624	2.14%	5.70%	-1.26%	2.19%	0.966	-3.36%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1072	0.080209	4.391	104.9355	111.3563	200.7568	207.7726	4.94%	11.36%	0.36%	3.89%	0.942	-5.77%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10648	0.080209	4.3615	109.0965	111.0185	200.5927	207.6031	9.10%	11.02%	0.30%	3.80%	0.983	-1.73%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10275	0.080209	4.2088	109.6653	129.2684	203.5974	210.7092	9.67%	29.27%	1.80%	5.35%	0.848	-15.16%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.098896	0.080209	4.0509	118.4949	141.0424	209.3646	216.698	18.49%	41.04%	4.69%	8.35%	0.840	-15.99%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10836	0.080209	4.4394	88.4866	83.444	166.6921	170.5091	-11.51%	-16.56%	-16.65%	-14.75%	1.060	6.04%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10783	0.080209	4.4169	89.8031	93.7175	174.3623	178.355	-10.20%	-6.28%	-12.82%	-10.82%	0.958	-4.18%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.1072	0.080209	4.391	92.0112	99.0065	177.2746	181.3338	-7.99%	-0.99%	-11.36%	-9.33%	0.929	-7.07%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10648	0.080209	4.3615	96.1153	98.3178	177.1456	181.2021	-3.88%	-1.68%	-11.43%	-9.40%	0.978	-2.24%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.10275	0.080209	4.2088	94.8741	116.4238	179.8769	183.9939	-5.13%	16.42%	-10.06%	-8.00%	0.815	-18.51%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.098896	0.080209	4.0509	102.4931	127.3691	185.0934	189.3294	2.49%	27.37%	-7.45%	-5.34%	0.805	-19.53%	

Table B.84. ¹²³I β-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC-Total			100	50	150				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10824	0.080231	4.4336	98.5975	47.7637	138.8228	143.6362	-1.40%	-4.47%	-7.45%	-4.24%	2.064	3.21%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10769	0.080231	4.4111	99.3647	57.6968	145.2142	150.2664	-0.64%	15.39%	-3.19%	0.18%	1.722	-13.89%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10706	0.080231	4.3853	101.2879	63.4878	147.8679	153.0184	1.29%	26.98%	-1.42%	2.01%	1.595	-20.23%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10634	0.080231	4.3558	105.0308	63.2286	147.4324	152.5651	5.03%	26.46%	-1.71%	1.71%	1.661	-16.94%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10262	0.080231	4.2033	104.9534	81.8084	150.8622	156.1304	4.95%	63.62%	0.57%	4.09%	1.283	-35.85%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.098769	0.080231	4.0456	113.494	93.7431	156.5955	162.066	13.49%	87.49%	4.40%	8.04%	1.211	-39.47%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10824	0.080231	4.4336	92.0883	37.1239	122.5413	125.3252	-7.91%	-25.75%	-18.31%	-16.45%	2.481	24.03%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10769	0.080231	4.4111	91.8351	46.8527	128.2089	131.1314	-8.16%	-6.29%	-14.53%	-12.58%	1.960	-2.00%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10706	0.080231	4.3853	93.165	52.3685	130.5715	133.5514	-6.83%	4.74%	-12.95%	-10.97%	1.779	-11.05%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10634	0.080231	4.3558	96.8474	51.7965	130.2011	133.171	-3.15%	3.59%	-13.20%	-11.22%	1.870	-6.51%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10262	0.080231	4.2033	94.9345	70.2915	133.3161	136.3665	-5.07%	40.58%	-11.12%	-9.09%	1.351	-32.47%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.098769	0.080231	4.0456	102.2465	81.425	138.4834	141.6532	2.25%	62.85%	-7.68%	-5.56%	1.256	-37.21%

Table B.85. ¹²³I β-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo6 hi6			
SC-Total			100	20	120				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10805	0.080336	4.4256	97.9565	19.4058	107.7067	111.3608	-2.04%	-2.97%	-10.24%	-7.20%	5.048	0.96%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1075	0.080336	4.4032	97.7337	28.9073	112.6062	116.4689	-2.27%	44.54%	-6.16%	-2.94%	3.381	-32.36%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10687	0.080336	4.3774	99.1358	34.7897	114.8699	118.8251	-0.86%	73.95%	-4.28%	-0.98%	2.850	-43.01%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10615	0.080336	4.348	102.6642	34.5724	114.5148	118.4546	2.66%	72.86%	-4.57%	-1.29%	2.970	-40.61%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10244	0.080336	4.1958	102.258	53.3105	117.5889	121.6747	2.26%	166.55%	-2.01%	1.40%	1.918	-61.64%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.098591	0.080336	4.0383	110.6737	65.4467	123.5409	127.843	10.67%	227.23%	2.95%	6.54%	1.691	-66.18%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10805	0.080336	4.4256	94.2683	9.3432	95.0502	97.1628	-5.73%	-53.28%	-20.79%	-19.03%	10.090	101.79%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1075	0.080336	4.4032	93.0835	18.7465	99.4067	101.6395	-6.92%	-6.27%	-17.16%	-15.30%	4.965	-0.69%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10687	0.080336	4.3774	93.8907	24.4042	101.4253	103.7133	-6.11%	22.02%	-15.48%	-13.57%	3.847	-23.05%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10615	0.080336	4.348	97.3587	23.893	101.127	103.4064	-2.64%	19.47%	-15.73%	-13.83%	4.075	-18.50%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.10244	0.080336	4.1958	95.1048	42.5743	103.9312	106.2973	-4.90%	112.87%	-13.39%	-11.42%	2.234	-55.32%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.098591	0.080336	4.0383	102.2726	53.9163	109.2921	111.7857	2.27%	169.58%	-8.92%	-6.85%	1.897	-62.06%

Table B.86. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo6 hi6			
SC-Total			40	40	80				139-144			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13								173-178			
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023962	0.017877	0.9815	39.7837	38.2759	76.7179	79.399	-0.54%	-4.31%	-4.10%	-0.75%	1.039	3.94%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023841	0.017877	0.97653	41.4143	41.4059	80.4794	83.2921	3.54%	3.51%	0.60%	4.12%	1.000	0.02%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023701	0.017877	0.9708	42.3299	42.9011	81.6349	84.488	5.82%	7.25%	2.04%	5.61%	0.987	-1.33%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023542	0.017877	0.96428	43.3772	42.9642	81.7551	84.6123	8.44%	7.41%	2.19%	5.77%	1.010	0.96%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022718	0.017877	0.93052	43.6113	47.2691	82.6952	85.5849	9.03%	18.17%	3.37%	6.96%	0.923	-7.74%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021865	0.017877	0.89559	45.4734	50.1448	83.8892	86.8205	13.68%	25.36%	4.86%	8.53%	0.907	-9.32%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023962	0.017877	0.9815	35.2727	33.6406	67.728	69.2789	-11.82%	-15.90%	-15.34%	-13.40%	1.049	4.85%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023841	0.017877	0.97653	36.5611	36.5651	71.0566	72.6837	-8.60%	-8.59%	-11.18%	-9.15%	1.000	-0.01%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023701	0.017877	0.9708	37.3172	37.9473	72.082	73.7326	-6.71%	-5.13%	-9.90%	-7.83%	0.983	-1.66%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.023542	0.017877	0.96428	38.3354	37.9183	72.1923	73.8455	-4.16%	-5.20%	-9.76%	-7.69%	1.011	1.10%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.022718	0.017877	0.93052	38.1417	42.1787	73.0418	74.7142	-4.66%	5.45%	-8.70%	-6.61%	0.904	-9.57%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.021865	0.017877	0.89559	39.7106	44.8791	74.1204	75.8174	-0.72%	12.20%	-7.35%	-5.23%	0.885	-11.52%

Table B.87. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	20	60					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024034	0.017744	0.96442	39.3211	19.3413	57.0751	59.0547	-1.70%	-3.29%	-4.87%	-1.58%	2.033	1.65%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023912	0.017744	0.97944	40.2785	22.1745	59.8704	61.9514	0.70%	10.87%	-0.22%	3.25%	1.816	-9.18%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023772	0.017744	0.97369	40.8339	23.7174	60.771	62.8852	2.08%	18.59%	1.29%	4.81%	1.722	-13.92%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023612	0.017744	0.96716	41.7174	23.7911	60.7896	62.904	4.29%	18.96%	1.32%	4.84%	1.753	-12.33%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022786	0.017744	0.9333	41.7011	28.1147	61.5147	63.6597	4.25%	40.57%	2.52%	6.10%	1.483	-25.84%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02193	0.017744	0.89827	43.4659	31.0139	62.7574	64.9473	8.66%	65.07%	4.60%	8.25%	1.401	-29.93%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.024034	0.017744	0.98442	36.6947	15.0939	50.3824	51.5274	-8.26%	-24.53%	-16.03%	-14.12%	2.431	21.55%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023912	0.017744	0.97944	37.3497	17.794	52.8576	54.0614	-6.63%	-11.03%	-11.90%	-9.90%	2.099	4.95%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023772	0.017744	0.97369	37.7459	19.2592	53.6581	54.8813	-5.64%	-3.70%	-10.57%	-8.53%	1.960	-2.01%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023612	0.017744	0.96716	38.6017	19.2567	53.6785	54.9018	-3.50%	-3.72%	-10.54%	-8.50%	2.005	0.23%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022786	0.017744	0.9333	38.1574	23.6608	54.339	55.5805	-4.61%	17.80%	-9.44%	-7.37%	1.620	-19.02%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02193	0.017744	0.89827	39.6254	26.2949	55.46	56.7279	-0.94%	31.47%	-7.57%	-5.45%	1.507	-24.65%	

Table B.88. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			40	8	48					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023993	0.017749	0.96275	39.0574	7.9948	44.7253	46.244	-2.36%	-0.07%	-6.82%	-3.66%	4.885	-2.29%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023671	0.017749	0.97777	39.6151	10.6548	46.9155	48.5199	-0.96%	33.19%	-2.26%	1.08%	3.718	-25.64%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023731	0.017749	0.97203	39.9625	12.2296	47.6343	49.2688	-0.09%	52.87%	-0.76%	2.64%	3.268	-34.65%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023572	0.017749	0.96551	40.7579	12.3132	47.5941	49.2268	1.89%	53.92%	-0.85%	2.56%	3.310	-33.80%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022747	0.017749	0.93171	40.6084	16.6905	48.2977	49.9676	1.52%	108.63%	0.62%	4.10%	2.433	-51.34%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021893	0.017749	0.89674	42.3032	19.6356	49.5385	51.2561	5.76%	145.45%	3.21%	6.78%	2.154	-56.91%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023993	0.017749	0.98275	37.5602	3.9784	39.4714	40.3494	-6.10%	-50.27%	-17.77%	-15.94%	9.441	88.82%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023671	0.017749	0.97777	37.8389	6.5482	41.413	42.3409	-5.40%	-18.15%	-13.72%	-11.79%	5.779	77.85%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023731	0.017749	0.97203	38.027	8.0655	42.0535	42.9989	-4.93%	0.82%	-12.39%	-10.42%	4.715	-28.52%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.023572	0.017749	0.96551	38.795	8.0814	42.0217	42.9662	-3.01%	1.02%	-12.45%	-10.49%	4.801	-19.95%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022747	0.017749	0.93171	38.2135	12.4514	42.6655	43.632	-4.47%	55.64%	-11.11%	-9.10%	3.069	-38.62%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021893	0.017749	0.89674	39.6068	15.2378	43.785	44.7796	-0.96%	90.47%	-8.78%	-6.71%	2.599	-48.02%	

Table B.89. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data, Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC-Total			84	84	168					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022806	0.016692	0.93414	83.502	80.8644	162.2958	167.9677	-0.59%	-3.73%	-3.40%	-0.02%	1.033	3.26%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022691	0.016692	0.92941	87.8468	85.7919	170.2698	176.2205	4.58%	2.13%	1.35%	4.89%	1.024	2.40%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022557	0.016692	0.92395	89.4397	87.7959	172.3545	178.3781	6.48%	4.52%	2.59%	6.18%	1.019	1.87%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022406	0.016692	0.91776	90.7409	88.1347	172.7217	178.7581	8.02%	4.92%	2.81%	6.40%	1.030	2.96%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021622	0.016692	0.88562	91.0605	92.7255	173.4616	179.5238	8.41%	10.39%	3.25%	6.86%	0.982	-1.80%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02081	0.016692	0.85238	92.5338	96.034	174.6787	180.7833	10.16%	14.33%	3.98%	7.61%	0.964	-3.64%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022806	0.016692	0.93414	73.9801	71.1269	143.278	146.559	-11.93%	-15.33%	-14.72%	-12.76%	1.040	4.01%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022691	0.016692	0.92941	77.7591	75.5517	150.3322	153.7747	-7.43%	-10.06%	-10.52%	-8.47%	1.029	2.92%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022557	0.016692	0.92395	79.1341	77.366	152.1799	155.6647	-5.79%	-7.90%	-9.42%	-7.34%	1.023	2.29%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022406	0.016692	0.91776	80.3743	77.5842	152.5085	156.0008	-4.32%	-7.64%	-9.22%	-7.14%	1.036	3.60%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021622	0.016692	0.88562	80.2388	82.1181	153.1816	156.6893	-4.48%	-2.24%	-8.82%	-6.73%	0.977	-2.29%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02081	0.016692	0.85238	81.3792	85.2851	154.2795	157.8123	-3.12%	1.53%	-8.17%	-6.06%	0.954	-4.58%	

Table B.90. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC-Total			84	42	126					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022739	0.016556	0.93138	82.5649	41.0436	121.2591	125.4657	-1.71%	-2.28%	-3.76%	-0.42%	2.012	0.58%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022624	0.016556	0.92667	85.5053	45.3595	127.2534	131.6736	1.79%	8.00%	0.99%	4.50%	1.885	-5.75%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022491	0.016556	0.92123	86.3539	47.4739	128.8982	133.3795	2.80%	13.03%	2.30%	5.86%	1.819	-9.05%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02234	0.016556	0.91505	87.3296	47.8423	129.1372	133.6256	3.96%	13.91%	2.49%	6.05%	1.825	-8.73%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021558	0.016556	0.88301	87.1643	52.5287	129.796	134.3044	3.77%	25.07%	3.00%	6.59%	1.659	-17.03%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020749	0.016556	0.84987	88.4262	55.926	130.9442	135.5064	5.27%	33.16%	3.92%	7.54%	1.581	-20.94%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022739	0.016556	0.93138	77.0059	32.1184	107.0405	109.4735	-8.33%	-23.53%	-15.05%	-13.12%	2.398	19.88%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022624	0.016556	0.92667	79.4626	36.0814	112.345	114.9018	-5.40%	-14.09%	-10.84%	-8.81%	2.202	10.12%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022491	0.016556	0.92123	80.0927	38.0799	113.8049	116.3967	-4.65%	-9.33%	-9.68%	-7.62%	2.103	5.16%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.02234	0.016556	0.91505	81.0088	38.3697	114.0203	116.6169	-3.56%	-8.67%	-9.51%	-7.45%	2.112	5.59%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021558	0.016556	0.88301	80.3782	43.0398	114.6136	117.228	-4.31%	2.48%	-9.04%	-6.96%	1.868	-6.62%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020749	0.016556	0.84987	81.2952	46.3194	115.659	118.2993	-3.22%	10.28%	-8.21%	-6.11%	1.755	-12.24%	

Table B.91. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 139-144 keV, 142-175 keV, 173-178 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo6 hi6			
SC-Total			84	16.8	100.8					139-144			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13									173-178			
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022746	0.016451	0.93168	81.9953	17.1425	95.5555	98.8018	-2.39%	2.04%	-5.20%	-1.98%	4.783	-4.34%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022631	0.016451	0.92696	84.0887	21.0884	100.3462	103.7715	0.11%	25.53%	-0.45%	2.95%	3.987	-20.25%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022498	0.016451	0.92152	84.4891	23.2672	101.6121	105.089	0.58%	38.50%	0.81%	4.25%	3.631	-27.37%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022347	0.016451	0.91534	85.2719	23.6504	101.741	105.2232	1.51%	40.78%	0.93%	4.39%	3.606	-27.89%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021565	0.016451	0.88329	84.8373	28.3728	102.2657	105.7835	1.00%	68.89%	1.45%	4.94%	2.990	-40.20%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020755	0.016451	0.85013	85.9792	31.6285	103.4498	107.0179	2.36%	89.46%	2.63%	6.17%	2.701	-45.97%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022746	0.016451	0.93168	78.8148	8.7057	84.3314	86.2083	-6.17%	-48.18%	-16.34%	-14.48%	9.053	81.06%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022631	0.016451	0.92696	80.4744	12.3903	88.5736	90.5544	-4.20%	-26.25%	-12.13%	-10.16%	6.495	29.90%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022498	0.016451	0.92152	80.6565	14.4961	89.7	91.7108	-3.98%	-13.71%	-11.01%	-9.02%	5.564	11.28%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.022347	0.016451	0.91534	81.381	14.8094	89.8181	91.8321	-3.12%	-11.65%	-10.89%	-8.90%	5.495	9.90%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.021565	0.016451	0.88329	80.479	19.5525	90.3034	92.3386	-4.19%	16.38%	-10.41%	-8.39%	4.116	-17.68%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.020755	0.016451	0.85013	81.2689	22.9029	91.3725	93.4373	-3.25%	36.33%	-9.35%	-7.30%	3.548	-29.03%	

Table B.92. ¹²³I β-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13_h10			
SC-Total			100	100	200					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13363	0.05932	5.4734	85.4401	92.0137	173.777	179.8498	-14.56%	-7.99%	-13.11%	-10.08%	0.929	-7.14%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13295	0.05932	5.4457	88.0216	99.956	181.7591	188.1097	-11.98%	-0.04%	-9.12%	-5.95%	0.881	-11.94%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13217	0.05932	5.4138	90.6672	102.7235	184.2469	190.6845	-9.33%	2.72%	-7.88%	-4.66%	0.883	-11.74%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13129	0.05932	5.3775	91.2102	105.9974	185.276	191.7488	-8.79%	6.00%	-7.36%	-4.13%	0.860	-13.95%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12669	0.05932	5.1892	96.3169	112.1593	184.8707	191.3291	-3.68%	12.16%	-7.56%	-4.34%	0.859	-14.12%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12193	0.05932	4.9944	101.253	122.365	189.635	196.2583	1.25%	22.37%	-5.18%	-1.87%	0.827	-17.25%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13363	0.05932	5.4734	74.777	81.874	153.4006	156.9131	-25.22%	-18.13%	-23.30%	-21.54%	0.913	-8.67%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13295	0.05932	5.4457	76.5209	89.4446	160.4628	164.1363	-23.48%	-10.56%	-19.77%	-17.93%	0.856	-14.45%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13217	0.05932	5.4138	78.8569	91.9069	162.6678	166.3919	-21.14%	-8.09%	-18.67%	-16.80%	0.858	-14.20%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.13129	0.05932	5.3775	79.0601	95.1025	163.5896	167.3343	-20.94%	-4.90%	-18.21%	-16.33%	0.831	-16.87%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12669	0.05932	5.1892	83.4986	100.767	163.2848	167.0224	-16.50%	0.77%	-18.36%	-16.49%	0.829	-17.14%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12193	0.05932	4.9944	87.3982	110.468	167.5509	171.3852	-12.60%	10.47%	-16.22%	-14.31%	0.791	-20.86%	

Table B.93. ^{123}I β -CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Io13 h0		
SC-Total			100	50	150					130-142		
Atten Coeff	0.1481									142-175		
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.13026	0.059651	5.3353	85.7867	43.5601	124.4859	128.8062	-14.21%	-12.88%	-17.01%	-14.13%	1.969	-1.53%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1296	0.059651	5.3083	87.4491	49.8642	129.9025	134.4204	-12.55%	-0.27%	-13.40%	-10.39%	1.754	-12.31%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12884	0.059651	5.2771	89.6748	52.0198	131.6085	136.187	-10.33%	4.04%	-12.26%	-9.21%	1.724	-13.81%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12797	0.059651	5.2417	90.398	54.4762	132.17	136.771	-9.60%	8.95%	-11.89%	-8.82%	1.659	-17.03%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12349	0.059651	5.0582	95.8005	59.3223	131.3233	135.8972	-4.20%	18.64%	-12.45%	-9.40%	1.615	-19.25%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11886	0.059651	4.8683	100.7006	69.1698	136.042	140.7889	0.70%	38.34%	-9.31%	-6.14%	1.456	-27.21%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.13026	0.059651	5.3353	79.9193	34.2625	109.8762	112.3747	-20.08%	-31.48%	-26.75%	-25.08%	2.333	16.63%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.1296	0.059651	5.3083	80.9213	40.313	114.673	117.2862	-19.08%	-19.37%	-23.55%	-21.81%	2.007	0.37%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12884	0.059651	5.2771	82.9036	42.2149	116.1864	118.8347	-17.10%	-15.57%	-22.54%	-20.78%	1.964	-1.81%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12797	0.059651	5.2417	83.3661	44.5888	116.6954	119.357	-16.63%	-10.83%	-22.20%	-20.43%	1.870	-6.51%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.12349	0.059651	5.0582	88.2254	48.9198	115.9972	118.6444	-11.77%	-2.16%	-22.67%	-20.90%	1.803	-9.83%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.11886	0.059651	4.8683	92.125	58.2641	120.222	122.9696	-7.88%	16.53%	-19.85%	-18.02%	1.581	-20.94%

Table B.94. ^{123}I β -CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW to13 h0			
SC-Total			100	20	120					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12623	0.059836	5.2523	85.9954	14.4855	93.5976	96.7605	-14.00%	-27.57%	-22.00%	-19.37%	5.937	18.73%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12758	0.059836	5.2257	87.1035	19.8078	97.5028	100.8227	-12.90%	-0.96%	-18.75%	-15.98%	4.397	-12.05%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12683	0.059836	5.195	89.0507	21.6209	98.5421	101.9031	-10.95%	8.10%	-17.88%	-15.08%	4.119	-17.63%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12598	0.059836	5.1602	89.7996	23.6477	98.7457	102.1212	-10.20%	18.24%	-17.71%	-14.90%	3.797	-24.05%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12157	0.059836	4.9795	95.2313	27.8727	97.445	100.7865	-4.77%	39.36%	-18.80%	-16.01%	3.417	-31.67%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11701	0.059836	4.7926	100.1157	37.4978	102.2005	105.7275	0.12%	87.49%	-14.83%	-11.89%	2.670	-46.60%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12623	0.059836	5.2523	83.0055	5.6936	62.5899	84.4141	-16.99%	-71.53%	-31.18%	-29.65%	14.579	191.57%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12758	0.059836	5.2257	83.5596	10.834	86.0502	87.9698	-16.44%	-45.83%	-28.29%	-26.69%	7.713	54.25%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12683	0.059836	5.195	85.301	12.4267	86.9746	88.9181	-14.70%	-37.87%	-27.52%	-25.90%	6.864	37.29%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12598	0.059836	5.1602	85.8316	14.3725	87.1689	89.1211	-14.17%	-28.14%	-27.36%	-25.73%	5.972	19.44%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.12157	0.059836	4.9795	90.7805	18.0856	86.0673	88.0009	-9.22%	-9.57%	-28.28%	-26.67%	5.019	0.39%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.11701	0.059836	4.7926	94.6868	27.2069	90.3262	92.3685	-5.31%	36.03%	-24.73%	-23.03%	3.480	-30.40%	

Table B.95. ¹²³I FP-CIT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW Iof13_h10			
SC-Total			40	40	80					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028787	0.013703	1.1791	35.2466	35.1511	69.5525	71.9833	-11.88%	-12.12%	-13.06%	-10.02%	1.003	0.27%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028641	0.013703	1.1731	36.766	37.4866	72.7867	75.3305	-8.09%	-6.28%	-9.02%	-5.84%	0.981	-1.92%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028473	0.013703	1.1663	37.706	38.194	73.6522	76.2262	-5.73%	-4.51%	-7.93%	-4.72%	0.987	-1.28%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028282	0.013703	1.1584	37.8873	39.0837	74.0484	76.6363	-5.28%	-2.29%	-7.44%	-4.20%	0.969	-3.06%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027292	0.013703	1.1179	39.3138	40.4165	73.9711	76.5562	-1.72%	1.04%	-7.54%	-4.30%	0.973	-2.73%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026267	0.013703	1.0759	40.3976	42.8773	75.0009	77.6219	0.99%	7.19%	-6.25%	-2.97%	0.942	-5.78%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028787	0.013703	1.1791	31.1258	31.0201	61.3966	62.8045	-22.19%	-22.45%	-23.25%	-21.49%	1.003	0.34%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028641	0.013703	1.1731	32.3872	33.1697	64.2597	65.7312	-19.03%	-17.08%	-19.68%	-17.84%	0.976	-2.36%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028473	0.013703	1.1663	33.2449	33.7721	65.0262	66.5153	-16.89%	-15.57%	-18.72%	-16.86%	0.984	-1.56%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.028282	0.013703	1.1584	33.3348	34.6351	65.3792	66.8763	-16.66%	-13.41%	-18.28%	-16.40%	0.962	-3.75%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027292	0.013703	1.1179	34.614	35.8298	65.3234	66.8192	-13.47%	-10.43%	-18.35%	-16.48%	0.966	-3.39%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026267	0.013703	1.0759	35.4496	38.1801	66.2471	67.764	-11.38%	-4.55%	-17.19%	-15.30%	0.928	-7.15%	

Table B.96. ¹²³I FP-CIT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW fo13_h10			
SC-Total			40	20	60				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028307	0.013741	1.1594	35.1587	17.3156	51.4777	53.2632	-12.10%	-13.42%	-14.20%	-11.23%	2.030	1.52%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028163	0.013741	1.1536	36.2835	19.2048	53.8565	55.7272	-9.29%	-3.98%	-10.24%	-7.12%	1.889	-5.54%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027998	0.013741	1.1468	37.04	19.7969	54.4909	56.3839	-7.40%	-1.02%	-9.18%	-6.03%	1.871	-6.45%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02781	0.013741	1.1391	37.2214	20.525	54.7582	56.6616	-6.95%	2.62%	-8.74%	-5.56%	1.813	-9.33%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026837	0.013741	1.0992	38.63	21.713	54.6824	56.584	-3.42%	8.57%	-8.86%	-5.69%	1.779	-11.04%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025829	0.013741	1.058	39.7306	24.0833	55.729	57.6695	-0.67%	20.42%	-7.12%	-3.88%	1.650	-17.51%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028307	0.013741	1.1594	32.8072	13.5158	45.4384	46.471	-17.98%	-32.42%	-24.27%	-22.55%	2.427	21.37%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028163	0.013741	1.1536	33.7239	15.2662	47.5433	48.6252	-15.69%	-23.67%	-20.76%	-18.96%	2.209	10.45%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027998	0.013741	1.1468	34.4118	15.7736	48.1054	49.2003	-13.97%	-21.13%	-19.82%	-18.00%	2.182	9.08%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.02781	0.013741	1.1391	34.5173	16.4779	48.3447	49.4458	-13.71%	-17.61%	-19.43%	-17.59%	2.095	4.74%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026837	0.013741	1.0992	35.7931	17.5306	48.2901	49.3904	-10.52%	-12.35%	-19.52%	-17.68%	2.042	2.08%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025829	0.013741	1.058	36.6537	19.7888	49.2288	50.3519	-8.37%	-1.06%	-17.95%	-16.08%	1.852	-7.39%

Table B.97. ¹²³I FP-CIT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW Iof13 hi6			
SC-Total			40	8	48				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027978	0.013821	1.146	35.113	6.6211	40.1723	41.5335	-12.22%	-17.24%	-16.31%	-13.47%	5.303	6.06%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027836	0.013821	1.1402	36.0022	8.2453	42.0432	43.4749	-9.99%	3.07%	-12.41%	-9.43%	4.366	-12.67%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027673	0.013821	1.1335	36.6534	8.7696	42.5333	43.9633	-8.37%	9.61%	-11.39%	-8.37%	4.180	-16.40%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027487	0.013821	1.1259	36.8432	9.3988	42.7281	44.1874	-7.89%	17.49%	-10.98%	-7.94%	3.920	-21.60%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026525	0.013821	1.0864	38.2688	10.5005	42.6309	44.09	-4.33%	31.26%	-11.19%	-8.15%	3.644	-27.11%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025529	0.013821	1.0457	39.4029	12.8258	43.6899	45.1921	-1.49%	60.32%	-8.98%	-5.85%	3.072	-38.56%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027978	0.013821	1.146	33.8222	3.0192	35.4494	36.2363	-15.44%	-62.26%	-26.15%	-24.51%	11.202	124.05%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027836	0.013821	1.1402	34.5332	4.5325	37.106	37.9338	-13.67%	-43.34%	-22.70%	-20.97%	7.619	261.90%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027673	0.013821	1.1335	35.1239	4.9828	37.5408	38.3793	-12.19%	-37.72%	-21.79%	-20.04%	7.049	204.90%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027487	0.013821	1.1259	35.247	5.5904	37.7166	38.5605	-11.88%	-30.12%	-21.42%	-19.67%	6.305	130.49%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026525	0.013821	1.0864	36.5483	6.5554	37.6441	38.4882	-8.63%	-18.06%	-21.57%	-19.82%	5.575	11.51%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025529	0.013821	1.0457	37.4465	8.765	38.5944	39.4638	-6.38%	9.56%	-19.60%	-17.78%	4.272	-14.55%

Table B.98. ¹²³I IBZM 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC-Total			84	84	168				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028854	0.013172	1.1818	75.1144	74.9078	149.0658	154.2754	-10.58%	-10.82%	-11.27%	-8.17%	1.003	0.28%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028707	0.013172	1.1759	78.912	79.0864	156.2332	161.6934	-6.06%	-5.85%	-7.00%	-3.75%	0.998	-0.22%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028539	0.013172	1.169	80.6221	80.3008	158.1358	163.6625	-4.02%	-4.40%	-5.87%	-2.58%	1.004	0.40%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028348	0.013172	1.1611	80.9015	81.7549	158.9479	164.5029	-3.69%	-2.67%	-5.39%	-2.08%	0.990	-1.04%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027355	0.013172	1.1205	82.6346	83.6633	159.52	165.0951	-1.63%	-0.40%	-5.05%	-1.73%	0.988	-1.23%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026328	0.013172	1.0784	83.5664	86.5535	160.5855	166.1976	-0.52%	3.04%	-4.41%	-1.07%	0.965	-3.45%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028854	0.013172	1.1818	66.3322	66.107	131.594	134.6074	-21.03%	-21.30%	-21.67%	-19.88%	1.003	0.34%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028707	0.013172	1.1759	69.6528	69.8433	137.9344	141.093	-17.08%	-16.85%	-17.90%	-16.02%	0.997	-0.27%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028539	0.013172	1.169	71.221	70.8674	139.6199	142.8171	-15.21%	-15.63%	-16.89%	-14.99%	1.005	0.50%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028348	0.013172	1.1611	71.3513	72.2768	140.3412	143.5549	-15.06%	-13.96%	-16.46%	-14.55%	0.987	-1.28%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027355	0.013172	1.1205	72.8761	74.0136	140.8635	144.0892	-13.24%	-11.89%	-16.15%	-14.23%	0.985	-1.54%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026328	0.013172	1.0784	73.519	76.8057	141.8205	145.068	-12.48%	-8.56%	-15.58%	-13.65%	0.957	-4.28%

Table B.99. ¹²³I IBZM 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)				TEW lo13 hi6			
SC-Total			84	42	126				130-142			
Atten Coeff	0.1481								142-175			
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028053	0.012927	1.1491	74.8946	37.3293	110.976	114.8261	-10.84%	-11.12%	-11.92%	-8.87%	2.006	0.32%
Atten Coeff	0.1481											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027911	0.012927	1.1432	77.8523	40.5558	116.4443	120.4876	-7.32%	-3.44%	-7.58%	-4.37%	1.920	-4.02%
Atten Coeff	0.1481											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027748	0.012927	1.1365	79.1628	41.516	117.8418	121.9341	-5.76%	-1.15%	-6.47%	-3.23%	1.907	-4.66%
Atten Coeff	0.1481											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027561	0.012927	1.1289	79.4277	42.6178	118.4214	122.5356	-5.44%	1.47%	-6.01%	-2.75%	1.864	-6.81%
Atten Coeff	0.1481											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026596	0.012927	1.0894	81.0885	44.1304	118.8429	122.9731	-3.47%	5.07%	-5.68%	-2.40%	1.837	-8.13%
Atten Coeff	0.1481											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025598	0.012927	1.0485	82.022	46.713	119.795	123.9617	-2.35%	11.22%	-4.92%	-1.62%	1.756	-12.21%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.028053	0.012927	1.1491	69.8407	29.2296	97.9601	100.1868	-16.86%	-30.41%	-22.25%	-20.49%	2.389	19.47%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027911	0.012927	1.1432	72.4273	32.1144	102.7981	105.1368	-13.78%	-23.54%	-18.41%	-16.56%	2.255	12.76%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027748	0.012927	1.1365	73.6268	32.9278	104.0366	106.4038	-12.35%	-21.60%	-17.43%	-15.55%	2.236	11.80%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.027561	0.012927	1.1289	73.7769	33.9924	104.5523	106.9322	-12.17%	-19.07%	-17.02%	-15.13%	2.170	8.52%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.026596	0.012927	1.0894	75.2688	35.3433	104.9405	107.3301	-10.39%	-15.85%	-16.71%	-14.82%	2.130	6.48%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.025598	0.012927	1.0485	75.9424	37.8294	105.7964	108.2075	-9.59%	-9.93%	-16.03%	-14.12%	2.007	0.37%

Table B.100. ¹²³I IBZM 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 130-142 keV, 142-175 keV, 000-000 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW lo13 hi6			
SC-Total			84	16.8	100.8					130-142			
Atten Coeff	0.1481									142-175			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027427	0.013008	1.1234	74.787	14.8068	87.3836	90.3482	-10.97%	-11.86%	-13.31%	-10.37%	5.051	1.02%	
Atten Coeff	0.1481												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027288	0.013008	1.1177	77.2477	17.4704	91.7956	94.9206	-8.04%	3.99%	-8.93%	-5.83%	4.422	-11.57%	
Atten Coeff	0.1481												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027128	0.013008	1.1111	78.328	18.286	92.8476	96.011	-6.75%	8.85%	-7.89%	-4.75%	4.283	-14.33%	
Atten Coeff	0.1481												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026946	0.013008	1.1037	78.5962	19.185	93.2859	96.468	-6.43%	14.20%	-7.45%	-4.30%	4.097	-18.06%	
Atten Coeff	0.1481												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026002	0.013008	1.065	80.2464	20.5162	93.6339	96.8322	-4.47%	22.12%	-7.11%	-3.94%	3.911	-21.77%	
Atten Coeff	0.1481												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025026	0.013008	1.0251	81.2167	22.9481	94.6742	97.8136	-3.31%	36.60%	-6.18%	-2.96%	3.539	-29.22%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027427	0.013008	1.1234	71.9673	7.1248	77.1152	78.8291	-14.32%	-57.59%	-23.50%	-21.80%	10.101	102.02%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027288	0.013008	1.1177	74.1197	9.5062	81.02	82.8269	-11.76%	-43.42%	-19.62%	-17.83%	7.797	55.94%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.027128	0.013008	1.1111	75.1063	10.2004	81.9533	83.7825	-10.59%	-39.28%	-18.70%	-16.88%	7.363	47.26%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026946	0.013008	1.1037	75.2791	11.0655	82.3445	84.1847	-10.38%	-34.13%	-18.31%	-16.48%	6.803	36.06%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.026002	0.013008	1.065	76.7785	12.237	82.6679	84.5178	-8.60%	-27.16%	-17.99%	-16.15%	6.274	25.49%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.025026	0.013008	1.0251	77.5032	14.5693	83.5138	85.3879	-7.73%	-13.28%	-17.15%	-15.29%	5.320	6.39%	

Table B.101. ^{99m}I TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	126.4	252.8					114-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025554	0.004482	0.10467	128.4074	128.2244	254.9849	264.4495	1.59%	1.44%	0.86%	4.61%	1.001	0.14%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025424	0.004482	0.10414	135.6278	133.4698	266.6251	276.5217	7.30%	5.59%	5.47%	9.38%	1.016	1.62%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025275	0.004482	0.10353	136.5808	137.0725	270.0693	280.0938	8.05%	8.44%	6.83%	10.80%	0.996	-0.36%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025106	0.004482	0.10283	136.9464	138.9665	271.3133	281.3839	8.34%	9.94%	7.32%	11.31%	0.995	-1.45%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024227	0.004482	0.099232	138.2019	141.1666	271.7951	281.8836	9.34%	11.68%	7.51%	11.50%	0.979	-2.10%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0023317	0.004482	0.095508	139.5987	141.7156	271.6774	281.7616	10.44%	12.12%	7.47%	11.46%	0.985	-1.49%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025554	0.004482	0.10467	111.1646	110.9725	220.7099	225.764	-12.05%	-12.21%	-12.69%	-10.69%	1.002	0.17%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025424	0.004482	0.10414	117.6491	115.3054	230.8094	236.0948	-6.92%	-8.78%	-8.70%	-6.61%	1.020	2.03%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025275	0.004482	0.10353	118.1778	118.7385	233.8053	239.1593	-6.50%	-6.06%	-7.51%	-5.40%	0.995	-0.47%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025106	0.004482	0.10283	118.3252	120.5595	234.8903	240.2691	-6.39%	-4.62%	-7.08%	-4.96%	0.981	-1.85%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024227	0.004482	0.099232	119.3205	122.6058	235.3219	240.7106	-5.60%	-3.00%	-6.91%	-4.78%	0.973	-2.68%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0023317	0.004482	0.095508	120.6565	123.002	235.2288	240.6153	-4.54%	-2.69%	-6.95%	-4.62%	0.981	-1.91%	

Table B.102. ^{99m}I TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	63.2	189.6					114-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027246	0.00404	0.1116	129.3225	63.4475	190.5369	197.5525	2.31%	0.39%	0.49%	4.19%	2.038	1.91%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027108	0.00404	0.11104	135.0164	67.1632	199.3044	206.645	6.82%	6.27%	5.12%	8.99%	2.010	0.51%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026949	0.00404	0.11038	135.384	70.3077	202.0091	209.4559	7.11%	11.25%	6.54%	10.47%	1.926	-3.72%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026768	0.00404	0.10964	135.0345	72.4289	202.8536	210.3364	6.83%	14.60%	6.99%	10.94%	1.864	-6.78%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025831	0.00404	0.1058	134.8709	75.4878	203.1532	210.6533	6.70%	19.44%	7.15%	11.10%	1.787	-10.67%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024862	0.00404	0.10183	135.5446	76.6416	203.1303	210.6311	7.23%	21.27%	7.14%	11.09%	1.769	-11.57%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027246	0.00404	0.1116	119.3432	47.5157	164.9079	168.6536	-5.58%	-24.82%	-13.02%	-11.05%	2.512	25.58%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027108	0.00404	0.11104	124.502	50.5224	172.5152	176.4349	-1.50%	-20.06%	-9.01%	-6.94%	2.464	23.21%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026949	0.00404	0.11038	124.5068	53.572	174.8693	178.8461	-1.50%	-15.23%	-7.77%	-5.67%	2.324	16.21%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026768	0.00404	0.10964	123.9242	55.699	175.6086	179.6048	-1.96%	-11.67%	-7.38%	-5.27%	2.225	11.24%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0025831	0.00404	0.1058	123.4159	58.7561	175.8812	179.8869	-2.36%	-7.03%	-7.24%	-5.12%	2.100	5.02%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0024862	0.00404	0.10183	123.9521	59.8447	175.869	179.8754	-1.94%	-5.31%	-7.24%	-5.13%	2.071	3.56%	

Table B.103. ^{99m}I TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 114-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	25.28	151.68					114-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027704	0.003759	0.11348	129.881	24.5909	150.3077	155.7083	2.75%	-2.73%	-0.90%	2.66%	5.282	5.63%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027564	0.003759	0.1129	134.6619	27.3917	157.2514	162.9124	6.54%	8.35%	3.67%	7.41%	4.916	-1.68%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027402	0.003759	0.11224	134.6795	30.2658	159.3976	165.1515	6.55%	19.72%	5.09%	8.88%	4.450	-11.00%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027218	0.003759	0.11149	133.9048	32.5259	159.9616	165.7483	5.94%	28.66%	5.46%	9.27%	4.117	-17.66%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026265	0.003759	0.10758	132.8879	36.1301	160.0757	165.8855	5.13%	42.92%	5.54%	9.37%	3.678	-26.44%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.002528	0.003759	0.10354	133.1891	37.6457	160.1759	165.9963	5.37%	48.91%	5.60%	9.44%	3.538	-29.24%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027704	0.003759	0.11348	124.2585	9.4498	130.0493	132.9315	-1.69%	-62.62%	-14.26%	-12.36%	13.149	162.99%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027564	0.003759	0.1129	128.6245	11.6634	136.0755	139.0971	1.76%	-53.86%	-10.29%	-8.30%	11.028	120.56%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027402	0.003759	0.11224	128.3161	14.4868	137.9465	141.018	1.52%	-42.69%	-9.05%	-7.03%	8.857	77.15%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0027218	0.003759	0.11149	127.299	16.7998	138.4444	141.5337	0.71%	-33.55%	-8.73%	-6.69%	7.577	51.55%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0026265	0.003759	0.10758	125.8839	20.4937	138.5589	141.6609	-0.41%	-18.93%	-8.65%	-6.61%	6.143	22.85%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.002528	0.003759	0.10354	126.0004	21.9912	138.6554	141.7634	-0.32%	-13.01%	-8.59%	-6.54%	5.730	14.59%	

Table B.104. ^{99m}I TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	126.4	252.8					121-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038806	0.005219	0.15895	130.304	126.4417	254.5439	263.992	3.09%	0.03%	0.69%	4.43%	1.031	3.05%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003961	0.005219	0.15815	136.933	132.3449	265.8976	275.767	8.33%	4.70%	5.18%	9.09%	1.035	3.47%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038383	0.005219	0.15722	138.1701	135.692	269.2438	279.2376	9.31%	7.35%	6.50%	10.46%	1.018	1.83%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038126	0.005219	0.15616	138.6465	137.4481	270.4815	280.5213	9.69%	8.74%	6.99%	10.97%	1.009	0.87%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036791	0.005219	0.1507	141.2105	138.2072	270.2986	280.3315	11.72%	9.34%	6.92%	10.89%	1.022	2.17%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003541	0.005219	0.14504	142.7504	138.6364	270.611	280.6555	12.94%	9.64%	7.05%	11.02%	1.028	2.82%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038806	0.005219	0.15895	113.2263	109.0094	220.3273	225.3726	-10.42%	-13.76%	-12.85%	-10.85%	1.039	3.87%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003961	0.005219	0.15815	119.0592	114.0513	230.1767	235.4475	-5.81%	-9.77%	-8.95%	-6.86%	1.044	4.39%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038383	0.005219	0.15722	119.8918	117.2043	233.0862	238.4237	-5.15%	-7.28%	-7.80%	-5.69%	1.023	2.29%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038126	0.005219	0.15616	120.1651	118.877	234.1653	239.5275	-4.93%	-5.95%	-7.37%	-5.25%	1.011	1.08%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036791	0.005219	0.1507	122.6165	119.3543	234.0218	239.3807	-2.99%	-5.57%	-7.43%	-5.31%	1.027	2.73%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003541	0.005219	0.14504	124.0822	119.8124	234.3047	239.6701	-1.83%	-5.21%	-7.32%	-5.19%	1.036	3.56%	

Table B.105. ^{99m}I TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	63.2	189.6					121-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040901	0.004672	0.16753	130.4417	62.2612	189.9669	196.9573	3.20%	-1.49%	0.19%	3.88%	2.095	4.75%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040694	0.004672	0.16668	135.7624	66.4266	198.4909	205.7969	7.41%	5.11%	4.69%	8.54%	2.044	2.19%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040455	0.004672	0.1657	136.3428	69.3452	200.9658	208.3917	7.87%	9.72%	6.01%	9.91%	1.966	-1.69%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040184	0.004672	0.16459	136.0965	71.2773	201.8395	209.2814	7.67%	12.78%	6.46%	10.38%	1.909	-4.53%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038777	0.004672	0.15883	137.0446	72.9019	201.6484	209.0858	8.42%	15.35%	6.35%	10.28%	1.880	-6.01%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037321	0.004672	0.15287	137.9057	73.8562	201.8773	209.3243	9.10%	16.86%	6.48%	10.40%	1.867	-6.64%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040901	0.004672	0.16753	120.5746	46.2259	164.412	168.1443	-4.61%	-26.86%	-13.28%	-11.32%	2.608	30.42%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040694	0.004672	0.16668	125.3177	49.7141	171.8069	175.7092	-0.86%	-21.34%	-9.38%	-7.33%	2.521	26.04%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040455	0.004672	0.1657	125.5546	52.5197	173.978	177.9328	-0.67%	-16.90%	-8.24%	-6.15%	2.391	19.53%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040184	0.004672	0.16459	125.0962	54.4491	174.7258	178.7	-1.03%	-13.85%	-7.85%	-5.75%	2.297	14.87%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038777	0.004672	0.15883	125.8457	55.9668	174.5716	178.5437	-0.44%	-11.44%	-7.93%	-5.83%	2.249	12.43%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037321	0.004672	0.15287	126.5885	56.8337	174.7793	178.7567	0.15%	-10.07%	-7.82%	-5.72%	2.227	11.37%	

Table B.106. ^{99m}I TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 121-126 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	25.28	151.68					121-126			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040325	0.004355	0.16517	130.5555	23.7838	149.6434	155.0147	3.29%	-5.92%	-1.34%	2.20%	5.489	9.79%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.004012	0.004355	0.16433	135.0993	26.9175	156.4601	162.0895	6.88%	6.48%	3.15%	6.86%	5.019	0.36%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0039885	0.004355	0.16337	135.2961	29.5898	158.4804	164.1969	7.04%	17.05%	4.48%	8.25%	4.572	-8.55%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0039617	0.004355	0.16227	134.6318	31.639	159.0984	164.8487	6.51%	25.15%	4.89%	8.68%	4.255	-14.90%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003823	0.004355	0.15659	134.688	33.8324	158.8393	164.5916	6.56%	33.83%	4.72%	8.51%	3.981	-20.36%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036795	0.004355	0.15071	135.201	35.0584	159.0909	164.8581	6.96%	38.68%	4.89%	8.69%	3.856	-22.87%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040325	0.004355	0.16517	125.0106	8.5826	129.4716	132.338	-1.10%	-66.05%	-14.64%	-12.75%	14.566	191.31%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.004012	0.004355	0.16433	129.1071	11.1482	135.3873	138.3919	2.14%	-55.90%	-10.74%	-8.76%	11.581	131.62%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0039885	0.004355	0.16337	128.9955	13.7549	137.1482	140.1996	2.05%	-45.59%	-9.58%	-7.57%	9.378	87.56%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0039617	0.004355	0.16227	128.112	15.8484	137.6935	140.7632	1.35%	-37.31%	-9.22%	-7.20%	8.084	61.67%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.003823	0.004355	0.15659	127.9117	18.031	137.4798	140.5508	1.20%	-28.67%	-9.36%	-7.34%	7.094	41.88%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0036795	0.004355	0.15071	128.2718	19.2125	137.7077	140.7869	1.48%	-24.00%	-9.21%	-7.18%	6.676	33.53%	

Table B.107. ^{99m}I TRODAT 1:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	126.4	252.8					123-128			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055911	0.005665	0.22901	127.5656	147.0167	272.6883	282.8065	0.92%	16.31%	7.87%	11.87%	0.8677	-13.23%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055628	0.005665	0.22785	134.1034	154.48	285.7694	296.3729	6.09%	22.22%	13.04%	17.24%	0.8681	-13.19%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055301	0.005665	0.22651	135.7102	158.2825	289.9702	300.729	7.37%	25.22%	14.70%	18.96%	0.8574	-14.26%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.005493	0.005665	0.22499	135.9614	160.8757	291.7859	302.6112	7.56%	27.28%	15.42%	19.70%	0.8451	-15.49%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0053007	0.005665	0.21712	138.836	164.2017	294.7733	305.7094	9.84%	29.91%	16.60%	20.93%	0.8455	-15.45%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0051017	0.005665	0.20897	142.0657	164.2348	295.7007	306.6726	12.39%	29.93%	16.97%	21.31%	0.865	-13.50%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055911	0.005665	0.22901	108.2197	129.461	236.0383	241.4415	-14.38%	2.42%	-6.63%	-4.49%	0.8359	-16.41%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055628	0.005665	0.22785	113.7918	136.0392	247.3895	253.0526	-9.97%	7.63%	-2.14%	0.10%	0.8365	-16.35%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0055301	0.005665	0.22651	114.9406	139.5939	251.0427	256.789	-9.07%	10.44%	-0.70%	1.58%	0.8234	-17.66%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.005493	0.005665	0.22499	114.8956	142.1169	252.6239	258.4059	-9.10%	12.43%	-0.07%	2.22%	0.8085	-19.15%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0053007	0.005665	0.21712	117.3608	145.0922	255.2442	261.0863	-7.15%	14.79%	0.97%	3.28%	0.8089	-19.11%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0051017	0.005665	0.20897	120.5665	144.7751	256.0668	261.9285	-4.62%	14.54%	1.29%	3.61%	0.8328	-16.72%	

Table B.108. ^{99m}I TRODAT 2:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW		
SC-Total			126.4	63.2	189.6					123-128		
Atten Coeff	0.1526									126-154		
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047866	0.004863	0.19606	128.1156	70.3272	195.9827	203.2149	1.36%	11.28%	3.37%	7.18%	1.822	-8.91%
Atten Coeff	0.1526											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047624	0.004863	0.19507	132.8049	75.5576	205.1196	212.6939	5.07%	19.55%	8.19%	12.18%	1.758	-12.12%
Atten Coeff	0.1526											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047344	0.004863	0.19392	133.6368	78.4115	207.8665	215.5464	5.73%	24.07%	9.63%	13.68%	1.704	-14.78%
Atten Coeff	0.1526											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047027	0.004863	0.19262	133.2597	80.6731	208.9	216.6221	5.43%	27.65%	10.18%	14.25%	1.652	-17.41%
Atten Coeff	0.1526											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.004538	0.004863	0.18588	134.3442	83.4811	210.2042	217.978	6.28%	32.09%	10.87%	14.97%	1.609	-19.54%
Atten Coeff	0.1526											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0043677	0.004863	0.1789	136.1225	83.9035	210.6595	218.4494	7.69%	32.76%	11.11%	15.22%	1.622	-18.88%
Atten Coeff	0.12											
Striata ROI Size	13 13											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047866	0.004863	0.19606	117.381	54.3921	169.6299	173.4916	-7.14%	-13.94%	-10.53%	-8.50%	2.158	7.90%
Atten Coeff	0.12											
Striata ROI Size	15 15											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047624	0.004863	0.19507	121.3897	58.992	177.5594	181.6043	-3.96%	-6.66%	-6.35%	-4.22%	2.058	2.89%
Atten Coeff	0.12											
Striata ROI Size	17 17											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047344	0.004863	0.19392	121.8835	61.7034	179.9434	184.0509	-3.57%	-2.37%	-5.09%	-2.93%	1.975	-1.23%
Atten Coeff	0.12											
Striata ROI Size	19 19											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0047027	0.004863	0.19262	121.258	63.972	180.8521	184.9763	-4.07%	1.22%	-4.61%	-2.44%	1.895	-5.23%
Atten Coeff	0.12											
Striata ROI Size	26 26											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.004538	0.004863	0.18588	122.0024	66.6458	182.0005	186.1528	-3.48%	5.45%	-4.01%	-1.82%	1.831	-8.47%
Atten Coeff	0.12											
Striata ROI Size	32 32											
Bkgd ROI Size	13 13											
Units:	microCi											
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err
0.0043677	0.004863	0.1789	123.7105	66.8886	182.4076	186.5688	-2.13%	5.84%	-3.79%	-1.60%	1.850	-7.52%

Table B.109. ^{99m}I TRODAT 5:1 Activity Uptake Ratio, Scatter-Corrected Emission and Camera Sensitivity Data; Energy Windows Used: 123-128 keV, 126-154 keV

Data Type	whole		LStr (microCi)	RStr (microCi)	Sum (microCi)					TEW			
SC-Total			126.4	25.28	151.68					123-128			
Atten Coeff	0.1526									126-154			
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0041021	0.004423	0.16802	128.4764	24.349	147.9875	153.3051	1.64%	-3.68%	-2.43%	1.07%	5.277	5.53%	
Atten Coeff	0.1526												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040613	0.004423	0.16717	132.0692	28.2506	154.6367	160.2113	4.49%	11.75%	1.95%	5.62%	4.675	-6.50%	
Atten Coeff	0.1526												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040574	0.004423	0.16619	132.4475	30.5477	156.3478	161.9957	4.78%	20.84%	3.08%	6.80%	4.336	-13.28%	
Atten Coeff	0.1526												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040302	0.004423	0.16508	131.7135	32.62	156.8727	162.5508	4.20%	29.03%	3.42%	7.17%	4.038	-19.24%	
Atten Coeff	0.1526												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038891	0.004423	0.1593	131.8097	35.1771	157.1475	162.8482	4.26%	39.15%	3.60%	7.36%	3.747	-25.06%	
Atten Coeff	0.1526												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037431	0.004423	0.15332	132.8006	35.8852	157.3438	163.0542	5.06%	41.95%	3.73%	7.50%	3.701	-25.99%	
Atten Coeff	0.12												
Striata ROI Size	13 13												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0041021	0.004423	0.16802	122.9063	9.3815	128.0455	130.8834	-2.76%	-62.89%	-15.58%	-13.71%	13.101	162.02%	
Atten Coeff	0.12												
Striata ROI Size	15 15												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040813	0.004423	0.16717	125.9862	12.8041	133.8183	136.7939	-0.33%	-49.35%	-11.78%	-9.81%	9.840	96.79%	
Atten Coeff	0.12												
Striata ROI Size	17 17												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040574	0.004423	0.16619	126.097	15.0202	135.3097	138.3246	-0.24%	-40.58%	-10.79%	-8.80%	8.395	67.90%	
Atten Coeff	0.12												
Striata ROI Size	19 19												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0040302	0.004423	0.16508	125.1416	17.1438	135.7725	138.8038	-1.00%	-32.18%	-10.49%	-8.49%	7.300	45.99%	
Atten Coeff	0.12												
Striata ROI Size	26 26												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0038891	0.004423	0.1593	124.935	19.6837	136.0242	139.0679	-1.16%	-22.14%	-10.32%	-8.31%	6.347	26.94%	
Atten Coeff	0.12												
Striata ROI Size	32 32												
Bkgd ROI Size	13 13												
Units:	microCi												
Ave bkgd Act	bkgd STD	bkgd Act per Area	LStr Act	RStr Act	Sum w/ Cor	Sum w/o Cor	LStr %err	RStr %err	Sum w Corr %err	Sum w/o Corr %err	Ratio	Ratio %err	
0.0037431	0.004423	0.15332	125.8248	20.3021	136.2041	139.2533	-0.46%	-19.69%	-10.20%	-8.19%	6.198	23.95%	

REFERENCES

- Abella, H. (2003). SPECT and PET Compete in Diagnosing Parkinson's. Diagnostic Imaging: 43-47, 71.
- Abi-Dargham, A., J. Rodenhiser, et al. (2000). "Increased Baseline Occupancy of D₂ Receptors By Dopamine in Schizophrenia." Proceedings of the National Academy of Sciences **97**(14): 8104-8109.
- Albanese, A. (2003). "Diagnostic Criteria for Parkinson's Disease." Neurological Sciences **24**: S23-S26.
- Albin, R., A. Young, et al. (1989). "The Functional Anatomy of Basal Ganglia Disorders." Trends in Neurosciences **12**: 366-375.
- Andreasen, N. C. (1995). "Symptoms, Signs, and Diagnosis of Schizophrenia." The Lancet **346**(8973): 477-481.
- Andreasen, N. C. (1996). Schizophrenia. Retrieved 5/03/2004, from <http://gtel.gatech.edu:2176/universe/printdoc>
- Anger, H. O. (1958). "Scintillation Camera." The Review of Scientific Instruments **29**(1): 27-33.
- Barnes, W. E. (1996). In Vivo Quantitation of Activity By Planar Imaging. Nuclear Medicine. M. A. B. Robert E. Henkin, Gary, L. Dillehay, James R. Halama, Stephen M. Karesh, Robert H. Wagner, A. Michael Zimmer, Mosby-Year Book, Inc.: 216-221.
- Benamer, H. T. S., J. Patterson, et al. (2000). "Accurate Differentiation of Parkinsonism and Essential Tremor Using Visual ASsessment of [123-I]-FP-CIT SPECT Imaging: The [123-I]-FP-CIT Study Group." Movement Disorders **15**(3): 503-510.
- Booij, J., J. T. G. M. Hemelaar, et al. (1999). "One-Day Protocol for Imaging of the Nigrostriatal Dopaminergic Pathway in Parkinson's Disease by [123-I]FP-CIT SPECT." The Journal of Nuclear Medicine **40**(5): 753-761.
- Booij, J., E. Sokole, et al. (1998). "Human Biodistribution and Dosimetry of [123-I]FP-CIT: A Potent Radioligand for Imaging of Dopamine Transporters." The European Journal of Nuclear Medicine **25**(1): 24-30.

Booij, J., J. D. Speelman, et al. (2001). "The Clinical Benefit of Imaging Striatal Dopamine Transporters with [123-I]FP-CIT SPECT in Differentiating Patients with Presynaptic Parkinsonism from Those with Other Forms of Parkinsonism." The European Journal of Nuclear Medicine **28**(3): 266-272.

Booij, J., G. Tissingh, et al. (1999). "Imaging of the Dopaminergic Neurotransmission System Using Single-Photon Emission Tomography and Positron Emission Tomography in Patients with Parkinsonism." The European Journal of Nuclear Medicine **26**(2): 171-182.

Bourguignon, M., C. Pauwel, et al. (1997). "Iodine-123 Labeled Radiopharmaceuticals and Single Photon Emission Tomography: A Natural Liaison." The European Journal of Nuclear Medicine **24**: 331-344.

Bower, J. H., D. M. Maraganore, et al. (1999). "Incidence and Distribution of Parkinsonism in Olmsted County, Minnesota, 1976-1990." Neurology **52**: 1214.

Briesmeister, J. F. (2000). MCNP--A General Monte Carlo N-Particle Transport Code, Version 4. Los Alamos, NM, Los Alamos National Laboratories.

Burn, D. J. and J. T. O'Brien. (2003). "Use of Functional Imaging in Parkinsonism and Dementia." Movement Disorders **18**(Suppl 6): S88-S95.

Buvat, I., H. Benali, et al. (1994). "Scatter Correction in Scintigraphy: The State of the Art." The European Journal of Nuclear Medicine **21**(7): 675-694.

Buvat, I., M. Rodriguea-Villafuerte, et al. (1995). "Comparative Assessment of Nine Scatter Correction Methods Based on Spectral Analysis Using Monte Carlo Simulations." The Journal of Nuclear Medicine **36**(8): 1476-1488.

Buvat, I., M. Soret, et al. (2001). "Respective Importance of Scatter, Attenuation, Collimator Response and Partial Volume Effect Corrections for Accurate Quantification in 123I Dopamine Receptor Imaging." 13-15 - 13-19.

Cabezas, G., Z. Israili, et al. (2003). "The Actions of Dopamine on the Airways." American Journal of Therapy **10**(6): 477-486.

Cabezas, G., Y. Lezama, et al. (2001). "Dopaminergic Modulation of Human Bronchial Tone." Archives of Medical Research **32**: 143-147.

Carlsson, A. and M. Lindqvist. (1963). "Effect of Chlorpromazine or Haloperidol on Formation of 3-Methoxytyramine and Normetanephrine in Mouse Brain." Acta Pharmacol Toxicol **20**.

Catafau, A. M. (2003). SPECT Imaging of Cerebral Neurotransmission. Diagnostic Nuclear Medicine. R. E. C. Martin P. Sandler, James A. Patton, Frans J. Th. Wackers, Alexander Gottschalk, Lippincott Williams & Wilkins: 821-833.

Costa, D. (2001). Movement Disorders - A New Role for Nuclear Medicine and SPECT. 48th SNM Annual Meeting, Society of Nuclear Medicine.

Cote, L. and M. D. Crutcher (1991). The Basal Ganglia, Chapter 42. Principles in Neural Science. E. e. a. Kandel, Appleton and Lange: 647-659.

Creese, I., D. R. Burt, et al. (1976). "Dopamine Receptor Binding Predicts Clinical and Pharmacological Potencies of Antischizophrenic Drugs." Science **192**: 481-483.

Davis, K. L., R. S. Kahn, et al. (1991). "Dopamine in Schizophrenia: A Review and Reconceptualization." American Journal of Psychiatry **148**(11): 1474-1486.

Depatie, L. and S. Lal. (2001). "Apomorphine and the Dopamine Hypothesis of Schizophrenia: A Dilemma?" Journal of Psychiatry and Neuroscience **26**(3): 203-220.

Deplon, G., L. Ferrer, et al. (2003). "Impact of Scatter and Attenuation Corrections for Iodine-131 Two Dimensional Quantitative Imaging in Patients." Cancer Biotherapy and Radiopharmaceuticals **18**(2): 191-199.

Dobbeleir, A., A. Hambye, et al. (1999). "Influence of High-Energy Photons on the Spectrum of Iodine-123 with Low- and Medium-Energy Collimators: Consequences for Imaging with 123-I Labelled Compounds in Clinical Practice." The European Journal of Nuclear Medicine **26**: 655-658.

Eary, J. F., F. L. Appelbaum, et al. (1989). "Preliminary Validation of the Opposing View Method for Quantitative Gamma Camera Imaging." Medical Physics **16**(3): 382-387.

Eary, J. F., O. W. Press, et al. (1990). "Imaging and Treatment of B-Cell Lymphoma." Journal of Nuclear Medicine **31**(8): 1257-1268.

Eidelberg. (2003). Lecture 28 - Basal Ganglia. Retrieved June 2003, from <http://www.columbia.edu/itc/hs/medical/neuralsci/2003/lectures/basal%20ganglia.doc>

Evans, R. (1937). "Radium Poisoning: The Quantitative Determination of the Radium Content and Radium Elimination Rate of Living Persons." American Journal of Roentgenology and Radium Therapy **37**: 368-378.

Fahn, S. (2003). "Description of Parkinson's Disease as a Clinical Syndrome." Annals of the New York Academy of Sciences **991**: 1-14.

- Fang, P., C. Wu, et al. (2000). "The Preclinical Pharmacologic Study of Dopamine Transporter Imaging Agent [^{99m}Tc] TRODAT-1." Nuclear Medicine & Biology **27**: 69-75.
- Frederickson, A. (1998). The Dopamine Hypothesis of Schizophrenia. Retrieved 5/18/2004, from <http://serendip.brynmawr.edu/bb/neuro/neuro98/202s98-paper2/Frederickson2.html>
- Frey, K. A. (2002). "Can SPET Imaging of Dopamine Uptake Sites Replace PET Imaging in Parkinson's Disease." European Journal of Nuclear Medicine **29**: 711-714.
- Gilland, D. R., R. L. Jaszczak, et al. (1994). "Volume and Activity Quantitation with Iodine-123 SPECT." Journal of Nuclear Medicine **35**(10): 1707-1713.
- Grosset, D. G. (2001). Movement Disorders - A New Role for Nuclear Medicine and SPECT. 48th Annual Meeting of the Society of Nuclear Medicine, Society of Nuclear Medicine.
- Halldin, C., B. Gulyas, et al. (2001). "Brain Radioligands - State of the Art and New Trends." The Quarterly Journal of Nuclear Medicine **45**(2): 139-152.
- Hammond, N. D., P. J. Moldofsky, et al. (1984). "External Imaging Techniques for Quantitation of Distribution of I-131 F(ab')₂ Fragments of Monoclonal Antibody in Humans." Medical Physics **11**(6): 778-783.
- Harris, C. C., K. L. Greer, et al. (1984). "Tc-99m Attenuation Coefficients in Water-Filled Phantoms Determined With Gamma Cameras." Medical Physics **11**(5): 681-685.
- Harrison, P. J. (2000). "Dopamine and Schizophrenia - Proof At Last." The Lancet **356**: 958-959.
- Hashimoto, J., T. Sasaki, et al. (1999). "Effects of Scatter and Attenuation Correction on Quantitative Analysis of β -CIT Brain SPET." Nuclear Medicine Communications **20**: 159-165.
- Haynor, D. R., M. S. Kaplan, et al. (1995). "Multiwindow Scatter Correction Techniques in Single-Photon Imaging." Medical Physics **22**(12): 2015-2024.
- Hokama, H., M. E. Shenton, et al. (1995). "Caudate, Putamen, and Globus Pallidus Volume in Schizophrenia: A Quantitative MRI Study." Psychiatry Research: Neuroimaging **61**(4): 209-229.
- Huang, W.-S., S.-Z. Lin, et al. (2001). "Evaluation of Early-Stage Parkinson's Disease with ^{99m}Tc-TRODAT-1 Imaging." The Journal of Nuclear Medicine **42**(9): 1303-1308.

Ichihara, T., K. Ogawa, et al. (1993). "Compton Scatter Compensation Using the Triple-Energy Window for Single- and Dual- Isotope SPECT." The Journal of Nuclear Medicine **34**(12): 2216-2221.

Ichise, M. and J. R. Ballinger. (1996). "SPECT Imaging of Dopamine Receptors." Journal of Nuclear Medicine **37**(10): 1591-1595.

Jangha, D. N., R. A. Mintzer, et al. (2001). MCNP-4c Simulations to Investigate Shielding of Scintillation Cameras for Conjugate Imaging of I-123 Labeled Brain Agents. IEEE Nuclear Science Symposium and Imaging Conference, San Diego CA.

Jaszczak, R. L., K. L. Greer, et al. (1984). "Improved SPECT Quantification Using Compensation for Scattered Photons." Journal of Nuclear Medicine **25**(8): 893-900.

Jersild, S. and M. Trevino. (2003). "Agent Improves Diagnosis of Parkinson's." Diagnostic Imaging **25**(2): 19-20.

Karesh, S. M. (1996). Principles of Radiopharmacy. Nuclear Medicine. M. A. B. Robert E. Henkin, Gary, L. Dillehay, James R. Halama, Stephen M. Karesh, Robert H. Wagner, A. Michael Zimmer, Mosby-Year Book Inc. **1**: 334-349.

Kim, K. M., A. Varrone, et al. (2003). "Contribution of Scatter and Attenuation Compensation to SPECT Images of Nonuniformly Distributed Brain Activities." The Journal of Nuclear Medicine **44**(4): 512-519.

King, M. and T. Farncombe. (2003). "An Overview of Attenuation and Scatter Correction of Planar and SPECT Data for Dosimetry Studies." Cancer Biotherapy and Radiopharmaceuticals **18**(2): 181-190.

King, M. A., D. T. Long, et al. (1991). "SPECT Volume Quantitation: Influence of Spatial Resolution, Source Size and Shape, and Voxel Size." Medical Physics **18**(5): 1016-1024.

Knable, M. B., D. W. Jones, et al. (1995). "Lateralized Differences in the Iodine-123-IBZM Uptake in the Basal Ganglia in Asymmetric Parkinson's Disease." The Journal of Nuclear Medicine **36**(7): 1216-1225.

Kojima, A., Y. Ohyama, et al. (2000). "Quantitative Planar Imaging Method for Measurement of Renal Activity by Using a Conjugate-Emission Image and Transmission Data." Medical Physics **27**(3): 608-615.

Koral, K. F., R. S. Adler, et al. (1982). "Two-orthogonal-view Method for Quantification of Rad Doase to Neck Lesions in Thyroid Cancer Patients." Medical Physics **9**(4): 497-503.

Kung, H., A. Alvai, et al. (1990). "In Vivo SPECT Imaging of CNS D-2 Dopamine Receptors: Initial Studies with Iodine-123 IBZM in Humans." The Journal of Nuclear Medicine **31**: 573-579.

Kung, H., H. J. Kim, et al. (1996). "Imaging of Dopamine Transporters in Humans with Technetium-99m TRODAT-1." The European Journal of Nuclear Medicine **31**: 1527-1530.

Kung, H. F., M. P. Kung, et al. (2003). "Radiopharmaceuticals for Single-Photon Emission Computed Tomography Brain Imaging." Seminars in Nuclear Medicine **33**(1): 2-13.

Laruelle, M. (2000). "The Role of Model-Based Methods in the Development of Single Scan Techniques." Nuclear Medicine and Biology **27**: 637-642.

Leenders, K. L. (2003). "Significance of Non-Presynaptic SPECT Tracer Methods in Parkinson's Disease." Movement Disorders **18**(Suppl. 7): S39-S42.

Leong, L. K., M. K. O' Connor, et al. (1999). "Quantitation of Iodine-123-B-CIT Dopamine Receptor Uptake in Phantom Model." The Journal of Nuclear Medicine Technology **27**(2): 117-122.

Leslie, W. D., D. N. Abrams, et al. (1996). "Comparison of Iodine-123-Epidopride and Iodine-123-IBZM for Dopamine D2 Receptor Imaging." Journal of Nuclear Medicine **37**(10): 1589-1591.

Lomena, F., A. M. Catafau, et al. (2004). "Striatal Dopamine D₂ Receptor Density in Neuroleptic-Naive and in Neuroleptic-Free Schizophrenic Patients: an ¹²³I-IBZM-SPECT Study." Psychopharmacology **172**: 165-169.

Mann, R. (1996). The Role of Dopamine Receptors in Schizophrenia. Retrieved from <http://www.chem.csustan.edu/chem44x0/SJBR/Mann.htm>

Marek, K. L., J. P. Seibly, et al. (1996). "[¹²³I]beta-CIT/SPECT Imaging Demonstrates Bilateral Loss of Dopamine Transporters in Hemi-Parkinson's Disease." Neurology **46**: 231-237.

Matthews, K., R. Mintzer, et al. (1997). Investigation of Conjugate Imaging Using Small Gamma Cameras. IEEE NSS/MIC, IEEE.

Mayo Clinic Staff. (2004). Parkinson's Disease. Retrieved June 2003, from Mayo Clinic Website: <http://www.mayoclinic.com/invoke.cfm?objectid=67352253-4C13-41A1-A6F790554DDA8474>

Meegalla, S. K., K. Plossl, et al. (1997). "Synthesis and Characterization of Technitium-99m-Labeled Tropanes as Dopamine Transporter-Imaging Agents." The Journal of Medicinal Chemistry **40**(1): 9-17.

Mental Health Channel. (2004). Schizophrenia. Retrieved February 16, 2002, from Mental Health Channel website:<http://www.mentalhealthchannel.net/schizophrenia/course.shtml>

Messa, C., M. A. Volonte, et al. (1998). "Differential Distribution of Striatal [I-123]B-CIT in Parkinson's Disease and Progressive Supranuclear Palsy, Evaluated with Single-Photon Emission Technology." The European Journal of Nuclear Medicine **25**(9): 1270-1276.

Mezey, E., G. Eisenhofer, et al. (1996). "A Novel Nonneuronal Catecholaminergic System: Exocrine Pancreas Synthesizes and Releases Dopamine." Proceedings of the National Academy of Sciences **93**: 10377-10382.

Mink, J. (1996). "The Basal Ganglia: Focused Selection and Inhibition of Competing Motor Programs." Progress in Neurobiology **50**: 381-425.

Mintzer, R., J. N. Aarsvold, et al. (2000). "Conjugate Imaging for Regional-Uptake Estimation: Some Analytical Considerations." The Journal of Nuclear Medicine **41**(5): 190P.

Moore, S. C., M. F. Kijewski, et al. (2001). "Evaluation of Scatter Compensation Methods By Their Effects On Parameter Estimation From SPECT Projections." Medical Physics **28**(2): 278-287.

Morrish, P. K. (2003). "How Valid Is Dopamine Transporter Imaging as a Surrogate Marker in Research Trials in Parkinson's Disease." Movement Disorders **18**(Suppl. 7): S63-S70.

Mozley, P. D., J. S. Schneider, et al. (2000). "Binding of [99m-Tc]TRODAT-1 to Dopamine Transporters in Patients with Parkinson's Disease and in Healthy Volunteers." The Journal of Nuclear Medicine **41**: 584-589.

Mozley, P. D., J. B. Stubbs, et al. (1998). "Biodistribution and Dosimetry of TRODAT-1: A Technetium-99m Tropane for Imaging Dopamine Transporters." The Journal of Nuclear Medicine **39**(12): 2069-2076.

Narita, Y., S. Eberl, et al. (1996). "Monte Carlo and Experimental Evaluation of Accuracy and Noise Properties of Two Scatter Correction Methods for SPECT." Physics in Medicine and Biology **41**: 2481-2496.

O'Donnell, P. and A. A. Grace. (1998). "Dysfunctions in Multiple Interrelated Systems as the Neurobiological Bases of Schizophrenic Symptoms Clusters." Schizophrenia Bulletin **24**: 267-283.

Oertel, W. H., A. Gerstner, et al. (2003). "Role of Dopamine Transporter SPECT for the Practitioner and the General Neurologist." Movement Disorders **18**(Suppl. 7): S9-S15.

Ogawa, K., Y. Harata, et al. (1991). "A Practical Method for Position-Dependent Compton Scatter Compensation in Single Photon Emission Computed Tomography." IEEE Transactions on Medical Imaging **10**: 408-412.

Okada, J., Y. Imal, et al. (1999). "⁶⁷Ga Planar Imaging With a Low-Energy Collimator and Scatter Correction Using the Triple Energy Window Method." Nuclear Medicine Communications **20**: 647-657.

Pestana, M., M. Faria, et al. (1997). "Assessment of Renal Dopaminergic System Activity During the Recovery of Renal Function in Human Kidney Transplant Recipients." Nephrology Dialysis Transplantation **12**: 2667-2672.

Piccini, P. P. (2003). "Dopamine Transporter: Basic Aspects and Neuroimaging." Movement Disorders **18**(Suppl. 7): S3-S8.

Pinborg, L., C. Videbaek, et al. (2002). "Quantification of [123-I]PE2I Binding to Dopamine Transporters with SPECT." The European Journal of Nuclear Medicine **29**(5): 623-631.

Pirker, W. (2003). "Correlation of Dopamine Transporter Imaging With Parkinsonian Motor Handicap: How Close Is It?" Movement Disorders **18**(Suppl. 7): S43-S51.

Poewe, W. and G. Wenning. (2002). "The Differential Diagnosis of Parkinson's Disease." European Journal of Neurology **9**(Suppl. 3): 23-30.

Pollard, K. R., T. K. Lewellen, et al. (1996). "Energy-Based Scatter Corrections for Scintillation Camera Images of Iodine-131." The Journal of Nuclear Medicine **37**(12): 2030-2037.

Quinn, N. (1994). "Multiple System Atrophy." Movement Disorders **3**: 262-281.

Rajput, A. H., B. Rozdilsky, et al. (1991). "Accuracy of Clinical Diagnosis in Parkinsonism - A Prospective Study." Canadian Journal of Neurological Sciences **18**(3): 275-278.

Rosenthal, M. S., J. Cullom, et al. (1995). "Quantitative SPECT Imaging: A Review and Recommendations by the Focus Committee of the Society of Nuclear Medicine Computer and Instrumentation Council." Journal of Nuclear Medicine **36**(8): 1489-1513.

Schrag, A., Y. Ben-Shlomo, et al. (2002). "How Valid is the Clinical Diagnosis of Parkinson's Disease in the Community?" Journal of Neurology, Neurosurgery, and Psychiatry **73**: 529-534.

Schwarz, J., W. H. Oertel, et al. (1996). "Iodine-123-Iodobenzamide Binding in Parkinsonism: Reduction by Dopamine Agonists but Not L-Dopa." Journal of Nuclear Medicine **37**(7): 1112-1115.

Sedvall, G. and L. Farde. (1995). "Chemical Brain Anatomy in Schizophrenia." The Lancet **346**: 743-749.

Seibyl, J. P. (2003). "Imaging Studies in Movement Disorders." Seminars in Nuclear Medicine **33**(2): 105-113.

Seibyl, J. P., K. Marek, et al. (1998). "Iodine-123-B-CIT and Iodine-123-FPCIT SPECT Measurement of Dopamine Transporters in Healthy Subjects and Parkinson's Patients." The Journal of Nuclear Medicine **39**(9): 1500-1508.

Seibyl, J. P., E. Wallace, et al. (1994). "Whole Body Distribution Radiation Absorbed Dose and Brain SPECT with Iodine-123-B-CIT in Healthy Human Subjects." The Journal of Nuclear Medicine **35**(7): 764-770.

Seibyl, J. P., S. W. Woods, et al. (1992). "Dynamic SPECT Imaging of Dopamine D2 Receptors in Human Subjects with Iodine-123 IBZM." The Journal of Nuclear Medicine **33**(11): 1964-1971.

Shulkin, B., J. Sisson, et al. (1988). "Conjugate View Gamma Camera Method for Estimating Tumor Uptake of Iodine-131 Metaiodobenzylguanidine." The Journal of Nuclear Medicine **29**(4).

Shung, K. K. (1992). Principles of Medical Imaging. New York, Academic Press, Inc.

Sjogreen, K., M. Ljungberg, et al. (2001). "Registration of Emission and Transmission Whole-Body Scintillation-Camera Images." Journal of Nuclear Medicine **42**(10): 1563-1570.

Sorenson, J. (1974). Quantitative Measurement of Radioactivity in vivo by Whole Body Counting. Instrumentation in Nuclear Medicine. F. J. Hine and J. A. Sorenson. New York, Academic Press. **II**.

Soret, M., P. M. Koulibaly, et al. (2003). "Quantitative Accuracy of Dopaminergic Neurotransmission Imaging with ¹²³I SPECT." Journal of Nuclear Medicine **44**(7): 1184-1193.

Stacy, M. and M. Jankovic. (1992). "Differential diagnosis of Parkinson's disease and the parkinsonism plus syndromes." Neurologic clinics **10**(2): 341-359.

Strange, P. (2001). "Antipsychotic Drugs: Importance of Dopamine Receptors for Mechanisms of Therapeutic Actions and Side Effects." Pharmacological Review **53**(1): 119-133.

Tatsch, K. and K. Frey. (2002). "Controversies: Can SPECT Imaging of Dopamine Uptake Sites Replace PET Imaging in Parkinson's Disease." The European Journal of Nuclear Medicine **29**(5): 711-714.

Tatsch, K. C., S. Asenbaum, et al. (2002). "European Association of Nuclear Medicine Procedure Guidelines for Brain Neurotransmission SPECT Using 123-I-Labelled Dopamine D2 Receptor Ligands." The European Journal of Nuclear Medicine **29**(10): BP23-BP29.

Tatsch, K. C., S. Asenbaum, et al. (2002). "European Association of Nuclear Medicine Procedure Guidelines for Brain Neurotransmission SPECT Using 123-I-Labelled Dopamine Transporter Ligands." The European Journal of Nuclear Medicine **29**(10): BP30-BP35.

Tatsch, K. M. D. (2001). Investigations of the Dopaminergic System In-Vivo: Available Radioligands for SPECT and Clinical Applications. The 48th SNM Annual Meeting 2001 Handout Book Tuesday and Wednesday Courses, SNM.

Testa, D., D. Monza, et al. (2001). "Comparison of Natural Histories of Progressive Supranuclear Palsy and Multiple System Atrophy." Neurological Sciences **22**: 247-251.

The HOPES Brain Tutorial. (2003). Basal Ganglia. Retrieved February 2006, from HOPES: The Huntington's Outreach Program for Education, at Stanford: <http://www.stanford.edu/group/hopes/basics/braintut/ab6.html>

Thomas, S., H. Maxon, et al. (1976). "In Vivo Quantitation of Lesion Radioactivity Using External Counting Methods." Medical Physics **3**: 253-255.

Tissingh, G., J. Booij, et al. (1998). "Iodine-123-N-w-Fluoropropyl-2-B-Carbomethoxy-3B-(4-Iodophenyl) Topane SPECT in Healthy Controls and Early-Stage, Drug-Naive Parkinson's Disease." The Journal of Nuclear Medicine **39**(7): 1143-1148.

Tonini, M., L. Cipollina, et al. (2004). "Review Article: Clinical Implications of Enteric and Central D2 Receptor Blockade By Antidopaminergic Gastrointestinal Prokinetics." Ailment Pharmacology Terapeutics **19**(4): 379-390.

Toru, M., A. Kurumaji, et al. (1994). "Excitatory amino acids: implications for psychiatric disorders research." Life Sciences **55**(22): 1683-1699.

Tsui, B. M. W. (1996). Quantitative SPECT. Nuclear Medicine. M. A. B. Robert E. Henkin, Gary, L. Dillehay, James R. Halama, Stephen M. Karesh, Robert H. Wagner, A. Michael Zimmer, Mosby-Year Book Inc.: 260-278.

Tsui, B. M. W., C.-T. Chen, et al. (1981). A Whole Body Scanning System for Collection of Quantitative In Vivo Distribution Data in Humans. Third International Radiopharmaceutical Dosimetry Symposium, HHS Publication FDA 81-8166.

Tzen, K.-Y., L. Chin-Song, et al. (2001). "Differential Diagnosis of Parkinson's Disease and Vascular Parkinsonism by 99mTc-TRODAT-1." The Journal of Nuclear Medicine **42**(3): 408-413.

Uhl, G. R. (2003). "Dopamine Transporter: Basic Science and Human Variation of a Key Molecule for Dopaminergic." Movement Disorders **18**(Suppl. 7): S71-S80.

Van_Laere, K., M. Koole, et al. (2000). "Nonuniform Transmission in Brain SPECT Using 201-Tl, 153-Gd, and 99m-Tc Static Line Sources: Anthropomorphic Dosimetry Studies and Influence on Brain Quantification." The Journal of Nuclear Medicine **41**: 2051-2062.

Van_Laere, K., M. Koole, et al. (2001). "Non-Uniform Versus Uniform Attenuation Coefficient in Brain Perfusion SPECT of Healthy Volunteers." The European Journal of Nuclear Medicine **28**: 90-98.

Vandenberghe, S., Y. D'Asseler, et al. (2001). "Iterative Reconstruction Algorithms in Nuclear Medicine." Computerized Medical Imaging and Graphics **25**: 105-111.

Vines, D. C., M. Ichise, et al. (2003). "Evaluation of 2 Scatter Correction Methods Using a Striatal Phantom for Quantitative Brain SPECT." Journal of Nuclear Medicine Technology **31**(3): 157-160.

Weinberger, D. R., K. F. Berman, et al. (1988). "Physiological Dysfunction of Dorsolateral Prefrontal Cortex in Schizophrenia: III A New Cohort and Evidence for a Monoaminergic Mechanism." Archives of General Psychiatry **45**: 609-615.

Wu, R. K. and J. A. Siegel. (1984). "Absolute Quantitation of Radioactivity Using the Buildup Factor." Medical Physics **11**(2): 189-192.

Yen, C.-K., A. D. Lim, et al. (1994). "A Count-Based Radionuclide Method for Volume Quantitation Using Conjugate Imaging and an External Reference Source: Theoretical Consideration and Phantom Study." Journal of Nuclear Medicine **35**(4): 644-651.

Zaidi, H. (1996). "Quantitative SPECT: Recent Developments in Detector Response, Attenuation, and Scatter Correction Techniques." Physica Medica **39**: 1719-1726.

Zaidi, H. and M.-L. Montandon. (2002). "Which Attenuation Coefficient to Use in Combined Attenuation and Scatter Correction for Quantitative Brain SPECT?" The European Journal of Nuclear Medicine and Molecular Imaging.

Zeng, G. L. (2001). "Image Reconstruction - A Tutorial." Computerized Medical Imaging and Graphics **25**: 97-103.

Zeng, L. (2003). Personal Communication in form of figures. University of Utah, Salt Lake City, Utah.