

Upper Extremity Robot Therapy for Individuals with Spinal Cord Injury

Bui KT, Sprigle S, Backus D

April 25, 2007

Background : Spinal Cord Injury (SCI)

- SCI annual incidence in the U.S.
 - 11,000 new cases per year
- SCI prevalence in the U.S. 2006
 - 253,000 people
- Statistics
 - Average age at injury: 28 yrs
 - Gender: 77.8% males
 - Etiology: 46.9% from motor vehicle accident
- Impairment or loss of motor/sensory function
 - Tetraplegia – impaired function of arms, trunk, legs, pelvic organs
 - Paraplegia – impaired function of trunk, legs, pelvic organs

Background: Recovery After SCI

- Functional recovery of some muscles can occur for several years after incomplete SCI
 - Recovery depends on the reorganization of preserved nerve connections and formation of new connections [Beekhuizen 2005]
- Cortical reorganization
 - Reorganization is associated with activity [Beekhuizen 2005, Fasioli 2003, Raintenteau 2001]
 - Cortical reorganization occurs after SCI as it does after stroke [Beekhuizen 2005]
- SCI recovery mechanism is similar to that in stroke subjects
 - Interventions used after stroke have demonstrated potential to improve strength after SCI [Beekhuizen 2005]

Background: Robotic therapy

- Current use: rehabilitation for motor impairment and limb weakness after stroke [Fasioli 2003]
- Therapy modes: passive & active-assisted, active-resisted [Stein 2004]
- Results
 - improved motor coordination and strength in early stroke subjects
 - reduced motor impairment in moderate to chronic stroke subjects [Fasioli 2003]
- Advantages: precise repetitive movements, quantitative performance monitoring, staffing/budget, compliments conventional therapy [Dikers 1991, Colombo 2005, Reinkensmeyer 2000]
- Disadvantages: generic treatment protocols, unable to identify specific benefits of each therapy modality [Stein 2004]

Robotic Device: Muscle Tech

■ Background


- ❑ Designed by Muscle Tech company in Israel

■ Functions

- ❑ Elbow flexion/extension
- ❑ Forearm pronation/supination

■ Setup

- ❑ Device stands independently to side of wheelchair
- ❑ Controlled through computer program



Objectives

■ Phase 1

- **Assess usability:** safety, comfort, ease of use, fit and exercise protocol in order to optimize the device


■ Phase 2

- Investigate effectiveness of robotic exercises for **strengthening** after SCI
-

Inclusion Criteria

Phase 1

- SCI level C5, C6, or C7
- ASIA impairment B, C, or D
- Male or female
- 18 - 59 years old
- Tolerate sitting upright at 90 degrees for at least one hour
- Demonstrate at least 50% of normal range of motion (ROM) for each joint involved in this therapy
- Palpable contraction (grade 1-3) in the biceps and triceps
- Chronic impairment





Study Protocol



■ Phase 1: 3 hr session

- Pre-treatment
 - Clinical Assessment
 - Robotic Evaluation
- Treatment
 - 1 hr robot exercises
- Post-study Questionnaire
 - Subjects



- Pre-treatment
- Robotic Exercises
- Post-study Questionnaire
- Modalities: Positional treatment
 - Design optimization
 - Passive assisted movement
 - Clinical Assessment
 - Comfort in use, motion, muscle
 - Robotic Evaluation
 - Extending
 - Treatment
 - Isometric contraction
 - 1 hr robot exercises
 - Therapist usefulness
 - Strengthening
 - 3 days/wk, 4 wks
 - Active resisted movement
 - Robotic Evaluation
 - Post-study Questionnaire
 - range of motion: active
 - Subjects
 - passive
 - Clinicians
 - muscle strength: kgF




Results – Phase 1: Pre-treatment

- Subjects

- ❑ 3 males, 2 females
- ❑ Ages: 19 – 46 yrs
- ❑ Injury year: 1980-2006

- IRB approval

- ❑ Georgia Institute of Technology
- ❑ Shepherd Center
- ❑ Emory University



- PROM average


Treatment Arm	Clinical	Robot
Elbow Flexion	135°	110°
Elbow Extension	4°	27°

- Device limitations

- ❑ Flexion limit: 110 degrees
- ❑ Extension: subject selected end point

Results – Phase 1: Treatment

- Subject #3
 - ❑ Passive-assisted
 - ❑ Isometric contraction
 - ❑ Active-resisted



Results – Phase 1: Questionnaire

- Safety (9/10)
 - ❑ “felt safe with padding and splint”
 - ❑ “didn’t feel it squeeze”
 - ❑ “ROM limits were good”
- Comfort (8.6/10)
 - ❑ “snug fit”
 - ❑ “exercising is uncomfortable, not the machine b/c it was padded well”
 - ❑ “arm cuff not wide enough for larger arm”
 - ❑ “pressure at wrist but with breaks between exercises it was not an issue”
- Design ideas
 - ❑ Flexible plastic arm section to expand and conform
 - ❑ Modify splint to redistribute pressure at wrist



Results – Phase 1: Questionnaire

- Ease of Use
 - Don (4.2/10)
 - “tetraplegics need assistance or initial set up”
 - “rings on velcro straps”
 - Use (7.6/10)
 - “need instructions for exercises”
 - Doff (6.4/10)
 - “rings on velcro straps”
 - “no buttons”
- Fit
 - Of arm in device (8.4/10)
 - “straps and pads helped the fit”
 - “one-size arm section was tight on larger arm”
 - Security during exercise (9.2/10)
- Design ideas
 - Rings on velcro straps



Results – Phase 1:

Questionnaire

■ Therapy

- Challenging (8.4/10)

- Motivation (7.4/10) "exercises were challenging but not biceps"

- "colored graphs made me try harder on the next set"



- Strengthening (7.6/10) "how many more repetitions?"

- "signals for stop/go not always clear"
- "enough exercise but not exhausting"

- "may be good for recently injured person"

- Design Ideas:

- "maintaining strength is more important than building up"
- "color code all graphs, repetition counter, stop/go signal on computer screen"
- "don't know how much I can build up the muscles I still use"
- "greater range of resistance"




Discussion



■ Future Design Considerations

- Robot exercises vs. hand weights
 - Safety – risk of stretching muscle
 - Independence
 - Measures of improvement
 - Interest in technology
- Home device vs. therapist tool
 - Level of assistance to set-up and operate
 - Portable or permanently fixed



Acknowledgements

- Advisors: Dr. Sprigle, Dr. Backus & Dr. Kong
- Kevin Grogg
- Therapists: Beth Pharo, Amy Brake, Teresa Foy
- Jon Jowers
- Study participants

