

# Power Wheelchairs 101



# Objectives

- Identify components of a power obility evaluation
- Recognize the functional benefits/limitations of power seating systems
- Identify specialty controls and discuss potential benefits and limitations of each.
- Recognize risks and discuss power mobility safety.
- Discuss driver training programs
- Discuss power mobility options

# Medicare Hierarchy of Mobility Aids

1. Cane or crutches
2. Walker
3. Rollator
4. Manual Wheelchair
5. Power assist
6. POV (Scooter)
7. Power wheelchair



# Power Wheelchairs

- POV/scooter
- Transportable power
- Basic power (non-adjustable)
- Power assist
- Power Base
- Off Road



# Manual vs. Power

- Client Considerations

- Function
- Cognitive Status
- Safety/Seizures
- Vision
- Pain
- Strength/Endurance
- Sensation

- Equipment Considerations

- Environment
- Range/speed
- Size
- Accessibility
- Transportation
- Cost/Funding



# Functional Status

- Upper Extremity Involvement/
- Strength / Symmetry
- Scapular Stability
- Endurance
- Pre-existing impairments
- Orthopedic limitations
- Age
- Mobility related activities of daily living

# Cognitive Status

- Understands cause and effect
- Free of seizures
- Psychological Stability

# Vision

- Vertical sitting tolerance
- Visual acuity
- Restrictions



# Transportation / Accessibility

- Safe environment?
- Does the person have the means to transport a power wheelchair/scooter?
- If not how do they plan to get around the community?
- How will they secure the wheelchair in the vehicle?
- Where do they plan to ride in the vehicle?
- Can they transfer independently?
- What is the overall width and length of the wheelchair?
- What is the turning radius?
- Knee and head clearance?
- Foot rest width?

# Funding of Power

- Prior authorization
- Medicare-must be needed for in the home mobility.
- Medicare just released new reimbursements and groupings of power wheelchairs.
- Assigned versus non-assigned
- Discuss the percentage that the individual will have to pay?

# Sequence of Evaluation

- Postural support and seating needs?
  - Will the person require power seat functions?
- Selection of Input Device
  - What is the most reliable and functional access point?
- Essential Controller Functions
  - What programming functions are needed for safe operation of the wheelchair (if any)
- Selection of Powered Mobility Base
  - Where does the person plan to use the chair and what performance characteristics are needed?

# Seating System

- Symmetrical sitting posture?
- Level of Balance?
- Vertical Tolerance?
- Stability for operation of Input Device?
- Need for Positional Changes?

# Interfacing Powered Seating Systems

- Tilt
- Recline
- Standing
- Seat Elevation
- Standing
- Power Elevating Leg Rest

# Medical Considerations

- Pressure Distribution
- Gravity Assisted positioning
- Accommodation of Contractures
- Spasticity
- Compensate for Postural Hypotension
- Intermittent Catheterization

# Electronics/Input Devices

- Proportional
- Non-proportional



# Questions When Choosing an Electronics System cont.

- Can I plug more than one input device into the system?
- Advanced programming, is it global or profile specific?
- Does the system allow me to grow with the user as they progress or digress?
- Is it an aesthetically pleasing system? Less wires and boxes.
- Does it use current technology and lend itself to easy interface with current specialty controls in the future?
- Know your options.

# Choosing the Appropriate Electronics Operating System



*This could be a day in itself!!*

# Basic Joystick Considerations

## Hand Function

- Joystick handles
- Toggles and Switch Jacks
- Ability to interface with joystick

## Lights and Actuators

## EADL and AAC Control

- Mouse Emulation



# Proportional Controls



Speed modulation is proportional to the input given to the control system. Like a gas pedal in a car.

- Chin Controls
- Mini-Joysticks
- Proportional Head Arrays (Rim)
- Touch pads
- Finger Steering
- Magitek



# Non-Proportional Controls

Even Better Than Proportional?

A Switch is Either Open or Closed

- No Matter What Type of Switch

- Proximity, fiber optic, infra red and mechanical switch systems
- Head Arrays
- Sip n Puff
- Single Switch Scanners
- Sip n Puff Head Arrays



A switch is a switch is a switch is a switch.....

# Positioning of Input Devices

- Can the person get to the on/off switch?
- Can person access the input device independently and move in all planes of motion required to operate the device?
- Can person maintain access throughout the range of the power seat functions?

# Basic Controller Functions

- Preprogrammed functions
- forward speed
- turn speed
- reverse speed
- acceleration
- deceleration
- turn acceleration
- turn deceleration
- Sensitivity
- Power level



# Additional Controller Functions

- momentary/latched
- Tremor dampening
- RIM control
- Switch direction swap
- Standby
- Reset/mode selection
- Joystick Throw
- Drive Selection

# Powered Mobility Technology

# Basic Considerations

**Growth** - Will the base and seating platform grow with the needs, present and future, of the user?

**Performance** - Will the base I choose accommodate the seating without sacrificing the performance of the chair?

**Modularity** - Can I add to my seating as needs change, can my seating adapt to changes made to the base or power positioning system?

**Accessories** - Will the seating and seating platform accept current and future accessories?

# Basic Considerations cont.

**Repair** - If my base goes down, will it be easy to move my seating to a temporary base?

**Ease of Adjustment** - How easy is it to make basic adjustments to COG, depth, STF, seat angle and width? Will I need to justify additional parts?

**Durability** - Will the seating platform and seating hold-up to the rigors of the individual user(s)?

**Caregiver** - What can I do to the seating platform and interface to make it easier on the caregiver?

**Aesthetics** - Can I match colors? Hide wires? Use minimal amounts of hardware?

# Rear Wheel Drive

- directional stability (tracks in a straight line)
- Largest turning radius
- Pivot point behind the user
- Excellent traction up inclines, poor down inclines
- Fastest speeds

# Front Wheel Drive

- Slightly slower
- May fish tail at higher speeds
- Good traction down hill
- Best climbing ability
- Smoothest on rough terrain because of long wheel base.
- Shorter turning radius than rear wheel drive

# Midwheel Drive

- Smallest turning radius
- Slower than rear wheel drive for safe control
- Traction about equal for up and down inclines
- Pivot point under the driver
- Additional casters for stability.
- May experience a more significant rocking motion over obstacles