

# Cushion use and performance in everyday life

Stephen Sprigle







- Surveying used cushions
- Documenting degradation
- Temperature and humidity
  - Controlled tests
  - Within everyday use


# Surveying used cushions

- Survey developed to document cushion status
  - Descriptions of cushion and cover;
  - Reasons for replacement
- Sent to Robert Bingham in Australia
- 209 surveys completed
- Flat and contoured foam

# Age of foam cushions




Cleanliness of foam




Percent within levels of agecode2.

Cover cleanliness




Panel variable: agecode2

# Physical and clinical signs of fatigue



After 6 months, 70% showed physical signs of fatigue



40% of cushions 7-12 mo of age showed clinical sign of fatigue

>12 months, >60% showed a clinical sign

# Why we should care

- Nice size data set on foam cushions
- Insight into a different delivery model
- Compression set occurs before clinical indicators of fatigue
  - High pressures or discomfort noted for 33% of cushions  $\leq 12$  mo
  - Compression set noted in 70% of cushions  $\leq 12$  mo
- Foam is in pretty good shape after 12 months of use
- Certain temporary wheelchair users may benefit from a foam cushion
  - i.e., stroke survivors are often d/c'd with orders for only a wheelchair; minimal cost might meet needs

# Documenting degradation

a collaborative project between



SHEPHERD CENTER  
A Catastrophic Care Hospital

- Objectives
  - Identify the expected lifespan of cushions and the significant predictors of cushion failure
  - Develop and validate a clinical measure of seat cushion degradation
- 138 different cushions studied (24 measured >1x)
  - Most common: 32 Jay2; 26 Roho HP; 14 Evolution
- Client eval, visual inspection & performance measures
- Mean age: 24 months (range: 1 day to 168 months)



# Testing cushions over time


Interview & physical exam



IPM with user




IPM using model



Visual inspection & dimensioning






Loaded contour depth



Impact dampening




# Model and Human IPM



Metrics cover:  
 magnitude  
 asymmetry  
 dispersion


### Dispersion Index vs cushion age



DI= ratio of IT pressures to total pressure



PPI= measure of pressure magnitude

### Model PPI vs cushion age




# Pressure magnitudes- ALL 162 cushions

Both model and subject pressures  
indicate NO relationship over time



# A tale of 2 predictors




← Cushion age

Subject weight



PPI<sub>asym</sub>=  
asymmetry of R & L peak values


normPPI= (normalized PPI)  
PPI divided by total pressure



# In fact....

- **Cushion age** has *not* been able to predict any IPM-related variable
  - For all cushions
  - Combining the 3 most tested cushions (Roho, Jay2, Evolution)
- Cushion age may predict certain performance in *Evolution*
  - *May be indicative of foam*

# Looking only at FOAM-based cushions



Pressure magnitudes tend to rise over time  
Huge variance in model testing

# Why we should care

- Tracking performance changes over time is needed to better understand “useful life”
- Extensive data on 138 cushions (and 162 measurements) is overwhelming
- Evidence suggests that Roho and Jay 2 cushion performance appears independent of age
  - For the cohort studied



# Temperature and humidity

- Humidity represents moisture
- Temperature represents temperature

# Friction and Moisture

- As moisture increases, friction increases
  - $\uparrow$  softness  $\rightarrow$   $\uparrow$  contact between surfaces
  - Want to learn more?- see cosmetics literature
- Excessive moisture weakens skin's ability to withstand load

# Temperature and it's impact on tissue viability


- ↑ tissue temperature ↑ metabolic demand
  - Added demand coupled with reduced nutrient delivery leaves tissues vulnerable
- Evidence suggests that reduced temperature has protective influence
  - Kokate (1995)
  - Patel (1999)
- Kokate: “At a given pressure, ... lower temperatures exert a significant protective influence with respect to the development of pressure ulcers”

# Temperature and pressure

- Lachenbruch (2005)
  - 2<sup>nd</sup> analysis of published data
  - 8°C decrease in skin temperature is equivalent to a 29% reduction in interface pressure
  - Rightly advocates attention to skin temperature

# Controlled testing- Ferrarin & Ludwig, 2000

- Sequence of images taken
  - Before sitting (T0)
  - After 15 of sitting (T15)
  - 5 & 15 minutes after transfer (T20 & T35)



# Controlled testing- Ferrarin & Ludwig, 2000

Roho heats the most and cools the quickest (*steepest slope*)  
R. Medica gel retains heat the most (*lowest slopes*)





Fig. 2 Average surface temperature ( $T_{sk}$ ) against time for each cushion


△ cushion A  
× cushion B  
■ cushion C  
◆ cushion D

What's one limitation of the study and conclusion

# Logging temperature & humidity

Controlled testing  
Monitoring daily life

- Accuracy
  - $\pm 0.1\text{ }^{\circ}\text{C}$
  - $\pm 2\%\text{ RH}$
- Inserted temperature and humidity sensors at cushion interface under buttocks




## MSR data logger


External temperature sensor



External humidity sensor



MSR 145WS



20 x 15 x 61 mm




# Controlled 45 min test- 4 cushions


*Same subject; same clothes, same room*

## Predictions?

Action Twister



Hi Profile Roho



Polyurethane foam




Silicone-impregnated foam



# Controlled interface temperature measurements


*Same subject; same clothes, same room*



Predictions of RH responses?


# Controlled interface humidity measurements

*Same subject; same clothes, same room*




# Monitoring in everyday life

- Attached logger and sensors
  - Everything fit within cover, on the side
- Monitor for 1 week
  - Occupancy switch and debriefing help contextualize data



One long bout  
Up @ 10:30- down at 11:30

Same person  
Two different days




Long day, 3 bouts  
No PRs  
Up @ 8 am – down @ 1am

In both instances  
Humidity hit 90%  
Temperature peaked <30° C  
BUT  
Sitting bouts were very long

# Data including many off-loading episodes

A 5-hour  
block



day	average bout length (min)	total occupancy (min)
1	34.51	483.17
2	39.60	712.75
3	34.53	552.42
4	33.03	495.50



## Temperature variations by bouts of sitting

Roho Harmony


Tighter temperature range

(Of bouts >30 min, range 27.5 to 31.3 [*outlier*])

Rapid rise to max

This subject sat  $\approx 12$  hr/day

Sat for >6 hours in a row 1+ times daily



← Note: different time scales


Jay 2

Wide temperature range

(Of bouts >30 min, range 29.5 to 35.2)

Slower rise to max

This subject sat  $\approx 11$  hr/day & got up frequently





# Why should we care?

- Tissue microclimate is important
- Cushions vary widely in microclimate management just like they vary widely in pressure management
- Moving is a good thing
  - unweights tissue so dissipates heat & alters normal and shear loading
  - Facilitate movement via education, proper positioning, bribes, threats
- If client reports sweating, we should seek other solutions
  - Shear, friction and temperature implications
- Pressure reliefs have at least 2 purposes:
  - Alleviate pressure and dissipate heat