

Cushion use and performance in everyday life

Stephen Sprigle

Or

Stephen Herbertson






Understanding the State-of-Use and Performance

- Surveying used cushions
- Documenting degradation
- Temperature and humidity
 - Controlled tests
 - Within everyday use

Surveying used cushions


- Survey developed to document cushion status
 - Cushion construction
 - Reasons for replacement
 - Cushion & cover inspection
 - Damage, fatigue, failure
- Sent to Robert Bingham in Australia
- 209 surveys completed
- Flat and contoured foam

Age of foam cushions (n=209)




90% were
<12 mo old

Majority of cushion
replaced because of
their age




All- but 1- was deemed
'clean' or 'slightly soiled'

>60% of covers deemed
'clean' or 'slightly soiled'




Physical and clinical signs of fatigue



After 6 months, 70% showed physical signs of fatigue

Compression set: when a material fails to return to its unloaded state



40% of cushions 7-12 mo of age showed clinical sign of fatigue

>12 months, >60% showed a clinical sign

Why we should care

- Nice size data set on foam cushions
 - *Opportunity to use peer pressure to collect more data*
 - Hopefully, on different types of cushions
- Insight into a different delivery model
- Compression set occurs before clinical indicators of fatigue
 - High pressures or discomfort noted for 33% of cushions ≤ 12 mo
 - Compression set noted in 70% of cushions ≤ 12 mo

Conclusion

- Foam is in pretty good shape after 12 months of use
 - Australia should be commended for this model
- Certain temporary wheelchair users may benefit from a foam cushion
 - i.e., stroke survivors are often d/c'd with orders for only a wheelchair; minimal cost might meet needs

Documenting degradation

a collaborative project between



SHEPHERD CENTER
A Catastrophic Care Hospital

- Objectives
 - Identify the expected lifespan of cushions and the significant predictors of cushion failure
 - Develop and validate a clinical measure of seat cushion degradation
- 138 different cushions studied (24 measured >1x)
 - Most common: 32 Jay2; 26 Roho HP; 14 Evolution
- Client eval, visual inspection & performance measures
- Mean age: 24 months (range: 1 day to 168 months)

Testing cushions over time


Interview & physical exam



IPM with user




IPM using model



Visual inspection & dimensioning





Loaded contour depth



Impact dampening




Model and Human IPM




Metrics cover:
magnitude
asymmetry
dispersion



Look only at ROHO and JAY2



DI= ratio of IT pressures to total pressure





PPI= measure of pressure magnitude




Pressure magnitudes- ALL 162 cushions

Both model and subject pressures
indicate NO relationship over time



Black: IPM using buttock model
Red: IPM using cushion user

Look at variability of **red** model data
compared to variability of **black** subject data





A tale of 2 predictors

← Cushion age

Subject weight
↓



Dispersion Index (DI)



In fact....


- **Cushion age** has *not* been able to predict any IPM-related variable
 - For all cushions
 - Combining the 3 most tested cushions (Roho, Jay2, Evolution)
- Regression reported cushion age has a weak performance relationship within *Evolution*
 - *May be indicative of foam*

Looking only at FOAM-based cushions




No relationship with age
Huge variance in model testing

Cushion Cleanliness



About $\frac{3}{4}$ of cushion deemed 'clean'
'Unclean' cushions ranged in age



Clinician's perception of cushion adequacy

“is the cushion good enough for the client to use?”

- Clinician's are 5 ½ times more likely to judge a foam cushion as *inadequate* compared to a non-foam cushion
- Clinician's were 3 times more likely to judge a 'clean' cushion as 'adequate' than one judged 'moderately clean' or 'unclean'

Why we should care

- Tracking performance changes over time is needed to better understand “useful life”
- We need a means to identify cushions that need to be replaced
- Extensive data on 138 cushions is overwhelming
- Evidence suggests that Roho and Jay 2 cushion performance appears independent of age
 - For the cohort studied
- Cushion usage is individualistic so identifying a global measure of fatigue is very difficult
 - As always, individualistic evaluation is indicated

Temperature and humidity

- Humidity represents moisture
- Temperature represents temperature

Both have impacts on tissue integrity

Friction and Moisture

- As moisture increases, friction increases
 - ↑ softness → ↑ contact between surfaces
 - Want to learn more?- see cosmetics literature
- Excessive moisture weakens skin's ability to withstand load

Temperature and it's impact on tissue viability


- ↑ tissue temperature ↑ metabolic demand
 - Added demand coupled with reduced nutrient delivery leaves tissues vulnerable
- Evidence suggests that reduced temperature has protective influence
 - Kokate (1995)
 - Patel (1999)
- Kokate: “At a given pressure, ... lower temperatures exert a significant protective influence with respect to the development of pressure ulcers”

Temperature and pressure

- Lachenbruch (2005)
 - 2nd analysis of published data
 - 8°C decrease in skin temperature is equivalent to a 29% reduction in interface pressure
 - Rightly advocates attention to skin temperature

Controlled testing- Ferrarin & Ludwig, 2000

- Sequence of images taken
 - Before sitting (T0)
 - After 15 of sitting (T15)
 - 5 & 15 minutes after transfer (T20 & T35)



Controlled testing- Ferrarin & Ludwig, 2000

Roho heats the most and cools the quickest (*steepest slope*)
R. Medica gel retains heat the most (*lowest slopes*)





Fig. 2 Average surface temperature (T_{sk}) against time for each cushion


△ cushion A
× cushion B
■ cushion C
◆ cushion D

What's one limitation of the study and conclusion

Logging temperature & humidity


Controlled testing
Monitoring daily life

- Accuracy
 - $\pm 0.1\text{ }^{\circ}\text{C}$
 - $\pm 2\%\text{ RH}$
- Inserted temperature and humidity sensors at cushion interface under buttocks



MSR data logger

External temperature sensor



External humidity sensor




MSR 145WS




1st question?:

Do skin interface temperatures equal cushion interface temperatures?



Does RH at skin interface equal cushion interface RH?




Controlled 45 min test- 4 cushions

Same subject; same clothes, same room

Predictions?

Action Twister



Hi Profile Roho



Polyurethane foam




Silicone-impregnated foam




Controlled interface temperature measurements

Same subject; same clothes, same room




Monitoring in everyday life

- Attached logger and sensors
 - Everything fit within cover, on the side
- Monitor for 1 week
 - Occupancy switch and debriefing help contextualize data



One long bout
Up @ 10:30- down at 11:30

Same person
Two different days




Long day, 3 bouts
No PRs
Up @ 8 am – down @ 1am

In both instances
Humidity hit 90%
Temperature peaked <30° C
BUT
Sitting bouts were very long


Data including many off-loading episodes

A 5-hour
block



day	average bout length (min)	total occupancy (min)
1	34.51	483.17
2	39.60	712.75
3	34.53	552.42
4	33.03	495.50

Single day: Long bout sandwiched between PR activities





Does movement matter?

Controlled testing- Action Exact cushion



Skin-mounted sensor



Why should we care?

- Tissue microclimate is important
- Cushions vary widely in microclimate management just like they vary widely in pressure management
- Moving is a good thing
 - unweights tissue so dissipates heat & alters normal and shear loading
 - Facilitate movement via education, proper positioning, bribes, threats
- If client reports sweating, we should seek other solutions
 - Shear, friction and temperature implications
- Pressure reliefs have at least 2 purposes:
 - Alleviate pressure and dissipate heat

Thanks & questions

